



CCD and CMOS Cameras

DCU223x, DCU224x

DCC1240x

DCC1545M, DCC1645C

DCC3240X

DCC3260X

**DCx Camera Functional
Description and SDK Manual**



2018

Version: 4.82
Date: 7/30/2018

Contents

Foreword	9
1 General Information	10
1.1 Safety	10
1.2 Ordering Codes and Accessories	12
1.3 Requirements	12
1.4 DCx Camera Family	13
1.5 Contents	15
1.6 What's New in this Version?	17
2 Camera Basics	18
2.1 Operating Modes	18
2.1.1 Freerun Mode	18
2.1.2 Trigger Mode	20
2.1.3 Standby	21
2.2 Image Display Modes	22
2.3 Sensor	25
2.3.1 Sensor Sizes	25
2.3.2 Micro Lenses	26
2.3.3 Color Filter (Bayer filter)	29
2.3.4 Hot Pixels	31
2.3.5 Shutter Methods	33
2.3.6 Line Scan Mode	37
2.4 Reading Out Partial Images	38
2.4.1 Area of Interest (AOI)	38
2.4.2 Subsampling	42
2.4.3 Binning	43
2.5 Digitizing Images	44
2.5.1 Characteristics and LUT	44
2.5.2 Bit Depth and Digital Contrast Adjustment	47
2.6 Camera Parameters	50
2.6.1 Pixel Clock, Frame Rate, Exposure Time	50
2.6.2 Gain and Offset	51
2.6.3 Automatic Image Control	51
2.6.4 Applying New Parameters	53
2.7 Firmware and Camera Start	53
2.8 Digital Inputs / Outputs	53
2.8.1 Using Digital Inputs/Outputs	54
2.8.2 Flash Timing (Trigger Mode)	54
2.8.3 Flash Timing (Freerun Mode)	56
2.8.4 Serial Interface RS-232 (DC3240x only)	57

2.9	USB Interface	58
2.9.1	History and Development	58
2.9.2	Structure and Topology	58
2.9.3	USB 2.0 Cabling and Connectors	59
2.9.4	USB 3.0 Cabling and Connectors	60
2.9.5	Data Transmission and Bandwidth	60
3	Operation	62
3.1	uc480 Quick Start	62
3.2	Installation and Connection	64
3.2.1	System Requirements	64
3.2.2	DCx Driver Compatibility	65
3.2.3	Connecting a DCx Camera	65
3.3	Application Notes by Camera Model	67
3.3.1	DCC3260x Application Notes	67
3.3.2	DCC1240x / DCC3240x Application Notes	68
3.3.3	DCC1545M Application Notes	71
3.3.4	DCC1645C Application Notes	72
3.3.5	DCU223x Application Notes	72
3.3.6	DCU224x Application Notes	72
3.4	Installed uc480 Programs	72
3.4.1	uc480 Camera Manager	73
3.4.1.1	Camera List	74
3.4.1.2	Control Center	74
3.4.1.3	General Information	76
3.4.1.4	Camera Information	77
3.4.1.5	Creating a Support File	77
3.4.1.6	Additional Functions	77
3.4.1.7	Parameters	81
3.4.2	uc480 Hotpixel Editor	82
4	Programming (SDK)	85
4.1	First Steps to uc480 Programming	86
4.2	How to Proceed	88
4.2.1	Preparing Image Capture	88
4.2.1.1	Querying Information	88
4.2.1.2	Opening and Closing the Camera	89
4.2.1.3	Allocating Image Memory	89
4.2.1.4	Image Memory Sequences	91
4.2.2	Selecting the Display Mode	92
4.2.3	Capturing Images	93
4.2.3.1	Image Capture Modes	93
4.2.3.2	Event / Message Handling	101

4.2.4	Setting Camera Parameters	105
4.2.4.1	Setting and Getting Parameters	105
4.2.4.2	Automatic Image Control	108
4.2.4.3	Image Pre-processing	109
4.2.4.4	Get Camera Status	109
4.2.4.5	Using the Camera EEPROM	109
4.2.5	Saving Images and Videos	110
4.2.5.1	Saving and Loading Single Frames	110
4.2.5.2	Capturing AVIs	110
4.2.6	Using Inputs and Outputs	113
4.2.6.1	Input/Output Control	113
4.3	Function Descriptions	116
4.3.1	is_AddToSequence	117
4.3.2	is_AllocImageMem	118
4.3.3	is_AOI	120
4.3.4	is_AutoParameter	130
4.3.5	is_Blacklevel	133
4.3.6	is_CameraStatus	136
4.3.7	is_CaptureStatus	138
4.3.8	is_CaptureVideo	141
4.3.9	is_ClearSequence	143
4.3.10	is_ColorTemperature	144
4.3.11	is_Configuration	148
4.3.12	is_Convert	155
4.3.13	is_CopyImageMem	157
4.3.14	is_CopyImageMemLines	158
4.3.15	is_DeviceFeature	159
4.3.15.1	Configuring the AOI Merge Mode	161
4.3.15.2	Using the Log Mode	165
4.3.15.3	Using Level Controlled Trigger	168
4.3.15.4	Switching the Shutter Mode	169
4.3.15.5	Using the Internal Image Memory	170
4.3.15.6	Using the Line Scan Mode	172
4.3.15.7	Configuring the Timestamp	173
4.3.16	is_DeviceInfo	174
4.3.17	is_DirectRenderer	177
4.3.18	is_DisableEvent	187
4.3.19	is_EdgeEnhancement	188
4.3.20	is_EnableAutoExit	191
4.3.21	is_EnableEvent	192
4.3.22	is_EnableMessage	195
4.3.23	is_ExitCamera	197

4.3.24	is_ExitEvent	198
4.3.25	is_ExitImageQueue	199
4.3.26	is_Exposure	200
4.3.26.1	Setting the Exposure Time	202
4.3.26.2	Exposure Time with Fine Increments	205
4.3.26.3	Setting the Long Exposure	207
4.3.26.4	Setting the Dual Exposure	208
4.3.27	is_ForceTrigger	210
4.3.28	is_FreelImageMem	212
4.3.29	is_FreezeVideo	214
4.3.30	is_Gamma	216
4.3.31	is_GetActiveImageMem	217
4.3.32	is_GetActSeqBuf	219
4.3.33	is_GetAutoInfo	221
4.3.34	is_GetBusSpeed	225
4.3.35	is_GetCameraInfo	226
4.3.36	is_GetCameraList	228
4.3.37	is_GetCameraLUT	230
4.3.38	is_GetColorConverter	231
4.3.39	is_GetColorDepth	232
4.3.40	is_GetDLLVersion	233
4.3.41	is_GetError	234
4.3.42	is_GetFramesPerSecond	235
4.3.43	is_GetFrameTimeRange	236
4.3.44	is_GetImageHistogram	238
4.3.45	is_GetImageInfo	241
4.3.46	is_GetImageMem	245
4.3.47	is_GetImageMemPitch	246
4.3.48	is_GetNumberOfCameras	248
4.3.49	is_GetOsVersion	249
4.3.50	is_GetSensorInfo	250
4.3.51	is_GetSensorScalerInfo	252
4.3.52	is_GetSupportedTestImages	253
4.3.53	is_GetTestImageValueRange	255
4.3.54	is_GetTimeout	256
4.3.55	is_GetUsedBandwidth	257
4.3.56	is_GetVsyncCount	258
4.3.57	is_HasVideoStarted	259
4.3.58	is_HotPixel	260
4.3.59	is_ImageFile	265
4.3.60	is_ImageFormat	269
4.3.61	is_InitCamera	276

4.3.62	is_InitEvent	280
4.3.63	is_InitImageQueue	282
4.3.64	is_InquireImageMem	284
4.3.65	is_IO	286
4.3.66	is_IsVideoFinish	299
4.3.67	is_LockSeqBuf	301
4.3.68	is_LUT	302
4.3.69	is_Measure	305
4.3.70	is_ParameterSet	309
4.3.71	is_PixelClock	312
4.3.72	is_ReadEEPROM	315
4.3.73	is_RenderBitmap	317
4.3.74	is_ResetToDefault	319
4.3.75	is_SetAllocatedImageMem	321
4.3.76	is_SetAutoParameter	324
4.3.77	is_SetBinning	332
4.3.78	is_SetCameraID	335
4.3.79	is_SetColorConverter	337
4.3.80	is_SetColorCorrection	339
4.3.81	is_SetColorMode	341
4.3.82	is_SetDisplayMode	345
4.3.83	is_SetDisplayPos	349
4.3.84	is_SetErrorReport	350
4.3.85	is_SetExternalTrigger	351
4.3.86	is_SetFrameRate	353
4.3.87	is_SetGainBoost	355
4.3.88	is_SetGamma	356
4.3.89	is_SetHardwareGain	358
4.3.90	is_SetHWGainFactor	361
4.3.91	is_SetImageMem	363
4.3.92	is_SetOptimalCameraTiming	364
4.3.93	is_SetRopEffect	366
4.3.94	is_SetSaturation	368
4.3.95	is_SetSensorScaler	369
4.3.96	is_SetSensorTestImage	372
4.3.97	is_SetSubSampling	374
4.3.98	is_SetTimeout	378
4.3.99	is_SetTriggerCounter	380
4.3.100	is_SetTriggerDelay	381
4.3.101	is_StopLiveVideo	382
4.3.102	is_Trigger	383
4.3.103	is_TriggerDebounce	385

4.3.104	is_UnlockSeqBuf	389
4.3.105	is_WaitEvent	390
4.3.106	is_WaitForNextImage	391
4.3.107	is_WriteEEPROM	393
4.4	AVI Function Descriptions	395
4.4.1	isavi_AddFrame	395
4.4.2	isavi_CloseAVI	396
4.4.3	isavi_DisableEvent	396
4.4.4	isavi_EnableEvent	397
4.4.5	isavi_ExitAVI	398
4.4.6	isavi_ExitEvent	399
4.4.7	isavi_GetAVIFileName	400
4.4.8	isavi_GetAVIFileNameW	401
4.4.9	isavi_GetAVISize	402
4.4.10	isavi_GetnCompressedFrames	403
4.4.11	isavi_GetnLostFrames	404
4.4.12	isavi_InitAVI	405
4.4.13	isavi_InitEvent	406
4.4.14	isavi_OpenAVI	407
4.4.15	isavi_OpenAVIW	408
4.4.16	isavi_ResetFrameCounters	409
4.4.17	isavi_SetFrameRate	410
4.4.18	isavi_SetImageQuality	411
4.4.19	isavi_SetImageSize	412
4.4.20	isavi_StartAVI	414
4.4.21	isavi_StopAVI	415
4.5	RAW function descriptions	416
4.5.1	israw_AddFrame	416
4.5.2	israw_CloseFile	416
4.5.3	israw_ExitFile	417
4.5.4	israw_GetFrame	417
4.5.5	israw_GetImageInfo	418
4.5.6	israw_GetSize	419
4.5.7	israw_InitFile	419
4.5.8	israw_OpenFile	420
4.5.9	israw_SeekFrame	421
4.5.10	israw_SetImageInfo	421
4.6	Obsolete functions	422
4.7	Programming Notes	425
4.7.1	Programming in C/C++	425
4.7.2	Programming in C#	426
4.7.3	Programming in VB.NET	427

4.7.4	Programming in Delphi	427
4.7.5	Programming with ActiveX	427
4.7.6	Thread Programming	427
4.8	Lists	429
4.8.1	Complete List of All Return Values	429
4.8.2	Error Codes of AVI Functions	433
4.8.3	Linux: Not Supported Functions	433
5	Specifications	435
5.1	Model Comparison	436
5.2	Model Naming Conventions	438
5.3	Camera and Sensor Data	438
5.3.1	DCC3260x	439
5.3.2	DCC1240x / DCC3240x	441
5.3.3	DCC1545M	445
5.3.4	DCC1645C	447
5.3.5	DCU223x	449
5.3.6	DCU224x	452
5.4	Mechanical Specifications	455
5.4.1	DCU223x, DCU224x	456
5.4.2	DCC1240x	457
5.4.3	DCC3240x	458
5.4.4	DCC1545M, DCC1645C	459
5.4.5	Flange Back Distance	460
5.4.5.1	Calculating the Flange Back Distance	460
5.4.5.2	Maximum Immersion Depth for Lenses	462
5.4.6	Position Accuracy of the Sensor	463
5.4.7	Filter Glasses	464
5.4.7.1	Filter Types	464
5.4.7.2	Mounting the Filter	468
5.4.7.3	Cleaning the Filter Glasses	469
5.4.8	Ambient Conditions	471
5.5	Camera Interface	472
5.5.1	DCU223x, DCU224x, DCC1240x	472
5.5.1.1	I/O Connector - Pin Assignment	472
5.5.1.2	Digital Input (Trigger) Circuit	474
5.5.1.3	Digital Output (Flash) Circuit	475
5.5.2	DCC3240x	477
5.5.2.1	I/O Connector Pin Assignment	477
5.5.2.2	GPIO Interface	478
5.5.2.3	Digital Input (Trigger) Circuit	480
5.5.2.4	Digital Output (Flash) Circuit	481
5.5.2.5	RS-232 Serial Interface	482

5.5.3	Camera EEPROM Specification _____	483
5.6	Accessories for DCx cameras _____	484
5.6.1	Accessories for DCU22xX / DCC1240X _____	484
5.6.2	Accessories for DCC1x45X _____	485
5.6.3	Accessories for DCC3240x / DCC3260x _____	485
6	Appendix _____	486
6.1	Troubleshooting/FAQ _____	486
6.1.1	PCs with Energy Saving CPU Technology _____	487
6.2	Status LED on USB DCx Cameras _____	488
6.3	Color and Memory Formats _____	490
6.4	uc480 Parameter File (ini file) _____	493
6.5	Definition of IP Protection Classes _____	500
6.6	History of API functions _____	500
6.7	Certifications and Compliances _____	502
6.8	Thorlabs 'End of Life' Policy (WEEE) _____	504
6.9	Exclusion of Liability and Copyright _____	505
6.10	Thorlabs Worldwide Contacts _____	505

Warning

Sections marked by this symbol explain dangers that might result in personal injury or death. Always read the associated information carefully, before performing the indicated procedure.

Attention

Paragraphs preceded by this symbol explain hazards that could damage the instrument and the connected equipment or may cause loss of data.

1 General Information

Thank you for purchasing a DCx camera!

You should first read the following chapters to get a quick overview on what is new in this software version and on getting started with your new camera.

Getting started

- [DCx quick-start](#)
- [First steps to DCx Camera programming](#)

Further important information

- [What is new in this version?](#)
- [Contents](#) of this Manual
- [The DCx camera family](#)
- [Specifications](#)

Enjoy your new DCx camera!

1.1 Safety

Attention

All statements regarding safety of operation and technical data in this instruction manual will only apply when the unit is operated correctly as it was designed for.

Prior to applying power to the DCx Camera Functional Description and SDK Manual, make sure that the protective conductor of the 3 conductor mains power cord is correctly connected to the protective earth ground contact of the socket outlet! Improper grounding can cause electric shock with damages to your health or even death!

The DCx Camera Functional Description and SDK Manual must not be operated in explosion endangered environments!

Do not remove covers! Do not obstruct the air ventilation slots in the housing!

Do not open the cabinet, there are no parts serviceable by the operator inside!

Refer servicing to qualified personnel!

Only with written consent from *Thorlabs Scientific Imaging* may changes to single components be made or components not supplied by *Thorlabs Scientific Imaging* be used.

This precision device is only serviceable if properly packed into the complete original packaging. If necessary, ask for a replacement package prior to return.

Attention

The following statement applies to the products covered in this manual, unless otherwise specified herein. The statement for other products will appear in the accompanying documentation.

This equipment has been tested and found to comply with the limits for a Class **A** digital device, pursuant to part 15 of the FCC Rules and meets all requirements of the Canadian Interference-Causing Equipment Standard ICES-003 for digital apparatus. These limits are

designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

Thorlabs Scientific Imaging is not responsible for any radio television interference caused by modifications of this equipment or the substitution or attachment of connecting cables and equipment other than those specified by Thorlabs Scientific Imaging. The correction of interference caused by such unauthorized modification, substitution or attachment will be the responsibility of the user.

The use of shielded I/O cables is required when connecting this equipment to any and all optional peripheral or host devices. Failure to do so may violate FCC and ICES rules.

Attention

Mobile telephones, cellular phones or other radio transmitters are not to be used within the range of three meters of this unit since the electromagnetic field intensity may then exceed the maximum allowed disturbance values according to IEC 61326-1.

This product has been tested and found to comply with the limits according to IEC 61326-1 for using connection cables shorter than 3 meters (9.8 feet).

1.2 Ordering Codes and Accessories

DCU223M	CCD camera, monochrome, 1024x768 pixel, C mount
DCU223C	CCD camera, color, 1280x1024 pixel, C mount
DCU224M	CCD camera, monochrome, 1280x1024 pixel, C mount
DCU224C	CCD camera, color, 1280x1024 pixel, C mount
DCC1545M	CMOS camera, monochrome, 1280x1024 pixel, CS mount
DCC1645C	CMOS camera, color, 1280x1024 pixel, CS mount
DCC1240M	CMOS camera, monochrome, 1280x1024 pixel, C mount
DCC1240C	CMOS camera, color, 1280x1024 pixel, C mount
DCC3240M	CMOS camera, monochrome, 1280x1024 pixel, C mount, USB 3.0
DCC3240C	CMOS camera, color, 1280x1024 pixel, C mount, USB 3.0
DCC3240N	CMOS camera, NIR enhanced, 1280x1024 pixel, C mount, USB 3.0
DCC3260M	CMOS camera, monochrome, 1936x1216 pixel, C mount, USB 3.0
DCC3260C	CMOS camera, color, 1936x1216 pixel, C mount, USB 3.0
CAB-DCU-T1	Trigger cable for DCU22xX and DCC1240X cameras (Trigger In/Out)
CAB-DCU-T2	Trigger cable for DCU22xX and DCC1240X cameras (Trigger In only)
CAB-DCU-T3	I/O cable for DC3240 CMOS USB 3.0 cameras

Thorlabs C Mount Camera Lenses (objectives): [See Thorlabs' website](#)

1.3 Requirements

For operating the DCx cameras, the following system requirements must be met:

	Recommended
CPU speed	>2.0 GHz Intel Core i5 or Core i7
Memory (RAM)	8 GByte
For USB DCx cameras: USB host controller	USB 3.0 Super Speed Intel® motherboard chipset
Graphics card	Dedicated AGP/PCIe graphics card Latest version of Microsoft DirectX Runtime 9.0c
Operating system	Windows 8.1 32 or 64 bit Windows 7 32 or 64 bit

Drivers for network cards

To ensure optimum performance of the network connection, you need to install the latest drivers for your network card. We recommend using the drivers of the following versions:

- Intel® chipsets: version 8.8 or higher
- Realtek chipsets: version 5.7 or higher

USB interface

- Onboard USB 2.0 ports usually provide significantly better performance than PCI and PCMCIA USB

adapters.

- Current generation CPUs with energy saving technologies can cause bandwidth problems on the USB bus. See section on PCs With Energy Saving CPU Technology.

Large multi-camera systems

Connecting a large number of cameras to a single PC may require a large working memory (RAM). This is especially the case when many cameras with high sensor resolution are used.

If you want to set up such a system we recommend to use PCs with 64 bit operating systems and more than 4 GB of RAM.



Note on color cameras with high frame rates

For uc480 color cameras, the color conversion is done by software in the PC. When you use a color camera with a high frame rate, the conversion might lead to a high CPU load. Depending on the PC hardware used you might not be able to reach the camera's maximum frame rate.

Direct3D graphics functions

The uc480 driver can use Direct3D to display the camera image with overlay information (Microsoft DirectX Runtime had to be installed). On Windows systems, you can use the supplied "DXDiag" diagnostic tool to check whether your graphics card supports Direct3D functions. To start the diagnostic tool, click "Run..." on the Windows start menu (shortcut: Windows+R) and enter "DXDiag" in the input box.

On the "Display" page of the diagnostic tool, click the button for testing the Direct3D functions.

OpenGL graphics functions

For OpenGL version 1.4 or higher must be installed. The OpenGL graphics functions do not work with QT under Linux.

1.4 DCx Camera Family

DCx cameras stand for a range of compact and cost-effective cameras for professional use in industrial, security and non-industrial applications. Equipped with the widely used USB 2.0 and particularly USB 3.0 ports, they can easily be interfaced with a vast variety of systems. The images are digitized in the camera and transmitted digitally to the PC. An additional frame grabber is not required.

DCU cameras have state-of-the-art CCD sensors while the DCC models are CMOS based. The CMOS models use either the [global or the rolling shutter](#) method; the CCD models use only the global shutter method.

The DCx cameras are available as monochrome and color versions, DC3240 series has a NIR version as well. The [Model Comparison](#) chapter shows the most important features of every series at a glance.

USB 3.0 DCC3260x and DCC3240x CMOS Cameras



Compact, fast and lightweight. The new **DCC3260x** and **DCC3240x**. The 29 x 29 x 29 mm small camera housing is not only ultra-compact, but due to its magnesium casing and a total camera weight of 43 g, it is also ultra-lightweight and robust. The powerful camera offers a bandwidth of 400 MByte/s via USB 3.0. Power is supplied via the USB bus, hence an extra power cable is obsolete.

With its lockable Micro USB connector the camera is perfectly suited even for rough environments. Offering trigger and flash as well as two GPIOs (General Purpose I/O), which can also be changed into a serial interface (RS232). Hence, peripheral devices can easily be triggered or controlled.

But also the camera's inner values are outstanding: brightness corrections are easily realized by a comfortable 12 bit lookup table and hardware gamma. 12 bit color depth offers a by factor 16 increased level of detail compared to the usual 8 bit. Hardware based data preprocessing saves additional CPU resources.

USB 2.0 DCC1240x (CMOS) and DCC22xX (CCD) Cameras

The **DCC1240X** and **DCC22xX** series feature a robust metal housing with a standard mini-B USB 2.0 connector. Connection is additionally possible via a lockable micro D-sub connector which also carries the opto-isolated I/O signals.

The USB 2.0 interface is meanwhile available in every standard PC and notebook/laptop and provides a gross bandwidth of 60 MByte/s. The camera is connected and powered through the USB port by just a single cable.



USB2.0 DCC1545M and DCC1645C Cameras



The **DCC1x45X** series features extremely compact cameras with high-speed CMOS sensors. The LE models are designed for professional use in non-industrial applications. Through the use of the widespread USB 2.0 technology, the cameras can easily be interfaced with a vast variety of systems. These cameras are available with a plastic housing with CS-mount lens adapter.

1.5 Contents

The DCx Camera Manual contains all the information you need for operating your DCx camera. It comprises the following parts:

Section A: Camera basics

- In this section you will find a lot of important information on the technical background of your USB camera. This section contains explanations on the DCx's [operating modes](#), on [sensor technology](#), important [camera parameters](#), and the [USB](#) interfaces. We recommend to read this chapter to become familiar with the general functionality of the DCx Cameras.

Section B: Operation

- [Quick start](#) to using your DCxCamera
- [Installing](#) and [Using DCx Camera software](#)
These sections show how to connect cameras and start operation using the software tools .
- [Application notes by camera model](#)
This section explains special features and limitations of some camera models.

Section C: Programming

- [First steps](#) to programming with your DCxCamera
- [How to proceed](#)
If you are not yet familiar with DCxCamera programming, we suggest that you first explore the basic functional flows in this chapter. The function blocks contain almost all the functions available for the uc480 API ordered by topics. The flowcharts help to easily find the appropriate API function for a certain task.
- [Description of functions/Description of AVI functions](#)
These chapters cover all the functions of the uc480 API in alphabetic order.
The AVI functions for video recording are implemented by the `uc480_tools.dll` which is also included in the DCxCamera software package.
- [Obsolete functions](#)
This chapter lists obsolete API functions and recommended alternatives.
- [Lists and programming notes](#)
In this chapter, you will find useful information on how to use the DCxCamera programming API. Programming environments, modes for DCxCamera color and image display as well as the automatic image control functions are discussed here.

Section D: Specifications

- [Specifications](#)
All information on the camera's [sensor and performance](#), [mechanical](#) as well as [electrical specifications](#) are contained in this section.
- [Accessories](#)
Here you will find a list of accessories for DCx Cameras sorted by model.

Appendix

- Information on [Troubleshooting](#)
- Status LEDs on [USB DCx](#) cameras
- [Color and memory formats](#)
- [uc480 parameter file \(ini file\)](#)

- [Definition of IP protection classes](#)

1.6 What's New in this Version?

Version 4.82 of the DCxCamera software package includes many new features and enhancements. The following table gives you an overview of the major new functions.

New in Version 4.82.00

Cameras & functions	Described in chapter
New models in the USB 3 DCC camera family: <ul style="list-style-type: none">• DCC3260	DCC3260

Older versions

See the History of uc480 Software Versions and [History of uc480 API functions](#) chapters.

2 Camera Basics

This chapter explains the basics of DCx Camera technology.

- [Operating modes](#)
- [Image display modes](#)
- [Sensor](#)
- [Reading out partial images](#)
- [Digitizing images](#)
- [Camera parameters](#)
- [Firmware and camera start-up](#)
- [Digital inputs/outputs](#)
- [USB interface](#)

2.1 Operating Modes

DCx Cameras support the following operating modes:

- [Freerun mode](#)
- [Trigger mode](#)
- [Standby](#)

2.1.1 Freerun Mode

In freerun mode, the camera sensor captures one image after another at the set frame rate. Exposure and readout/transfer of the image data are performed in parallel. This allows the maximum camera frame rate to be achieved. The frame rate and the exposure time can be set separately. The captured images can be transferred one by one or continuously to the PC. If trigger mode is active, you need to disable it before activating freerun mode.


Note

Note on the schematic diagrams: These illustrations show a schematic view of the image capture sequence. The sensor exposure and readout times and the transmission times depend on the camera model and settings. The pre-processing time depends on the API functions you are using (e.g. color conversion, edge enhancement).

For more information on flash timing see the [Digital In-/Output \(Trigger/Flash\)](#) chapter.

Continuous mode (live mode)

Images are captured and transferred continuously. You can use the DCxCamera flash outputs.




* Flash function optional. See also [Digital in/output \(trigger/flash\)](#).

Note

In freerun mode the flash function starts with the second image as the setting of the flash timing depends on the finish of the first image. If you change the flash timing during operation, the freerun mode will restart. Therefore the first image after the change is black.

Single frame mode (snap mode)

The next image exposed by the sensor will be transferred. In this mode, flash is not making sense (only manually).



See also:

- Basics: [Shutter methods](#)
- Basics: [Trigger mode](#)
- Basics: [Applying new parameters](#)

Programming:

- [Capture modes](#)

2.1.2 Trigger Mode

In trigger mode, the sensor is on standby and starts exposing on receipt of a trigger signal. A trigger event can be initiated by a software command (software trigger) or by an electrical signal via the camera's digital input (hardware trigger).

This chapter describes the different trigger modes you can use with the DCx Cameras. To choose a mode, go to the Settings icon in the ThorCam application or use the [API command](#).


Note

Note on the schematic diagrams: These illustrations show a schematic view of the image capture sequence. The sensor exposure and readout times and the transmission times depend on the camera model and settings. The pre-processing time depends on the API functions you are using (e.g. color conversion, edge enhancement).

For more information on flash timing see the [Digital In-/Output \(Trigger/Flash\)](#) chapter.

Software trigger mode

When this mode is enabled, calling the "Snap" function triggers the capture of an image, which is then transferred to the PC. If you call the "Live" function in this mode, the image capture is triggered continuously and images are transferred continuously.



* Optional flash function. See also [Digital input/output \(trigger/flash\)](#)


Hardware trigger mode

When this mode is enabled, calling the [is_FreezeVideo\(\)](#) (Snap) function makes the camera ready for triggering just once. When the camera receives an electrical trigger signal, one image is captured and transferred.

If you call the [is_CaptureVideo\(\)](#) (Live) function, the camera is made ready for triggering continuously. An image is captured and transferred each time an electrical trigger signal is received; the camera is then ready for triggering again (recommended procedure).

Attention

When you use triggered image capture, the camera is only ready to process the next trigger signal after completion of the data transfer to the PC. Trigger events that occur during image exposure or data transfer are ignored. An internal counter records the number of ignored trigger events and can be read out from the PC.



* Optional flash function. See also [Digital input/output \(trigger/flash\)](#)

Frame rate in trigger mode

With many sensors, the maximum frame rate is lower in trigger mode than in freerun mode because these sensors expose and transfer sequentially. Which frame rate is possible in trigger mode therefore depends on the exposure time. The time required for capturing a frame in trigger mode can be approximated with the following formula:

$$t_{\text{capture}} = \text{Current exposure time} + \left(\frac{1}{\text{max.frame rate}} \right)$$

Example: At the maximum exposure time, the frame rate is about half as high as in freerun mode; at the minimum exposure time, the frame rate is about the same.

Freerun synchronization

This mode is currently not supported by DCx Cameras.

See also:

- Basics: [Freerun mode](#)
- Basics: [Digital input/output \(trigger/flash\)](#)
- ThorCam: Settings > Trigger

Programming:

- [Image capture modes: Trigger](#)

2.1.3 Standby

DCx Cameras can be set to a power-saving standby mode. Standby mode switches off the sensor of CMOS cameras and the timing board of CCD cameras. The camera remains open in the software.

In standby mode, the camera cools down and the number of hot pixels visible when longer exposure times are used is reduced.

Standby is the default state when the camera is not open in the software. When you open the camera or switch to a different mode (freerun or trigger mode), the camera wakes up from standby mode.

Note

In standby mode, you can continue to use the camera's digital inputs or outputs.

2.2 Image Display Modes

The uc480 driver provides different modes for displaying the captured images on Windows systems. We recommend using the Bitmap mode or the Direct3D functions, depending on your specific application.

For a list of API functions for image display see [How to proceed: Image display](#).


Attention

The "DirectDraw BackBuffer" and "DirectDraw Overlay Surface" display modes are obsolete. Please use the Direct3D functions instead (see also [Obsolete functions](#)).

1. Bitmap mode (Device Independent Bitmap, DIB)

In Bitmap mode, images captured by the DCxCamera are written to the random access memory of the PC. Programming the image display is up to the user. The application software uses the [is_RenderBitmap\(\)](#) function to initiate the image display by the graphics card. This may result in a slightly higher CPU load as compared to the Direct3D display.


The advantage of Bitmap mode is that it is compatible with all graphics cards and that image data in the memory is directly accessible. Programming of overlay functions is up to the user. Since the operating system controls the image display, the image may be completely or partly overlapped by other windows and dialog boxes.



2. Direct3D mode (only under Windows with DirectX)


In this mode, the uc480 driver writes the image data to the invisible area of the graphics card. This process runs automatically and does not have to be controlled by the application software. It requires an installed Direct3D driver, sufficient memory on the graphics card and Direct3D function support by the graphics card (see [System requirements](#)). For this purpose, graphics cards generally provide better performance than graphics chips integrated on the mainboard. In Direct3D mode, the CPU load may be lower than in Bitmap mode. You can display overlay data and also scale the video image.

The Direct3D mode and the overlay functions can be configured using the [is_DirectRenderer\(\)](#) API function.



OpenGL mode

OpenGL stands for Open Graphics Library and it is an interface specification for graphics hardware. Unlike Direct3D OpenGL is not only available under Windows but also under Linux and Mac OS, if OpenGL is supported by the graphics hardware. There are several implementations of OpenGL, such as e.g. NVIDIA® or AMD/ATI. The implementations are always dependent on the graphics card manufacturer.



Comparison of the display modes

The following table illustrates the major differences between the display modes:

	Bitmap mode	Direct3D mode	OpenGL mode
Graphics card requirements	Low. No special graphics hardware required. Runs on all systems.	High. Graphics card has to support Direct3D. Does not run on all systems.	High. Graphics card has to support OpenGL.
Operating system	Windows, Linux	Only Windows with DirectX	Cross-platform
Programming effort	Greater. Memory management, event handling and display performed by the application.	Low. Memory management, event handling and display performed by DirectX.	High. OpenGL itself does not provide functions for opening windows or reading files. However, there are related libraries, e.g. GLUT.
CPU load	Slightly increased by copying of data.	Low. Display performed by graphics card.	Low. Display performed by graphics card.
Overlay functions	Not available. A simple overlay can be	Integrated. Complex overlays can be	Integrated.

	programmed by the user.	displayed without flicker.	
Access to image memory	Direct access possible. Image data already provided in user memory.	Possible using Steal Mode. Single images can be copied to the user memory.	Direct access to graphics card and image memory.


2.3 Sensor

- [Sensor sizes](#)
- [Micro lenses](#)
- [Color filter \(Bayer filter\)](#)
- [Hot pixels](#)
- [Shutter methods](#)
- [Line scan mode](#)

2.3.1 Sensor Sizes

The size of a digital camera sensor is usually specified in inches. However, the specified value does not indicate the actual size of the active sensor area. The sensor size specifications date back to the formerly used tube systems: The curvature of the imaging surface of the camera tube caused distortions to the display, reducing the usable capture area of a 1" tube to a rectangle with a diagonal of 16 mm.

With the introduction of the semiconductor sensor technology, the dimensional specifications were taken over from tube systems. For this reason, a sensor whose active area diagonal measures 16 mm is specified as a 1-inch sensor. The following illustrations show the most common sensor sizes. The diameter in inch multiplied with $2/3$ equals approximately the actual sensitiv area in millimeters.



Common sensor sizes (in inch)



1" Sensor



2/3" Sensor



1/2" Sensor



1/3" Sensor


Comparison of common sensor sizes and examples for different fields of view

The size of each single sensor cell (pixel) depends on the size of the active sensor area and the resolution. In general, less pixels over the same sensor area (or a larger sensor area with the same resolution) will result in greater photo sensitivity of the sensor.


2.3.2 Micro Lenses

Micro lenses improve the fill factor

The fill factor is the percentage of the pixel area that is exposed to light during exposure. Ideally this would be 100 %. Since other elements are located on the sensor surface besides the light-sensitive photodiodes, this value may be reduced to approx. 30–50 %, depending on the sensor technology. The use of micro lenses compensates for this and increases the fill factor to 90 % or more. Micro lenses collect the light that falls onto a photocell, thus increasing the useable sensor area.






Using micro lenses to increase the effective fill factor



CMOS pixel design with Bayer filter (red) and micro lens

Micro lenses with CRA correction

Some sensors have micro lenses offset to the sensor edge. They compensate for shading created by obliquely incident light. The angle of incident light is called Chief Ray Angle (CRA), the micro lens offset is thus called CRA correction. The amount of micro lens shift is specified in degrees and refers to the micro lenses in the corners of the sensor.




Note

Using parallel light on sensors with CRA correction may cause slight color variations. These may occur, for example, if telecentric lenses are used. The following models are equipped with sensors with offset micro lenses:

- [DCC1240x/DCC3240x](#)
- [DCC1645C](#)

2.3.3 Color Filter (Bayer filter)

For technical reasons, digital image sensors can only detect brightness information, but no color information. To produce color sensors, a color filter is applied to each photocell (pixel). The arrangement of the color filters is illustrated in the following figure. Two out of every four pixels have a green filter, one pixel has a red filter and one has a blue filter. This color distribution corresponds to the color sensitivity of the human eye, and is called the Bayer filter pattern. With the help of the Bayer pattern the correct brightness and color information can be calculated for each pixel. Full sensor resolution is retained.



Bayer RGB filter pattern

Bayer conversion

A Bayer conversion, also referred to as de-Bayering, is carried out to determine the color information from the raw sensor data (raw Bayer). By default all DCx Cameras transmit the image data to the PC in raw Bayer format. The PC then uses the functions of the uc480 API to convert the image data to the color format you need for displaying or further processing the data.

To convert the colors, a filter mask moves over the image and calculates a color value for each pixel from the surrounding pixels. The uc480 API provides two filter masks that differ in image quality and CPU load.


- Normal Quality (Mode `IS_CONV_MODE_SOFTWARE_3x3`/`IS_CONV_MODE_HARDWARE_3x3`)
A smaller filter mask is used for conversion. This algorithm has a low load on the CPU. The filter's averaging function may cause a slight blur. Noise is reduced. This filter is recommended for image processing tasks.
- High Quality (Mode `IS_CONV_MODE_SOFTWARE_5x5`)
A large filter mask is used for conversion. This algorithm offers very accurate color positioning and an increased level of detail. The CPU load is higher than with the normal filter. This filter is recommended for visualization applications.

Note

Software conversion with high quality should only be used for sensors whose green pixels have the same sensitivity. This applies to the following sensors:

- DCU223C / DCU224C
- DCC1240C, DCC3240C

For all other sensors, we recommend using the standard filter mask.



Bayer conversion using the standard mask

See also:


- Color conversion: [is_SetColorConverter\(\)](#)
- ThorCam: Settings > Camera > Pixel Data Format

2.3.4 Hot Pixels

Definition

Hot pixels (or in a broader sense, defective pixels) are pixels that do not react linearly to incident light – or do not react at all. They occur for various reasons, such as contamination during sensor production or sensor age, and with both CCD and CMOS sensors. CCD sensors generally have fewer hot pixels than CMOS sensors under the same operating conditions. With darkened sensors and prolonged exposure times, hot pixels are visible as individual bright dots in the image. The following factors promote the occurrence of hot pixels:

- Long exposure times
- High gain settings
- High sensor operating temperature



Hot pixels detected in a monochrome camera



Hot pixels detected in a color camera

Hot pixel correction

During the manufacture of our cameras, all sensors that will be used in DCx Cameras are checked for hot pixels. In the process, images are taken with a darkened sensor and long exposure times. Pixels with a brightness higher than a specific value are classified as hot pixels. A list of the

coordinates of each hot pixel is stored in the camera EEPROM. The hotpixel correction is done in the uc480 driver. However, some sensors also provide an internal hotpixel correction.

The maximum number of hot pixels stored in a DCx camera is:

DCx model	max. hot pixels stored
DCC1240x, DCC1545M, DCC1645C, DCC3240x (CMOS)	768
DCU223x, DCU224x (CCD)	20

How many hot pixels are on the camera's internal list depends above all on the defined threshold values. It is not an indication of the quality of the sensors used.

When you enable the "Hotpixel correction" function in the DCx software, the software automatically corrects the hot pixels in the captured image by calculating the average from the brightness value of two neighboring pixels. When using color sensors, the hot pixel is corrected with the appropriate color in raw Bayer format, i.e. before color conversion. The correction does not work with activated subsampling and binning factors greater than 2x.

Note

The sensors are tested during manufacturing also for cold pixel and dead pixels. Sensors with dead pixel clusters (more than two neighboring defective pixels of the same color) are rejected by our quality control. When the camera is operated in very warm ambient conditions, other defective pixels can occur, however.

Defining additional hot pixels

If additional hot pixels occur during use of the camera, you can add them to the camera's internal hot pixel list. To do this, use the API function given below.

See also:

- ThorCam: Hot pixel correction
- [uc480 Hotpixel Editor](#)
- Programming: [is_HotPixel\(\)](#)

2.3.5 Shutter Methods

- [Global shutter](#)
- [Rolling shutter](#)
- [Rolling shutter with global start](#)

Note

Note on the schematic diagrams: These illustrations show a schematic view of the image capture sequence. The sensor exposure and readout times and the transmission times depend on the camera model and settings.

For more information on flash timing see the [Digital In-/Output \(Trigger/Flash\)](#) chapter.

General

The image is recorded in the sensor in four phases:

- Reset pixels of the rows to be exposed
- Exposure of pixel rows
- Charge transfer to sensor
- Data readout


The sensor cells must not be exposed during the readout process. The sensors of the DCx Cameras have no mechanical shutters, but work with electronic shutter methods instead. Depending on the sensor type, either the rolling shutter method or the global shutter method is used.

Global shutter


On a global shutter sensor, all pixel rows are reset and then exposed simultaneously. At the end of the exposure, all rows are simultaneously moved to a darkened area of the sensor. The pixels are then read out row by row.

Exposing all pixels simultaneously has the advantage that fast-moving objects can be captured without geometric distortions. Sensors that use the global shutter system are more complex in design than rolling shutter sensors.

All CCD sensors as well as some CMOS sensors use the global shutter method.



Global shutter sensor in live mode



Global shutter sensor in trigger mode

* Optional flash function. The start time and duration are defined by the flash delay and duration parameters (see also Camera settings: I/O).

Rolling shutter

With the rolling shutter method, the pixel rows are reset and exposed one row after another. At the end of the exposure, the lines are read out sequentially. As this results in a time delay between the exposure of the first and the last sensor rows, captured images of moving objects are distorted.




Example for the rolling shutter effect with a moving car

To counteract this effect, the DCx Camera software provides a global flash window where you set the time by which flash activation is delayed. You can also specify the flash duration. This allows implementing a global flash functionality which exposes all rows of a rolling shutter sensor simultaneously.


Rolling shutter sensors offer a higher pixel density compared to global shutter CMOS sensors. The rolling shutter system is used in DCC Cameras with high-resolution CMOS sensors.

Note


Some CMOS sensors with global shutter can be operated also with rolling shutter. The operation in the rolling shutter mode is used to reduce the image noise. This function is only supported from the camera models [DCC1240x/DCC3240x](#).



Rolling shutter sensor in live mode



Rolling shutter sensor in live mode with global flash window




Rolling shutter sensor in triggered mode with global flash window

* Optional flash function. The start time and duration are defined by the flash delay and duration parameters (see also Camera settings: I/O).

Rolling shutter with global start

Some rolling shutter sensors also provide a global start mode, which starts exposure of all rows simultaneously (see illustration). For best results, use a flash for this mode. No light is allowed to fall on the sensor outside the flash period because otherwise the image brightness will be distributed unevenly.



Rolling shutter sensor in trigger mode with global start function

* Optional flash function. The start time and duration are defined by the flash delay and duration parameters (see also Camera settings: I/O).

2.3.6 Line Scan Mode

Area scan sensor (matrix)

The sensors of area scan cameras have a matrix of many (usually between several hundred and several thousand) rows and columns of pixels. State-of-the-art area scan sensors use only square pixels with a consistent pixel pitch.

Area scan cameras are suitable for applications in which stationary or moving objects should be captured as completely as possible in one image capture.

Line scan mode

In some applications, however, it is necessary to read out and transfer only one sensor line instead of the entire sensor area. This applies, for example, to endless web inspection systems. These systems often use line scan cameras for this reason. Their sensors have only one pixel row, which they can read out at very high speeds in the kilohertz (kHz) range. Some DCxCamera models have area scan sensors that optionally also offer a line scan mode. This mode can read out any pixel row of the sensor at high speed.

There are two line scan modes to distinguish:

- **Fast line scan**

In this mode, the sensor achieves very high line scan rates. Several hundred to thousand lines are combined and transferred in one frame. The camera can be triggered on the beginning of a frame, but not on each individual line. You can choose any line of the area scan sensor for the line scan mode. Color images are not supported in this mode because [Bayer color sensors](#) need at least two neighboring lines for color calculation.

- **Triggered line scan**

In this mode, the sensor achieves lower line scan rates than in fast line scan mode. The camera can be triggered on each individual line. Several hundred to thousand lines are combined and transferred in one frame.

Color images are possible in this mode because Bayer color sensors can use two lines.

Note

The line scan mode is currently only supported by the monochrome [DCC1240M_and DCC3240M,N](#) models in form of the fast line scan mode. The triggered line scan mode is not supported by any camera model yet.

See also:

- ThorCam: Settings > Shutter

Programming:

- Function: [is_DeviceFeature\(\)](#)

2.4 Reading Out Partial Images

The camera sensors have defined resolutions which are given as the number of pixels (width x height). However, for some applications it may be necessary to read out only a selected part of the sensor area or to reduce the local resolution. For this purpose, the DCx Cameras provide various functions:



- [Area of interest \(AOI\)](#)
- [Subsampling](#) (skipping) pixels
- [Binning](#) (combining) pixels

These functions reduce the amount of data to be transferred and thus allow you to increase the frame rate considerably, depending on the camera model.

2.4.1 Area of Interest (AOI)

Using this function, you can set the size and position of an area of interest (AOI) within an image. In this case, only data included in this AOI will be read out and transferred to the computer. The smaller partial image enables the camera to use a higher frame rate.

For information on the AOI position grid and the frame rates that your camera model can achieve with AOI, see the model specifications in the [Camera and sensor data](#) chapter.



Note

Step widths for AOI definition (position grid): The available step widths for the position and size of image AOIs depend on the sensor. The values defining the position and size of an AOI have to be integer multiples of the allowed step widths.

For details on the AOI grids of the individual camera models, please see [Camera and sensor data](#)

and click a camera model.

Please note that, after defining an AOI, the resulting image may be darker if the camera cannot maintain the originally set exposure time due to the increased frame rate.

Multi AOI

The Multi AOI function allows defining more than one AOI in an image and transferring these AOIs all at the same time. Only DCC1240x models support this feature. In the Multi AOI mode you can define two or four AOIs in one image and transfer them simultaneously. The AOIs are positioned side by side or one below the other, and share the same X or Y axis. This feature is not supported in ThorCam at this time.

Sequence AOI mode

Apart from the multi AOI mode, DCC1240x and DCC3240x also support the sequence AOI mode. This mode allows to define up to four AOIs, which need to have the same size but may differ in position, exposure time or gain settings. This feature is not supported in ThorCam at this time.

See also:

- ThorCam: Crop to Region of Interest
- Programming: [is_AOI\(\)](#)

2.4.2 Subsampling



Subsampling is a technique that skips multiple sensor pixels when reading out image data. This reduces the amount of data to be transferred and enables higher camera frame rates. The captured image has a lower resolution but still the same field of view compared to the full-resolution image. This mode can be used as a fast preview mode for high-resolution cameras.

Color subsampling as performed by most color sensors skips pixels while maintaining colors (see illustration). For some monochrome sensors, the camera also performs color subsampling, resulting in slight artifacts.

Monochrome sensors and some color sensors ignore the Bayer pattern and the color information gets lost (mono subsampling).

Depending on the model, DCx Cameras support different subsampling factors. Subsampling of horizontal and vertical pixels can be enabled independently.

The [Camera and sensor datas](#) chapter lists the subsampling methods and factors supported by each camera model.



2.4.3 Binning

Binning is a function that averages or adds multiple sensor pixels to obtain a single value. This reduces the amount of data to be transferred and enables higher camera frame rates. The captured image has a lower resolution but still the same field of view compared to the full-resolution image. This mode can be used as a fast preview mode for high-resolution cameras.



Color binning, as performed by most color sensors, combines only pixels of the same color (see also the [Color filter \(Bayer filter\)](#) chapter). For some monochrome sensors, the camera also performs color binning, resulting in slight artifacts.

Most monochrome sensors and some color sensors combine neighboring Bayer pattern pixels; in this case, the color information gets lost (mono binning).

With CCD sensors, binning makes the images brighter because the pixel values are added up. With CMOS sensors, pixel values are usually averaged; this reduces image noise.

Depending on the model, DCx Cameras support different binning factors. Binning of horizontal and vertical pixels can be enabled independently.

The [Camera and sensor data](#) chapter lists the binning methods and factors the individual camera models support.



2.5 Digitizing Images


- [Characteristics and LUT](#)
- [Bit depth and digital contrast adjustment](#)

2.5.1 Characteristics and LUT

When perceiving or imaging a scene, the form of the imaging characteristic is crucial for displaying the differences in brightness. With image processing (e.g. applications such as edge detection and character recognition), linear characteristics are generally required. The human eye, on the other hand, perceives differences in brightness based on a logarithmic characteristic, which often approximates a gamma characteristic in practice. All three forms will be shown in the following.

Linear characteristic

If a system (e.g. a camera with a conventional CCD sensor) yields double the output value for double the brightness, the system features a linear characteristic:




Gamma characteristic

Gamma characteristics (or gamma curves) are named after the Greek formula symbol γ . Gamma curves are power functions of the form

$$y = x^{\frac{1}{\gamma}}$$

and are often used in photography or image display on computer screens. A gamma value of 1 generates a linear characteristic again. A curve with the value $\gamma = 2.2$ used for computer screens is shown in the figure below.



Such a gamma characteristic brightens dark areas of an image, which corresponds more to the perception of the human eye. In light areas of an image, the differences in brightness are condensed for this.



Image with linear characteristic





Image with gamma characteristic

Logarithmic characteristic

The effect of the logarithmic characteristic is even stronger. Here, the characteristic follows the function

$$y = \lg(x)$$

The following diagram illustrates how very large jumps in brightness in light areas of a scene only cause small changes in image brightness. This explains why image sensors with a logarithmic characteristic, in particular, are ideal for imaging scenes with very high dynamic range.




Imaging with logarithmic characteristic

Lookup table (LUT)

With a lookup table (LUT) it is easy to apply characteristic curves to digital images. A LUT is a table which assigns an output value to every possible input value. The figure below shows a LUT which would binarize an image: For an 8 bit image, for example, this LUT would replace all pixels with gray values 0...127 with value 0 and all pixels with gray values 128...255 with value 255.


Using LUTs has the advantage that calculations can be done very fast. Typical applications of LUTs are enhancing image contrast, or gamma characteristics.



2.5.2 Bit Depth and Digital Contrast Adjustment

Digitizing

Image sensor pixels first generate an analog voltage signal proportional to the amount of light that strikes them. The image is digitized for further processing, i.e. the stepless signal is converted to a digital numerical value. The following figure shows this using a gray gradient as an example



Various bit depths using a gray-scale gradient as an example

If the stepless gradient is imaged in a digital range in 2 bits, for example, the result is $2^2 = 4$ levels; for 4 bits, it is $2^4 = 16$ levels, and so on. The intermediate brightness values of the original gradient are irreversibly lost after digitization.

With around 200 levels or more, the jumps in brightness can no longer be discerned with the eye, which is why current monitors and digital cameras use 8 bits (256 levels) per color channel (fully adequate for visualization).

Bit depth in image processing

If digital image data undergoes further image processing, a bit depth greater than 8 may be necessary. The computer is able to differentiate between these very fine differences in brightness (no longer discernable by the eye) and process them. This is why industrial cameras often use 12 bits.

Bit depth	Brightness levels
8	$2^8 = 256$
10	$2^{10} = 1024$
12	$2^{12} = 4096$
14	$2^{14} = 16.384$

Note

Greater bit depths require extremely low-noise image sensors, however. As soon as the differences in brightness created by noise are greater than the digitization levels, no further data is gained.

Bit depth by sensors

Platform	USB 2.0	USB 3.0
CMOS sensors	8 bit	10 bit
CCD sensors	8 bit	-

Note

Color formats with a bit depth of more than 8 bits per channel are only supported by USB 3 DCC3240x camera models. Using color formats with higher bit depth increases the bandwidth used by a camera.

Histogram and contrast

The brightness distribution of digital images is represented in a histogram. If an image has optimum contrast, the histogram includes practically all brightness values between 0 and the highest value (255 in 8-bit images). If an image has low contrast, the histogram only includes a small number of the values; the image appears dull:




Image capture and histogram with minimal contrast (l.e.) and with optimum contrast after a contrast adjustment (r.i.)


For improved display on the screen or when printed, the histogram can be spread to optimally utilize the possible brightness levels. For this purpose, the dark parts of the image are further darkened via an LUT characteristic and the light parts of the image are brightened. Thus the human eye can better differentiate between the different brightness levels; the image has more contrast.

It must be noted, however, that subsequent processing with a computer will not yield more data. Therefore, subsequent contrast adjustment via software is not necessary for use in image processing. The computer can differentiate between the differences in brightness without contrast adjustment.

Advantage of greater bit depth with contrast adjustment


The bit depth in the output image is crucial for contrast adjustment. The following figures illustrate this. In the first example the 8 bit output image contains fewer than 100 brightness levels, as there are no dark or very bright parts. The image is low-contrast.

With a contrast adjustment, the values of the histogram are spread in such a way as to create a contrast-rich image. The fewer than 100 brightness values are now distributed across levels 0 to 255; gaps arise in the histogram and are visible as jumps in brightness in the resulting image.



The second example shows the same output image with a 10 bit bit depth right at the time of capture. This image also has low contrast, as it features only average brightness values. The greater bit depth allows the brightness values of the image to be imaged over 500 different digital levels, however. The entire histogram includes 1,024 values in the 10 bit image (in contrast to 256 values with 8 bits).

This means that a contrast adjustment can now be made for screen display without a reduction in quality. The 500 values of the output image are distributed over the 256 values of the 8-bit target image in such a way that optimum contrast is the result. The large number of output values means that there are no gaps in the histogram.



Note

This type of contrast adjustment can already be done in the camera when an image is digitized in 10 bits and transferred in 8 bit. In this case, optimum utilization of the 8 bit data is also important for image processing.

2.6 Camera Parameters

- [Pixel clock, frame rate, exposure time](#)
- [Gain and offset](#)
- [Automatic image control](#)
- [Applying new parameters](#)

2.6.1 Pixel Clock, Frame Rate, Exposure Time

Pixel clock

The basic parameter for camera timing is the pixel clock. It determines the speed at which the sensor cells can be read out.

Attention

We recommend not setting the pixel clock any higher than necessary to achieve the desired frame rate.

An excessive pixel clock can cause delays or transmission errors. If the data is read from the sensor at a higher speed (high pixel clock), you will also need a faster transmission over the data connection. Thus, by controlling the pixel clock, you can also influence the bandwidth required for a camera.

The pixel clock influences the connected load and consequently the temperature inside the camera.

Frame rate

The possible range of settings for the frame rate depends on the currently selected pixel clock. You can select a lower frame rate without changing the pixel clock. To set a higher frame rate, however, you need to increase the pixel clock.

Exposure time

The exposure time depends on the currently selected frame rate and is preset to its reciprocal value. You can select a shorter exposure time without changing the frame rate. To set a longer exposure time, however, you need to reduce the frame rate.

Note

The increments for setting the exposure time depend on the sensor's current timing settings (pixel clock, frame rate). The exposure time values are rounded down to the nearest valid value, if required. For this reason, the actual exposure time can deviate slightly from the exposure time you have selected.

See also:

- ThorCam: Settings > Camera
- [is_PixelClock\(\)](#)
- [is_SetFrameRate\(\)](#)
- [is_Exposure\(\)](#)

2.6.2 Gain and Offset

Gain

In digital imaging, a voltage proportional to the amount of incident light is output by the sensor. To increase image brightness and contrast, this signal can be amplified by an analog gain and offset before the digitizing process. The results of analog signal processing are usually better than the results of digital post-processing.

Analog amplification of the read-out pixel values increases overall image brightness and contrast. Depending on the sensor type, a global gain value for all pixels (master gain) or a separate gain value for each color (RGB gain) can be set.

Note

Using Sensor Gain: A signal gain will also result in a noise gain. High gain settings are therefore not recommended.

We suggest the following gain settings:

1. Enable the Gain boost function ([is_SetGainBoost\(\)](#)).
2. If required, adjust the gain setting with the master gain control.

Note

Linearity of sensor gain: You can set the gain factor in increments from 0 to 100.

- For CCD sensors the gain increases usually not linear but disproportionate.
- For CMOS sensors the gain increases linear. Some sensors have only 32 or fewer levels, so not each step is assigned to a level.

The maximum gain factor settings also vary from sensor to sensor (see [Camera- and sensor data](#)).

Offset

Every digital image sensor has light-insensitive cells next to the active image area. These dark pixels are used to measure a reference voltage (black level) which is subtracted from the image signal. This compensates thermally generated voltages on the sensor which would otherwise falsify the signals.

Normally, the sensor adjusts the black level automatically. If the environment is very bright or if exposure times are very long, it may be necessary to adjust the black level manually.

2.6.3 Automatic Image Control

The uc480 driver provides various options to automatically adjust the image capture parameters to the lighting situation. These include:

- Auto exposure shutter (AES)
- Auto gain control (AGC)
- Auto white balance (AWB)
- Auto frame rate (AFR)

The auto functions are used to adjust the average brightness and color rendering of the camera image to their setpoint values, while trying to keep the frame rate at the highest possible value.

All controls are configured using the [is_SetAutoParameter\(\)](#) SDK function.

Auto exposure shutter (AES)

The control of the average brightness is preferably achieved by adjusting the exposure, i.e. you set the highest possible exposure time before gain is controlled. The auto exposure feature always

uses the current exposure range which results from the selected pixel clock frequency and the frame rate. You can set separate control range limits for exposure and gain.

Auto gain control (AGC)

The auto gain feature controls the camera master gain in a range from 0-100 %. You can set separate control range limits for exposure and gain.

Auto frame rate (AFR)

With the exposure control function enabled, you can still change the frame rate manually or automatically to maintain a dynamic exposure control range. A lower frame rate allows for longer exposure times, but then the live image display may exhibit jitter. The objective of the automatic frame rate control is to set the frame rate to an optimum value. This way, in all situations, the automatic exposure control can use the required control range at the highest possible frame rate.

Auto white balance (AWB)

Depending on the lighting source, light can have different color temperatures so that the images may have a color cast. At low color temperatures (e.g. light from incandescent lamps), the white content is offset towards a red hue. At high color temperatures (e.g. light from fluorescent lamps), the white content is offset towards a blue hue.

The white balance control feature uses the RGB gain settings of the camera to correct the white level. This is achieved by adjusting the gain controls within the 0-100 % range until the red or blue channel matches the average brightness of the green channel. In order to manually influence the color rendering, you can adjust the setpoint values for the red and blue channels relative to the green channel by using an offset value (see also ThorCam > Histogram).

Automatically disabling the control function

You can disable the control functionality automatically once the target value has been reached (approximately) and after 3 regulations no improvement has been reached (API parameters `IS_SET_AUTO_WB_ONCE` and `IS_SET_AUTO_BRIGHTNESS_ONCE`). An event/a message notifies the system of this (see also [is_InitEvent\(\)](#)). Alternatively, you can keep the control feature enabled so that it responds to deviations from the target value.

Control speed

You can set the auto function speeds in a 0–100 % range. This influences the control increments. High speed (100 %) causes a little attenuation of a fast-responding control and vice versa. The control functions for average brightness and for color rendering use separate speeds.

In trigger mode, every frame is evaluated for automatic control. The freerun mode skips a number of frames by default because in that mode, changes to the image parameters only become effective after one or more image captures (see also [Applying new parameters](#)). With the "Skip Frames" parameter (API parameter `IS_SET_AUTO_SKIPFRAMES`), you can select how many frames should be skipped in freerun mode (default: 4). This parameter strongly influences the control speed. Choosing small values can destabilize the automatic control.

Note

For higher frame rates select for the "Skip frames" parameter a bigger value. This reduces the number of automatic adjustments that must be done by the camera.

Hysteresis

The automatic control feature uses a hysteresis function for stabilization. Automatic control is stopped when the actual value lies in a range between (setpoint - hysteresis value) and (setpoint + hysteresis value). It is resumed when the actual value drops below (setpoint - hysteresis value) or exceeds (setpoint + hysteresis value). If the hysteresis value is increased, the control function will stop sooner. This can be useful in some situations.

See also:

- [is_SetAutoParameter\(\)](#)

2.6.4 Applying New Parameters

New capture parameters (such as exposure time or gain settings) can be transferred to the camera via software at any time. Depending on the operating mode, these settings will not always be immediately effective for next image, however.

- Freerun mode

In freerun mode, the camera is internally busy with capturing the next image while new parameters are transmitted to the camera. Depending on the exact time of transmission, new parameters might only come into effect two or even three images later.

- Trigger mode

In this mode, the camera reverts to idle state between two images. When you change the camera parameters, the new settings will be applied immediately to the next image.

2.7 Firmware and Camera Start

Every DCx camera has its own firmware that handles internal processes in the camera. The camera firmware varies from model to model.

USB DCx Cameras have a two-tier firmware that is uploaded to the camera each to you connect it to a PC:

1. Common firmware (uc480 boot)

The general firmware identifies what camera model you have connected, and uploads the corresponding firmware.

2. Model-specific firmware (e.g.: uc480 DC1240x series)

The model-specific firmware is named after the camera type and provides the functions of the relevant model.

Note

When you connect a USB DCxCamera with a Windows PC or a new USB port for the first time, it is detected as a new device. This is normal standard behavior of the operating system.

The USB DCx Cameras firmware is part of the driver. The automatic upload always loads the firmware that matches the driver installed in the camera.

2.8 Digital Inputs / Outputs

All DCx cameras (see [Model comparison](#)), except DCC1545M and DCC1645C, come with opto-isolated inputs/outputs that can be used for triggering the camera and for flash control. DC3240x cameras have in addition general purpose I/Os (GPIO). Use of the GPIOs for flash control is possible to a certain degree. External triggering via the GPIO is not supported.

See also:

- Basics: [Trigger mode](#)
- ThorCam: Settings > Input/output
- Specification: [Electrical specifications](#)

Programming:

- [is_IO\(\)](#)

2.8.1 Using Digital Inputs/Outputs

Digital input (trigger)

Models with optocoupler input can use the digital input for externally triggering the image capture, or query the applied signal level.

In [trigger mode](#), a digital signal is applied to the camera's input. You can determine whether the camera will respond to the rising or falling edge of the digital signal. After an internal delay, the sensor is exposed for the defined exposure time. The captured image is then transferred to the PC.

On models with general purpose I/Os (GPIO), you can query a voltage level at these inputs (TTL compatible).

Digital output (flash)

The digital outputs can be used in both freerun mode and trigger mode. You can synchronize the output level to the exposure time or set it statically.

Models with optocoupler output allow control of a DC voltage applied to the output. This allows controlling a flash, either directly or via a separate flash controller unit. Models with general purpose I/Os (GPIO) can output a voltage at these outputs (TTL compatible).

Note

Please read the notes on I/O wiring for your camera model in the [Electrical specifications](#) chapter.

Note

The settings specified for the digital output will be reset when the camera is disconnected from the PC or the PC is powered down.

2.8.2 Flash Timing (Trigger Mode)

When using the digital output for flash control, you can set the delay and the duration of the flash. The flash timing can be adjusted manually or automatically by the camera driver.

Note



Sensor latency and delay times The sensor latency is due to a number of technical factors, including sensor type, image geometry, pixels clock and, with CCD sensors, the exposure time. The latency is constant for a specific combination of parameters.

[Trigger delay](#), flash delay and flash duration are optional and can be set by software.

The following illustrations show a schematic view of the image capture sequence. The sensor exposure and readout times and the transmission times depend on the camera model and the current parameter settings.

Automatic flash

If flash delay = 0 and flash duration = 0, the flash signal is automatically synchronized to the exposure time. The automatic flash feature has the advantage that the flash is synchronized automatically if the settings for image geometry or camera timing are changed. The disadvantage is that the flash signal is active slightly longer than the exposure time. The flash duration with automatic flash is longer for rolling shutter sensors than for global shutter sensors.



Manual flash synchronization

If one of the flash delay or flash duration parameters is set to a value greater than 0, you can shift the flash signal to any point in the exposure time or change its duration. In this case, the flash delay will be calculated exactly from the start of the exposure time (after the sensor latency time). When manually synchronizing the flash signal to the exposure time, you can use the [is_IO\(\)](#) function to query the data you need.

The advantage of manual flash synchronization is that the flash can be precisely controlled based on the start of exposure. This applies to both rolling and global shutter sensors. You can thus achieve a higher accuracy with the manual flash synchronization than with the automatic flash feature.



The disadvantage is that the flash signal has to be resynchronized whenever any settings for image geometry or camera timing change.

Note

With rolling shutter sensors, you can avoid the [rolling shutter effect](#) by selecting suitable delay and duration settings (global flash function). Using [is_IO\(\)](#), you can query the appropriate values.

Note


The flash output is reset with the start of the next image capture. This also applies if you have set a longer flash duration.



2.8.3 Flash Timing (Freerun Mode)

Automatic flash

In freerun mode, the automatic flash feature works in the same way as in trigger mode. As a result, the flash output is continuously or almost continuously active (see illustration below).




Manual flash synchronization

Note

It is recommended to synchronize the flash manually in freerun mode. This applies to both rolling and global shutter sensors.

In freerun mode, the manual flash synchronization works in the same way as in trigger mode.



2.8.4 Serial Interface RS-232 (DC3240x only)

DC3240x cameras are equipped with a serial interface (RS-232). It provides functionality for communication with peripheral devices (e.g. lighting controller, lens controller or the serial port of a PC). Before you can send data through the camera's serial interface, one or more virtual COM ports have to be defined on the PC. Once defined, they can be used for data communication with appropriate software just like any physical COM port.

To set up and use the serial interface, the ["Additional functions"](#) dialog box is provided in the uc480 Camera Manager. For the serial interface specifications, please refer to the [DC3240x Serial Interface Wiring \(RS-232\)](#) chapters.

2.9 USB Interface

- [History and development](#)
- [Structure and topology](#)
- [USB 2.0 cabling and connection](#)
- [USB 3.0 cabling](#)
- [Data transmission and bandwidth](#)

2.9.1 History and Development

The **Universal Serial Bus** (USB) is an interface which enables you to easily connect various devices to a PC. As all data exchange is controlled by the PC, no additional interface controller is needed. Further advantages of USB are:

- The PC does not have to be shut down when connecting USB devices (hot plugging).
- USB devices can be supplied with power from the PC.
- High bandwidth for data transmission.

The USB standard was developed by a group of companies including Compaq, IBM, Intel, and Microsoft. Version 1.0 was presented in 1995. The slightly faster USB 1.1 standard followed in 1998.


At first, the USB interface was designed to connect peripheral devices such as printers, mice, or keyboards. With the introduction of USB 2.0 in 2000, the transfer rate increased to 480 Mbit/s, making USB 2.0 suitable for connecting devices with higher data volumes (such as mass storage devices, scanners, or cameras).

In 2008, with USB 3.0 a new version of the interface has been published, which is significantly faster than USB 2.0 (400 MByte/s).

2.9.2 Structure and Topology

USB uses a tree topology and is host-controlled. That means that a PC with host functionality is mandatory for using USB. Therefore, it is not possible to directly connect two USB devices (with the exception of USB on-the-go compliant devices). Neither is it possible to connect a camera to a PDA device.

Theoretically, 127 devices can be connected to a host controller. Using external hubs or repeaters, even more devices can be connected, and from a greater distance. Provided that a maximum of 5 hubs/repeaters may be daisy-chained, USB devices can be connected in up to seven levels.



Note

The maximum bandwidth of 480 Mbit/s per USB 2.0 host or 400 MByte/s per USB 3.0 host cannot be exceeded. Therefore, the maximum possible frame rate will be reduced if image data from multiple USB cameras is transferred simultaneously.

The available bandwidth might also be decreased when you use hubs or repeaters. You can reduce the bandwidth required for each camera by lowering the frame rate or the image size.


2.9.3 USB 2.0 Cabling and Connectors

In order to comply with the specifications, the maximum length of USB 2.0 cables is limited to 5 m. Longer cables may be connected if you use high-quality material. The USB bus provides power supply with 5 V and 500 mA max. Many USB devices use the bus power and do not need external power supply (bus-powered devices).

Cable design

The following illustration shows the basic design of a shielded USB cable:

- D+/D-: data transfer
- +5 V/GND: power supply



Connector types

On the PC side, USB 2.0 cables are equipped with a standard A type plug (four pins) and on the device side either with a standard B plug (four pins) or a mini-B plug (five pins).


2.9.4 USB 3.0 Cabling and Connectors

In order to comply with the specifications, the maximum length of USB 3.0 cables is limited to 3-8 m. With the use of repeaters cable lengths up to 20 m are possible. With signal conversion into optical signals cable lengths up to 100 m are possible. The USB bus provides power supply with 5 V and 900 mA max.

Cable design

The following illustration shows the basic design of a shielded USB cable:


- SSTX+/-: SuperSpeed Transmit (data transfer from host to the device)
- SSRX+/-: SuperSpeed Receive (data transfer from device to the host)
- D+/D-: data transfer (USB 2.0)
- +5 V/GND: power supply



Basic design of a USB 3.0 cable

Connector types

On the PC side, USB 3.0 cables are equipped with a standard A type plug (8 pins) and on the device side either with a standard B plug or a micro-B plug.



While a USB 3.0 standard A plug and a USB 2.0 port can be used together (with the restriction that only USB 2.0 speed is possible), the standard B and micro B plug are no longer compatible with USB 2.0 ports.

2.9.5 Data Transmission and Bandwidth

USB 2.0

The USB 2.0 standard specifies an overall bandwidth of 480 Mbit/s shared between different transmission modes. DCx Cameras use the USB 2.0 bulk mode for transmitting images. This mode uses error correction to ensure correct delivery of the image data, but does not guarantee a

fixed bandwidth. To ensure error-free communication with all connected devices at all times, the maximum bandwidth for payload data is limited to 416 Mbit/s.

Theoretically, up to 50 MByte/s of data can be transmitted in this mode, but in practice, this value is hardly ever reached. A high-performance desktop PC can transmit about 40 MByte/s, most notebooks or embedded PC systems even less than that.

The overall bandwidth can be increased by the use of USB 2.0 expansion cards. These cards are available for the PCI and PCIe buses and have their own host controller chip.

Note on hardware selection

To achieve optimum USB bandwidth, it is important to use a powerful mainboard chipset. The mainboard chipsets from e.g. Intel® or NVIDIA® provide very good results.

Note on image content and bandwidth usage

For USB cameras, you can use a white test image to check the camera's maximum load on the USB bus. Due to the transmission process, completely white camera images require a somewhat more bandwidth on the USB bus than completely black images.

USB 3.0

Compared to USB 2.0, USB 3.0 offers a tenfold increased bandwidth of 5 Gbit/s, i.e. images can be transmitted with a bandwidth of 400 Mbytes/s.

3 Operation

This chapter explains how to connect the DCx camera and how to use the applications contained in the uc480 software package.

- [Quick start](#)
- [Installation and connection](#)
- Installed uc480 programs
 - ThorCam: A comprehensive viewer for exploring the camera functionality (See separate ThorCam User's Manual)
 - [uc480 Camera Manager](#): The central tool for managing all connected DCx Cameras.
 - [uc480 Hotpixel Editor](#): A tool to edit the sensor hot pixel list stored in the camera.

3.1 uc480 Quick Start

This chapter show how to quickly get started with your DCx camera. You will learn how easy it is to connect the camera and explore important functions. For further steps of integrating the DCx camera into your own applications please also see the [First steps to uc480 programming](#) chapter.



Connect the camera

Install the latest version of the DCx software. Then connect the DCx camera with the PC. USB cameras are automatically detected as new hardware under Windows. Check the status LEDs on your camera to see if the camera has been correctly identified.

See also:

- [Installation and connection](#)
- [Troubleshooting](#)
- [Connection - Status LED](#)



Configure the camera

USB DCx Cameras are ready for use right out of the box. You can assign a unique ID to your camera with the uc480 Camera Manager.

See also:

- [uc480 Camera Manager](#)
- [Assigning a camera ID in the camera manager](#)
- [Firmware and camera start](#)



Capture images

The uc480 software package includes many sample programs that you can use to try out the extensive functionality of your DCx camera. We recommend starting off with the ThorCam application. To run the application, simply double-click the corresponding icon on your Windows desktop.

See also:

- ThorCam Software User Manual (Separate Document)
- [Camera basics: Operating modes](#)



Customize the key camera properties

Select the Settings icon on the menu bar to open the dialog box for modifying the camera properties.

The "Camera" tab provides all the parameters for adjusting the camera's speed. You can increase the pixel clock to run the camera at a higher frame rate. Reduce the pixel clock if transmission errors occur too often.

On the "Image" tab, you find various sensor gain controllers. Use the "Master gain" controller to increase image brightness if no longer exposure time setting is possible. Switch to the "AES/AGC" tab to enable the Auto Exposure Shutter (AES) and Auto Gain Control (AGC) features.

Tip: Select a low sensor gain to minimize visible noise.

If you are using a color camera, you should activate sensor color correction on the "Color" tab in order to achieve rich vibrant colors for on-screen display. To adapt a color camera to the ambient light conditions, it is essential to carry out Auto White Balance (AWB). Aim the camera at a surface of a uniform gray color, then enable the "Image white balance: Enable" and "Run once" check boxes on the "AWB" tab.

See also:

- ThorCam Software User Manual (Separate Document)
- [Camera basics: Camera parameters](#)



Activate trigger and flash modes

DCx Cameras provide the possibility to trigger the image capture and to have the flash controlled by the camera. To switch the camera to trigger mode, go to the camera properties as described above, select the "Trigger" tab and enable the desired mode. To trigger on "falling edges" or "rising edges", a digital signal has to be applied to the camera. When you are finished with the trigger settings, select "Start Continuous Trigger" or "Triggered Snapshot" icons on the menu bar to start the triggered image capture.

If you have connected the digital output on your DCx camera to a flash controller, you can configure the flash function on "Input/Output" tab. Enable "Flash high active" and "Global exposure window". This way, the DCx camera automatically activates the flash during the exposure time.

See also:

- ThorCam Software User Manual (Separate Document)
- [Camera basics: Digital input/output](#)
- [Specifications: Electrical specifications](#)



Save the camera settings and images

With ThorCam, saving single frames or videos is very easy to do. Just choose the relevant option on the "File" menu. If you have recorded AVI videos, you can play them through Windows Media Player, ImageJ, or other AVI-compliant viewer with Motion-JPEG decoding capability.

See also:

- ThorCam Software User Manual (Separate Document)

3.2 Installation and Connection

- [System requirements](#)
- Installing uc480 software under Linux
- [Connecting a USB DCx camera](#)

3.2.1 System Requirements

For operating the DCx cameras, the following system requirements must be met:

	Recommended
CPU speed	>2.0 GHz Intel Core i5 or Core i7
Memory (RAM)	8 GByte
For USB DCx cameras: USB host controller	USB 3.0 Super Speed Intel® motherboard chipset
Graphics card	Dedicated AGP/PCIe graphics card Latest version of Microsoft DirectX Runtime 9.0c
Operating system	Windows 8.1 32 or 64 bit Windows 7 32 or 64 bit

Drivers for network cards

To ensure optimum performance of the network connection, you need to install the latest drivers for your network card. We recommend using the drivers of the following versions:

- Intel® chipsets: version 8.8 or higher
- Realtek chipsets: version 5.7 or higher

USB interface

- Onboard USB 2.0 ports usually provide significantly better performance than PCI and PCMCIA USB adapters.
- Current generation CPUs with energy saving technologies can cause bandwidth problems on the USB bus. See section on PCs With Energy Saving CPU Technology.

Large multi-camera systems

Connecting a large number of cameras to a single PC may require a large working memory (RAM). This is especially the case when many cameras with high sensor resolution are used.

If you want to set up such a system we recommend to use PCs with 64 bit operating systems and more than 4 GB of RAM.



Note on color cameras with high frame rates

For uc480 color cameras, the color conversion is done by software in the PC. When you use a color camera with a high frame rate, the conversion might lead to a high CPU load. Depending on the PC hardware used you might not be able to reach the camera's maximum frame rate.

Direct3D graphics functions

The uc480 driver can use Direct3D to display the camera image with overlay information (Microsoft DirectX Runtime had to be installed). On Windows systems, you can use the supplied "DXDiag" diagnostic tool to check whether your graphics card supports Direct3D functions. To start the diagnostic tool, click "Run..." on the Windows start menu (shortcut: Windows+R) and enter "DXDiag" in the input box.

On the "Display" page of the diagnostic tool, click the button for testing the Direct3D functions.

OpenGL graphics functions

For OpenGL version 1.4 or higher must be installed. The OpenGL graphics functions do not work with QT under Linux.

3.2.2 DCx Driver Compatibility

Attention


Support of older DCC1545M cameras by driver versions 3.10 and higher

From driver version 3.10 onwards, only USB board revision 2.0 or higher are supported. To operate a camera with an earlier USB board revision, you will need the *uc480* driver version 2.40. Please contact [Thorlabs](#)

The LED on the back of the camera housing also indicates the USB board version (see [DCx_Status_LED](#)). In addition, the *uc480 Camera Manager* version 3.10 or higher provides information about the compatibility (see [Camera Manager](#)). An incompatible camera will be displayed as *free* and *not available*.

3.2.3 Connecting a DCx Camera

Please install the software first as described in the Quick Start Guide. Connect the DCx camera to the PC, using the USB cable. The camera will be recognized automatically and the necessary driver software is being installed:




When the camera has been correctly installed, the LED on the back of the camera lights up green.

Note

The first time you connect a USB DCx camera to a USB port under Windows, two driver files will be registered. The first file (uc480 boot) contains the generic driver, the second file the model-specific driver.

The model will be immediately recognized whenever you connect the camera to this port again. If you use a different port, the registration will be repeated. Under Windows the camera will show up in the [uc480 Camera Manager's](#) camera list.



The DCx Cameras can be connected to a USB port either directly or via hubs and repeaters. A wide range of different hubs and repeaters are available commercially. The USB 2.0 hubs being used must be "full powered" hubs that are able to provide 500 mA per USB port. "Low Powered" hubs, in comparison, only supply 100 mA per port, which is not sufficient for DCx Cameras.

Note

To use maximum bandwidth, we recommend connecting the cameras directly to the USB ports on the mainboard. Many USB ports on PC/PCIe cards and the USB ports on the front of the PC often supply lower bandwidth.

Attention

USB cables with non-standard connectors must be connected to the camera first and then to the PC. Otherwise the camera might not be recognized correctly.

3.3 Application Notes by Camera Model

Cameras with CMOS sensors

- [DCC1240x / DCC3240x Application Notes](#)
- [DCC1545M Application Note](#)
- [DCC1645C Application Notes](#)

Cameras with CCD sensors

- [DC223x Application Notes](#)
- [DC224x Application Notes](#)

3.3.1 DCC3260x Application Notes



For the technical specifications of this model go to: [Camera and sensor data > DCC3260](#).

Sensor, pixel clock

- Changes on the capturing parameters (such as exposure or gain) only have an effect on the image after the next which is caused by the sensor. This is independent of the operating mode.
- The pixel clock can only be adjusted in steps. The steps get bigger with higher values.
- The adjustable pixel clock changes between 8/10 bit and 12 bit.
- Note the bandwidth requirement for bit depth greater than 8 bits per pixel.
- It is recommended to use packed color formats in 10 bit mode with high pixel clocks.
- If the camera is operated without image memory there are additional pixel clocks in 10 MHz steps. The image quality may be reduced in the lower image area because of the slow reading-out. It is recommended to operate the camera with image memory in case of bandwidth problems or in multi-camera operation.

AOI

- The sensor speed does not increase by the use of a horizontal AOI.

Frame rate

- In 10 bit mode, the sensor will not get faster than in 8 bit mode.

Binning

- The sensor does not support binning.

Subsampling

- Subsampling does not increase the frame rate.

Trigger

- The internal sensor delay is about 2-3 lines when triggering. The line period depends on the selected pixel clock. The higher the pixel clock, the smaller the line period is.

Vertical AOI merge mode

- For color sensors, the AOI merge mode can only be used with at least 2 lines.

White balance

- The white level may not be reached for very long exposure times and minimum gain. Increase the gain by one level if necessary.

3.3.2 DCC1240x / DCC3240x Application Notes

For the technical specifications of this model go to: Camera and sensor data > [DCC1240x / DCC3240x](#).

Shutter modes

The following table displays the four shutter modes and their advantages and disadvantages in different situations:

	Global shutter (default)	Global shutter (alternative timing)	Rolling shutter	Rolling shutter (global start)
Black level constancy	-	+	++	++
Capturing of moving objects	++	+	-	+
Hotpixels	-	-	+	+
Image quality with high gain	-	+	++	+

- The rolling shutter mode offers a better signal/noise ratio and a more consistent black level compared to the global shutter mode.
- If the sensor is used in global shutter mode at a low pixel clock frequency and a high gain, the bottom pixel rows might become brighter for technical reasons. Color distortion will occur for the color sensor with active white balance. In this case, use a higher pixel clock frequency, less gain or the rolling shutter mode.
- The "Global shutter (alternative timing)" mode offers a more consistent black level compared to the global shutter mode. This mode should not be used with a frame rate below 2 fps. This mode is also not suitable for bright, moving image contents in combination with long exposure times.
- The rolling shutter mode with global start is suitable for capturing moving objects with flash.
- When using flash in rolling shutter mode make sure to set the flash duration accordingly ($(1/\text{maxFramerate}) + \text{exposure}$) or that the global time window is available by a long exposure time ($2 * (1/\text{MaxFramerate}) + \text{FlashDuration}$). For flashing into this time window use the flash delay ($1/\text{MaxFramerate}$).
- The hardware sensor gamma curve is piecewise linear with three sections. This allows evaluating four times more details per pixel for lower gray level values and in 8 bit per pixel mode in dark image areas.
- In global shutter mode the shutter efficiency of 1:3000 have a negative impact with bright conditions and the usage of exposure times under 100 μs . In this case, set the pixel clock to the maximum possible value and close the aperture a little bit. Also enabling the Log mode with low values achieve huge improvements.

Black level

- The black level can also be set to negative values. Therefore, the factory setting of the offset control is nearly in the middle of the range.
- Use of the gain functions can lead to slight fluctuations of the black level. In global shutter mode the black level can also vary slightly.
- In global shutter mode the black level can also vary slightly between two image captures.

-
- When enabling the rolling or global shutter mode the black level is set to a fixed factory-provided value. Therefore, the black level can individually adjusted after switching the shutter mode.
 - Depending on the internal black level and shutter mode the offset control shows no additional cumulative effect at the top.
 - The factory setting of the offset control are so selected that the black level is always slightly increased to avoid losing image information by cutting underneath the origin. For linearity measurements the black level must be adjusted to the origin with the offset control before the measurement is done.

Color sensor

- The color sensor's black level cannot be adjusted manually, as the RGB gains are downstreamed and an adjustment would cause color errors.
- Automatic black level correction is always enabled.
- The RGB gains work analog.
- The [fast line scan](#) mode is disabled.

NIR sensor

- It is recommended to use a IR-coated and IR-corrected high-quality lens, especially for non-monochromatic light.
- In the high IR wavelength range picture blur can occur with strong contrasts. This reduced MTF (modulation transfer function) is a characteristic of the sensor pixels.
- The master gain of the NIR sensor is adapted in comparison to the monochrome sensor. If both sensors are compared directly, the gain of the NIR sensor must be set to the double factor of the monochrome sensor. This can be done via the master gain or the gain boost.
- In the Log mode the guaranteed dynamic range of the NIR sensor is reached with a gain value of at least 1. For the monochrome sensor you need at least a value of 3. With very short exposure times, lower values can result in even higher dynamics.

Gain, pixel clock

- Master gain uses a combination of coarsely scaled analog gain factors and finer digital scaling. To achieve optimum homogeneity of the gray level, use only the gain factors 0, 33, 66, and 100.
- The gain boost has the factor 2. When using the master gain a maximum factor of 8 is possible caused by the sensor.
- For global shutter mode the pixel clock should be set to the maximum possible value to increase the image quality.
- In the 10 bit mode the usage of the digital gain intermediate level produces missing pixel values as the sensor works internally with maximally 10 bit.

Hot pixel

- In the rolling shutter mode, there are less hot pixels, as the pixel charges are not buffered in the sensor.
- The sensor corrects hot pixels dynamically. Neighbourhood hot pixels in diagonal direction cannot be corrected effectively. These positions are covered by the factory-made hot pixel correction and are eliminated by the software hot pixel correction. Therefore, the hardware hot pixel correction is a prerequisite and should not be deactivated.
- On the color sensor the hot pixel correction works with the appropriate color neighbours.
- In the global shutter mode, increased hot pixels can become visible in the lower image area with

log exposure time, a high gain and disabled hot pixel correction.

- Both hot pixel corrections should be disabled if extremely fine structures are captured with a high-quality lenses.
- For the measurements of noise characteristics both hot pixel corrections should be disabled.
- The activation of the factory-made hot pixel correction reduces the frame rate slightly. Here, the desired frame rate has to set **after** the activation or deactivation of the factory-made hot pixel correction.
- Increased hot pixels can occur in the entire image border area.

Binning

- 2x binning makes the image brighter by a factor of about two. It also reduces image noise.
- Binning does not result in a higher frame rate. Using binning allows higher pixel clock frequencies for USB DCx cameras with the USB 2.0 and GigE interface. To achieve the maximum frame rate, activate first the binning and then change the maximum pixel clock frequency.
- When using binning the frame is slightly shifted horizontally.
- Binning can only be enabled for both horizontal and vertical pixels. For this reason, the parameters of the [is_SetBinning\(\)](#) function have to be passed together (`IS_BINNING_2X_VERTICAL` | `IS_BINNING_2X_HORIZONTAL`) to enable binning.

Scaling, AOI

- The digital scaling functions result in a higher possible frame rate. The maximum frame rate is increased approximately proportionally to the scaling factor. When using the scaling functions with USB 2.0 and GigE cameras, you can set higher pixel clock frequencies. To achieve the maximum frame rate, select first the scaling factor and then the maximum pixel clock frequency.
- For improved image quality without undersampling artifacts, the digital scaling feature permanently uses an anti-aliasing filter.
- On color sensors the scaler works in consideration the color information.
- Reducing the horizontal resolution does not result in a higher frame rate.
- The available step widths for the position and size of image AOIs depend on the sensor. The values defining the position and size of an AOI have to be integer multiples of the allowed step widths. For detailed information on the AOI grid see the [DCC1240x / DCC3240x](#) chapter.

Multi AOI

- When the Multi AOI function is enabled, no changes can be made to the image size settings (e.g. binning, subsampling, scaling).
- If sensor hotpixel correction and Multi AOI are enabled, the sensor displays a four pixel wide black line between the AOIs.

Line scan mode

- [Fast line scan](#) mode: The exposure time is fixed to the readout time of one sensor line. Exposure time cannot be changed in this line scan mode.
- Fast line scan mode: The time stamp is generated for the complete image.
- Fast line scan mode: There is a time gap between two frames with the line scan information. This corresponds to the duration of 15 lines at a frame rate set to maximum.
- In the fast line scan mode color images are not possible as [Bayer color sensors](#) needs at least two neighboring lines for color calculating. Therefore, only monochrome models support the line scan mode.

Log mode

- The Log mode shows visible effect only for short exposure times (< 5 ms)
- In global shutter mode the Log mode can help to increase the shutter efficiency for extremely short exposure time (< 100 μ s).
- To find the right operating point of the Log mode use the following procedure:
 1. Set the Log mode gain to the minimum value.
 2. Find the operating point via the Log mode value. The higher the value, the more bright image areas are damped and more details become visible.
 3. Set the image as bright as possible via the Log mode gain. A typical value is 2 or 3 for monochrome sensors and 0 or 1 for NIR sensors.
- The master gain is disabled in Log mode.
- See ThorCam: Settings > Shutter > Log Mode

Anti blooming

- Activation: With exposure times over 10 ms and no use of gain or gain boost it can occur that bright image areas do not reach saturation and so no white level which is caused by the sensor. A visible, firm pattern is formed in bright image areas. For color sensors with enabled white balancing bright image areas gets purple. In this case, disable the anti blooming mode or increase the master gain from 1x to 1.5x.
- The anti blooming mode should not be enabled for short exposure times (< 5 ms). Depending on the shutter mode a "Black Sun" effect or overexposure occurs.

Micro lenses

- The sensor has non-removable micro lenses on each pixel. These micro lenses focus the incoming light for the subjacent smaller photodiode. This lens has a directive efficiency.
- To the corners of the active image area the micro lenses are slightly shifted to the photodiode. So the unavoidable non-vertical light incidence of C mount lenses is compensated. The shift is constant from the center to the corner and has a maximum of 12 degrees. When using a telecentric lens or parallel light incidence the shift must be considered as little shading effects may occur.

3.3.3 DCC1545M Application Notes

For the technical specifications of this model go to: Camera and sensor data > [DCC1545M](#).

Sensor

- Sensor speed does not increase for AOI width <240 pixels.
- Extreme overexposure shifts the black level. Please deactivate the Auto offset function in this case.
- At very long exposure times and minimum gain, the white level may not be reached. The gain should be increased by one step in this case.
- Monochrome version only: The sensor internally works like the color version. This might lead to artifacts when subsampling is used.
- The brightness of the first and last line might deviate due to the sensor.
- Gain values between 59 and 99 may lead to image inhomogeneity.
- When using very narrow AOIs, the sensor may not be able to calculate the correct black level. Use manual black level offset when problems with the black level occur.

Calibration

- Cameras with a [date of manufacture](#) after Dec. 9, 2008: The offset control has been calibrated internally. The calibration corrects offset errors when gain is used. In calibrated cameras, automatic black level correction is disabled by default. The calibration can only be used with uc480 driver version 3.31 or higher.
- Cameras with a [date of manufacture](#) before Dec. 9, 2008: If manual offset control is used, fixed pattern noise and horizontal lines may become visible. High gain values may shift the black level and therefore should be avoided. Offset increases the black level every 7th step. The steps in-between change the appearance of fixed pattern noise.

3.3.4 DCC1645C Application Notes

For the technical specifications of this model go to: Camera and sensor data > [DCC1645C](#).

Sensor

- At very long exposure times and minimum gain, the white level may not be reached. The gain should be increased by one step in this case.
- The RGB gain controls have no effect for values >90.

3.3.5 DCU223x Application Notes

For the technical specifications of this model go to: Camera and sensor data > [DCU223x](#).

Sensor

- Long exposure times will increase the number of hotpixels.
- High temperatures will increase the black level of individual pixels.

3.3.6 DCU224x Application Notes

For the technical specifications of this model go to: Camera and sensor data > [DCU224x](#).

Sensor

- Long exposure times will increase the number of hotpixels.
- High temperatures will increase the black level of individual pixels.
- When vertical 4x binning is activated, the minimum image width increases to 640 pixels.

3.4 Installed uc480 Programs

- [uc480 Camera Manager](#): The central tool for managing all connected DCx cameras.
- ThorCam: A comprehensive viewer for exploring the camera functionality.
- [uc480 Hotpixel Editor](#): A tool to edit the sensor hot pixel list stored in the camera.

3.4.1 uc480 Camera Manager

The DCx Camera Manager is the central tool for managing all DCx cameras. It displays information on the connected USB DCx Cameras and provides options for configuring them.

On Windows systems the DCx Camera Manager can be accessed as follows:

- Start > Programs > Thorlabs > Scientific Imaging > DCx Camera Manager


Click in the figure to get help on the functions.




- Camera list
The [camera list](#) displays information on the connected DCx Cameras.

Attention

Under Windows XP (64 bit) and Windows Vista DC3240x Cameras are displayed in the DCx Camera Manager, but they cannot be opened, because they are supported by Windows 7 and Windows XP (32 bit) only

- Control center
In the [control center](#), you can access the configuration and display detailed information on the connected DCx Cameras.
- In the  drop down box, you can choose the language for the DCx Camera Manager. This setting is saved and remains effective even after you close and reopen the program.
For proper display of Asian languages on Windows systems, the support package for East Asian languages has to be installed on your system (in "Control Panel > Regional and Language Options").

- Click  to close the application; any settings you have made are saved.
- The status box at the bottom indicates the current status of the selected camera. If it is available, the status message is shown in black. Otherwise, the status message is shown in red. If an error has occurred in a camera, a black exclamation mark on a yellow background is shown next the camera. The status box then indicates the cause of the error and suggests remedies.

3.4.1.1 Camera List


When a camera is activated (switched on or connected to the PC), it appears in the camera list of the uc480 Camera Manager after a few seconds.

Free	Avail.	Type	Cam.ID	Dev.ID	Model	SerNo.
------	--------	------	--------	--------	-------	--------


The data shown in the camera list can be sorted in ascending or descending order by left-clicking on the respective column header.

- Free/Avail.
Free: indicates whether a camera is currently in use.
Avail. (Available): indicates whether a camera can be opened by this computer with the current setup (computer and camera).
Cameras shown with a red x are currently in use (Free = No) and are not available (Avail. = No).
Cameras shown with an exclamation mark are not in use, but are currently unavailable for various reasons, such as:
 - The camera is not compatible with the driver. Please update the uc480 driver.
 - The driver has not properly detected (initialized) the camera. Please disconnect the camera from the PC and then reconnect it.
 - The camera is currently being removed from the Manager.
 - The camera reports that it is "Not operational".
- Type
This column indicates the USB camera type.
- Cam.ID
The [camera ID](#) assigned by the user.
- Dev.ID
Unique device identifier sequentially assigned by the system. DCx cameras are assigned device IDs from 1 upwards. After deactivating a DCx camera (switching it off or disconnecting it from the network), the device ID is no longer valid and can be assigned again by the system.
- Model
Model name of the camera
- SerNo.
Serial number of the camera.

3.4.1.2 Control Center

- Expert mode
When you select the  check box, the uc480 Camera Manager additionally displays the [Parameters](#) box on the right. There you will find detailed information on the DCx camera selected in the [camera list](#).

Click in the figure to get help on the functions.



- Automatic parameter refresh


If you select the Automatic parameter refresh check box, the data shown in the tree structure is updated periodically. If the option is disabled, the data in the tree structure is only updated when a different camera is selected.

All other "Control Center" buttons are described in detail in the following sections:

- [General information](#)
- [Camera information](#)
- [Creating a support file](#)
- [Additional functions \(COM port\)](#)

3.4.1.3 General Information


This dialog box provides information on the installed uc480 drivers and the available USB controllers and network adapters.



- **uc480 drivers**
This list shows the location and version of the uc480 driver files installed on your system.
- **3rd party drivers**
This list shows the location and version of the uc480 interface files that have been installed on your system for third-party software.
- **USB controller and network adapters**
All USB controllers and network adapters that are available in your system are shown in a tree structure.

3.4.1.4 Camera Information

In the "Camera information" dialog box, you can assign a unique ID to the selected camera and write to the user area of the EEPROM. The data you enter is retained in the camera memory even when the camera is disconnected from the PC or power supply.



- **Camera ID**
The camera ID identifies a camera in multi-camera operation. The ID can range from 1 to 254. The default value for the camera ID is 1. The same ID can be assigned to multiple cameras. You do not have to assign sequential ID numbers to all connected cameras.
- **User EEPROM (max. 64 characters)**
Every uc480 has a 64-byte user area in its EEPROM (Electrically Erasable and Programmable Read Only Memory) to which you can write text of your choice.

The "Camera information" dialog box displays two additional boxes that are for your information only and cannot be edited:


- Manufacturer (e.g.OEMINC)
- Date of QC (date of final camera quality test)

Notes

1. Setting a camera ID and writing to the EEPROM is possible only, if the camera is marked "Free" and "Available" in the Camera Manager (see also [Camera list](#)).
2. If software accesses cameras through the uc480 DirectShow interface, the camera IDs must be in a range from 1 to 24.
3. If software accesses cameras through the uc480 Cognex VisionPro interface, the camera IDs must be assigned consecutively beginning with 1.

3.4.1.5 Creating a Support File

A uc480 support file is a binary file with the extension `.bin`. The file contains camera and driver details that are required for diagnostics by our [Technical Support](#) team. No personal computer data or user data is stored in this file.

The  button opens the "Save as" dialog box, where you can save the displayed camera information and additional driver information to a file.

3.4.1.6 Additional Functions

CPU idle states

Windows only: Processor operating states (idle states/C-states)

Modern processors have various operating states, so-called C-states, that are characterized by

different power requirements. When the operating system selects an operating state with low power consumption (unequal C0), the USB transmission efficiency may be affected (see also [is_Configuration\(\)](#) and section [troubleshooting](#)).

Camera parameters when camera is opened

Here, you can set whether to apply the parameters stored on the camera automatically when opening the camera. You must first store the camera parameters on the camera using the [is_ParameterSet\(\)](#) function.

This setting applies to all connected cameras. If no parameters are stored on the camera, the standard parameters of this camera model are applied (see also [is_Configuration\(\)](#)).

Boot boost

This mode is not related to DCx Cameras.

Bulk Transfer Size

Via "Bulk Transfer Size" the behavior of the USB sub-system can be set.

Warning

Contact our [technical support](#) before changing the value under "Bulk Transfer Size".

COM Ports


The "Additional functions" dialog box allows installing virtual COM ports for communication through the serial interface of a DCC3240x camera. The following sections show you how to set up and use the [serial interface](#).

Note

This feature is only available for DCC3240x cameras.




You need administrator privileges to install a virtual COM port.

The selected in the uc480 Camera Manager DCC3240x camera has to be marked "Free" and "Available".



Setting up the serial interface on the DCC3240x camera

Before using the serial interface on the camera, one or more virtual COM ports have to be installed on the PC. Most systems support up to 255 COM ports; COM1 to COM8 are often assigned operating system functions by default. You can check the current port assignment in the Device Manager on your computer. Some older systems may not have more than eight ports; in that case you will need to assign the camera to one of these ports.

- **COM port**
In the drop down box, select the number of the port you want to install (default: 100). COM ports in use are marked "(used)" in the list.
- 
Click this button to install the selected virtual COM port.
During the first installation of a virtual COM port, an additional broadcast port with number 255 is installed. Data sent to this port will be forwarded to all paired cameras.
You can install any number of virtual COM ports on a single system.
- 
With this button, you can release a COM port that is marked "used." If the port number has been saved in that camera, it will be deleted in the camera, too. To release a COM port, select it in the drop down box and then click this button.
- 
Click this button to assign the selected port number to the camera. The port number is saved in the camera's non-volatile memory and retained even when the camera is switched off. You can look up the assigned port number in the Camera Manager's expert mode. A COM port number can also be saved in a camera without a virtual COM port installed on the PC.

Note

If you want to control more than one DCC3240x camera from a PC, each camera should be assigned a unique port number. If multiple cameras are assigned the same port number, only the port of the first camera that is opened will be used.

To send data via the serial interfaces of multiple cameras, you can use the broadcast port with number 255. Before connecting to the broadcast port, ensure that all the cameras that are to receive the broadcast have been opened.

Testing the serial interface on the DCC3240x camera


Note

To avoid transmission errors, please ensure that both the camera and the receiving end use the same communication parameters (baud rate, data bits, stop bits, parity). Further information on the communication parameters is provided in the [Serial interface DCC3240x](#) chapter.

- 

Clicking this button opens a dialog box for transferring data through the COM port. The dialog box is provided as the `uc480ComportDemo.exe` sample program together with the C++ source code and is included in the uc480 SDK.

This program allows sending ASCII characters to the COM port assigned to a camera. The characters are output unchanged on the camera's serial port. To check the proper functionality, you can connect a PC to the camera's serial port and read the transmitted characters on the PC's COM port.



- Baud
In this drop down box, you can change the data transfer rate of the serial interface.
- Append
This drop down box allows appending the special characters "CR" (Carriage Return) and "LF" (Line Feed) to the ASCII text you want to transmit. Some devices with serial interface require ASCII strings to be terminated with CR/LF.
- Send file
Using these functions, you can send a file in either direction (output on the camera's virtual COM port or output on the PC's COM port).


Note

Since the sample program has to open the camera, please make sure the selected camera is not used by other applications at the same time.

3.4.1.7 Parameters

This box displays the parameters of the camera you have selected in the camera list. The parameters box is only shown when [Expert mode](#) is active.

The parameters are organized in a tree structure. Only the information that applies to the selected camera is shown. The data displayed in the camera list is not repeated in the "Parameters" box. The data shown in the tree structure cannot be changed.



The screenshot shows a window titled "Parameters" with a table of camera parameters. The table has two columns: "Parameter" and "Value". The parameters are organized into a tree structure with expandable sections: "Device", "USB", and "Local driver".

Parameter	Value
Device	
Sensor ID	81
USB	
1. Hub	USB Root Hub (Port 2)
2. Controller	Intel(R) ICH9 Family US...
Local driver	
USB version	4.20.06.01

- Device
 - Sensor ID
- USB
 - Hub
 - Indicates which hub and port a USB camera is connected to. In addition, the full path through all hubs to the USB controller on the computer is displayed.
 - Controller
 - Indicates the USB controller to which the camera is connected.
- Local driver
 - Indicates the USB version of the camera driver


3.4.2 uc480 Hotpixel Editor

During manufacturing, each DCx camera is tested and calibrated for hot pixels (aka bad pixels), which are caused by technical reasons (see also [Camera basics: Sensor: Hot pixels](#)). In some cases, users may wish to extend this factory calibration. Using the uc480 Hotpixel Editor, you can now edit the sensor hot pixel list stored in the camera.

The uc480 Hotpixel Editor can be accessed as follows:

- C:\Program Files\Thorlabs\Scientific Imaging\DCx Camera Support\Program\HotpixelEditor\uc480HotpixelEditor.exe

After program start, the uc480 Hotpixel Editor shows the following window:



In the toolbar on the top the following functions are provided:

- Open camera
Opens the camera that is marked in the camera list.
- Close camera
Closes the camera that is marked in the camera list.
- Badpixel correction
Turns the hot pixel correction on/off.
- [Detect Badpixel](#)
Opens the dialog for the automatic hot pixel detection.
- Write Badpixel
Writes the list of the hot pixel to the EEPROM.

- Add Badpixel
Writes the marked hot pixel in the hot pixel list of the program.
- Delete Badpixel
Deletes a hot pixel from the hot pixel list of the program.
- Highlight
Mark bad pixel with a circle in order to improve localization. Yellow marker stands for a factory calibrated bad pixel, green marker - for user defined bad pixel.
- 1:1
Zoom factor of the shown image.

Camera and bad pixel list

Tab Camera List:

Recognized cameras are shown with status in the camera list. Via the context menu (right click to the appropriate camera), a camera can be opened and closed. Cameras with the status "in use" cannot be opened.


Tab Badpixel List:

In the bad pixel list the coordinates and the class of bad pixels are displayed.

Class "user": added by the user


Class "ignore": factory calibrated hot pixels

Magnifier function




The magnifier allows to exactly mark hot pixels. "Value" displays the relative pixel intensity.

Camera properties



Here, the camera settings can be changed for a better hot pixels visualization. By default, exposure time is set to max. and master gain to 50. Note that the frame rate is limited to lower values.

Detecting hot pixel automatically



Detect hotpixel opens the "Badpixel detection" dialog. Bad pixels are detected and can be added to the "user" class list.

- "Threshold": Any pixel with an intensity exceeding that threshold will be recognized as a bad (hot) pixel, where the threshold is the minimum intensity difference to the reference intensity. The reference is the averaged over the most recent 10 frames total pixel intensities. Threshold can be set between 5 and 255.
- The "Clear old badpixel list" option removes the old user list before the search.

Note

The number of bad pixels, stored to the EEPROM, is limited. If the number of detected bad pixels exceeds the memory capacity, please increase the bad pixel detection threshold and repeat the detection procedure.

4 Programming (SDK)

In addition to the drivers, the uc480 software Development Kit (SDK) includes the [uc480 Camera Manager](#), the uc480 Viewer and the uc480 API programming interface for creating your own DCx programs under Windows and Linux. Numerous demo applications make it easy for you to get started with programming. The uc480 API offers you over 150 commands with which you can access all the parameters and functions of your DCx camera. This chapter contains all the information you need to integrate the DCx camera in your own applications using the uc480 API.

Note

Older Functions

We are continuously extending and enhancing the uc480 API. The resulting product upgrades sometimes require replacing obsolete functions with new ones. If it is necessary to continue working with the older functions, it is possible to add the `uc480_deprecated.h` header file additionally to the `ueye.h` header file. The `uc480_deprecated.h` header file contains all obsolete function definitions and constants which are no longer part of the `ueye.h` header file.

See also:

- [First steps to uc480 programming](#)
- [How to proceed](#)
- [Function descriptions](#)
- [AVI function descriptions](#)
- [Obsolete functions](#)
- [Programming notes](#)
- [Lists](#)

4.1 First Steps to uc480 Programming

This chapter shows the most important functions of the uc480 API for integrating your camera into your own applications. You will find comprehensive lists of the API functions, sorted by task, in the [How to proceed](#) chapter.

The **uc480 SimpleLive** and **uc480 SimpleAcquire C++** programming samples included in the SDK illustrate the steps described below.

For information on required include files (uc480 API and header) see [Programming notes](#) chapter.



▣ Select a display mode

The uc480 API provides different modes you can use to display the camera's images on the PC. To quickly show a live image under Windows, it is easiest to use the [Direct3D mode](#). Under Linux the [OpenGL mode](#) can be used.

This mode has the advantage that no image memory has to be allocated, and that image capture is handled by the driver. Call [is_SetDisplayMode\(\)](#) to select the display mode. You can then customize the Direct3D mode by using [is_DirectRenderer\(\)](#).

For advanced users:

You can also access the image data directly by selecting the Bitmap (DIB) mode. To use DIB mode, you first have to allocate one or more memories by using [is_AllocImageMem\(\)](#), add them to a [memory sequence](#), if required, and then activate a memory with [is_SetImageMem\(\)](#) before each image capture. To show the image on-screen, call the [is_RenderBitmap\(\)](#) function after each completed image capture. From the [events or messages](#) you can see when an image is available for display.

See also:

- How to proceed: [Display mode selection](#)



▣ Capture images

Recording live images with the DCxCamera is very simple. Just call the [is_CaptureVideo\(\)](#) function and the camera captures the live images at the default frame rate. To capture single frames, use the [is_FreezeVideo\(\)](#) function. Every DCx camera of course also provides different trigger modes for image capture. Use [is_SetExternalTrigger\(\)](#) to activate the desired mode before starting the image capture.

See also:

- How to proceed: [Image capture](#)



▣ Adjust the frame rate, brightness and colors

To change the frame rate, for example, you call [is_SetFrameRate\(\)](#). With [is_SetColorMode\(\)](#) you set the color mode. Image brightness is adjusted through the exposure time set with [is_Exposure\(\)](#). You can also implement [automatic control](#) of image brightness and other parameters by using [is_SetAutoParameter\(\)](#).

If you are using a color camera, you should activate color correction in order to achieve rich vibrant colors for on-screen display ([is_SetColorCorrection\(\)](#)). To adapt a color camera to the ambient light conditions, it is essential to carry out white balancing. This is also done using the [is_SetAutoParameter\(\)](#) function.

See also:

- How to proceed: [Setting camera parameters](#)



▣ Save an image

Use the [is_ImageFile\(\)](#) function to save the current image as a BMP or JPEG file. To save a specific image, it is better to use the Snap function (single frame mode) than the Live function (continuous mode).

See also:

- How to proceed: [Saving images and videos](#)




▣ Close the camera

When you want to exit your application, close the camera with [is_ExitCamera\(\)](#). The camera and the allocated memory are automatically released. All previously set camera parameters will be lost, however. So, if you want to save specific settings, use the [is_ParameterSet\(\)](#) function before closing the camera. The next time you start the application, you can simply load the settings again by using the same function.

4.2 How to Proceed

This chapter shows function blocks and workflows for important camera functions. The charts are structured as follows:



Hint

Click the function names in the flowcharts to open the corresponding function description!

4.2.1 Preparing Image Capture

- [Querying information](#): Before you open one or more DCx Cameras, we recommend querying some key information.
- [Opening and closing a camera](#)
- [Allocating an image memory](#): This is necessary if you want to access image memory contents or if you are not using Direct3D or OpenGL for image display.
- Creating an [image memory sequence](#) is required when capturing live images.

4.2.1.1 Querying Information

It is recommended to query the following important information before opening one or more DCx Cameras.

is_GetNumberOfCameras()	Determines the number of cameras connected to the system.
is_GetCameraList()	Returns information on all connected cameras.

It is also very useful to have the message boxes for error output enabled during the programming process

is_SetErrorReport()	Enables/disables dialog messages for error output.
-------------------------------------	--

With the following functions, you can read out additional information on cameras and software.

is_CameraStatus()	Returns the event counters and other information. Enables standby mode.
-----------------------------------	--

is_GetCameraInfo()	Returns the camera information of an opened camera.
is_GetDLLVersion()	Returns the version of the <code>uc480.dll</code> .
is_GetOsVersion()	Returns the operating system version.

4.2.1.2 Opening and Closing the Camera

The following functions are required to open and close a DCx camera.

is_InitCamera()	Hardware initialization
is_ExitCamera()	Closes the camera and releases the created image memory.

When multiple cameras are used on one system you should assign every camera a unique camera ID.

is_SetCameraID()	Sets a new camera ID.
----------------------------------	-----------------------

4.2.1.3 Allocating Image Memory

When you are programming an application that

- requires direct access to the image data in stored in memory, or
- uses Bitmap mode (DIB) for display

use the following functions to allocate and manage image memories (see also [Quick start: Image display](#)).

is_AllocImageMem()	Allocates an image memory.
is_SetAllocatedImageMem()	The user provides pre-allocated memory for image capturing.
is_FreeImageMem()	Releases an allocated image memory.

An image memory has to be activated before each image capture:

is_SetImageMem()	Makes an image memory active.
----------------------------------	-------------------------------

To query image memory information and access the data in the image memories, you can use these functions:


is_CopyImageMem()	Copies the image to the user-defined memory.
is_CopyImageMemLines()	Copies selected image lines to the user-defined memory.
is_GetActiveImageMem()	Returns the number and address of the active image memory.
is_GetImageMem()	Returns the pointer to the starting address of the image memory.
is_GetImageMemPitch()	Returns the line offset used in the image memory.
is_InquireImageMem()	Returns the properties of an image memory.

Note

[Image memory sequences](#) should be used for frame sequence capture.

Flowchart: Allocating memory

Click in the figure to get help on the functions.




Calculate image buffer size

Allocate user memory

Get pointer to user memory

Activate user memory for the uc480 SDK. User memory pointer must be passed.



4.2.1.4 Image Memory Sequences

When you are capturing and displaying frame sequences (e.g. live display), it is advisable to provide appropriate image memory sequences. The uc480 driver offers a set of easy-to-use features for this purpose. For example, the system automatically cycles through the specified sequence of image memories and can generate an [event](#) when it reaches the end of a sequence cycle.

Before you can use a memory sequence, you have to allocate the relevant image memories (see [Allocating image memory](#)).

is_AddToSequence()	Adds image memory to the sequence list.
is_ClearSequence()	Deletes the entire sequence list.
is_GetActSeqBuf()	Determines the image memory currently used for the sequence.
is_SetImageMem()	Makes the indicated image memory the active memory.
is_LockSeqBuf()	Protects the sequence image memory from being overwritten.
is_UnlockSeqBuf()	Releases the sequence image memory for overwriting.

4.2.2 Selecting the Display Mode

Note

The uc480 driver provides different modes for displaying the captured images. We recommend using the Bitmap mode or the Direct3D functions (only Windows) or OpenGL functions, depending on your specific application.

For further information on the different display modes, see [Basics: Image display modes](#).

Select the desired mode. The display mode has to be set before you start image capture.


is_SetDisplayMode()	Selection of the display mode.
-------------------------------------	--------------------------------

When Bitmap mode (DIB) is active, image display has to be called explicitly for each image.

is_RenderBitmap()	Outputs the contents of the active image memory to a window.
is_SetDisplayPos()	Enables offsetting the image output inside the window.

Flowchart: Image display in DIB mode

Click in the figure to get help on the functions.



Find the buffer ID according to the pointer.

Lock buffer

Perform processing / display.

Optional: Copy image buffer to the application memory for processing

Unlock the buffer

4.2.3 Capturing Images

- DCx Cameras support the capture of single frames (snap) and frame sequences (live) in trigger mode and untriggered (freerun) mode. Select the [image capture mode](#) that best meets your application requirements.
- Using [events or messages](#), the uc480 driver can provide information to an application, e.g. about the completion of image capture. You will need events and messages, for example, when you are using [image memory sequences](#).

4.2.3.1 Image Capture Modes

For more information on the capture modes of the DCx Cameras see also [Camera basics: Freerun](#) and [Camera basics: Trigger](#).

Freerun mode

In freerun mode, the camera sensor captures one image after another at the set frame rate. Exposure of the current image and readout/transfer of the previous image data are performed simultaneously. This allows the maximum camera frame rate to be achieved. The frame rate and the exposure time can be set separately. The captured images can be transferred one by one or continuously to the PC.

If trigger mode is active, you need to disable it with [is_SetExternalTrigger\(\)](#) before activating freerun mode.

- Single frame mode (snap mode)
When [is_FreezeVideo\(\)](#) is called, the next image exposed by the sensor is transferred. You cannot use the DCxCamera flash outputs in this mode.
- Continuous mode (live mode)
When [is_CaptureVideo\(\)](#) is called, images are captured and transferred continuously. You can use the DCxCamera flash outputs.

Trigger mode

In trigger mode, the sensor is on standby and starts exposing on receipt of a trigger signal. A trigger event can be initiated by a software command (software trigger) or by an electrical signal via the camera's digital input (hardware trigger). For the specifications of the electrical trigger signals, see the [Specifications: Electrical specifications](#) chapter.

The trigger mode is selected using [is_SetExternalTrigger\(\)](#).

- Software trigger mode
When this mode is enabled, calling [is_FreezeVideo\(\)](#) immediately triggers the capture of an image and then transfers the image to the PC. If [is_CaptureVideo\(\)](#) is called, the triggering of image capture and the transfer of images are performed continuously.
- Hardware trigger mode
When this mode is enabled, calling [is_FreezeVideo\(\)](#) makes the camera ready for triggering just once. When the camera receives an electrical trigger signal, one image is captured and transferred.
If you call [is_CaptureVideo\(\)](#), the camera is made ready for triggering continuously. An image is captured and transferred each time an electrical trigger signal is received; the camera is then ready for triggering again (recommended procedure).
- Freerun synchronization
In this mode, cameras running in freerun mode (live mode, see above) can be synchronized with an external trigger signal. The cameras still remain in freerun mode. The trigger signal stops and restarts the current image capture process. You can use this mode to synchronize multiple

cameras that you are operating in the fast live mode. Not all camera models support this mode (see [is_SetExternalTrigger\(\)](#)).

Notes

- The freerun synchronization mode is currently not supported by DCx Cameras.
- In trigger mode, the maximum frame rate is lower than in freerun mode because the sensors expose and transfer sequentially. The possible frame rate in trigger mode depends on the exposure time.

Example: At the maximum exposure time, the frame rate is about half as high as in freerun mode; at the minimum exposure time, the frame rate is about the same.

Overview on image capture modes

Image capture	Trigger	Function calls	Allowed flash modes		Frame rate
			Standard	Global Start	
Continuous	Off	is_SetExternalTrigger (OFF) is_CaptureVideo()	X		Freely selectable
	Software	is_SetExternalTrigger (SOFTWARE) is_CaptureVideo()	X	X	Depending on exposure time and trigger delay
	Hardware	is_SetExternalTrigger (e.g. HI_LO) is_CaptureVideo()	X	X	Depending on exposure time and trigger delay
	Freerun sync.	is_SetExternalTrigger (e.g. HI_LO_SYNC) is_CaptureVideo()	X		Freely selectable
Single frame	Off	is_SetExternalTrigger (OFF) is_FreezeVideo()			Freely selectable
	Software	is_SetExternalTrigger (SOFTWARE) is_FreezeVideo()	X	X	Depending on exposure time and trigger delay
	Hardware	is_SetExternalTrigger (e.g. HI_LO) is_FreezeVideo()	X	X	Depending on exposure time and trigger delay

Timeout values for image capture

When you call [is_FreezeVideo\(\)](#) or [is_CaptureVideo\(\)](#), the timeout value for the image capture is determined from the `Wait` parameter. If no image arrives within this timeout period, a timeout error message is issued. Under Windows, a dialog box is displayed if you have enabled error reports (see [is_SetErrorReport\(\)](#)). Information on the error cause can be queried using `is_CaptureStatus()`.

The following table shows the effect of the Wait parameter depending on the image capture mode:

Parameter Wait	Image capture mode	Function returns	Timeout for 1st image	Timeout for subsequent images* ¹
IS_DONT_WAIT	HW trigger	Immediately	API default or user-defined value* ³	API default or user-defined value* ³
IS_WAIT	HW trigger	When 1st image in memory	API default or user-defined value* ³	API default or user-defined value* ³
Time t Value range [4... 429496729]	HW trigger	When 1st image in memory	Time t in steps of 10 ms (40 ms to approx. 1193 h)	API default or user-defined value* ³
IS_DONT_WAIT	Freerun/SW trigger	Immediately	Calculated internally by API* ²	Calculated internally by API* ²
IS_WAIT	Freerun/SW trigger	When 1st image in memory	Calculated internally by API* ²	Calculated internally by API* ²
Time t Value range [4... 429496729]	Freerun/SW trigger	When 1st image in memory	Time t in steps of 10 ms (40 ms to approx. 1193 h)	Calculated internally by API* ²

*1 Only with continuous image capture using `is_CaptureVideo()`

*2 The timeout is calculated from the exposure time setting, the image transfer time (depending on the pixel clock) and the optional trigger delay (see `is_SetTriggerDelay()`); it is at least 40 ms.


*3 The default value of the uc480 API is 60 s. User-defined values can be set using the `is_SetTimeout()` function.

Function list

<code>is_CaptureVideo()</code>	Captures a live video.
<code>is_FreezeVideo()</code>	Captures an image and writes it to the active image memory.
<code>is_ForceTrigger()</code>	Forces image capture in hardware trigger mode.
<code>is_HasVideoStarted()</code>	Returns whether the capture process has been started or not.
<code>is_IsVideoFinish()</code>	Returns whether the capture process has been terminated or not.
<code>is_SetSensorTestImage()</code>	Enables test image output from sensor (all cameras).
<code>is_StopLiveVideo()</code>	Terminates the capturing process (live video or single frame).

Flowchart: Single Capture

Click in the figure to get help on the functions.



Open Camera an get device handle

Get useful information about the sensor and the camera

- Sensor resolution
- Colour format (monochrome, Bayer)
- Rolling/global shutter
- Camera model


Set color mode for the acquisition
(e.g. IS_SET_CM_Y8, IS_SET_CM_RGB24)

Set image size or AOI

Memory allocation by driver: continue with 6
Memory allocation by user: continue with 7


Set Display mode to DIB

Flowchart - Single Capture (1 of 2)




Flowchart - Single Capture (2 of 2)

Flowchart: Sequence Capture



Flowchart - Capturing a frame sequence (1 of 2)



Flowchart - Capturing a frame sequence (2 of 2)


4.2.3.2 Event / Message Handling

Triggering events for single frame capture

The following figure shows the time sequence when triggering the `IS_SET_EVENT_EXTTRIG` and `IS_SET_EVENT_FRAME` events. The camera is prepared for triggered capture using the `is_SetExternalTrigger()` command. An incoming trigger signal at the camera starts the exposure and the subsequent image transfer. Upon completion of the data transfer, the `IS_SET_EVENT_EXTTRIG` event signals that the camera is ready for the next capture. The `IS_SET_EVENT_FRAME` event is set once pre-processing (e.g. color conversion) is complete and the finished image is available in the user memory.

Note


The following illustrations show a schematic view of the image capture sequence. The sensor exposure and readout times and the transmission times depend on the camera settings. The pre-processing time depends on the API functions you are using (e.g. color conversion, edge enhancement).



* Optional function. The start time and duration of the flash signal are defined by the "Flash delay" and "Duration" parameters (see `is_IO()`).

Events in live mode (image sequence)

The following figure shows the time sequence when triggering the `IS_SET_EVENT_FRAME` and `IS_SET_EVENT_SEQ` events. The camera is set to live mode using `is_CaptureVideo()` so that it continuously captures frames. The `IS_SET_EVENT_FRAME` event is set once pre-processing (e.g. color conversion) is complete and a finished image is available in the user memory. The `IS_SET_EVENT_SEQ` event is set after one cycle of a storing sequence has been completed (see also `is_AddToSequence()`).




* Optional function. The start time and duration of the flash signal are defined by the "Flash delay" and "Duration" parameters (see [is_IO\(\)](#)).

Function List

is_DisableEvent()	Disables a single event object.
is_EnableEvent()	Enables a single event object.
is_EnableMessage()	Turns the Windows messages on/off.
is_ExitEvent()	Closes the event handler (Windows only)
is_InitEvent()	Initializes the event handler (Windows only)
is_EnableAutoExit()	Automatically releases the camera resources when the camera is disconnected from the PC.
is_WaitEvent()	Waits for DCxCamera events (Linux only)

Flowchart: Enable events

Click in the figure to get help on the functions.



Create an event and receive the handle of the event object



Initialize the uc480 event handling. The event handle from CreateEvent is required.

Enable a specific event.

Using a thread (recommended) continue with **5** else continue with **6**.

Windows functions:
WaitForSingleObject()
WaitForMultipleObjects()

Flowchart: Enabling Messages



MS Windows only

Enable the message handling.

This may be different in each programming language. Hook on the message map.

When the messages are sent to the applications main window and the GUI is locked for some reasons the uc480 messages cannot be processed.

Disable the message handling.

4.2.4 Setting Camera Parameters

- [Setting and getting parameters](#): Using these functions, you can make settings for the camera and for image capture and preprocessing.
- The DCxCamera's [automatic image control](#) features allow automatically adjusting image brightness and image color to changing ambient conditions.
- [Image preprocessing](#): These functions specify e.g. how color images are processed after image capture.
- [Querying the camera status](#): With these functions, you can query additional useful information on the camera status.
- [Using the camera EEPROM](#): All DCx Cameras have a non-volatile EEPROM where you can save the camera settings or any other information.

4.2.4.1 Setting and Getting Parameters

Capture parameters

This set of functions specifies the camera's image capture parameters, such as exposure, pixel clock and frame rate:

is_ColorTemperature()	Sets the color temperature
is_Exposure()	Returns the adjustable exposure range.
is_GetFramesPerSecond()	Returns the current frame rate in live mode.
is_GetFrameTimeRange()	Returns the adjustable frame rate range.
is_PixelClock()	Returns the adjustable pixel clock range.
is_SetAutoParameter()	Enables/disables automatic imaging functions.
is_Blacklevel()	Turns black level correction on / off.
is_Exposure()	Sets the exposure time.
is_SetFrameRate()	Sets the frame rate.
is_SetGainBoost()	Sets additional sensor hardware gain boost.
is_SetGamma()	Sets the gamma value (digital post-processing).
is_SetHardwareGain()	Enables the sensor hardware gain.
is_SetHWGainFactor()	Sets the sensor hardware gain factor.
is_PixelClock()	Sets the pixel clock frequency.
is_ResetToDefault()	Resets the camera parameters to its default values.

Image geometry

This set of functions lets you influence the image geometry for image capture, e.g. the area of interest:

is_ImageFormat()	Sets a predefined image size
is_AOI()	Sets the size and position of an area of interest (AOI) or of a reference AOI for auto imaging functions.
is_SetBinning()	Sets the binning modes.
is_SetRopEffect()	Makes real-time geometry changes to an image

	(Rop = raster operation)
<u>is_SetSensorScaler()</u>	Scales the image in the camera
<u>is_SetSubSampling()</u>	Sets the subsampling modes.


Processing image data

The following set of functions refers to the further processing of image data in the PC:

<u>is_GetColorDepth()</u>	Determines the desktop color mode set in the graphics card.
<u>is_GetTimeout()</u>	Returns the user-defined timeout values.
<u>is_HotPixel()</u>	Enables and configures the hot pixel correction.
<u>is_SetColorConverter()</u>	Selects Bayer conversion mode.
<u>is_SetColorCorrection()</u>	Sets color correction.
<u>is_SetColorMode()</u>	Selects a color mode.
<u>is_Convert()</u>	Conversion parameters for raw Bayer conversion.
<u>is_EdgeEnhancement()</u>	Sets edge enhancement.
<u>is_SetSaturation()</u>	Sets the image saturation (digital post-processing).
<u>is_SetSensorTestImage()</u>	Enables test image output from sensor.
<u>is_SetTimeout()</u>	Sets user-defined timeout values.

Flowchart: Changing camera timing

Click in the figure to get help on the functions.



Recommended order of function calls for timing changes.

Return the possible range for the pixel clock.

The pixel clock determines the range of frame rate settings.

Return the possible range of the frame rate.

The frame rate defines the range for the exposure time.

Return the possible exposure range and set the exposure time.


4.2.4.2 Automatic Image Control

The uc480 driver provides various options to automatically adjust the image capture parameters to the lighting situation. All controls are configured using the [is_SetAutoParameter\(\)](#) SDK function.

For more information on the automatic image control see [Camera basics: Automatic image control](#).

Flowchart: Enable Auto Brightness

Click in the figure to get help on the functions.



4.2.4.3 Image Pre-processing

Bayer conversion

The following functions enable and adjust the Bayer conversion (see [Color filter \(Bayer filter\)](#)).

is_Convert()	Converts a Bayer raw image into the desired output format
is_GetColorConverter()	Returns the currently set Bayer conversion mode
is_SetBayerConversion()	Sets the algorithm for Bayer conversion
is_SetColorConverter()	Sets the algorithm for Bayer conversion in the camera (not applicable for DCx Cameras)

Lookup table

Using lookup table (LUT) functions, you can e. g. adjust brightness or contrast after the acquisition.

is_GetCameraLUT()	Read out current hardware LUT
is_GetImageHistogram()	Computes a histogram for the image buffer passed to the function

4.2.4.4 Get Camera Status

Using these functions, you can read out additional useful information on the camera status.

is_CameraStatus()	Returns the event counters and other information. Enables standby mode.
is_GetAutoInfo()	Returns status information on the auto features.
is_GetCameraList()	Returns information on all connected cameras.
is_GetCameraType()	Returns the camera type.
is_CaptureStatus()	Displays information on errors that have occurred.
is_GetError()	Displays errors that have occurred.
is_GetUsedbandwidth()	Returns the bus bandwidth (in Mbyte/s) currently used by all initialized or selected cameras.
is_GetVsyncCount()	Returns the VSYNC counter. It will be incremented by 1 each time the sensor starts capturing an image.
is_SetErrorReport()	Enables/disables dialog messages for error output.

4.2.4.5 Using the Camera EEPROM

The non-volatile EEPROM of every DCx camera can hold user data or camera settings.

is_GetCameraInfo()	Returns the factory-set information (e.g. revision information for the individual DCxCamera components).
is_GetSensorInfo()	Returns the sensor information.
is_ReadEEPROM()	Reads out the writable data area of the EEPROM.
is_WriteEEPROM()	Writes user data to the EEPROM.

4.2.5 Saving Images and Videos

Using the uc480 API, you can

- [Save and load single frames](#)
- [Capture an AVI frame sequence](#)

4.2.5.1 Saving and Loading Single Frames

With the [is_ImageFile\(\)](#) function you can save the image data of the current image memory to a BMP, PNG or JPG file, and load saved image data into an image memory.

4.2.5.2 Capturing AVIs

The functions of the `uc480_tools.dll` enable you to save images captured with the DCxCamera as sequences to an AVI file. In order to reduce the file size, the single frames are stored in the AVI container using an adjustable JPEG compression. It is possible to extract single frames from the AVI file.

AVI Capture Workflow

First initialize the AVI interface and then create a empty AVI file.

isavi_InitAVI()	Initializes the AVI interface.
isavi_ExitAVI()	Terminates and closes the AVI interface.
isavi_OpenAVI()	Opens an AVI file for capturing.
isavi_CloseAVI()	Closes an AVI file.
isavi_GetAVIFileName()	Returns the name of the current AVI file.

The following settings should also be done prior to starting the recording.

isavi_SetFrameRate()	Sets the frame rate of the AVI video.
isavi_SetImageQuality()	Sets the compression level/image quality of the AVI video.
isavi_SetImageSize()	Sets the size and offset of the input image memory.

Once the AVI file has been created, captured images are placed in a buffer. Then, the images are compressed and added to the AVI file which is stored on the hard disk. These operations are not performed in the same thread as the capturing process. If you capture more images while a compression or write operation is in progress, the new images will be discarded.

isavi_StartAVI()	Starts AVI recording.
isavi_AddFrame()	Adds a compressed image to the AVI file.
isavi_StopAVI()	Stops AVI recording.

With these functions, you can query additional information on the ongoing recording.

isavi_GetAVISize()	Returns the size of the current AVI file.
isavi_GetnCompressedFrames()	Returns the number of frames in the current AVI file.
isavi_GetnLostFrames()	Returns the number of frames that have been discarded so far.
isavi_ResetFrameCounters()	Resets the counters for discarded and saved frames to 0.

Events can be used to get signalled when a frame was added.

isavi_DisableEvent()	Disables a AVI event.
isavi_EnableEvent()	Enables a AVI-Event.

isavi_ExitEvent()	Turns off AVI event handling.
isavi_InitEvent()	Turns on AVI event handling.

Supported color formats

The supported input color formats are RGB32, RGB24, Y8 and raw Bayer. The output file will always be in RGB24 format, regardless of the input data format. You can adjust the size of the images to be stored by defining a freely selectable area of interest (AOI).

Capture speed

The possible speed of capture depends on the selected color format, the image size and the compression level of the AVI file as well as the PC performance.

Playback in external applications


AVI files you have captured using the `uc480_tools.dll` can also be played back in external applications, such as Windows Media Player. To do this, you need to install the uc480 MJPEG codec on your system:

- Open the uc480 installation directory (default: `C:\Program Files\Thorlabs\DCx Cameras\Tools32` or `C:\Program Files\Thorlabs\DCx Cameras\Tools64`).
- Right-click the `uc480Mjpeg.inf` (`uc480Mjpeg_64.inf`) file.
- Select "Install". The codec is installed automatically.

In player or recording software, the codec will show up as "Intermedia-X MJPEG Codec".

Flowchart: AVI capture

Click in the figure to get help on the functions.



Init the AVI file engine

Set AVI format options:

- Image size and position
- Colour mode

Optional:

- Compression / img. quality

Create new or open existing AVI file

Set AVI frame rate

Start recording

Compress current image and write to AVI file

Stop recording

Close AVI file

Release AVI file engine

Release uc480_tools.dll

4.2.6 Using Inputs and Outputs

Depending on the model, DCx Cameras have one or more digital inputs and outputs designed for different purposes.

- [Input/output control](#): Here, you will find functions for setting the DCxCamera's I/Os and for using the trigger and flash modes.

4.2.6.1 Input/Output Control

With these functions you can use the camera's digital in-/outputs for trigger and flash control.


is_SetExternalTrigger()	Enables the digital input for trigger operation or returns the applied signal level.
is_IO()	Sets the digital output for flash control or a static output level.
is_IO()	Sets the delay and power-on time of the flash output.
is_SetTriggerDelay()	Sets the trigger signal delay time.
is_IO()	Determines the delay and power-on times of the flash output to obtain a global shutter effect when using rolling shutter sensors.
is_ForceTrigger()	Simulates a trigger signal in hardware trigger mode.

With these commands you can activate additional functions or use GPIOs on some DCx Cameras.


is_IO()	Sets the additional digital outputs (GPIO).
is_IO()	Defines each port as a digital input or output (GPIO).
is_IO()	Toggles the color of the status LED for DCU22x and DCC1240x cameras.

Flowchart: Digital input

Click in the figure to get help on the functions.





Flowchart: Digital output



4.3 Function Descriptions

To integrate the DCx Cameras into your own programs, you can use the functions and parameters provided by the uc480 SDK. These are described in this chapter. The descriptions are listed alphabetically by function and are structured as follows:

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

This table shows the availability of the function. For both Windows and Linux the table shows which DCx camera series supports the function.

Syntax

Prototype of the function from the [uc480.h](#) header file.

Description

Description of the function with cross-references to related functions.

Input parameters

Description of the function parameters including their value ranges.

Return value

Description and value range of the return value. If a function returns the `IS_NO_SUCCESS (-1)` value, you can get information on the error from the [is_GetError\(\)](#) function.

Related functions

List with similar or related SDK functions.

Example



For some functions, C++ programming samples are have been added.

Sample programs

Some descriptions include references to uc480 SDK sample programs. When you install the uc480 software, the demo applications are copied to the `C:\Program Files\Thorlabs\DCx Cameras\Samples` directory. The associated source code can be found under `C:\Program Files\Thorlabs\DCx Cameras\Develop\Source`.

All sample programs are described in the **uc480 Samples Manual**.

4.3.1 is_AddToSequence

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_AddToSequence (HCAM hCam, char* pcImgMem, INT nID)
```

Description

`is_AddToSequence()` adds an image memory to the list of image memories used for ring buffering. The image memory must have been previously requested using [is_AllocImageMem\(\)](#). Using the [is_SetAllocatedImageMem\(\)](#) function, you can set a memory that has been allocated before as image memory. Image memories that are used for ring buffering must all have been allocated with the same color depth (bits per pixel).

Input parameters

hCam	Camera handle
pcMem	Pointer to image memory
nID	Image memory ID

Return values

IS_INVALID_CAMERA_HANDLE	Invalid camera handle
IS_NO_SUCCESS	General error message
IS_SUCCESS	Function executed successfully



Related functions

- [is_AllocImageMem\(\)](#)
- [is_InitImageQueue\(\)](#)
- [is_SetImageMem\(\)](#)
- [is_SetAllocatedImageMem\(\)](#)

Sample programs

- uc480Sequence (C++)

4.3.2 is_AllocImageMem

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_AllocImageMem (HCAM hCam, INT width, INT height, INT bitspixel, char** ppcImgMem, INT* pid)
```

Description

`is_AllocImageMem()` allocates an image memory for an image having its dimensions defined by `width` and `height` and its color depth defined by `bitspixel`. The memory size is at least:

`size = [width * int((bitspixel + 7) / 8) + adjust] * height` (for details on `adjust`, see below)

The line increment is calculated as:

```
line = width * int[(bitspixel + 7) / 8]
```

```
lineinc = line + adjust
```

`adjust = 0`, if line can be divided by 4 without remainder

`adjust = 4 - rest(line / 4)`, if line cannot be divided by 4 without remainder

To read out the line increment, you can use the [is_GetImgMemPitch\(\)](#) function.

The starting address of the memory area is returned in `ppcImgMem`.

`pid` returns an ID for the allocated memory. A newly allocated memory is not directly active, i.e. digitised images will not be stored immediately in this new memory. It must first be made active using [is_SetImageMem\(\)](#).

The returned pointer must be write-protected and may not be altered because it will be used for all further `ImageMem` functions. To release the memory, you can use [is_FreeImageMem\(\)](#).

Notes

- In the Direct3D or OpenGL modes, image memory allocation is not necessary.
- RGB16 and RGB15 require the same amount of memory, but can be distinguished by the `bitspixel` parameter. For information on the bit depths of different color formats please refer to the [Appendix: Color and memory formats](#) chapter.
- In case the operating system is short of physical memory, today's OS versions swap individual areas of the RAM that have not been used for some time out to the slower hard disk. This can slow down image capture if more image memory has been allocated than can be provided by the RAM at a time.

Input parameters

<code>hCam</code>	Camera handle
<code>width</code>	Image width
<code>height</code>	Image height
<code>bitspixel</code>	Image bit depth (bits per pixel).
<code>ppcImgMem</code>	Returns the pointer to the memory starting address
<code>pid</code>	Returns the ID of this memory



Return values

IS_CANT_ADD_TO_SEQUENCE	The image memory is already included in the sequence and cannot be added again.
IS_INVALID_CAMERA_HANDLE	Invalid camera handle
IS_INVALID_MEMORY_POINTER	Invalid pointer or invalid memory ID
IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
IS_NO_SUCCESS	General error message
IS_OUT_OF_MEMORY	No memory could be allocated.
IS_SEQUENCE_BUF_ALREADY_LOCKED	The memory could not be locked. The pointer to the buffer is invalid.
IS_SUCCESS	Function executed successfully

Related functions

- [is_FreeImageMem\(\)](#)
- [is_AddToSequence\(\)](#)
- [is_SetImageMem\(\)](#)
- [is_SetAllocatedImageMem\(\)](#)
- [is_GetColorDepth\(\)](#)
- [is_GetImgMemPitch\(\)](#)

4.3.3 is_AOI

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_AOI (HCAM hCam, UINT nCommand, void* pParam, UINT nSizeOfParam)
```

Description

`is_AOI()` can be used to set the size and position of an [area of interest \(AOI\)](#) within an image. The following AOIs can be defined:

- Image AOI – display of an image portion
- Auto Brightness AOI – reference area of interest for automatic brightness control
- Auto Whitebalance AOI – reference area of interest for automatic white balance control

The `nCommand` input parameter is used to select the function mode. The `pParam` input parameter depends on the selected function mode. If you select functions for setting or returning a value, `pParam` contains a pointer to a variable of the `UINT` type. The size of the memory area to which `pParam` refers is specified in the `nSizeOfParam` input parameter.

Note

Previous AOI functions

The `is_AOI()` function comprises all the functions for setting and positioning an AOI. The following uc480 API commands are therefore obsolete:

- `is_SetAOI()`
- `is_SetImageAOI()`
- `is_SetImageSize()`
- `is_SetImagePos()`

See also [Obsolete functions](#)

AOI for automatic image control

The AOI for automatic brightness control (AES/AGC) and automatic white balance (AWB) defaults to the same size as the current image (i.e. the image AOI).

After changes to the image geometry (by resetting an image AOI, by binning or subsampling), the AOIs for automatic image control will always be reset to the image AOI value. This means that it might be necessary to set the AOIs for auto brightness/auto white balance again manually.

Fast changes of AOI position

Using the `IS_AOI_IMAGE_SET_POS_FAST` command, you can change the positions of AOIs very quickly. Executing this command takes just a few milliseconds. When using this command, a few special requirements have to be met:

- The command is currently not supported by all DCx Cameras. With the `IS_AOI_IMAGE_SET_POS_FAST_SUPPORTED` command, you can check whether your sensor supports fast position changes.
- Hot pixel correction has to be disabled (see [is_HotPixel\(\)](#)).
- Image capture is not suspended for fast AOI position changes. As a result, when you call the command, a number of images might still be transferred with the old AOI position if they were in the driver buffer at that moment.

Notes

1. Changing the image size

- When changing the size of the AOI, please make sure that the selected image memory is large enough. If it isn't, allocate a new image memory (see [is_AllocImageMem\(\)](#)).
- Changes to the image size affect the value ranges of the frame rate and exposure time. After executing `is_AOI()`, calling the following functions is recommended in order to keep the defined camera settings:
 - [is_SetFrameRate\(\)](#)
 - [is_Exposure\(\)](#)
 - If you are using the flash function: [is_IO\(\)](#)

2. Step widths for AOI definition (position grid)

The available step widths for the position and size of image AOIs depend on the sensor. The values defining the position and size of an AOI have to be integer multiples of the allowed step widths.

For details on the AOI grids of the individual camera models, please see [Camera and sensor data](#).

3. AOI in combination with high frame rates

With very small AOI and therefore high frame rate and maximum possible frame rate set, it is possible that the USB camera transfers in freerun mode only half frame rates. This is a signal for a camera-internal overload. In this case it is recommended to set the frame rate to maximum of 98 %

Multi AOI function of the DCC1240x and DCC3240x models

The sensor of a.m. cameras supports multiple AOIs in one image capture. The AOIs are transferred together as one image. In this mode you can create 2 or 4 AOIs, which have either the same X axis or the same Y axis (see also uc480 Viewer: Multi AOI). The sensor is faster in this mode. It is possible to switch the AOI in the horizontal direction.

Sequence AOI mode for DCC1240x and DCC3240x camera models

A.m. camera models have a special AOI mode. In this mode you can define besides the normal AOI (AOI 1) up to 3 further AOI on the sensor (see uc480 Viewer: Sequence AOI). When activating the sequence mode, note that only the following combinations are possible:

1. All additional AOIs are off. AOI 1 is always active.
2. AOI 2 (+ AOI 1)
3. AOI 2 and 3 (+ AOI 1)
4. AOI 2, 3 and 4 (+ AOI 1)

It is not possible to have a combination e.g. of AOI 2 and AOI 4.

The parameters of AOI 2, 3 and 4 are defined by the [AOI_SEQUENCE_PARAMS](#) structure. In the version 4.82 binning, subsampling and scaler are not supported.

Input parameters

hCam	Camera handle
<input type="checkbox"/> nCommand IS_AOI_IMAGE_SET_AOI	Sets the position and size of the image by using an object of the <code>IS_RECT</code> type. Sample 1 for AOI You can define the start position of the AOI in the memory by ORing <code>IS_AOI_IMAGE_POS_ABSOLUTE</code> with the X or Y position. Sample for setting the AOI position

hCam	Camera handle
IS_AOI_IMAGE_GET_AOI	Returns the AOI in an IS_RECT object. Sample 2 for AOI
IS_AOI_IMAGE_SET_POS	Sets the AOI position by using an object of the IS_POINT_2D type. IS_AOI_IMAGE_POS_ABSOLUTE can be ORed here, as well. Sample for setting the AOI position
IS_AOI_IMAGE_GET_POS	Returns the position in an IS_POINT_2D object.
IS_AOI_IMAGE_SET_SIZE	Sets the AOI size by using an object of the IS_SIZE_2D type.
IS_AOI_IMAGE_GET_SIZE	Returns the size in an IS_SIZE_2D object.
IS_AOI_IMAGE_GET_POS_MIN	Returns the minimum possible position in an IS_POINT_2D object.
IS_AOI_IMAGE_GET_SIZE_MIN	Returns the smallest possible size in an IS_SIZE_2D object.
IS_AOI_IMAGE_GET_POS_MAX	Returns the maximum possible position in an IS_POINT_2D object.
IS_AOI_IMAGE_GET_SIZE_MAX	Returns the largest possible size in an IS_SIZE_2D object.
IS_AOI_IMAGE_GET_POS_INC	Returns the increment for the position in an IS_POINT_2D object.
IS_AOI_IMAGE_GET_SIZE_INC	Returns the increment for the size in an IS_SIZE_2D object.
IS_AOI_IMAGE_GET_POS_X_ABS	Returns an UINT object indicating whether IS_AOI_IMAGE_POS_ABSOLUTE is set for the X position. Sample 3 for AOI
IS_AOI_IMAGE_GET_POS_Y_ABS	Returns an UINT object indicating whether IS_AOI_IMAGE_POS_ABSOLUTE is set for the Y position.
IS_AOI_IMAGE_GET_ORIGINAL_AOI	Returns the AOI in an IS_RECT object without binning, subsampling or scaling.
IS_AOI_IMAGE_SET_POS_FAST	Allows changing the AOI position very quickly by using an IS_POINT_2D object. Hot pixel correction has to be disabled (see information above).
IS_AOI_IMAGE_SET_POS_FAST_SUPPORT	Returns an UINT object indicating whether fast AOI position changes are supported. The passed variable returns 0 if the function is not supported by the sensor.
IS_AOI_AUTO_BRIGHTNESS_SET_AOI	Sets the AOI for automatic brightness control (similar to IS_AOI_IMAGE_SET_AOI).
IS_AOI_AUTO_BRIGHTNESS_GET_AOI	Returns the AOI for automatic brightness control (similar to IS_AOI_IMAGE_GET_AOI).
IS_AOI_AUTO_WHITEBALANCE_SET_AOI	Sets the AOI for automatic white balance (similar to IS_AOI_IMAGE_SET_AOI).
IS_AOI_AUTO_WHITEBALANCE_GET_AOI	Returns the AOI for automatic white balance (similar to IS_AOI_IMAGE_GET_AOI).
IS_AOI_MULTI_GET_SUPPORTED_MODES	Returns the supported multi AOI modes in an UINT object.
IS_AOI_MULTI_SET_AOI	Sets the multi AOI mode. The mode you want to use has to be ORed with IS_AOI_MULTI_SET_AOI. The axes are passed in an UINT array: <ul style="list-style-type: none"> • Array[0] - Array[3] = X1...X4 • Array[4] - Array[8] = Y1...Y4 Sample 1 for multi AOI
IS_AOI_MULTI_GET_AOI	Returns the set multi AOI mode. The mode that is used has to be ORed with IS_AOI_MULTI_SET_AOI.

hCam	Camera handle
	Sample 2 fo multi AOI
IS_AOI_MULTI_MODE_X_Y_AXES	Multi AOI mode of the camera models DCC1240x/DCC3240x with up to AOIs (up to 4 x and y axes). The axes are passed by a UINT array: <ul style="list-style-type: none"> • array[0] - array[3] = X1...X4 • array[4] - array[8] = Y1...Y4 Attention: This parameter has been renamed in version 4.20. In formerly versions this parameter was named IS_AOI_MULTI_MODE_AXES.
IS_AOI_MULTI_DISABLE_AOI	Disables Multi AOI. The mode that is used has to be ORed with IS_AOI_MULTI_SET_AOI. Sample 3 for multi AOI
IS_AOI_SEQUENCE_GET_SUPPORTED	Returns a bitmask with the supported AOIs (only DCC1240x/DCC3240x camera models) Sample 1 for AOI sequence mode
IS_AOI_SEQUENCE_SET_PARAMS	Sets the parameters of AOI 2, 3 or 4 (only DCC1240x/DCC3240x camera models) Sample 2 for sequence AOI mode
IS_AOI_SEQUENCE_GET_PARAMS	Returns the parameters of AOI 2, 3 or 4 (only DCC1240x/DCC3240x camera models) Sample 2 for sequence AOI mode
IS_AOI_SEQUENCE_SET_ENABLE	Set a bitmask defining which AOIs should be active (only DCC1240x/DCC3240x camera models). Note: IS_AOI_SEQUENCE_SET_PARAMS must be called after IS_AOI_SEQUENCE_SET_ENABLE, with enabling the sequence AOI mode all AOIs are set to the same value and therefore the parameters are lost. Sample 3 for sequence AOI mode
IS_AOI_SEQUENCE_GET_ENABLE	Returns the bitmask (only DCC1240x/DCC3240x camera models) Sample 3 for sequence AOI mode
pParam	Pointer to a function parameter, whose function depends on nCommand.
nSizeOfParam	Size (in bytes) of the memory area to which pParam refers.

Contents of the IS_RECT structure

INT	s32X	X position of the AOI
INT	s32Y	Y position of the AOI
INT	s32Width	Width of the AOI
INT	s32Height	Height of the AOI

Contents of the IS_POINT_2D structure

INT	s32X	X position of the AOI
INT	s32Y	Y position of the AOI

Contents of the IS_SIZE_2D structure

INT	s32Width	Width of the AOI
INT	s32Height	Height of the AOI

Content of the AOI_SEQUENCE_PARAMS structure

INT	s32AOIIndex	Index of the AOI
INT	s32NumberOfCycleRepetitions	Number of readout cycles
INT	s32X	X position of the AOI
INT	x32Y	Y position of the AOI
Double	dblExposure	Exposure
INT	s32Gain	Gain
INT	s32BinningMode	Binning mode (not supported in version 4.82)
INT	s32SubsamplingMode	Subsampling mode (not supported in version 4.82)
INT	s32DetachImageParameters	<ul style="list-style-type: none"> • 0 = every change of the exposure time and the master gain is copied from AOI 1 to the additional AOIs (default). As a change of AOI 1 also reset the exposure time, this change is also transferred to AOI 2, 3 and 4. • 1 = a change of exposure time, gain or position of AOI 1 does not affect the parameters of AOI 2, 3 and 4.
Double	dblScalerFactor	Scaling factor (not supported in version 4.82)
BYTE	byReserved[64]	Reserved

Return values

IS_BAD_STRUCTURE_SIZE	An internal structure has an incorrect size.
IS_CANT_ADD_TO_SEQUENCE	The image memory is already included in the sequence and cannot be added again.
IS_CANT_COMMUNICATE_WITH_DRIVER	Communication with the driver failed because no driver has been loaded.
IS_CANT_OPEN_DEVICE	An attempt to initialize or select the camera failed (no camera connected or initialization error).
IS_CAPTURE_RUNNING	A capturing operation is in progress and must be terminated before you can start another one.
IS_DR_CANNOT_CREATE_SURFACE	The image surface or overlay surface could not be created.
IS_DR_CANNOT_CREATE_TEXTURE	The texture could not be created.
IS_DR_CANNOT_CREATE_VERTEX_BUFFER	The vertex buffer could not be created.
IS_DR_DEVICE_OUT_OF_MEMORY	Not enough graphics memory available.
IS_DR_LIBRARY_NOT_FOUND	The DirectRenderer library could not be found.
IS_INVALID_BUFFER_SIZE	The image memory has an inappropriate size to store the image in the desired format.
IS_INVALID_CAMERA_TYPE	The camera type defined in the .ini file does not match the current camera model.
IS_INVALID_CAPTURE_MODE	The function can not be executed in the current camera operating mode (free run, trigger or standby).
IS_INVALID_CAMERA_HANDLE	Invalid camera handle
IS_INVALID_MEMORY_POINTER	Invalid pointer or invalid memory ID
IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
IS_IO_REQUEST_FAILED	An IO request from the uc480 driver failed. Possibly the versions of the uc480.dll (API) and the driver file (uc480_usb.sys) do not match.
IS_NO_ACTIVE_IMG_MEM	No active image memory available. You must set the memory to active using the <code>is_SetImageMem()</code> function or create a sequence using the <code>is_AddToSequence()</code> function.
IS_NO_SUCCESS	General error message
IS_NOT_CALIBRATED	The camera does not contain any calibration data.
IS_NOT_SUPPORTED	The camera model used here does not support this function or setting.
IS_OUT_OF_MEMORY	No memory could be allocated.
IS_SEQUENCE_BUF_ALREADY_LOCKED	The memory could not be locked. The pointer to the buffer is invalid.
IS_SUCCESS	Function executed successfully
IS_TIMED_OUT	A timeout occurred. An image capturing process could not be terminated within the allowable period.

IS_TRIGGER_ACTIVATED	The function cannot be used because the camera is waiting for a trigger signal.
----------------------	---

Related functions

- [is_ImageFormat\(\)](#)
- [is_SetBinning\(\)](#)
- [is_SetSubSampling\(\)](#)

Sample 1 for AOI

```
// Sets the position and size of the image by using an object of the IS_RECT type.
IS_RECT rectAOI;

rectAOI.s32X      = 100;
rectAOI.s32Y      = 100;
rectAOI.s32Width  = 200;
rectAOI.s32Height = 100;

INT nRet = is_AOI( hCam, IS_AOI_IMAGE_SET_AOI, (void*)&rectAOI, sizeof(rectAOI));
```

Sample 2 for AOI

```
// Returns the AOI position and size by using an object of the IS_RECT type.
IS_RECT rectAOI;

INT nRet = is_AOI(hCam, IS_AOI_IMAGE_GET_AOI, (void*)&rectAOI, sizeof(rectAOI));
if (nRet == IS_SUCCESS)
{
    INT x      = rectAOI.s32X;
    INT y      = rectAOI.s32Y;
    INT width  = rectAOI.s32Width;
    INT height = rectAOI.s32Height;
}

```

Sample 3 for AOI


```
// Returns an UINT object indicating whether IS_AOI_IMAGE_POS_ABSOLUTE is set for the X position.
UINT nAbsPos = 0;

INT nRet = is_AOI(hCam, IS_AOI_IMAGE_GET_POS_X_ABS, (void*)&nAbsPos, sizeof(nAbsPos));
if (nRet == IS_SUCCESS)
{
    if (nAbsPos == IS_AOI_IMAGE_POS_ABSOLUTE)
    {
        // set
    }
    else if (nAbsPos == 0)
    {
        // not set
    }
}

```

Examples for setting absolute AOI positions in memory

```
x = 100
y = 100
```



AOI without absolute memory positioning

```
x = 100 | IS_AOI_IMAGE_POS_ABSOLUTE
y = 100
```


Sample 2 for sequence AOI mode

```

AOI_SEQUENCE_PARAMS Param;

// Set parameters of AOI 2
Param.s32AOIIndex = IS_AOI_SEQUENCE_INDEX_AOI_2;
Param.s32NumberOfCycleRepetitions = 1;
Param.s32X = 100;
Param.s32Y = 200;
...

INT nRet = is_AOI(m_hCam, IS_AOI_SEQUENCE_SET_PARAMS, (void*)&Param, sizeof(Param));

// Get parameters of AOI 2
Param.s32AOIIndex = IS_AOI_SEQUENCE_INDEX_AOI_2;

nRet = is_AOI(m_hCam, IS_AOI_SEQUENCE_GET_PARAMS, (void*)&Param, sizeof(Param));

```

Sample 3 for sequence AOI mode

```

INT nMask = 0;

// Enable AOI 1, Disable AOI 2, 3 and 4
nMask = IS_AOI_SEQUENCE_INDEX_AOI_1;

INT nRet = is_AOI(m_hCam, IS_AOI_SEQUENCE_SET_ENABLE, (void*)&nMask, sizeof(nMask));

// Enable AOI 1 and 2
nMask = IS_AOI_SEQUENCE_INDEX_AOI_1 |
        IS_AOI_SEQUENCE_INDEX_AOI_2;

nRet = is_AOI(m_hCam, IS_AOI_SEQUENCE_SET_ENABLE, (void*)&nMask, sizeof(nMask));

// Enable AOI 1, 2 and 3
nMask = IS_AOI_SEQUENCE_INDEX_AOI_1 |
        IS_AOI_SEQUENCE_INDEX_AOI_2 |
        IS_AOI_SEQUENCE_INDEX_AOI_3;

nRet = is_AOI(m_hCam, IS_AOI_SEQUENCE_SET_ENABLE, (void*)&nMask, sizeof(nMask));



// Enable AOI 1, 2, 3 and 4
nMask = IS_AOI_SEQUENCE_INDEX_AOI_1 |
        IS_AOI_SEQUENCE_INDEX_AOI_2 |
        IS_AOI_SEQUENCE_INDEX_AOI_3 |
        IS_AOI_SEQUENCE_INDEX_AOI_4;

nRet = is_AOI(m_hCam, IS_AOI_SEQUENCE_SET_ENABLE, (void*)&nMask, sizeof(nMask));

// Get current AOI mask
INT nRet = is_AOI(m_hCam, IS_AOI_SEQUENCE_GET_ENABLE, (void*)&nMask, sizeof(nMask));

```

4.3.4 is_AutoParameter

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_AutoParameter(HCAM hCam, UINT nCommand, void* pParam, UINT cbSizeOfParam)
```

Description

This functions enables/disables the auto white balance. With this function, you can require all supported types for white balance. In addition to the older white balance with the Gray-World algorithm, there is also a color temperature control according to Kelvin. In addition to the function the supported color spaces are queried and set.

The `nCommand` input parameter is used to select the function mode. The `pParam` input parameter depends on the selected function mode. If you select functions for setting or returning a value, `pParam` contains a pointer to a variable of the `UINT` type. The size of the memory area to which `pParam` refers is specified in the `cbSizeOfParam` input parameter.

Note

In a later version the `is_AutoParameter()` function will replace the [is_SetAutoParameter\(\)](#) function. In version 4.82, it is only partly replaced.

Input parameters

hCam	Camera handle
<input type="checkbox"/> nCommand	
IS_AWB_CMD_GET_SUPPORTED_TYPES	Returns the supported types for auto white balance (Example 1) <ul style="list-style-type: none"> IS_AWB_GREYWORLD: 0x0001 IS_AWB_COLOR_TEMPERATURE: 0x0002
IS_AWB_CMD_GET_TYPE	Returns the current set type of the auto white balance (Example 2) <ul style="list-style-type: none"> IS_AWB_GREYWORLD: 0x0001 IS_AWB_COLOR_TEMPERATURE: 0x0002
IS_AWB_CMD_SET_TYPE	Sets the type of the auto white balance (Example 2) <ul style="list-style-type: none"> IS_AWB_GREYWORLD: 0x0001 IS_AWB_COLOR_TEMPERATURE: 0x0002
IS_AWB_CMD_GET_ENABLE	Returns if the auto white balance is enabled (Example 3) <ul style="list-style-type: none"> IS_AUTOPARAMETER_DISABLE: 0 IS_AUTOPARAMETER_ENABLE: 1 IS_AUTOPARAMETER_ENABLE_RUNONCE: 2
IS_AWB_CMD_SET_ENABLE	Enables/Disables the auto white balance (Example 3) <ul style="list-style-type: none"> IS_AUTOPARAMETER_DISABLE: 0 IS_AUTOPARAMETER_ENABLE: 1 IS_AUTOPARAMETER_ENABLE_RUNONCE: 2
IS_AWB_CMD_GET_SUPPORTED_RGB_COLOR_M	Returns the supported color spaces for the auto white

hCam	Camera handle
ODELS	balance (Example 4) <ul style="list-style-type: none"> • RGB_COLOR_MODEL_SRGB_D50: 0x0001 • RGB_COLOR_MODEL_SRGB_D65: 0x0002 • RGB_COLOR_MODEL_CIE_RGB_E: 0x0004 • RGB_COLOR_MODEL_ECI_RGB_D50: 0x0008 • RGB_COLOR_MODEL_ADOBE_RGB_D65: 0x0010
IS_AWB_CMD_GET_RGB_COLOR_MODEL	Returns the current color space for the auto white balance (Example 5) <ul style="list-style-type: none"> • RGB_COLOR_MODEL_SRGB_D50: 0x0001 • RGB_COLOR_MODEL_SRGB_D65: 0x0002 • RGB_COLOR_MODEL_CIE_RGB_E: 0x0004 • RGB_COLOR_MODEL_ECI_RGB_D50: 0x0008 • RGB_COLOR_MODEL_ADOBE_RGB_D65: 0x0010
IS_AWB_CMD_SET_RGB_COLOR_MODEL	Sets the color space for the auto white balance (Example 5) <ul style="list-style-type: none"> • RGB_COLOR_MODEL_SRGB_D50: 0x0001 • RGB_COLOR_MODEL_SRGB_D65: 0x0002 • RGB_COLOR_MODEL_CIE_RGB_E: 0x0004 • RGB_COLOR_MODEL_ECI_RGB_D50: 0x0008 • RGB_COLOR_MODEL_ADOBE_RGB_D65: 0x0010
pParam	Pointer to a function parameter, whose function depends on nCommand.
cbSizeOfParam	Size (in bytes) of the memory area to which pParam refers.

Return values

IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
IS_NO_SUCCESS	General error message
IS_NOT_SUPPORTED	The camera model used here does not support this function or setting.
IS_SUCCESS	Function executed successfully

Related functions

- [is_SetAutoParameter\(\)](#)

Example 1

```

UINT nSupportedTypes = 0;
INT nRet = is_AutoParameter(m_hCam,
    IS_AWB_CMD_GET_SUPPORTED_TYPES,
    (void*)&nSupportedTypes,
    sizeof(nSupportedTypes)
);

if (nRet == IS_SUCCESS)
{
    if ((nSupportedTypes & IS_AWB_COLOR_TEMPERATURE) != 0)
    {
        // AWB type "Color Temperature" is supported
    }
}

```

```

    if ((nSupportedTypes & IS_AWB_GREYWORLD) != 0)
    {
        // AWB type "Greyworld" is supported
    }
}

```

Example 2

```

UINT nType = 0;

// Read current type
INT nRet = is_AutoParameter(m_hCam, IS_AWB_CMD_GET_TYPE, (void*)&nType, sizeof(nType));

// Write new type
nType = IS_AWB_GREYWORLD;
nRet = is_AutoParameter(m_hCam, IS_AWB_CMD_SET_TYPE, (void*)&nType, sizeof(nType));

```

Example 3

```

// Is AWB enabled?
UINT nEnable;
INT nRet = is_AutoParameter(m_hCam, IS_AWB_CMD_GET_ENABLE, (void*)&nEnable, sizeof(nEnable));

// Enable AWB (once)
nEnable = IS_AUTOPARAMETER_ENABLE_RUNONCE;
nRet = is_AutoParameter(m_hCam, IS_AWB_CMD_SET_ENABLE, (void*)&nEnable, sizeof(nEnable));

```

Example 4

```

UINT nSupportedRGBColorModels = 0;
nRet = is_AutoParameter(m_hCam,
    IS_AWB_CMD_GET_SUPPORTED_RGB_COLOR_MODELS,
    (void*)&nSupportedRGBColorModels,
    sizeof(nSupportedRGBColorModels)
);

if (nRet == IS_SUCCESS)
{
    if ((nSupportedRGBColorModels & RGB_COLOR_MODEL_SRGB_D50) != 0)
    {
        // Color model SRGB D50 is supported. See uc480.h for color model defines
    }
}

```

Example 5

```



UINT nRGBColorModel = 0;
INT nRet = is_AutoParameter(m_hCam,
    IS_AWB_CMD_GET_RGB_COLOR_MODEL,
    (void*)&nRGBColorModel,
    sizeof(nRGBColorModel)
);

nRGBColorModel = RGB_COLOR_MODEL_CIE_RGB_E;

nRet = is_AutoParameter(m_hCam,
    IS_AWB_CMD_SET_RGB_COLOR_MODEL,
    (void*)&nRGBColorModel,
    sizeof(nRGBColorModel)
);

```

4.3.5 is_Blacklevel


	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_Blacklevel(HCAM hCam, UINT nCommand, void* pParam, UINT cbSizeOfParam)
```

Beschreibung

`is_Blacklevel()` controls the black level correction of the camera which might improve the image quality under certain circumstances. By default, the sensor adjusts the black level value for each pixel automatically. If the environment is very bright, it can be necessary to adjust the black level manually.



The `nCommand` input parameter is used to select the function mode. The `pParam` input parameter depends on the selected function mode. If you select functions for setting or returning a value, `pParam` contains a pointer to a variable of the `UINT` type. The size of the memory area to which `pParam` refers is specified in the `cbSizeOfParam` input parameter.

Input parameters

hCam	Camera handle
<input type="checkbox"/> nCommand	
IS_BLACKLEVEL_CMD_GET_CAPS	Returns the black level feature of the camera (Example 1): <ul style="list-style-type: none"> IS_BLACKLEVEL_CAP_SET_AUTO_BLACKLEVEL: The state of the automatic black level can be changed. The flag does not indicate whether the camera is running with auto black level by default or not. For this purpose, use <code>IS_BLACKLEVEL_CMD_GET_MODE_DEFAULT</code>. IS_BLACKLEVEL_CAP_SET_OFFSET: The offset can be changed. The flag does not indicate whether the camera has set an offset by default or not. For this purpose, use <code>IS_BLACKLEVEL_CMD_GET_OFFSET_DEFAULT</code>.
IS_BLACKLEVEL_CMD_GET_MODE_DEFAULT	Returns the default black level mode (Example2)
IS_BLACKLEVEL_CMD_GET_MODE	Returns the current black level mode (Example2)
IS_BLACKLEVEL_CMD_SET_MODE	Sets the black level mode (Example 2) <ul style="list-style-type: none"> IS_AUTO_BLACKLEVEL_OFF: The automatic black level mode is switched off.

hCam	Camera handle
	<ul style="list-style-type: none"> IS_AUTO_BLACKLEVEL_ON: The automatic black level mode is switched on.
IS_BLACKLEVEL_CMD_GET_OFFSET_DEFAULT	Returns the default offset (Example 3)
IS_BLACKLEVEL_CMD_GET_OFFSET_RANGE	Returns the range of the offset (Example 3)
IS_BLACKLEVEL_CMD_GET_OFFSET	Returns the current offset (Example 3)
IS_BLACKLEVEL_CMD_SET_OFFSET	Sets the offset (Example 3)
pParam	Pointer to a function parameter, whose function depends on nCommand.
cbSizeOfParam	Size (in bytes) of the memory area to which pParam refers.

Return values

IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
IS_NO_SUCCESS	General error message
IS_NOT_SUPPORTED	The camera model used here does not support this function or setting.
IS_SUCCESS	Function executed successfully

Example 1

```

INT nBlacklevelCaps;

nRet = is_Blacklevel(hCam, IS_BLACKLEVEL_CMD_GET_CAPS,
                    (void*)&nBlacklevelCaps, sizeof(nBlacklevelCaps));
if (nRet == IS_SUCCESS) {

    // The user can changed the state of the auto blacklevel
    BOOL bSetAutoBlacklevel = (nBlacklevelCaps & IS_BLACKLEVEL_CAP_SET_AUTO_BLACKLEVEL) != 0;

    // The user can change the offset
    BOOL bSetBlacklevelOffset = (nBlacklevelCaps & IS_BLACKLEVEL_CAP_SET_OFFSET) != 0;
}

```

Example 2

```

INT nMode = IS_AUTO_BLACKLEVEL_OFF;

// Get default blacklevel mode
INT nRet = is_Blacklevel(hCam, IS_BLACKLEVEL_CMD_GET_MODE_DEFAULT, (void*)&nMode, sizeof(nMode));

// Get current blacklevel mode
nRet = is_Blacklevel(hCam, IS_BLACKLEVEL_CMD_GET_MODE, (void*)&nMode, sizeof(nMode));

// Set new mode (enable auto blacklevel)
nMode = IS_AUTO_BLACKLEVEL_ON;
nRet = is_Blacklevel(hCam, IS_BLACKLEVEL_CMD_SET_MODE, (void*)&nMode, sizeof(nMode));

```

Example 3

```

INT nOffset = 0;

// Get default blacklevel offset
INT nRet = is_Blacklevel(hCam, IS_BLACKLEVEL_CMD_GET_OFFSET_DEFAULT,
                        (void*)&nOffset, sizeof(nOffset));

// Get offset range
IS_RANGE_S32 nRange;
nRet = is_Blacklevel(hCam, IS_BLACKLEVEL_CMD_GET_OFFSET_RANGE, (void*)&nRange, sizeof(nRange));
INT nOffsetMin = nRange.s32Min;
INT nOffsetMax = nRange.s32Max;
INT nOffsetInc = nRange.s32Inc;



```



```
// Get current blacklevel offset
nRet = is_Blacklevel(hCam, IS_BLACKLEVEL_CMD_GET_OFFSET, (void*)&nOffset, sizeof(nOffset));

// Set new offset
nOffset = 100;
nRet = is_Blacklevel(hCam, IS_BLACKLEVEL_CMD_SET_OFFSET, (void*)&nOffset, sizeof(nOffset));
```

4.3.6 is_CameraStatus

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

ULONG is_CameraStatus (HCAM hCam, INT nInfo, ULONG ulValue)

Description

Using `is_CameraStatus()`, you can query and partly set various status information and settings.

Input parameters

hCam	Camera handle
<div style="border: 1px solid gray; padding: 2px;"> <input type="checkbox"/> nInfo </div>	
IS_FIFO_OVR_CNT	Number of FIFO overruns. Is increased if image data gets lost because the USB bus is congested.
IS_SEQUENCE_CNT	Returns the sequence count. For is_CaptureVideo() , this parameter is set to 0. Each time the sequence buffer (image counter) changes, the counter is increased by 1.
IS_SEQUENCE_SIZE	Returns the number of sequence buffers.
IS_EXT_TRIGGER_EVENT_CNT	Returns the camera's internal count of external trigger events.
IS_TRIGGER_MISSED	Returns the number of unprocessed trigger signals. Is reset to 0 after each call.
IS_LAST_CAPTURE_ERROR	Returns the last image capture error. For a list of all possible error events, see is_CaptureStatus() .
IS_PARAMETER_SET_1	Indicates whether parameter set 1 including camera settings is present on the camera (read-only). See also is_ParameterSet() . Return values: TRUE Parameter set 1 present FALSE Parameter set 1 not present
IS_PARAMETER_SET_2	Indicates whether parameter set 2 including camera settings is present on the camera (read-only). See also is_ParameterSet() . Return values: TRUE Parameter set 2 present FALSE Parameter set 2 not present
IS_STANDBY	Sets the camera to standby mode. Return values: TRUE Camera changes to standby mode FALSE The camera changes to freerun mode
IS_STANDBY_SUPPORTED	Queries whether the camera supports standby mode (read-only). Return values: TRUE The camera supports standby mode FALSE The camera does not support standby mode
<div style="border: 1px solid gray; padding: 2px;"> <input type="checkbox"/> ulValue </div>	
IS_GET_STATUS	Returns the information specified by <code>nInfo</code> .



Return values

Only if <code>ulValue = IS_GET_STATUS</code>	Returns the information specified by <code>nInfo</code>
When used with <code>IS_LAST_CAPTURE_ERROR</code>	Returns the last image capture error. For a list of all possible error events, see <code>is_GetCaptureErrorInfo()</code> .
<code>IS_BAD_STRUCTURE_SIZE</code>	An internal structure has an incorrect size.
<code>IS_CANT_COMMUNICATE_WITH_DRIVER</code>	Communication with the driver failed because no driver has been loaded.
<code>IS_CANT_OPEN_DEVICE</code>	An attempt to initialize or select the camera failed (no camera connected or initialization error).
<code>IS_INVALID_CAMERA_TYPE</code>	The camera type defined in the .ini file does not match the current camera model.
<code>IS_INVALID_EXPOSURE_TIME</code>	This setting is not available for the currently set exposure time.
<code>IS_INVALID_CAMERA_HANDLE</code>	Invalid camera handle
<code>IS_INVALID_PARAMETER</code>	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
<code>IS_IO_REQUEST_FAILED</code>	An IO request from the uc480 driver failed. Possibly the versions of the <code>uc480.dll</code> (API) and the driver file (<code>uc480_usb.sys</code>) do not match.
<code>IS_NO_SUCCESS</code>	General error message
<code>IS_NOT_CALIBRATED</code>	The camera does not contain any calibration data.
<code>IS_NOT_SUPPORTED</code>	The camera model used here does not support this function or setting.
<code>IS_OUT_OF_MEMORY</code>	No memory could be allocated.
<code>IS_SUCCESS</code>	Function executed successfully
<code>IS_TIMED_OUT</code>	A timeout occurred. An image capturing process could not be terminated within the allowable period.

Related functions

- [is_GetCameraInfo\(\)](#)
- [is_GetError\(\)](#)
- [is_SetErrorReport\(\)](#)
- [is_SetTriggerCounter\(\)](#)

4.3.7 is_CaptureStatus

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_CaptureStatus (HCAM hCam, UINT nCommand, void* pParam, UINT cbSizeOfParam)
```

Description

The function returns information on errors that occurred during an image capture. All errors are listed that occurred since the last reset of the function.

The `nCommand` input parameter is used to select the function mode. The `pParam` input parameter depends on the selected function mode. If you select functions for setting or returning a value, `pParam` contains a pointer to a variable of the `UINT` type. The size of the memory area to which `pParam` refers is specified in the `cbSizeOfParam` input parameter.

Note

The following functions are obsolete by the `is_CaptureStatus()` function:

- `is_GetCaptureErrorInfo()`
- `is_ResetCaptureErrorInfo()`

See also [Obsolete functions](#).

Input parameters

<code>hCam</code>	Camera handle
<input type="checkbox"/> <code>nCommand</code>	
<code>IS_CAPTURE_STATUS_INFO_CMD_GET</code>	Returns the CaptureStatus information (Example 1)
<code>IS_CAPTURE_STATUS_INFO_CMD_RESET</code>	Resets the CaptureStatus information (Example 2)
<code>pParam</code>	Pointer to a function parameter, whose function depends on <code>nCommand</code> .
<code>cbSizeOfParam</code>	Size (in bytes) of the memory area to which <code>pParam</code> refers.

Content of the `uc480_CAPTURE_STATUS` enumeration

<code>IS_CAP_STATUS_API_NO_DEST_MEM</code>	<p>There is no destination memory for copying the finished image.</p> <p><input type="checkbox"/> Possible cause/remedy</p> <p>Not enough destination memory allocated or all destination buffers locked by the application.</p> <ul style="list-style-type: none"> • Release locked destination memory • Allocate more destination memory • Reduce the frame rate so that there is more time to process the filled destination memory
<code>IS_CAP_STATUS_API_CONVERSION_FAILED</code>	The current image could not be processed correctly.

	<ul style="list-style-type: none"> ▣ Possible cause <p>Internal error during internal processing of the image</p>
IS_CAP_STATUS_API_IMAGE_LOCKED	<p>The destination buffers are locked and could not be written to.</p> <ul style="list-style-type: none"> ▣ Possible cause/remedy <p>All destination buffers locked by the application</p> <ul style="list-style-type: none"> • Release locked destination memory • Allocate more destination memory • Reduce the frame rate so that there is more time to process the filled destination memory
IS_CAP_STATUS_DRV_OUT_OF_BUFFERS	<p>No free internal image memory is available to the driver. The image was discarded.</p> <ul style="list-style-type: none"> ▣ Possible cause/remedy <p>The computer takes too long to process the images in the uc480 API (e.g. color conversion)</p> <ul style="list-style-type: none"> • Reduce the frame rate so that there is more time to process the filled image memory of the driver • Disable resource-intensive API image pre-processing functions (e.g. edge enhancement, color correction, choose smaller filter mask for software color conversion)
IS_CAP_STATUS_DRV_DEVICE_NOT_READY	<p>The camera is no longer available. It is not possible to access images that have already been transferred.</p> <ul style="list-style-type: none"> ▣ Possible cause <p>The camera has been disconnected or closed.</p>
IS_CAP_STATUS_USB_TRANSFER_FAILED	<p>The image was not transferred over the USB bus.</p> <ul style="list-style-type: none"> ▣ Possible cause/remedy <p>Not enough free bandwidth on the USB bus for transferring the image</p> <ul style="list-style-type: none"> • Reduce the pixel clock frequency • Operate fewer cameras simultaneously on a USB bus • Check the quality of the USB cabling and components
IS_CAP_STATUS_DEV_TIMEOUT	<p>The maximum allowable time for image capturing in the camera was exceeded.</p> <ul style="list-style-type: none"> ▣ Possible cause/remedy <p>The selected timeout value is too low for image capture</p> <ul style="list-style-type: none"> • Reduce the exposure time • Increase the timeout

Contents of the uc480_CAPTURE_STATUS_INFO structure

DWORD	dwCapStatusCnt_Total	Returns the total number of errors occurred since the last reset.
BYTE	reserved[60]	Reserved for an internal function
DWORD	adwCapStatusCnt_Detail[256]	This array returns the current count for each possible error. The possible errors are listed above.

Return values

IS_INVALID_CAMERA_HANDLE	Invalid camera handle
IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
IS_NO_SUCCESS	General error message
IS_SUCCESS	Function executed successfully

Related functions

- [is_GetError\(\)](#)
- [is_CameraStatus\(\)](#)
- [is_SetErrorReport\(\)](#)

Example 1



```
uc480_CAPTURE_STATUS_INFO CaptureStatusInfo;
INT nRet = is_CaptureStatus(m_hCam,
                           IS_CAPTURE_STATUS_INFO_CMD_GET,
                           (void*)&CaptureStatusInfo,
                           sizeof(CaptureStatusInfo));

if (nRet == IS_SUCCESS)
{
    UINT nConversionFailed = CaptureStatusInfo.adwCapStatusCnt_Detail[IS_CAP_STATUS_API_CONVERSION_FAILED];
    UINT nTotalInfos = CaptureStatusInfo.dwCapStatusCnt_Total;
}
```

Example 2

```
INT nRet = is_CaptureStatus(m_hCam,
                           IS_CAPTURE_STATUS_INFO_CMD_RESET,
                           NULL,
                           0);
```

4.3.8 is_CaptureVideo

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_CaptureVideo (HCAM hCam, INT Wait)
```

Description

`is_CaptureVideo()` activates the camera's live video mode (free run mode). The driver transfers the images to an allocated image memory or, if Direct3D/OpenGL is used, to the graphics card. The image data (DIB mode) is stored in the memory created using [is_AllocImageMem\(\)](#) and designated as active image memory using [is_SetImageMem\(\)](#). Using [is_GetImageMem\(\)](#), you can query the memory address.

If ring buffering is used, the image capturing function cycles through all image memories used for storing the images of a capture sequence in an endless loop. Sequence memories locked by [is_LockSeqBuf\(\)](#) will be skipped. If the last available sequence memory has been filled, the sequence event or message will be triggered. Capturing always starts with the first element of the sequence.

For further information on the image capture modes, see the [How to proceed: Image capture](#) section.

Input parameters

hCam	Camera handle
<input type="checkbox"/> Wait	
IS_DONT_WAIT	Timeout value for image capture (see also the How to proceed: Timeout values for image capture section)
IS_WAIT	
Time t	
IS_GET_LIVE	Returns if live capture is enabled.

Return values

When used with IS_GET_LIVE	TRUE if live capture is enabled
IS_BAD_STRUCTURE_SIZE	An internal structure has an incorrect size.
IS_CANT_COMMUNICATE_WITH_DRIVER	Communication with the driver failed because no driver has been loaded.
IS_CANT_OPEN_DEVICE	An attempt to initialize or select the camera failed (no camera connected or initialization error).
IS_CAPTURE_RUNNING	A capturing operation is in progress and must be terminated before you can start another one.
IS_INVALID_BUFFER_SIZE	The image memory has an inappropriate size to store the image in the desired format.
IS_INVALID_CAMERA_TYPE	The camera type defined in the .ini file does not match the current camera model.
IS_INVALID_EXPOSURE_TIME	This setting is not available for the currently set exposure time.
IS_INVALID_CAMERA_HANDLE	Invalid camera handle

IS_INVALID_MEMORY_POINTER	Invalid pointer or invalid memory ID
IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
IS_IO_REQUEST_FAILED	An IO request from the uc480 driver failed. Possibly the versions of the uc480.dll (API) and the driver file (uc480_usb.sys) do not match.
IS_NO_ACTIVE_IMG_MEM	No active image memory available. You must set the memory to active using the is_SetImageMem() function or create a sequence using the is_AddToSequence() function.
IS_NO_SUCCESS	General error message
IS_NOT_CALIBRATED	The camera does not contain any calibration data.
IS_NOT_SUPPORTED	The camera model used here does not support this function or setting.
IS_OUT_OF_MEMORY	No memory could be allocated.
IS_SUCCESS	Function executed successfully
IS_TIMED_OUT	A timeout occurred. An image capturing process could not be terminated within the allowable period.
IS_TRIGGER_ACTIVATED	The function cannot be used because the camera is waiting for a trigger signal.



Related functions

- [is_FreezeVideo\(\)](#)
- [is_StopLiveVideo\(\)](#)
- [is_SetExternalTrigger\(\)](#)
- [is_ForceTrigger\(\)](#)
- [is_SetTimeout\(\)](#)
- [is_CaptureStatus\(\)](#)

Sample programs

- SimpleLive (C++)

4.3.9 is_ClearSequence

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_ClearSequence (HCAM hCam)
```

Description

`is_ClearSequence()` removes all image memories from the sequence list that were added using [is_AddToSequence\(\)](#). After a call of `is_ClearSequence()`, there is no more active image memory. To make an image memory the active memory, call [is_SetImageMem\(\)](#).

Input parameters

hCam	Camera handle
------	---------------



Return values

IS_INVALID_CAMERA_HANDLE	Invalid camera handle
IS_NO_SUCCESS	General error message
IS_SEQUENCE_LIST_EMPTY	The sequence list is empty and cannot be deleted.
IS_SUCCESS	Function executed successfully

Related functions

- [is_AddToSequence\(\)](#)
- [is_FreeImageMem\(\)](#)
- [is_SetImageMem\(\)](#)

4.3.10 is_ColorTemperature

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_ColorTemperature (HCAM hCam, UINT nCommand,
                        void *pParam, UINT nSizeOfParam)
```

Description

Using `is_ColorTemperature()` you can fix a setting (in kelvins) for the color temperature of an image when you are using a color camera. The function will use the sensor's hardware gain controls for the setting, as far as possible. In addition, you can choose between different color spaces. A specific color temperature will result in slightly differing RGB values, depending on the selected color space.

The `nCommand` input parameter is used to select the function mode. The `pParam` input parameter depends on the selected function mode. If you select functions for setting or returning a value, `pParam` contains a pointer to a variable of the `UINT` type. The size of the memory area to which `pParam` refers is specified in the `nSizeOfParam` input parameter.

The color temperature is the temperature to which a black body radiator has to be heated to glow and give off light in the corresponding color. Warm light (reddish) has a lower color temperature than cold light (bluish). The following table lists a few example values:

Light source	Color temperature
Light bulb (100 W)	2800
Halogen lamp	3200
Fluorescent lamp (cold white)	4000
Morning and evening sunlight	5000
Noon sunlight	5500-5800
Flash strobe	6000
Overcast daylight	6500-7500
Fog	8000

Note

The `is_ColorTemperature()` function cannot be used simultaneously with the automatic white balance function in [is_SetAutoParameter\(\)](#)/[is_AutoParameter\(\)](#).

Input parameters

hCam	Camera handle
<input type="checkbox"/> nCommand	
Setting the color space	
COLOR_TEMPERATURE_CMD_GET_SUPPORTED_RGB_COLORS_MODELS	Returns the supported color spaces . <input type="checkbox"/> More details pParam: Pointer to a bit mask of type <code>UINT</code> The bit mask returns the supported modes,

hCam	Camera handle
	linked by logical ORs (see Color spaces table). nSizeOfParam: 4
COLOR_TEMPERATURE_CMD_SET_RGB_COLOR_MODEL	Sets a color space . <ul style="list-style-type: none"> ▣ More details pParam: Pointer to variable of type <code>UINT</code> that passes the value to be set. nSizeOfParam: 4
COLOR_TEMPERATURE_CMD_GET_RGB_COLOR_MODEL	Returns the set color space . <ul style="list-style-type: none"> ▣ More details pParam: Pointer to variable of type <code>UINT</code> returning the current value. nSizeOfParam: 4
COLOR_TEMPERATURE_CMD_GET_RGB_COLOR_MODEL_DEFAULT	Returns the default color space. <ul style="list-style-type: none"> ▣ More details pParam: Pointer to variable of type <code>UINT</code> returning the default value. nSizeOfParam: 4
Setting the color temperature	
COLOR_TEMPERATURE_CMD_SET_TEMPERATURE	Sets a color temperature. More details pParam: Pointer to variable of type <code>UINT</code> that passes the value to be set. nSizeOfParam: 4
COLOR_TEMPERATURE_CMD_GET_TEMPERATURE	Returns the set color temperature. More details pParam: Pointer to variable of type <code>UINT</code> returning the current value. nSizeOfParam: 4
COLOR_TEMPERATURE_CMD_GET_TEMPERATURE_MIN	Returns the minimum value for the color temperature. More details pParam: Pointer to variable of type <code>UINT</code> returning the minimum value. nSizeOfParam: 4
COLOR_TEMPERATURE_CMD_GET_TEMPERATURE_MAX	Returns the maximum value for the color temperature. More details pParam: Pointer to variable of type <code>UINT</code> returning the maximum value. nSizeOfParam: 4
COLOR_TEMPERATURE_CMD_GET_TEMPERATURE_INC	Returns the increment for setting the color temperature. More details pParam: Pointer to variable of type <code>UINT</code> returning the increment.

hCam	Camera handle
	nSizeOfParam: 4
COLOR_TEMPERATURE_CMD_GET_TEMPERATURE_DEFAULT	Returns the default value for the color temperature. More details pParam: Pointer to variable of type <code>UINT</code> returning the default value. nSizeOfParam: 4
pParam	Pointer to a function parameter, whose function depends on <code>nCommand</code> .
nSizeOfParam	Size (in bytes) of the memory area to which <code>pParam</code> refers.

Color Spaces

RGB_COLOR_MODEL_SRGB_D50	sRGB (standard RGB) color space with a white point of 5000 kelvins (warm light)
RGB_COLOR_MODEL_SRGB_D65	sRGB (standard RGB) color space with a white point of 6500 kelvins (mid daylight)
RGB_COLOR_MODEL_CIE_RGB_E	CIE-RGB color space with standard illumination E
RGB_COLOR_MODEL_ECI_RGB_D50	ECI-RGB color space with a white point of 5000 kelvins (warm light)
RGB_COLOR_MODEL_ADOBE_RGB_D65	Adobe RGB color space with a white point of 6500 kelvins (mid daylight). The Adobe RGB color space is larger than the sRGB color space, but not all devices can render it.



Return values

IS_INVALID_CAMERA_HANDLE	Invalid camera handle
IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
IS_NO_SUCCESS	General error message
IS_NOT_SUPPORTED	The camera model used here does not support this function or setting.
IS_SUCCESS	Function executed successfully

Related functions

- [is_SetHardwareGain\(\)](#)
- [is_AutoParameter\(\)](#)
- [is_SetAutoParameter\(\)](#)

4.3.11 is_Configuration

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_Configuration (UINT nCommand, void* pParam, UINT cbSizeOfParam)
```

Description

Use `is_Configuration()` to set various system-wide options:

- Windows only: Processor operating states (idle states/C-states)**
 Modern processors have various operating states, so-called C-states, that are characterized by different power requirements. When the operating system selects an operating state with low power consumption (unequal C0), the USB transmission efficiency may be affected. Use the function parameters `IS_CONFIG_CPU_IDLE_STATES_CMD...` to disable these low power consumption operating states and improve USB transmission efficiency. The uc480 driver changes the current energy settings of the operating system when the first USB DCx camera is opened. After the last USB DCx camera is closed, the uc480 driver restores the original settings. The settings are valid for the current user only.
- Windows only: Activate OpenMP (Open Multi-Processing)**
 OpenMP is a programming interface that supports distributed computing on multi-core processors. When you activate OpenMP support, intensive computing operations, such as the [Bayer conversion](#), are distributed across several processor cores to accelerate execution. The use of OpenMP, however, increases CPU load.
- Load camera parameters during installation**
 Use the function parameters `IS_CONFIG_INITIAL_PARAMETERSET...` to indicate whether to apply the parameters stored on the camera automatically when opening the camera. You must first store the camera parameters on the camera using the [is_ParameterSet\(\)](#) function or via the corresponding function in the uc480 demo. This setting applies to all connected cameras. If no parameters are stored on the camera, the standard parameters of this camera model are applied.

The `nCommand` input parameter is used to select the function mode. The `pParam` input parameter depends on the selected function mode. If you select functions for setting or returning a value, `pParam` contains a pointer to a variable of the `UINT` type. The size of the memory area to which `pParam` refers is specified in the `cbSizeOfParam` input parameter.

Note

Settings for processor operating states: The settings for processor operating states are available only on Windows operating systems.

Input parameters

hCam	Camera handle		
<div style="border: 1px solid gray; padding: 5px;"> <div style="border-bottom: 1px solid gray; padding-bottom: 5px;"> [-] nCommand </div> <div style="padding: 5px;"> <table border="0" style="width: 100%;"> <tr> <td style="width: 40%; vertical-align: top;"> <div style="border: 1px solid gray; padding: 5px;"> IS_CONFIG_CMD_GET_CAPABILITIES </div> </td> <td style="vertical-align: top;"> <p>Returns the configuration options supported by the system.</p> <div style="border: 1px solid gray; padding: 5px;"> [-] Additional information <ul style="list-style-type: none"> • <code>pParam</code>: Pointer to a <code>UINT</code> bitmask. The status flags from CONFIGURATION_CAPS are returned in the bitmask. </div> </td> </tr> </table> </div> </div>		<div style="border: 1px solid gray; padding: 5px;"> IS_CONFIG_CMD_GET_CAPABILITIES </div>	<p>Returns the configuration options supported by the system.</p> <div style="border: 1px solid gray; padding: 5px;"> [-] Additional information <ul style="list-style-type: none"> • <code>pParam</code>: Pointer to a <code>UINT</code> bitmask. The status flags from CONFIGURATION_CAPS are returned in the bitmask. </div>
<div style="border: 1px solid gray; padding: 5px;"> IS_CONFIG_CMD_GET_CAPABILITIES </div>	<p>Returns the configuration options supported by the system.</p> <div style="border: 1px solid gray; padding: 5px;"> [-] Additional information <ul style="list-style-type: none"> • <code>pParam</code>: Pointer to a <code>UINT</code> bitmask. The status flags from CONFIGURATION_CAPS are returned in the bitmask. </div>		

hCam	Camera handle
	<ul style="list-style-type: none"> • nSizeOfParam: 4 Example 1
IS_CONFIG_CPU_IDLE_STATES_CMD_GET_ENABLE	<p>Returns whether the current settings allow low power consumption operating states (unequal C0).</p> <ul style="list-style-type: none"> ▣ Additional information • pParam: Pointer to a <code>UINT</code> bitmask. The status flags from CONFIGURATION_SEL are returned in the bitmask. • nSizeOfParam: 4 Example 2
IS_CONFIG_CPU_IDLE_STATES_CMD_SET_DISABLE_ON_OPEN	<p>Changes the energy settings of the operating system so that low power consumption operating states (unequal C0) are disabled.</p> <ul style="list-style-type: none"> ▣ Additional information • pParam: Pointer to a <code>UINT</code> variable, see CONFIGURATION_SEL. • nSizeOfParam: 4 <p>Note: To apply a new setting, you must close all open DCx Cameras and then reopen at least one camera.</p> Example 3
IS_CONFIG_CPU_IDLE_STATES_CMD_GET_DISABLE_ON_OPEN	<p>Returns the current setting for <code>IS_CONFIG_CPU_IDLE_STATES_CMD_SET_DISABLE_ON_OPEN</code>.</p> <ul style="list-style-type: none"> ▣ Additional information • pParam: Pointer to a <code>UINT</code> bitmask. The status flags from CONFIGURATION_SEL are returned in the bitmask. • nSizeOfParam: 4 Example 3
IS_CONFIG_OPEN_MP_CMD_GET_ENABLE	<p>Returns whether OpenMP support is enabled.</p> <ul style="list-style-type: none"> ▣ Additional information • pParam: Pointer to a <code>UINT</code> bitmask. The status flags from CONFIGURATION_SEL are returned in the bitmask. • nSizeOfParam: 4 Example 4
IS_CONFIG_OPEN_MP_CMD_SET_ENABLE	<p>Enables OpenMP support.</p> <ul style="list-style-type: none"> ▣ Additional information • pParam: Pointer to a <code>UINT</code> variable, <ul style="list-style-type: none"> to enable: <code>IS_CONFIG_OPEN_MP_ENABLE</code> to disable: <code>IS_CONFIG_OPEN_MP_DISABLE</code> • nSizeOfParam: 4 <p>Note: The settings are lost after the application is closed and must be set again the next time the</p>

hCam	Camera handle
	camera is started. Example 4
IS_CONFIG_OPEN_MP_CMD_GET_ENABLE_DEFAULT	Returns the default setting for OpenMP support. <ul style="list-style-type: none"> ▣ Additional information • pParam: Pointer to a <code>UINT</code> bitmask. The status flags from CONFIGURATION_SEL are returned in the bitmask. • nSizeOfParam: 4 Example 4
IS_CONFIG_INITIAL_PARAMETERSET_CMD_SET	Sets the parameter set to read and apply from the camera EEPROM when the camera is opened. <ul style="list-style-type: none"> ▣ Additional information • pParam: Pointer to a <code>UINT</code> variable, <code>IS_CONFIG_INITIAL_PARAMETERSET_NONE</code>, <code>IS_CONFIG_INITIAL_PARAMETERSET_1</code>, <code>IS_CONFIG_INITIAL_PARAMETERSET_2</code> • nSizeOfParam: 4 Example 5
IS_CONFIG_INITIAL_PARAMETERSET_CMD_GET	Returns which parameter set will be read and applied from the camera EEPROM when the camera is opened. <ul style="list-style-type: none"> ▣ Additional information • pParam: Pointer to a <code>UINT</code> variable, <code>IS_CONFIG_INITIAL_PARAMETERSET_NONE</code>, <code>IS_CONFIG_INITIAL_PARAMETERSET_1</code>, <code>IS_CONFIG_INITIAL_PARAMETERSET_2</code> • nSizeOfParam: 4 Example 5
pParam	Pointer to a function parameter, whose function depends on <code>nCommand</code> .
cbSizeOfParam	Size (in bytes) of the memory area to which <code>pParam</code> refers.

Contents of the CONFIGURATION_CAPS Structure

INT	IS_CONFIG_CPU_IDLE_STATES_CAP_SUPPORTED	Function parameters for setting the processor operating states are supported.
INT	IS_CONFIG_OPEN_MP_CAP_SUPPORTED	Function parameters to configure OpenMP are supported.
INT	IS_CONFIG_INITIAL_PARAMETERSET_CAP_SUPPORTED	Function parameters to load camera parameters during initialization are supported.

Contents of the CONFIGURATION_SEL Structure

INT	IS_CONFIG_CPU_IDLE_STATES_BIT_AC_VALUE	Set/recover processor operating states for power supply unit operation
INT	IS_CONFIG_CPU_IDLE_STATES_BIT_DC_VALUE	Set/recover processor operating states for

	E	battery operation
INT	IS_CONFIG_OPEN_MP_DISABLE	OpenMP support disabled
INT	IS_CONFIG_OPEN_MP_ENABLE	OpenMP support enabled
INT	IS_CONFIG_INITIAL_PARAMETERSET_NONE	Load camera parameters during initialization disabled
INT	IS_CONFIG_INITIAL_PARAMETERSET_1	Load camera parameter set 1 during initialization
INT	IS_CONFIG_INITIAL_PARAMETERSET_2	Load camera parameter set 2 during initialization

Return values

IS_CANT_OPEN_REGISTRY	Error opening a Windows registry key
IS_CANT_READ_REGISTRY	Error reading settings from the Windows registry
IS_ERROR_CPU_IDLE_STATES_CONFIGURATION	The configuration of the CPU idle has failed.
IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
IS_NO_IMAGE_MEM_ALLOCATED	The driver could not allocate memory.
IS_NO_SUCCESS	General error message
IS_NOT_SUPPORTED	The camera model used here does not support this function or setting.
IS_OPERATING_SYSTEM_NOT_SUPPORTED	Operating system not supported
IS_SUCCESS	Function executed successfully

Related functions

- [is_ParameterSet\(\)](#)

Example 1

```
UINT nCaps = 0;
INT nRet = is_Configuration(IS_CONFIG_CMD_GET_CAPABILITIES, (void*)&nCaps, sizeof(UINT));
if (nRet == IS_SUCCESS)
{
    if (nCaps & IS_CONFIG_CPU_IDLE_STATES_CAP_SUPPORTED)
    {
        // CPU idle states supported
    }

    if (nCaps & IS_CONFIG_OPEN_MP_CAP_SUPPORTED)
    {
        // OpenMP supported
    }

    if (nCaps & IS_CONFIG_INITIAL_PARAMETERSSET_CAP_SUPPORTED)
    {
        // Initial parameter set supported
    }
}
```

Example 2

```
INT nCurrentCpuStates = 0;
INT nRet = is_Configuration(IS_CONFIG_CPU_IDLE_STATES_CMD_GET_ENABLE,
                          (void*)&nCurrentCpuStates,
                          sizeof(nCurrentCpuStates)
                          );

if (nRet == IS_SUCCESS)
{
    if ((nCurrentCpuStates & IS_CONFIG_CPU_IDLE_STATES_BIT_AC_VALUE) == 0)
    {
        // The CPU idle states for mains power is already deactivated
    }

    if ((nCurrentCpuStates & IS_CONFIG_CPU_IDLE_STATES_BIT_DC_VALUE) == 0)
    {
        // The CPU idle states for battery power is already deactivated
    }
}
```

Example 3

```
UINT nCpuStates = IS_CONFIG_CPU_IDLE_STATES_BIT_AC_VALUE | IS_CONFIG_CPU_IDLE_STATES_BIT_DC_VALUE;

INT nRet = is_Configuration(IS_CONFIG_CPU_IDLE_STATES_CMD_SET_DISABLE_ON_OPEN,
                          (void*)&nCpuStates,
                          sizeof(nCpuStates)
                          );

if (nRet == IS_SUCCESS)
{
    nCpuStates = 0;
    nRet = is_Configuration(IS_CONFIG_CPU_IDLE_STATES_CMD_GET_DISABLE_ON_OPEN,
                          (void*)&nCpuStates,
                          sizeof(nCpuStates)
                          );

    if (nRet == IS_SUCCESS)
    {
        if (nCpuStates & IS_CONFIG_CPU_IDLE_STATES_BIT_AC_VALUE)
        {
            // CPU idle states for mains power are deactivated when camera is opened
        }

        if (nCpuStates & IS_CONFIG_CPU_IDLE_STATES_BIT_DC_VALUE)
        {
            // CPU idle states for battery power are deactivated when camera is opened
        }
    }
}
```

```
}  
}
```

Example 4

```
UINT nEnabled = 0;
INT nRet = is_Configuration(IS_CONFIG_OPEN_MP_CMD_GET_ENABLE, (void*)&nEnabled, sizeof(nEnabled));
if (nRet == IS_SUCCESS)
{
    if (nEnabled == IS_CONFIG_OPEN_MP_ENABLE)
    {
        // OpenMP enabled
    }
}



nEnabled = 0;
nRet = is_Configuration(IS_CONFIG_OPEN_MP_CMD_GET_ENABLE_DEFAULT, (void*)&nEnabled, sizeof(nEnabled));
if (nRet == IS_SUCCESS)
{
    nRet = is_Configuration(IS_CONFIG_OPEN_MP_CMD_SET_ENABLE, (void*)&nEnabled, sizeof(nEnabled));
    if (nRet == IS_SUCCESS)
    {
        // Default value set
    }
}
```

Example 5

```
UINT nNumber = 0;
INT nRet = is_Configuration(IS_CONFIG_INITIAL_PARAMETERSET_CMD_GET, (void*)&nNumber, sizeof(nNumber));
if (nRet == IS_SUCCESS)
{
    if (nNumber == IS_CONFIG_INITIAL_PARAMETERSET_NONE)
    {
        // No parameter set specified
    }
    else if (nNumber == IS_CONFIG_INITIAL_PARAMETERSET_1)
    {
        // Parameter set 1
    }
    else if (nNumber == IS_CONFIG_INITIAL_PARAMETERSET_2)
    {
        // Parameter set 2
    }
}

nNumber = IS_CONFIG_INITIAL_PARAMETERSET_2;
is_Configuration(IS_CONFIG_INITIAL_PARAMETERSET_CMD_SET, (void*)&nNumber, sizeof(nNumber));
if (nRet == IS_SUCCESS)
{
    // Set to parameter set 2
}
```

4.3.12 is_Convert

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
IDSEXP is_Convert(HCAM hCam, UINT nCommand, void* pParam, UINT cbSizeOfParam)
```

Description

`is_Convert()` is a general function for conversions. In version 4.82 a raw Bayer image is converted to the desired format. You can set all parameters, which are important for software conversion:

- Pixel format
- Pixel converter (3×3, 5×5)
- Color correction
- Gamma
- Saturation
- Edge enhancement

The target buffer must be allocated with the [is_AllocImageMem\(\)](#) function and must have the right size.

The `nCommand` input parameter is used to select the function mode. The `pParam` input parameter depends on the selected function mode. If you select functions for setting or returning a value, `pParam` contains a pointer to a variable of the `UINT` type. The size of the memory area to which `pParam` refers is specified in the `cbSizeOfParam` input parameter.

Note

The following functions are obsolete by the `is_Convert()` function:

- `is_ConvertImage()`
- `is_SetConvertParam()`

See also [Obsolete functions](#)

Input parameters

<code>hCam</code>	Camera handle
<input type="checkbox"/> <code>nCommand</code>	
<code>IS_CONVERT_CMD_APPLY_PARAMS_AND_CONVERT_BUFFER</code>	Converts a raw Bayer buffer with the passed conversion parameters
<code>pParam</code>	Pointer to a function parameter, whose function depends on <code>nCommand</code> .
<code>cbSizeOfParam</code>	Size (in bytes) of the memory area to which <code>pParam</code> refers.

Contents of the BUFFER_CONVERSION_PARAMS structure

char*	pSourceBuffer	Pointer to the raw Bayer buffer which was created with the is_AllocImageMem() function.
char*	pDestBuffer	Pointer to the target buffer with the converted data which was created with the is_AllocImageMem() function.
INT	nDestPixelFormat	Color mode of the target image; see is_SetColorMode() for the possible modes
INT	nDestPixelConverter	Conversion mode of the target image; see is_SetColorConverter() for the possible modes
INT	nDestGamma	Sets the gamma correction, see is_SetGamma()
INT	nDestEdgeEnhancement	Sets the edge enhancement, see is_EdgeEnhancement()
INT	nDestColorCorrectionMode	Sets the color correction, see is_SetColorCorrection()
INT	nDestSaturationU	Sets the color saturation (saturation U), see is_SetSaturation()
INT	nDestSaturationV	Sets the color saturation (saturation V), see is_SetSaturation()

Return value

IS_INVALID_BUFFER_SIZE	The image memory has an inappropriate size to store the image in the desired format.
IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
IS_NO_SUCCESS	General error message
IS_NOT_SUPPORTED	The camera model used here does not support this function or setting.
IS_SUCCESS	Function executed successfully

Related functions

- [is_SetColorMode\(\)](#)
- [is_SetBayerConversion\(\)](#)

Example

```

BUFFER_CONVERSION_PARAMS conversionParams;



conversionParams.nDestPixelFormat           = IS_CM_BGRA8_PACKED;
conversionParams.nDestPixelConverter       = IS_CONV_MODE_SOFTWARE_3X3;
conversionParams.nDestColorCorrectionMode  = IS_CCOR_DISABLE;
conversionParams.nDestGamma                = 100;
conversionParams.nDestSaturationU          = 100;
conversionParams.nDestSaturationV         = 100;
conversionParams.nDestEdgeEnhancement     = 0;

conversionParams.pSourceBuffer             = pSourceBuffer;
conversionParams.pDestBuffer               = pDestBuffer;

INT nRet = is_Convert(m_hCam,
                    IS_CONVERT_CMD_APPLY_PARAMS_AND_CONVERT_BUFFER,
                    (void*)&conversionParams,
                    sizeof(conversionParams)
                    );

```

4.3.13 is_CopyImageMem

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_CopyImageMem (HCAM hCam, char* pcSource, INT nID, char* pcDest)
```

Description

`is_CopyImageMem()` copies the contents of the image memory described by `pcSource` and `nID` to the memory area to whose starting address `pcDest` points.

Attention

The allocated memory must be large enough to accommodate the entire image in its current format (bits per pixel).

Input parameters

<code>hCam</code>	Camera handle
<code>pcSource</code>	Pointer to the image memory
<code>nID</code>	ID of this image memory
<code>pcDest</code>	Pointer to the destination memory to copy the image to



Return values

<code>IS_CANT_COMMUNICATE_WITH_DRIVER</code>	Communication with the driver failed because no driver has been loaded.
<code>IS_CANT_OPEN_DEVICE</code>	An attempt to initialize or select the camera failed (no camera connected or initialization error).
<code>IS_INVALID_CAMERA_HANDLE</code>	Invalid camera handle
<code>IS_INVALID_MEMORY_POINTER</code>	Invalid pointer or invalid memory ID
<code>IS_INVALID_PARAMETER</code>	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
<code>IS_IO_REQUEST_FAILED</code>	An IO request from the <code>uc480</code> driver failed. Possibly the versions of the <code>uc480.dll</code> (API) and the driver file (<code>uc480_usb.sys</code>) do not match.
<code>IS_NO_SUCCESS</code>	General error message
<code>IS_SUCCESS</code>	Function executed successfully

Related functions

- [is_AllocImageMem\(\)](#)
- [is_SetAllocatedImageMem\(\)](#)

4.3.14 is_CopyImageMemLines

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_CopyImageMemLines (HCAM hCam, char* pcSource,
                          INT nID, INT nLines, char* pcDest)
```

Description

`is_CopyImageMemLines()` copies the contents of the image memory described by `pcSource` and `nID` to the memory area to whose starting address `pcDest` points. The function only copies the number of lines indicated by `nLines`.

Attention

The allocated memory must be large enough to accommodate the in `nLines` given number of image lines considering the image width and format (Bits per Pixel).

Input parameters

<code>hCam</code>	Camera handle
<code>pcSource</code>	Pointer to the image memory
<code>nID</code>	ID of this image memory
<code>nLines</code>	Number of lines to be copied
<code>pcDest</code>	Pointer to the destination memory to copy the image to



Return values

<code>IS_CANT_COMMUNICATE_WITH_DRIVER</code>	Communication with the driver failed because no driver has been loaded.
<code>IS_CANT_OPEN_DEVICE</code>	An attempt to initialize or select the camera failed (no camera connected or initialization error).
<code>IS_INVALID_CAMERA_HANDLE</code>	Invalid camera handle
<code>IS_INVALID_MEMORY_POINTER</code>	Invalid pointer or invalid memory ID
<code>IS_INVALID_PARAMETER</code>	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
<code>IS_IO_REQUEST_FAILED</code>	An IO request from the uc480 driver failed. Possibly the versions of the <code>uc480.dll</code> (API) and the driver file (<code>uc480_usb.sys</code>) do not match.
<code>IS_NO_SUCCESS</code>	General error message
<code>IS_SUCCESS</code>	Function executed successfully

Related functions

- [is_AllocImageMem\(\)](#)
- [is_SetAllocatedImageMem\(\)](#)

4.3.15 is_DeviceFeature

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_DeviceFeature (HCAM hCam, UINT nCommand, void* pParam, UINT cbSizeOfParam)
```

Description

Using `is_DeviceFeature()` you can configure special camera functions provided by specific uc480 models:

- Configuring the AOI merge mode (DCC3240 and DCC3260)
- Using the Log mode (DCC1240/DCC3240 and DCC3260)
- Using level controlled trigger (DCC3260)
- Switching the shutter mode (DCC1240/DCC3240)
- Using the internal image memory (DCC3260)
- Using the line scan mode (DCC1240/DCC3240)
- Configuring the timestamp (DCC3260)

The `nCommand` input parameter is used to select the function mode. The `pParam` input parameter depends on the selected function mode. If you select functions for setting or returning a value, `pParam` contains a pointer to a variable of the `UINT` type. The size of the memory area to which `pParam` refers is specified in the `cbSizeOfParam` input parameter.

Input parameters

<code>hCam</code>	Camera handle
<input type="checkbox"/> <code>nCommand</code>	
<code>IS_DEVICE_FEATURE_CMD_GET_SUPPORTED_FEATURES</code>	Returns the functions supported by the camera. See Status flags table (Example)
<code>pParam</code>	Pointer to a function parameter, whose function depends on <code>nCommand</code> .
<code>cbSizeOfParam</code>	Size (in bytes) of the memory area to which <code>pParam</code> refers.

Status flags from DEVICE_FEATURE_MODE_CAPS

<code>IS_DEVICE_FEATURE_CAP_BLACK_REFERENCE</code>	Displays the black level reference in the image which can be used to calibrate the black level.
<code>IS_DEVICE_FEATURE_CAP_EXTENDED_PIXELCLOCK_RANGE</code>	The use of the extended pixel clock range is supported.
<code>IS_DEVICE_FEATURE_CAP_FPN_CORRECTION</code>	Fixed pattern noise correction is supported.
<code>IS_DEVICE_FEATURE_CAP_IMAGE_EFFECT</code>	Image effects are supported.
<code>IS_DEVICE_FEATURE_CAP_JPEG_COMPRESSION</code>	JPEG compression is supported.
<code>IS_DEVICE_FEATURE_CAP_LEVEL_CONTROLLED_TRIGGER</code>	The use of the level controlled trigger is supported.

IS_DEVICE_FEATURE_CAP_LINESCAN_MODE_FAST	Fast line scan mode is supported/Set mode
IS_DEVICE_FEATURE_CAP_LINESCAN_NUMBER	Line number at fast line scan mode is supported/Set number
IS_DEVICE_FEATURE_CAP_LOG_MODE	Log mode is supported/Set mode
IS_DEVICE_FEATURE_CAP_MEMORY_MODE	Internal image memory is supported.
IS_DEVICE_FEATURE_CAP_MULTI_INTEGRATION	Multi integration mode is supported.
IS_DEVICE_FEATURE_CAP_NOISE_REDUCTION	Noise suppression is supported.
IS_DEVICE_FEATURE_CAP_REPEATED_START_CONDITION_I2C	Enabling/disabling the I ² C stop bit for the read command is supported.
IS_DEVICE_FEATURE_CAP_SEND_EXTERNAL_INTERFACE_DATA	Querying the timestamp via the I ² C bus is supported.
IS_DEVICE_FEATURE_CAP_SENSOR_BIT_DEPTH	Sets the bit depth of the sensor. With the IS_SENSOR_BIT_DEPTH_AUTO setting the software selects the appropriate sensor bit depth to the chosen image format. The bit depth can also be selected independent from the image format. Attention: As you can choose combinations that does not fit, this function should be used by experts only. In most cases the auto control are sufficient.
IS_DEVICE_FEATURE_CAP_SENSOR_SOURCE_GAIN	Analog master gain of the sensor is supported.
IS_DEVICE_FEATURE_CAP_SHUTTER_MODE_GLOBAL	Global shutter mode is supported/Set mode
IS_DEVICE_FEATURE_CAP_SHUTTER_MODE_GLOBAL_ALTERNATIVE_TIMING	Global shutter mode with different timing parameters is supported/Set mode
IS_DEVICE_FEATURE_CAP_SHUTTER_MODE_ROLLING	Rolling shutter mode is supported/Set mode
IS_DEVICE_FEATURE_CAP_SHUTTER_MODE_ROLLING_GLOBAL_START	Rolling shutter mode with global start is supported/Set mode
IS_DEVICE_FEATURE_CAP_TEMPERATURE	Display of the internal camera temperature is supported.
IS_DEVICE_FEATURE_CAP_TEMPERATURE_STATUS	Monitoring the temperature status is supported.
IS_DEVICE_FEATURE_CAP_TIMESTAMP_CONFIGURATION	Configuration of the timestamp (e.g. reset the timestamp to 0 at a signal on the trigger pin).
IS_DEVICE_FEATURE_CAP_VERTICAL_AOI_MERGE	Special AOI merge mode which combines the lines of an AOI to a new image. Example: The AOI merge mode height is 1 line for a monochrome sensor. That means that the image is composed of 4000 "merged" AOIs. If the AOI merge mode height is 100 lines, the image is composed of 400 "merged" AOIs.
IS_DEVICE_FEATURE_CAP_WIDE_DYNAMIC_RANGE	The use of the wide dynamic range is supported.

Return values

IS_CANT_COMMUNICATE_WITH_DRIVER	Communication with the driver failed because no driver has been loaded.
IS_CANT_OPEN_DEVICE	An attempt to initialize or select the camera

	failed (no camera connected or initialization error).
IS_INVALID_CAPTURE_MODE	The function can not be executed in the current camera operating mode (free run, trigger or standby).
IS_INVALID_CAMERA_HANDLE	Invalid camera handle
IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
IS_IO_REQUEST_FAILED	An IO request from the uc480 driver failed. Possibly the versions of the ueye_api.dll (API) and the driver file (ueye_usb.sys or ueye_eth.sys) do not match.
IS_NO_SUCCESS	General error message
IS_NOT_SUPPORTED	The camera model used here does not support this function or setting.
IS_OUT_OF_MEMORY	No memory could be allocated.
IS_SUCCESS	Function executed successfully



Example

```

INT nSupportedFeatures;
INT nRet = is_DeviceFeature(hCam, IS_DEVICE_FEATURE_CMD_GET_SUPPORTED_FEATURES,
                           (void*)&nSupportedFeatures, sizeof(nSupportedFeatures));
if (nRet == IS_SUCCESS)
{
    if (nSupportedFeatures & IS_DEVICE_FEATURE_CAP_LINESCAN_MODE_FAST)
    {
        // Enable line scan mode
        INT nMode = IS_DEVICE_FEATURE_CAP_LINESCAN_MODE_FAST;
        nRet = is_DeviceFeature(hCam, IS_DEVICE_FEATURE_CMD_SET_LINESCAN_MODE, (void*)&nMode,
                               sizeof(nMode));
        // Disable line scan mode
        nMode = 0;
        nRet = is_DeviceFeature(hCam, IS_DEVICE_FEATURE_CMD_SET_LINESCAN_MODE, (void*)&nMode,
                               sizeof(nMode));
        // Return line scan mode
        nRet = is_DeviceFeature(hCam, IS_DEVICE_FEATURE_CMD_GET_LINESCAN_MODE, (void*)&nMode,
                               sizeof(nMode));
    }
    if (nSupportedFeatures & IS_DEVICE_FEATURE_CAP_LINESCAN_NUMBER)
    {
        // Set line number
        INT nLineNumber = 512;
        nRet = is_DeviceFeature(hCam, IS_DEVICE_FEATURE_CMD_SET_LINESCAN_NUMBER,
                               (void*)&nLineNumber, sizeof(nLineNumber));
        nRet = is_DeviceFeature(hCam, IS_DEVICE_FEATURE_CMD_GET_LINESCAN_NUMBER,
                               (void*)&nLineNumber, sizeof(nLineNumber));
    }
}

```

4.3.15.1 Configuring the AOI Merge Mode

	
USB 3.0	USB 3.0

Syntax

```

INT is_DeviceFeature (HCAM hCam, UINT nCommand, void* pParam, UINT cbSizeOfParam)

```

Description

Using [is_DeviceFeature\(\)](#) you can configure special camera functions provided by specific uc480 models:



The AOI merge mode is currently supported by the following camera models:

- DCC3240
- DCC3260

Flashing per line

Note that depending on the set line rate, the optocoupler of the flash output may be too slow if the flash output is set to flash high active or flash low active. In this case, use a GPIO as flash output because the GPIO have no optocoupler.

The `nCommand` input parameter is used to select the function mode. The `pParam` input parameter depends on the selected function mode. If you select functions for setting or returning a value, `pParam` contains a pointer to a variable of the `UINT` type. The size of the memory area to which `pParam` refers is specified in the `cbSizeOfParam` input parameter.

Input parameters

hCam	Camera handle
<div style="border: 1px solid gray; padding: 5px;"> <input type="checkbox"/> nCommand </div>	

IS_DEVICE_FEATURE_CMD_GET_VERTICAL_AOI_MERGE_MODE_SUPPORTED_LINE_MODES	<p>Returns the supported modes of the AOI merge mode (freerun, software trigger or hardware trigger):</p> <ul style="list-style-type: none"> • 1 = IS_VERTICAL_AOI_MERGE_MODE_LINE_FREERUN: AOI merge mode in freerun mode is supported. • 2 = IS_VERTICAL_AOI_MERGE_MODE_LINE_SOFTWARE_TRIGGER: AOI merge mode with software trigger is supported. • 4 = IS_VERTICAL_AOI_MERGE_MODE_LINE_GPIO_TRIGGER: AOI merge mode with hardware trigger is supported. <p>(Example 3)</p>
IS_DEVICE_FEATURE_CMD_GET_VERTICAL_AOI_MERGE_MODE_DEFAULT	Returns the default value of the AOI merge mode
IS_DEVICE_FEATURE_CMD_GET_VERTICAL_AOI_MERGE_MODE	Returns the current set AOI merge mode (Example 1)
IS_DEVICE_FEATURE_CMD_SET_VERTICAL_AOI_MERGE_MODE	<p>Sets the AOI merge mode:</p> <ul style="list-style-type: none"> • 0 = IS_VERTICAL_AOI_MERGE_MODE_OFF: Disables the AOI merge mode • 1 = IS_VERTICAL_AOI_MERGE_MODE_FREERUN: The sensor runs with maximum speed. • 2 = IS_VERTICAL_AOI_MERGE_MODE_TRIGGERED_SOFTWARE: The sensor is triggered via software. • 3 = IS_VERTICAL_AOI_MERGE_MODE_TRIGGERED_FALLING_GPIO1: The sensor is triggered on GPIO 1 (falling edge). • 4 = IS_VERTICAL_AOI_MERGE_MODE_TRIGGERED_RISING_GPIO1: The sensor is triggered on GPIO 1 (rising edge). • 5 = IS_VERTICAL_AOI_MERGE_MODE_TRIGGERED_FALLING_GPIO2: The sensor is triggered on GPIO 2 (falling edge). • 6 = IS_VERTICAL_AOI_MERGE_MODE_TRIGGERED_RISING_GPIO2: The sensor is triggered on GPIO 2 (rising edge).
IS_DEVICE_FEATURE_CMD_GET_VERTICAL_AOI_MERGE_POSITION_DEFAULT	Returns the default value for the AOI merge position
IS_DEVICE_FEATURE_CMD_GET_VERTICAL_AOI_MERGE_POSITION_RANGE	Returns the range for the AOI merge position (IS_RANGE_S32)
IS_DEVICE_FEATURE_CMD_GET_VERTICAL_AOI_MERGE_POSITION	Returns the position of the two sensor lines (default = 0, i.e. the two top lines)
IS_DEVICE_FEATURE_CMD_SET_VERTICAL_AOI_MERGE_POSITION	Sets the position of the two sensor lines

pParam	Pointer to a function parameter, whose function depends on nCommand.
cbSizeOfParam	Size (in bytes) of the memory area to which pParam refers.

Status flags from DEVICE_FEATURE_MODE_CAPS

IS_DEVICE_FEATURE_CAP_VERTICAL_AOI_MERGE	Special AOI merge mode which combines the lines of an AOI to a new image. Example: The AOI merge mode height is 1 line for a monochrome sensor. That means that the image is composed of 4000 "merged" AOIs. If the AOI merge mode height is 100 lines, the image is composed of 400 "merged" AOIs.
--	--

Example 1

```

INT nVerticalAoiMergeMode = 0;
INT nVerticalAoiMergePosition = 0;

/* Read current vertical AOI merge mode */
nRet = is_DeviceFeature(hCam, IS_DEVICE_FEATURE_CMD_GET_VERTICAL_AOI_MERGE_MODE,
                       (void*)&nVerticalAoiMergeMode, sizeof(nVerticalAoiMergeMode));
if (nRet == IS_SUCCESS)
{
    if (nVerticalAoiMergeMode == IS_VERTICAL_AOI_MERGE_MODE_FREERUN)
    {
        /* Read current AOI merge position */
        nRet = is_DeviceFeature(hCam, IS_DEVICE_FEATURE_CMD_GET_VERTICAL_AOI_MERGE_POSITION,
                               (void*)&nVerticalAoiMergePosition,
                               sizeof(nVerticalAoiMergePosition));
    }
}

/* Set vertical AOI merge mode */
nVerticalAoiMergeMode = IS_VERTICAL_AOI_MERGE_MODE_FREERUN;
nRet = is_DeviceFeature(hCam, IS_DEVICE_FEATURE_CMD_SET_VERTICAL_AOI_MERGE_MODE,
                       (void*)&nVerticalAoiMergeMode, sizeof(nVerticalAoiMergeMode));

if (nRet == IS_SUCCESS)
{
    INT nVerticalAoiMergePosition = 100;
    nRet = is_DeviceFeature(hCam, IS_DEVICE_FEATURE_CMD_SET_VERTICAL_AOI_MERGE_POSITION,
                           (void*)&nVerticalAoiMergePosition, sizeof(nVerticalAoiMergePosition));
}

```

Example 2

```

/* Get the default value for the vertical AOI merge mode height */
INT nHeight = 0;
nRet = is_DeviceFeature(hCam, IS_DEVICE_FEATURE_CMD_GET_VERTICAL_AOI_MERGE_HEIGHT_DEFAULT,
                       (void*)&nHeight, sizeof(nHeight));

/* Get current value of the vertical AOI merge mode height */
nRet = is_DeviceFeature(hCam, IS_DEVICE_FEATURE_CMD_GET_VERTICAL_AOI_MERGE_HEIGHT,
                       (void*)&nHeight, sizeof(nHeight));

/* Get the number of elements in the vertical AOI merge mode height list */
INT nVerticalAoiMergeModeHeightNumber;
nRet = is_DeviceFeature(hCam, IS_DEVICE_FEATURE_CMD_GET_VERTICAL_AOI_MERGE_HEIGHT_NUMBER,
                       (void*)&nVerticalAoiMergeModeHeightNumber,

```

```

        sizeof(nVerticalAoiMergeModeHeightNumber));
if (nRet == IS_SUCCESS)
{
    UINT* pVerticalAoiMergeModeHeightList = new UINT[nVerticalAoiMergeModeHeightNumber];

    /* Get the vertical AOI merge mode height list */
    nRet = is_DeviceFeature(hCam, IS_DEVICE_FEATURE_CMD_GET_VERTICAL_AOI_MERGE_HEIGHT_LIST,
        (void*)pVerticalAoiMergeModeHeightList,
        nVerticalAoiMergeModeHeightNumber * sizeof(UINT));

    if (nRet == IS_SUCCESS)
    {
        /* Set the maximum possible vertical AOI merge mode height depending on the current image height */
        UINT nMaxHeight = pVerticalAoiMergeModeHeightList[nVerticalAoiMergeModeHeightNumber - 1];
        nRet = is_DeviceFeature(hCam, IS_DEVICE_FEATURE_CMD_SET_VERTICAL_AOI_MERGE_HEIGHT,
            (void*)&nMaxHeight, sizeof(nMaxHeight));
    }
}
}

```

Example 3

```

UINT nModes;

/* Get mask with all supported vertical AOI merge line modes */
INT nRet = is_DeviceFeature(hCam, IS_DEVICE_FEATURE_CMD_GET_VERTICAL_AOI_MERGE_MODE_SUPPORTED_LINE_MODES);



if (nRet == IS_SUCCESS)
{
    if (nModes & IS_VERTICAL_AOI_MERGE_MODE_LINE_FREERUN)
    {
        // freerun supported
    }

    if (nModes & IS_VERTICAL_AOI_MERGE_MODE_LINE_SOFTWARE_TRIGGER)
    {
        // software triggered supported
    }

    if (nModes & IS_VERTICAL_AOI_MERGE_MODE_LINE_GPIO_TRIGGER)
    {
        // GPIO triggered supported
    }
}

```

4.3.15.2 Using the Log Mode

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_DeviceFeature (HCAM hCam, UINT nCommand, void* pParam, UINT cbSizeOfParam)
```

Description

Using [is_DeviceFeature\(\)](#) you can configure special camera functions provided by specific uc480 models:

- DCC1240/DCC3240

The `nCommand` input parameter is used to select the function mode. The `pParam` input parameter depends on the selected function mode. If you select functions for setting or returning a value, `pParam` contains a pointer to a variable of the `UINT` type. The size of the memory area to which `pParam` refers is specified in the `cbSizeOfParam` input parameter.

Input parameters

hCam	Camera handle
------	---------------

nCommand	
IS_DEVICE_FEATURE_CMD_GET_LOG_MODE_DEFAULT	Returns the default settings for the Log mode (Example 1)
IS_DEVICE_FEATURE_CMD_GET_LOG_MODE	Returns the current Log mode (Example 2)
IS_DEVICE_FEATURE_CMD_SET_LOG_MODE	Sets the Log mode: <ul style="list-style-type: none"> • 0 = IS_LOG_MODE_FACTORY_DEFAULT: Factory setting for the Log mode • 1 = IS_LOG_MODE_OFF: Log mode off • 2 = IS_LOG_MODE_MANUAL: Manual Log mode. In this case the Log mode value and the Log mode gain are effective. • 3 = IS_LOG_MODE_AUTO: Automatic Log mode (default setting)
IS_DEVICE_FEATURE_CMD_GET_LOG_MODE_MANUAL_VALUE_DEFAULT	Returns the default settings for the manual value of the Log mode (Example 3)
IS_DEVICE_FEATURE_CMD_GET_LOG_MODE_MANUAL_VALUE_RANGE	Returns the range of the manual value of the Log mode.
IS_DEVICE_FEATURE_CMD_GET_LOG_MODE_MANUAL_VALUE	Returns the current manual value of the Log mode (Example 4)
IS_DEVICE_FEATURE_CMD_SET_LOG_MODE_MANUAL_VALUE	Sets the manual value of the Log mode.
IS_DEVICE_FEATURE_CMD_GET_LOG_MODE_MANUAL_GAIN_DEFAULT	Returns the default settings for the manual gain for the Log mode (Example 5)
IS_DEVICE_FEATURE_CMD_GET_LOG_MODE_MANUAL_GAIN_RANGE	Returns the range for the manual gain of the Log mode.
IS_DEVICE_FEATURE_CMD_GET_LOG_MODE_MANUAL_GAIN	Returns the current manual gain of the Log mode (Example 6)
IS_DEVICE_FEATURE_CMD_SET_LOG_MODE_MANUAL_GAIN	Sets the manual gain of the Log mode.
pParam	Pointer to a function parameter, whose function depends on nCommand.
cbSizeOfParam	Size (in bytes) of the memory area to which pParam refers.

Status flags from DEVICE_FEATURE_MODE_CAPS

IS_DEVICE_FEATURE_CAP_LOG_MODE	Log mode is supported/Set mode
--------------------------------	--------------------------------

Example 1

```

/* Read and set default Log mode */
UINT nDefaultLogMode = 0;
INT nRet = is_DeviceFeature(hCam, IS_DEVICE_FEATURE_CMD_GET_LOG_MODE_DEFAULT,
    (void*)&nDefaultLogMode, sizeof(nDefaultLogMode));
if (nRet == IS_SUCCESS) {
    nRet = is_DeviceFeature(hCam, IS_DEVICE_FEATURE_CMD_SET_LOG_MODE,
        (void*)&nDefaultLogMode, sizeof(nDefaultLogMode));
}
    
```


Example 2

```

/* Read current Log pixel mode */
UINT nLogMode = 0;
INT nRet = is_DeviceFeature(hCam, IS_DEVICE_FEATURE_CMD_GET_LOG_MODE,
                            (void*)&nLogMode, sizeof(nLogMode));

```

Example 3

```

/* Read default Log pixel mode manual value */
UINT nDefaultLogModeManualValue = 0;
INT nRet = is_DeviceFeature(hCam, IS_DEVICE_FEATURE_CMD_GET_LOG_MODE_MANUAL_VALUE_DEFAULT,
                            (void*)&nDefaultLogModeManualValue,
                            sizeof(nDefaultLogModeManualValue));

/* Get the range of the manual value */
IS_RANGE_S32 nLogModeManualValueRange;
nRet = is_DeviceFeature(hCam, IS_DEVICE_FEATURE_CMD_GET_LOG_MODE_MANUAL_VALUE_RANGE,
                        (void*)&nLogModeManualValueRange,
                        sizeof(nLogModeManualValueRange));

if (nRet == IS_SUCCESS) {
    INT nMin = nLogModeManualValueRange.s32Min;
    INT nMax = nLogModeManualValueRange.s32Max;
    INT nInc = nLogModeManualValueRange.s32Inc;
}

```

Example 4

```

UINT nLogModeValue = 0;

/* Read current Log mode manual value */
nRet = is_DeviceFeature(hCam, IS_DEVICE_FEATURE_CMD_GET_LOG_MODE_MANUAL_VALUE,
                        (void*)&nLogModeValue, sizeof(nLogModeValue));

/* Set log pixel mode value */
nRet = is_DeviceFeature(hCam, IS_DEVICE_FEATURE_CMD_SET_LOG_MODE_MANUAL_VALUE,
                        (void*)&nLogModeValue, sizeof(nLogModeValue));

```

Example 5

```

/* Read default Log mode manual gain */
UINT nDefaultLogModeManualGain = 0;
INT nRet = is_DeviceFeature(hCam, IS_DEVICE_FEATURE_CMD_GET_LOG_MODE_MANUAL_GAIN_DEFAULT,
                            (void*)&nDefaultLogModeManualGain,
                            sizeof(nDefaultLogModeManualGain));

/* Get the range of the manual value */
IS_RANGE_S32 nLogModeManualGainRange;
nRet = is_DeviceFeature(hCam, IS_DEVICE_FEATURE_CMD_GET_LOG_MODE_MANUAL_GAIN_RANGE,
                        (void*)&nLogModeManualGainRange,
                        sizeof(nLogModeManualGainRange));

if (nRet == IS_SUCCESS) {
    INT nMin = nLogModeManualGainRange.s32Min;
    INT nMax = nLogModeManualGainRange.s32Max;
    INT nInc = nLogModeManualGainRange.s32Inc;
}

```

Example 6

```



UINT nLogModeGain = 0;

/* Read current Log mode gain */
nRet = is_DeviceFeature(hCam, IS_DEVICE_FEATURE_CMD_GET_LOG_MODE_MANUAL_GAIN,
                        (void*)&nLogModeGain, sizeof(nLogModeGain));

/* Set Log mode gain*/
nRet = is_DeviceFeature(hCam, IS_DEVICE_FEATURE_CMD_SET_LOG_MODE_MANUAL_GAIN,
                        (void*)&nLogModeGain, sizeof(nLogModeGain));

```

4.3.15.3 Using Level Controlled Trigger

	
USB 3.0	USB 3.0

Syntax

```
INT is_DeviceFeature (HCAM hCam, UINT nCommand, void* pParam, UINT cbSizeOfParam)
```

Description

Using [is_DeviceFeature\(\)](#) you can configure special camera functions provided by specific uc480 models:

- DCC3260

By using the level-controlled trigger, you can control the exposure time directly via the duration of the trigger signal.



The level controlled trigger cannot be used in combination with the software trigger mode, trigger prescaler, or the AOI merge mode.

The `nCommand` input parameter is used to select the function mode. The `pParam` input parameter depends on the selected function mode. If you select functions for setting or returning a value, `pParam` contains a pointer to a variable of the `UINT` type. The size of the memory area to which `pParam` refers is specified in the `cbSizeOfParam` input parameter.

Input parameters

<code>hCam</code>	Camera handle
<div style="border: 1px solid gray; padding: 2px;"> <input type="checkbox"/> <code>nCommand</code> </div>	
<code>IS_DEVICE_FEATURE_CMD_SET_LEVEL_CONTROL_LED_TRIGGER_INPUT_MODE</code>	Enables/disables the use of the level controlled trigger (Example): <ul style="list-style-type: none"> • 1 = <code>IS_LEVEL_CONTROLLED_TRIGGER_INPUT_ON</code>: Enables the level controlled trigger • 0 = <code>IS_LEVEL_CONTROLLED_TRIGGER_INPUT_OFF</code>: Disables the level controlled trigger
<code>IS_DEVICE_FEATURE_CMD_GET_LEVEL_CONTROL_LED_TRIGGER_INPUT_MODE</code>	Returns the current state of the level controlled trigger (Example).
<code>IS_DEVICE_FEATURE_CMD_GET_LEVEL_CONTROL_LED_TRIGGER_INPUT_MODE_DEFAULT</code>	Returns the default settings for the level controlled trigger (Example).
<code>pParam</code>	Pointer to a function parameter, whose function depends on <code>nCommand</code> .
<code>cbSizeOfParam</code>	Size (in bytes) of the memory area to which <code>pParam</code> refers.

Status flags from `DEVICE_FEATURE_MODE_CAPS`

<code>IS_DEVICE_FEATURE_CAP_LEVEL_CONTROLLED_TRIGGER</code>	The use of the level controlled trigger is supported.
---	---

Example

```

/* Query the current settings */
INT nMode = IS_LEVEL_CONTROLLED_TRIGGER_INPUT_OFF;
INT nRet = is_DeviceFeature(hCam, IS_DEVICE_FEATURE_CMD_GET_LEVEL_CONTROLLED_TRIGGER_INPUT_MODE,
                           (void*)&nMode, sizeof(nMode));

if (IS_SUCCESS == nRet)
{
    if (IS_LEVEL_CONTROLLED_TRIGGER_INPUT_OFF == nMode)
    {
        /* Level controlled trigger is disabled */
    }
    else if (IS_LEVEL_CONTROLLED_TRIGGER_INPUT_ON == nMode)
    {
        /* Level controlled trigger is enabled */
    }
}

/* Enable level controlled trigger */
nMode = IS_LEVEL_CONTROLLED_TRIGGER_INPUT_ON;
nRet = is_DeviceFeature(hCam, IS_DEVICE_FEATURE_CMD_SET_LEVEL_CONTROLLED_TRIGGER_INPUT_MODE,
                       (void*)&nMode, sizeof(nMode));



if (IS_SUCCESS == nRet)
{
    /* Level controlled trigger is enabled */
}

/* Get default setting */
nRet = is_DeviceFeature(hCam,
                       IS_DEVICE_FEATURE_CMD_GET_LEVEL_CONTROLLED_TRIGGER_INPUT_MODE_DEFAULT,
                       (void*)&nMode, sizeof(nMode));

if (IS_SUCCESS == nRet)
{
    /* nMode returns the default setting */
}

```

4.3.15.4 Switching the Shutter Mode

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_DeviceFeature (HCAM hCam, UINT nCommand, void* pParam, UINT cbSizeOfParam)
```

Description

Using [is_DeviceFeature\(\)](#) you can configure special camera functions provided by specific uc480 models:

- DCC1240/DCC3240 models: Toggle between shutter modes, see [Basics: Shutter methods](#).

The `nCommand` input parameter is used to select the function mode. The `pParam` input parameter depends on the selected function mode. If you select functions for setting or returning a value, `pParam` contains a pointer to a variable of the `UINT` type. The size of the memory area to which `pParam` refers is specified in the `cbSizeOfParam` input parameter.

Input parameters



hCam	Camera handle
<input type="checkbox"/> nCommand	
IS_DEVICE_FEATURE_CMD_SET_SHUTTER_MODE	Sets the shutter mode, see <code>DEVICE_FEATURE_MODE_CAPS</code>

IS_DEVICE_FEATURE_CMD_GET_SHUTTER_MODE	Returns the shutter mode
pParam	Pointer to a function parameter, whose function depends on nCommand.
cbSizeOfParam	Size (in bytes) of the memory area to which pParam refers.

Status flags from DEVICE_FEATURE_MODE_CAPS

IS_DEVICE_FEATURE_CAP_SHUTTER_MODE_ROLLING	Rolling shutter mode is supported/Set mode
IS_DEVICE_FEATURE_CAP_SHUTTER_MODE_ROLLING_GLOBAL_START	Rolling shutter mode with global start is supported/Set mode
IS_DEVICE_FEATURE_CAP_SHUTTER_MODE_GLOBAL	Global shutter mode is supported/Set mode
IS_DEVICE_FEATURE_CAP_SHUTTER_MODE_GLOBAL_ALTERNATIVE_TIMING	Global shutter mode with different timing parameters is supported/Set mode

4.3.15.5 Using the Internal Image Memory

	
USB 3.0	USB 3.0

Syntax

INT is_DeviceFeature (HCAM hCam, UINT nCommand, void* pParam, UINT cbSizeOfParam)

Description

Using [is_DeviceFeature\(\)](#) you can configure special camera functions provided by specific uc480 models:

- DCC3260: Enable or disable the image memory of the camera. By default, the camera is operated without image memory. Note that the camera has to be closed if you enable or disable the image memory.



Hint: Note that from driver version 4.80 on the internal image memory can only be enabled or disabled when the device ID and not the camera handle is used for calling. Please refer to the [Application Note on our website](#) for more information on this topic.

The nCommand input parameter is used to select the function mode. The pParam input parameter depends on the selected function mode. If you select functions for setting or returning a value, pParam contains a pointer to a variable of the UINT type. The size of the memory area to which pParam refers is specified in the cbSizeOfParam input parameter.

Input parameters

hCam	Camera handle
nCommand	
IS_DEVICE_FEATURE_CMD_GET_MEMORY_MODE_ENABLED	Returns if the image memory of the camera is enabled (Example 3).
IS_DEVICE_FEATURE_CMD_SET_MEMORY_MODE_ENABLED	Enables/disables the image memory of the camera (Example 2): <ul style="list-style-type: none"> • 1 = IS_MEMORY_MODE_ON • 0 = IS_MEMORY_MODE_OFF

IS_DEVICE_FEATURE_CMD_GET_MEMORY_MODE_ENABLE_DEFAULT	Returns the default setting (Example 4).
IS_DEVICE_FEATURE_CMD_GET_MEMORY_MODE_ENABLE_SUPPORTED	Returns if the internal image memory is supported (Example 1).
pParam	Pointer to a function parameter, whose function depends on nCommand.
cbSizeOfParam	Size (in bytes) of the memory area to which pParam refers.

Status flags from DEVICE_FEATURE_MODE_CAPS

IS_DEVICE_FEATURE_CAP_MEMORY_MODE	Internal image memory is supported.
-----------------------------------	-------------------------------------

Related functions

- Querying the temperature state
- [is_EnableEvent\(\)](#)
- [is_EnableMessage\(\)](#)

Example 1

```
// Is using the image memory supported?
// Via camera handle
int nMemoryMode;
nRet = is_DeviceFeature(hCam, IS_DEVICE_FEATURE_CMD_GET_MEMORY_MODE_ENABLE_SUPPORTED, &nMemoryMode, sizeof(nMemoryMode));
if (nMemoryMode == IS_MEMORY_MODE_ON)
{
    // is supported...
}

// Is using the image memory supported?
// Via device ID
int nMemoryMode;
nRet = is_DeviceFeature(nDevID | IS_USE_DEVICE_ID, IS_DEVICE_FEATURE_CMD_GET_MEMORY_MODE_ENABLE_SUPPORTED, &nMemoryMode, sizeof(nMemoryMode));
if (nMemoryMode == IS_MEMORY_MODE_ON)
{
    // is supported...
}
```

Example 2

```
// Enable memory mode, device ID must be used for this purpose
UINT nEnable = IS_MEMORY_MODE_ON;
INT nRet = is_DeviceFeature(nDevID | IS_USE_DEVICE_ID, IS_DEVICE_FEATURE_CMD_SET_MEMORY_MODE_ENABLE, nEnable, sizeof(nEnable));
```

Example 3

```
// Query the status of the image memory via camera handle
int nMemoryMode;
nRet = is_DeviceFeature(hCam, IS_DEVICE_FEATURE_CMD_GET_MEMORY_MODE_ENABLE, &nMemoryMode, sizeof(nMemoryMode));
if (nMemoryMode == IS_MEMORY_MODE_ON)
{
    ...
}



// Query the status of the image memory via device ID
int nMemoryMode;
nRet = is_DeviceFeature(nDevID | IS_USE_DEVICE_ID, IS_DEVICE_FEATURE_CMD_GET_MEMORY_MODE_ENABLE, &nMemoryMode, sizeof(nMemoryMode));
if (nMemoryMode == IS_MEMORY_MODE_ON)
{
    ...
}
```

Example 4

```
// Query the default setting for the image memory via camera handle
int nMemoryMode;
nRet = is_DeviceFeature(hCam, IS_DEVICE_FEATURE_CMD_GET_MEMORY_MODE_ENABLE_DEFAULT, &nMemoryMode, sizeof(nMemoryMode));
if (nMemoryMode == IS_MEMORY_MODE_ON)
{
    ...
}

// Query the default setting for the image memory via device ID
int nMemoryMode;
nRet = is_DeviceFeature(nDevID | IS_USE_DEVICE_ID, IS_DEVICE_FEATURE_CMD_GET_MEMORY_MODE_ENABLE_DEFAULT, &nMemoryMode, sizeof(nMemoryMode));
if (nMemoryMode == IS_MEMORY_MODE_ON)
{
    ...
}
```

4.3.15.6 Using the Line Scan Mode

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_DeviceFeature (HCAM hCam, UINT nCommand, void* pParam, UINT cbSizeOfParam)
```

Description

Using [is_DeviceFeature\(\)](#) you can configure special camera functions provided by specific uc480 models:

- DCC1240/DCC3240



Only the monochrome camera models DCC1240/DCC3240 currently support the "fast line scan" mode.

The `nCommand` input parameter is used to select the function mode. The `pParam` input parameter depends on the selected function mode. If you select functions for setting or returning a value, `pParam` contains a pointer to a variable of the `UINT` type. The size of the memory area to which `pParam` refers is specified in the `cbSizeOfParam` input parameter.



Input parameters

<code>hCam</code>	Camera handle
<input type="checkbox"/> <code>nCommand</code>	
<code>IS_DEVICE_FEATURE_CMD_SET_LINESCAN_MODE</code>	Sets the line scan mode, see <code>DEVICE_FEATURE_MODE_CAPS</code>
<code>IS_DEVICE_FEATURE_CMD_GET_LINESCAN_MODE</code>	Returns the current set line scan mode
<code>IS_DEVICE_FEATURE_CMD_SET_LINESCAN_NUMBER</code>	Sets the scan line used for the line scan mode
<code>IS_DEVICE_FEATURE_CMD_GET_LINESCAN_NUMBER</code>	Returns the scan line used for the line scan mode
<code>pParam</code>	Pointer to a function parameter, whose function depends on <code>nCommand</code> .
<code>cbSizeOfParam</code>	Size (in bytes) of the memory area to which <code>pParam</code> refers.

Status flags from DEVICE_FEATURE_MODE_CAPS

IS_DEVICE_FEATURE_CAP_LINESCAN_MODE_FAST	Fast line scan mode is supported/Set mode
IS_DEVICE_FEATURE_CAP_LINESCAN_NUMBER	Line number at fast line scan mode is supported/Set number

4.3.15.7 Configuring the Timestamp

	
USB 3.0	USB 3.0

Syntax

```
INT is_DeviceFeature (HCAM hCam, UINT nCommand, void* pParam, UINT cbSizeOfParam)
```

Description

Using [is_DeviceFeature\(\)](#) you can configure special camera functions provided by specific uc480 models:

The `nCommand` input parameter is used to select the function mode. The `pParam` input parameter depends on the selected function mode. If you select functions for setting or returning a value, `pParam` contains a pointer to a variable of the `UINT` type. The size of the memory area to which `pParam` refers is specified in the `cbSizeOfParam` input parameter.

Input parameters

<code>hCam</code>	Camera handle
<input type="checkbox"/> <code>nCommand</code>	
IS_DEVICE_FEATURE_CMD_GET_TIMESTAMP_CONFIGURATION	Returns the current configuration of the timestamp (Example)
IS_DEVICE_FEATURE_CMD_SET_TIMESTAMP_CONFIGURATION	Sets the new configuration of the timestamp
<code>pParam</code>	Pointer to a function parameter, whose function depends on <code>nCommand</code> .
<code>cbSizeOfParam</code>	Size (in bytes) of the memory area to which <code>pParam</code> refers.

Contents of the IS_TIMESTAMP_CONFIGURATION structure

INT	<code>s32Mode</code>	Mode to be set: <ul style="list-style-type: none"> IS_RESET_TIMESTAMP_ONCE: Sets once the camera timestamp to 0.
INT	<code>s32Pin</code>	Pin for setting the timestamp: <ul style="list-style-type: none"> TIMESTAMP_CONFIGURATION_PIN_NONE: Timestamp is not reset. TIMESTAMP_CONFIGURATION_PIN_TRIGGER: Timestamp is set by signal on the trigger pin. TIMESTAMP_CONFIGURATION_PIN_GPIO_1: Timestamp is set by

		<p>signal on GPIO 1.</p> <ul style="list-style-type: none"> • <code>TIMESTAMP_CONFIGURATION_PIN_GPIO_2</code>: Timestamp is set by signal on GPIO 2.
INT	s32Edge	<p>The timestamp is set by rising or falling edge:</p> <ul style="list-style-type: none"> • <code>TIMESTAMP_CONFIGURATION_EDGE_FALLING</code> • <code>TIMESTAMP_CONFIGURATION_EDGE_RISING</code>

Status flags from `DEVICE_FEATURE_MODE_CAPS`

<code>IS_DEVICE_FEATURE_CAP_TIMESTAMP_CONFIGURATION</code>	Configuration of the timestamp (e.g. reset the timestamp to 0 at a signal on the trigger pin).
--	--

Example

```
IS_TIMESTAMP_CONFIGURATION timestampConfiguration;

/* Returns the current timestamp configuration */
nRet = is_DeviceFeature(hCam, IS_DEVICE_FEATURE_CMD_GET_TIMESTAMP_CONFIGURATION,
    (void*) &timestampConfiguration, sizeof(timestampConfiguration));



/* A falling edge on the trigger pin will reset the device timestamp once. All 3 parameters will be set at t
timestampConfiguration.s32Mode = IS_RESET_TIMESTAMP_ONCE;
timestampConfiguration.s32Edge = TIMESTAMP_CONFIGURATION_EDGE_FALLING; timestampConfiguration.s32Pin = TIMES
nRet = is_DeviceFeature(hCam, IS_DEVICE_FEATURE_CMD_SET_TIMESTAMP_CONFIGURATION,
    (void*) &timestampConfiguration, sizeof(timestampConfiguration));

/* Set only the edge */
timestampConfiguration.s32Mode = IS_IGNORE_PARAMETER;
timestampConfiguration.s32Edge = TIMESTAMP_CONFIGURATION_EDGE_RISING;
timestampConfiguration.s32Pin = IS_IGNORE_PARAMETER;
nRet = is_DeviceFeature(hCam, IS_DEVICE_FEATURE_CMD_SET_TIMESTAMP_CONFIGURATION,
    (void*)&timestampConfiguration, sizeof(timestampConfiguration));

/* Set only the pin */
timestampConfiguration.s32Mode = IS_IGNORE_PARAMETER;
timestampConfiguration.s32Edge = IS_IGNORE_PARAMETER;
timestampConfiguration.s32Pin = TIMESTAMP_CONFIGURATION_PIN_GPIO_1;
nRet = is_DeviceFeature(hCam, IS_DEVICE_FEATURE_CMD_SET_TIMESTAMP_CONFIGURATION,
    (void*)&timestampConfiguration, sizeof(timestampConfiguration));

/* Set only the mode */
timestampConfiguration.s32Mode = IS_RESET_TIMESTAMP_ONCE;
timestampConfiguration.s32Edge = IS_IGNORE_PARAMETER;
timestampConfiguration.s32Pin = IS_IGNORE_PARAMETER;
nRet = is_DeviceFeature(hCam, IS_DEVICE_FEATURE_CMD_SET_TIMESTAMP_CONFIGURATION,
    (void*)&timestampConfiguration, sizeof(timestampConfiguration));
```

4.3.16 is_DeviceInfo

	
USB 3	USB 3

Note

This command is supported by DCC3240x USB3 cameras only!

Syntax

```
INT is_DeviceInfo (HCAM hCam, UINT nCommand, void* pParam, UINT cbSizeOfParam)
```


Description

Using `is_DeviceInfo()`, you can query information about connected USB 3 Cameras. The resulting information is written to the `IS_DEVICE_INFO` structure. For this purpose, the cameras need not be opened.

The `nCommand` input parameter is used to select the function mode. The `pParam` input parameter depends on the selected function mode. If you select functions for setting or returning a value, `pParam` contains a pointer to a variable of the `UINT` type. The size of the memory area to which `pParam` refers is specified in the `cbSizeOfParam` input parameter.

Input parameters

<code>hCam</code>	Camera handle
<input type="checkbox"/> <code>nCommand</code>	
<code>IS_DEVICE_INFO_CMD_GET_DEVICE_INFO</code>	Returns a information structure about the specified device (Example 1)
<code>pParam</code>	Pointer to a function parameter, whose function depends on <code>nCommand</code> .
<code>cbSizeOfParam</code>	Size (in bytes) of the memory area to which <code>pParam</code> refers.

Contents of the `IS_DEVICE_INFO` structure

<code>IS_DEVICE_INFO_HEARTBEAT</code>	<code>infoDevHeartbeat</code>	Camera-related data retrieved from the camera (from the heartbeat telegram) See IS_DEVICE_INFO_HEARTBEAT
<code>IS_DEVICE_INFO_CONTROL</code>	<code>infoDevControl</code>	Camera-related driver data See IS_DEVICE_INFO_CONTROL
<code>BYTE</code>	<code>reserved[240]</code>	reserved

Contents of the `IS_DEVICE_INFO::IS_DEVICE_INFO_HEARTBEAT` structure

<code>BYTE</code>	<code>reserved_1[36]</code>	reserved
<code>DWORD</code>	<code>dwRuntimeFirmwareVersion</code>	Firmware version
<code>BYTE</code>	<code>reserved_2[8]</code>	reserved
<code>WORD</code>	<code>wTemperature</code>	Camera temperature in degree Celsius Bit 15: algebraic sign Bit 14...11: filled according to algebraic sign Bit 10...4: temperature (places before the decimal point) Bit 3...0: temperature (places after the decimal point) See the Ambient conditions chapter for permissible temperature range.
<code>WORD</code>	<code>wLinkSpeed_Mb</code>	Transfer rate: <ul style="list-style-type: none"> <code>IS_USB_HIGH_SPEED</code> = 480 <code>IS_USB_SUPER_SPEED</code> = 4000

BYTE	reserved[208]	reserved
------	---------------	----------

Contents of the IS_DEVICE_INFO::IS_DEVICE_INFO_CONTROL structure

DWORD	dwDeviceId	Internal device ID of the camera
BYTE	reserved[148]	reserved

Return values

IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
IS_IO_REQUEST_FAILED	An IO request from the uc480 driver failed. Possibly the versions of the uc480.dll (API) and the driver file (uc480_usb.sys) do not match.
IS_SUCCESS	Function executed successfully

Related functions

- [is_GetCameraList\(\)](#)



Example 1

```
// The camera has the device ID 1
UINT nDeviceId = 1;

IS_DEVICE_INFO deviceInfo;
memset(&deviceInfo, 0, sizeof(IS_DEVICE_INFO));

INT nRet = is_DeviceInfo((HCAM)(nDeviceId | IS_USE_DEVICE_ID),
                        IS_DEVICE_INFO_CMD_GET_DEVICE_INFO,
                        (void*)&deviceInfo, sizeof(deviceInfo));
if (nRet == IS_SUCCESS)
{
    WORD wTemperature = deviceInfo.infoDevHeartbeat.wTemperature;
}
```

4.3.17 is_DirectRenderer

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_DirectRenderer (HCAM hCam, UINT nMode, void* pParam, UINT nSize)
```

Description

Note

The `is_DirectRenderer()` functions works under Linux only in OpenGL mode.

`is_DirectRenderer()` provides a set of advanced rendering functions and allows inserting overlay data into the camera's live image without flicker. The graphics card functions of the Direct3D library are supported under Windows.

The second input parameter `nMode` specifies the effect of the `is_DirectRenderer()` call.

The value of the third parameter `pParam` depends on the mode selected with `nMode`: For example, when setting the overlay size (`nMode = DR_SET_OVERLAY_SIZE`), a pointer to an array of two values (x and y) is passed (see [code samples](#)). When you load a bitmap image (`nMode = DR_LOAD_OVERLAY_FROM_FILE`), `pParam` passes the path to the file (see [code samples](#)). The required parameters are illustrated in the sample codes at the end of this section.

Notes

1. System requirements

- To use the Direct3D functionality, the appropriate version of the Microsoft DirectX Runtime has to be installed in your PC.
- When you are using high-resolution cameras, the maximum texture size supported by the graphics card should be at least 4096 x 4096 pixels. You can check the maximum texture size by reading out the `D3D_GET_MAX_OVERLAY_SIZE` parameter.
- The Direct3D mode automatically uses the Windows Desktop color depth setting for the display.

Please also read the notes on graphics cards which are provided in the [System requirements](#) chapter.

2. Displaying monochrome or raw data formats

To display monochrome or Bayer raw data in Direct3D, please set the appropriate constants using the [is_SetDisplayMode\(\)](#) function.

The `nCommand` input parameter is used to select the function mode. The `pParam` input parameter depends on the selected function mode. If you select functions for setting or returning a value, `pParam` contains a pointer to a variable of the `UINT` type. The size of the memory area to which `pParam` refers is specified in the `nSizeOfParam` input parameter.

Input parameters

hCam	Camera handle
<input type="checkbox"/> nMode	
DR_GET_SUPPORTED	Returns either if Direct3D or OpenGL is

hCam	Camera handle
	<p>supported by the graphics card.</p> <ul style="list-style-type: none"> ▣ More details <ul style="list-style-type: none"> • IS_SET_DM_DIRECT3D: Tests if Direct3D is supported. • IS_SET_DM_OPENGL: Tests if OpenGL is supported. <p>Example</p>
DR_GET_OVERLAY_DC	<p>Direct3D only: Returns the device context (DC) handle to the overlay area of the graphics card.</p> <ul style="list-style-type: none"> ▣ More details <p>In Direct3D mode, the DR_GET_OVERLAY_DC mode returns the device context (DC) handle of the overlay area. Using this handle, it is possible to access the overlay using the Windows GDI functionality. Thus, all Windows graphics commands (e.g. Line, Circle, Rectangle, TextOut) are available. To transfer the drawn elements to the overlay, release the DC handle by calling DR_RELEASE_OVERLAY_DC.</p> <p>Example</p>
DR_RELEASE_OVERLAY_DC	<p>Direct3D only: Releases the device context (DC) handle.</p> <ul style="list-style-type: none"> ▣ More details <p>Using DR_RELEASE_OVERLAY_DC, you can release the DC handle and update the overlay data.</p> <p>Example</p>
DR_GET_MAX_OVERLAY_SIZE	<p>Returns the width x and height y of the maximum overlay area supported by the graphics card.</p> <p>Example</p>
DR_SET_OVERLAY_SIZE	<p>Defines the size of the overlay area (default: current camera image size). Example</p>
DR_GET_OVERLAY_SIZE	<p>Returns the size of the overlay area. (Sample: see DR_SET_OVERLAY_SIZE)</p>
DR_SET_OVERLAY_POSITION	<p>Defines the position of the overlay area. Example</p>
DR_GET_OVERLAY_KEY_COLOR	<p>Returns the RGB values of the current key color (default: black). Example</p>
DR_SET_OVERLAY_KEY_COLOR	<p>Defines the RGB values of the key color.</p> <ul style="list-style-type: none"> ▣ More details <p>The key color specifies where the camera image will be visible in the overlay area. For example: if you fill the complete overlay with the key color, the whole camera image will be visible. If you fill part of the overlay with a different color, the camera image will be covered by the overlay in those places.</p> <p>The key color has no effect in semi-transparent mode!</p>

hCam	Camera handle
	Example
DR_SHOW_OVERLAY	Enables overlay display on top of the current camera image. Example
DR_HIDE_OVERLAY	Disables overlay display. Example
DR_ENABLE_SCALING	Enables real-time scaling of the image to the size of the display window. The overlay is scaled together with the camera image. Example
DR_ENABLE_IMAGE_SCALING	Direct3D only: Enables real-time scaling of the image to the size of the display window. The overlay is not scaled. (Sample: see DR_ENABLE_SCALING)
DR_DISABLE_SCALING	Disables real-time scaling. Example
DR_ENABLE_SEMI_TRANSPARENT_OVERLAY	<p>Enables a semi-transparent display of the overlay area.</p> <ul style="list-style-type: none"> ☐ More details <p>In semi-transparent mode, the values of the camera image and the overlay data are added up for each pixel. Since black has the value 0, the complete camera image will be visible if the overlay is black; if the overlay is white, only the overlay will be visible. With all other colors, the camera image will be visible with the overlay superimposed.</p> <p>The key color has no effect in semi-transparent mode!</p> <p>Example</p>
DR_DISABLE_SEMI_TRANSPARENT_OVERLAY	Disables the semi-transparent display of the overlay area. Example
DR_SET_VSYNC_AUTO	Enables synchronization of the image display with the monitor's image rendering. The image is displayed upon the monitor's next VSYNC signal. Example
DR_SET_VSYNC_OFF	Disables image display synchronization. The image is displayed immediately. Example
DR_SET_USER_SYNC	<p>Direct3D only: Enables synchronization of the image display with a monitor pixel row specified by the user.</p> <ul style="list-style-type: none"> ☐ More details <p>When displaying very large camera images, the auto VSYNC function might not always optimally synchronize image rendering. In this case, you can eliminate flicker by manually setting a suitable position for synchronization. The position needs to be determined individually, based on the camera type and the graphics card.</p> <p>Example</p>
DR_GET_USER_SYNC_POSITION_RANGE	Direct3D only: Returns the minimum and

hCam	Camera handle
	maximum row position for <code>DR_SET_USER_SYNC</code> . Example
<code>DR_LOAD_OVERLAY_FROM_FILE</code>	Direct3D only: Loads a bitmap image (*.BMP file) into the overlay area. If the bitmap image is larger than the overlay area, the bitmap image is clipped. Example
<code>DR_CLEAR_OVERLAY</code>	Deletes the data of the overlay area by filling it with black color. Example
<code>DR_STEAL_NEXT_FRAME</code>	Copies the next image to the active user memory (Steal function). <ul style="list-style-type: none"> ▣ More details Using the <code>pParam</code> parameter, you specify when the function should return: <ul style="list-style-type: none"> • <code>IS_WAIT</code>: The function waits until the image save is complete. • <code>IS_DONT_WAIT</code>: The function returns immediately. Example
<code>DR_SET_STEAL_FORMAT</code>	Defines the color format for the Steal function. <ul style="list-style-type: none"> ▣ More details For a list of all available color formats, see the function description for is_SetColorMode() . The default is <code>IS_CM_BGRA8_PACKED</code> (RGB 32). Example
<code>DR_GET_STEAL_FORMAT</code>	Returns the color format setting for the Steal function. Example
<code>DR_SET_HWND</code>	Sets a new window handle for image output in Direct3D. Example
<code>DR_CHECK_COMPATIBILITY</code>	Returns whether the graphics card supports the uc480 Direct3D functions. Example
<code>DR_GET_OVERLAY_DATA</code>	OpenGL only: Returns a pointer to the overlay.
<code>DR_UPDATE_OVERLAY_DATA</code>	OpenGL only: Updates the overlay.
<code>pParam</code>	Pointer to a function parameter, whose function depends on <code>nCommand</code> .
<code>nSize</code>	Size (in bytes) of the data object or array.

Return values

When used with DR_CHECK_COMPATIBILITY	IS_DR_DEVICE_CAPS_INSUFFICIENT The graphics hardware does not fully support the uc480 Direct3D functions.
IS_DR_CANNOT_CREATE_SURFACE	The image surface or overlay surface could not be created.
IS_DR_CANNOT_CREATE_TEXTURE	The texture could not be created.
IS_DR_CANNOT_CREATE_VERTEX_BUFFER	The vertex buffer could not be created.
IS_DR_CANNOT_GET_OVERLAY_DC	Could not get the device context handle for the overlay.
IS_DR_CANNOT_LOCK_OVERLAY_SURFACE	The overlay surface could not be locked.
IS_DR_CANNOT_RELEASE_OVERLAY_DC	Could not release the device context handle for the overlay.
IS_DR_CANNOT_UNLOCK_OVERLAY_SURFACE	The overlay surface could not be unlocked.
IS_DR_DEVICE_CAPS_INSUFFICIENT	Function is not supported by the graphics hardware.
IS_DR_DEVICE_OUT_OF_MEMORY	Not enough graphics memory available.
IS_DR_NOT_ALLOWED_WHILE_DC_IS_ACTIVE	A device context handle is still open in the application.
IS_INVALID_CAMERA_HANDLE	Invalid camera handle
IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
IS_NO_SUCCESS	General error message
IS_NOT_SUPPORTED	The camera model used here does not support this function or setting.
IS_SUCCESS	Function executed successfully
IS_TIMED_OUT	A timeout occurred. An image capturing process could not be terminated within the allowable period.

Related functions

- [is_SetDisplayMode\(\)](#)
- [is_SetColorMode\(\)](#)
- [is_SetImageMem\(\)](#)
- [is_RenderBitmap\(\)](#)

Example Supported function

```

UINT nType = IS_SET_DM_DIRECT3D;
if (is_DirectRenderer(m_pMainView->GetCameraHandle(), DR_GET_SUPPORTED,
                    (void*)&nType, sizeof(nType)) == IS_SUCCESS)
{
    // Direct3D is supported
}
nType = IS_SET_DM_OPENGL;
if (is_DirectRenderer(m_pMainView->GetCameraHandle(), DR_GET_SUPPORTED,
                    (void*)&nType, sizeof(nType)) == IS_SUCCESS)
{
    // OpenGL is supported

```

```
}
```

Example DC handle

```
//-----  
//          DC-Handle  
//-----  
  
// Get DC handle for Overlay  
HDC hDC;  
is_DirectRenderer (hCam, DR_GET_OVERLAY_DC, (void*)&hDC, sizeof (hDC));  
  
// Release DC handle  
is_DirectRenderer (hCam, DR_RELEASE_OVERLAY_DC, NULL, NULL);
```


Example overlay size and position

```
//-----
//          Size of overlay
//-----

// Query maximum size of overlay area
UINT OverlaySize[2];
is_DirectRenderer (hCam, DR_GET_MAX_OVERLAY_SIZE,
                  (void*)OverlaySize, sizeof(OverlaySize));

INT nWidth = OverlaySize[0];
INT nHeight = OverlaySize[1];

// Set size of overlay area
UINT Size[2];
Size[0] = 100;
Size[1] = 120;
is_DirectRenderer (hCam, DR_SET_OVERLAY_SIZE,
                  (void*)Size, sizeof (Size));

// Set position of overlay area
UINT Position[2];
Position[0] = 20;
Position[1] = 0;
is_DirectRenderer (hCam, DR_SET_OVERLAY_POSITION,
                  void*)Position, sizeof (Position));
```

Example key color

```
//-----
//          Key color
//-----

// Get current key color
UINT OverlayKeyColor[3];
is_DirectRenderer (hCam, DR_GET_OVERLAY_KEY_COLOR,
                  (void*)OverlayKeyColor, sizeof(OverlayKeyColor));

INT nRed = OverlayKeyColor[0];
INT nGreen = OverlayKeyColor[1];
INT nBlue = OverlayKeyColor[2];

// Set new key color
OverlayKeyColor[0] = GetRValue(m_rgbKeyColor);
OverlayKeyColor[1] = GetGValue(m_rgbKeyColor);
OverlayKeyColor[2] = GetBValue(m_rgbKeyColor);
is_DirectRenderer (hCam, DR_SET_OVERLAY_KEY_COLOR,
                  (void*)OverlayKeyColor, sizeof(OverlayKeyColor));
```

Example display

```
//-----
//          Display
//-----

// Show overlay
is_DirectRenderer (hCam, DR_SHOW_OVERLAY, NULL, NULL);

// Hide overlay
is_DirectRenderer (hCam, DR_HIDE_OVERLAY, NULL, NULL);
```

Example scaling

```
//-----
//          Scaling
//-----

// Enable scaling
is_DirectRenderer (hCam, DR_ENABLE_SCALING, NULL, NULL);

// Disable scaling
is_DirectRenderer (hCam, DR_DISABLE_SCALING, NULL, NULL);
```

Example transparency

```
//-----  
//          Transparency  
//-----  
  
// Enable semi-transparent overlay  
is_DirectRendererer (hCam, DR_ENABLE_SEMI_TRANSPARENT_OVERLAY, NULL, NULL);  
  
// Disable semi-transparent overlay  
is_DirectRendererer (hCam, DR_DISABLE_SEMI_TRANSPARENT_OVERLAY, NULL, NULL);
```

Example synchronization

```
//-----  
//          Synchronization  
//-----  
  
// Enable auto-synchronization  
is_DirectRendererer (hCam, DR_SET_VSYNC_AUTO, NULL, NULL);  
  
// User defined synchronization: Query range and set position  
UINT UserSync[2];  
is_DirectRendererer (hCam, DR_GET_USER_SYNC_POSITION_RANGE,  
                    (void*)UserSync, sizeof (UserSync));  
INT Min = UserSync[0];  
INT Max = UserSync[1];  
INT SyncPosition = 400;  
is_DirectRendererer (hCam, DR_SET_USER_SYNC,  
                    (void*)&SyncPosition, sizeof (SyncPosition));  
  
// Disable synchronization  
is_DirectRendererer (hCam, DR_SET_VSYNC_OFF, NULL, NULL);
```

Example overlay with BMP

```
//-----  
//          BMP file  
//-----  
  
// Load overlay from BMP file  
is_DirectRendererer (hCam, DR_LOAD_OVERLAY_FROM_FILE,  
                    (void*)"c:\test.bmp", NULL);  
  
//-----  
//          Delete overlay  
//-----  
  
// Delete overlay area  
is_DirectRendererer (hCam, DR_CLEAR_OVERLAY, NULL, NULL);
```

Example steal mode

```
//-----  
//          Steal mode  
//-----  
  
// Get and set color mode for image to be copied  
INT nColorMode;  
is_DirectRendererer (hCam, DR_GET_STEAL_FORMAT,  
                    (void*)&nColorMode, sizeof (nColorMode));  
  
nColorMode = IS_CM_MONO8;  
is_DirectRendererer (hCam, DR_SET_STEAL_FORMAT,  
                    (void*)&nColorMode, sizeof (nColorMode));  
  
// Copy image with function returning immediately  
INT nwait = IS_DONT_WAIT;  
is_DirectRendererer(hCam, DR_STEAL_NEXT_FRAME,  
                    (void*)&nwait, sizeof (nwait));
```

Example window handle

```
//-----  
//          Handle to window  
//-----  
  
// Set new window handle for image display  
is_DirectRenderer (hCam, DR_SET_HWND,  
                  (void*)&hWnd, sizeof (hWnd));
```

Example compatibility

```
//-----  
//          Compatibility  
//-----  
  
// Check graphics card compatibility  
INT nRet = is_DirectRenderer (hCam, DR_CHECK_COMPATIBILITY, NULL, NULL);  
  
if (nRet == IS_DR_DEVICE_CAPS_INSUFFICIENT )  
// Graphics card does not support Direct3D
```

Example OpenGL under Linux

```
//OpenGL initialize
OPENGL_DISPLAY display;
display.pDisplay = NULL;
display.nWindowID = 0 /* window id */
```

```
is_InitCamera(&hCam, (void*)&display);
```

Example under Linux (with usage of the Cairo library)

```
UINT Size[2] = { 480, 480 };
is_DirectRenderer (hCam, DR_SET_OVERLAY_SIZE, (void*)Size, sizeof (Size));

char *pOverlayBuffer;
is_DirectRenderer(hCam, DR_GET_OVERLAY_DATA, (void*)&pOverlayBuffer, sizeof(pOverlayBuffer));



cairo_surface_t *surface = 0;
cairo_t *cr = 0;
int w, h;
w = Size[0];
h = Size[1];
surface = cairo_image_surface_create_for_data(buffer, CAIRO_FORMAT_ARGB32, w, h, w * 4);
cr = cairo_create(surface);
cairo_set_line_width (cr, 6);
cairo_rectangle (cr, 12, 12, 232, 70);
cairo_new_sub_path (cr); cairo_arc (cr, 64, 64, 40, 0, 2* 3.14);
cairo_new_sub_path (cr); cairo_arc_negative (cr, 192, 64, 40, 0, -2*3.14);
cairo_set_fill_rule (cr, CAIRO_FILL_RULE_EVEN_ODD);
cairo_set_source_rgb (cr, 0, 0.7, 0); cairo_fill_preserve (cr);
cairo_set_source_rgb (cr, 0, 0, 0); cairo_stroke (cr);
cairo_translate (cr, 0, 128);
cairo_rectangle (cr, 12, 12, 232, 70);
cairo_new_sub_path (cr); cairo_arc (cr, 64, 64, 40, 0, 2*3.14);
cairo_new_sub_path (cr); cairo_arc_negative (cr, 192, 64, 40, 0, -2*3.14);
cairo_set_fill_rule (cr, CAIRO_FILL_RULE_WINDING);
cairo_set_source_rgb (cr, 0, 0, 0.9); cairo_fill_preserve (cr);
cairo_set_source_rgb (cr, 0, 0, 0); cairo_stroke (cr);
cairo_select_font_face (cr, "Sans", CAIRO_FONT_SLANT_NORMAL,
                        CAIRO_FONT_WEIGHT_BOLD );
cairo_set_font_size (cr, 90.0);
cairo_move_to (cr, 10.0, 135.0);
cairo_show_text (cr, "Hello");
cairo_move_to (cr, 70.0, 165.0);
cairo_set_font_size (cr, 150.0);
cairo_text_path (cr, "uc480");
cairo_set_source_rgb (cr, 0.5, 0.5, 1);
cairo_fill_preserve (cr);
cairo_set_source_rgb (cr, 0, 0, 0);
cairo_set_line_width (cr, 2.56);
cairo_stroke (cr);
/* draw helping lines */
cairo_set_source_rgba (cr, 1, 0.2, 0.2, 0.6);
cairo_arc (cr, 10.0, 135.0, 5.12, 0, 2*3.14);
cairo_close_path (cr);
cairo_arc (cr, 70.0, 165.0, 5.12, 0, 2*3.14);
cairo_fill (cr);
cairo_destroy (cr);
cairo_surface_destroy (surface);

// update overlay
is_DirectRenderer(hCam, DR_UPDATE_OVERLAY_DATA, NULL, 0);
```

Sample programs

- uc480DirectRenderer
- uc480Steal

4.3.18 is_DisableEvent

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_DisableEvent (HCAM hCam, INT which)
```

Description

Using `is_DisableEvent()`, you disable the event indicated here. The event (e.g. image capture completed) will usually still occur, but will no longer trigger an event signal. Disabled events are no longer signaled to the application. You can re-enable the desired event using [is_EnableEvent\(\)](#). See also [is_InitEvent\(\)](#).

Input parameters

hCam	Camera handle
which	ID of the event to be disabled. See also is_EnableEvent() .



Return values

IS_INVALID_CAMERA_HANDLE	Invalid camera handle
IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
IS_NO_SUCCESS	General error message
IS_SUCCESS	Function executed successfully

Related functions

- [is_EnableEvent\(\)](#)
- Windows only: [is_InitEvent\(\)](#)
- Windows only: [is_ExitEvent\(\)](#)
- Linux only: [is_WaitEvent\(\)](#)

4.3.19 is_EdgeEnhancement

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_EdgeEnhancement(HCAM hCam, UINT nCommand, void* pParam, UINT cbSizeOfParam)
```

Description

This function enables/disables a software edge filter.

Due to Bayer format color conversion, the original edges of a color image may easily become blurred. By enabling the digital edge filter, you can optimize edge representation. This function causes a higher CPU load.

The `nCommand` input parameter is used to select the function mode. The `pParam` input parameter depends on the selected function mode. If you select functions for setting or returning a value, `pParam` contains a pointer to a variable of the `UINT` type. The size of the memory area to which `pParam` refers is specified in the `cbSizeOfParam` input parameter.

Note

The following function is obsolete by the `is_EdgeEnhancement()` function:

- `is_SetEdgeEnhancement()`

See also [Obsolete functions](#)

Input parameters

<code>hCam</code>	Camera handle
<div style="border: 1px solid gray; padding: 2px;"> <input type="checkbox"/> <code>nCommand</code> </div>	
<code>IS_EDGE_ENHANCEMENT_CMD_GET_RANGE</code>	Returns the range of the edge enhancement (Example 1)
<code>IS_EDGE_ENHANCEMENT_CMD_GET_DEFAULT</code>	Returns the standard value of the edge enhancement (Example 2)
<code>IS_EDGE_ENHANCEMENT_CMD_GET</code>	Returns the current set edge enhancement (Example 3)
<code>IS_EDGE_ENHANCEMENT_CMD_SET</code>	Sets the edge enhancement (Example 4) 0: no edge enhancement
<code>pParam</code>	Pointer to a function parameter, whose function depends on <code>nCommand</code> .
<code>cbSizeOfParam</code>	Size (in bytes) of the memory area to which <code>pParam</code> refers.

Return values

<code>IS_INVALID_PARAMETER</code>	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
<code>IS_NO_SUCCESS</code>	General error message
<code>IS_SUCCESS</code>	Function executed successfully

Related functions

- [is_SetColorMode\(\)](#)
- [is_SetColorConverter\(\)](#)

Example 1

```
UINT nRange[3];
ZeroMemory(nRange, sizeof(nRange));

INT nRet = is_EdgeEnhancement(m_hCam,
                             IS_EDGE_ENHANCEMENT_CMD_GET_RANGE,
                             (void*)nRange,
                             sizeof(nRange)
                             );

if (nRet == IS_SUCCESS)
{
    UINT nEdgeEnhancementMin = nRange[0];
    UINT nEdgeEnhancementMax = nRange[1];
    UINT nEdgeEnhancementInc = nRange[2];
}
```

Example 2

```
UINT nDefault;
INT nRet = is_EdgeEnhancement(m_hCam,
                             IS_EDGE_ENHANCEMENT_CMD_GET_DEFAULT,
                             (void*)&nDefault,
                             sizeof(nDefault)
                             );
```



Example 3

```
UINT nEdgeEnhancement;
INT nRet = is_EdgeEnhancement(m_hCam,
                             IS_EDGE_ENHANCEMENT_CMD_GET,
                             (void*)&nEdgeEnhancement,
                             sizeof(nEdgeEnhancement)
                             );
```

Example 4

```
UINT nEdgeEnhancement = 4;
INT nRet = is_EdgeEnhancement(m_hCam,
                             IS_EDGE_ENHANCEMENT_CMD_SET,
                             (void*)&nEdgeEnhancement,
                             sizeof(nEdgeEnhancement)
                             );
```


4.3.20 is_EnableAutoExit

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_EnableAutoExit (HCAM hCam, INT nMode)
```

Description

`is_EnableAutoExit()` enables automatic closing of the camera handle after a camera has been removed on-the-fly. Upon closing of the handle, the entire memory allocated by the driver will be released.

Input parameters

hCam	Camera handle
<input type="checkbox"/> nMode	
IS_ENABLE_AUTO_EXIT	Enables automatic closing
IS_DISABLE_AUTO_EXIT	Disables automatic closing
IS_GET_AUTO_EXIT_ENABLED	Returns the current setting



Return values

Current setting when used together with <code>IS_GET_AUTO_EXIT_ENABLED</code>	
IS_INVALID_CAMERA_HANDLE	Invalid camera handle
IS_NO_SUCCESS	General error message
IS_SUCCESS	Function executed successfully

Related functions

- [is_ExitCamera\(\)](#)

4.3.21 is_EnableEvent

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

INT is_EnableEvent (HCAM hCam, INT which)

Description

Using `is_EnableEvent()`, you release an event object. Following the release, the event messages for the created event object are enabled. Depending on the operating system different functions are to call.

Windows	<ul style="list-style-type: none"> • Event has to be provided by the application program • Event has to be declared by is_InitEvent() • Event has to be activated by <code>is_EnableEvent()</code> • You have to wait for the event in the application program by <code>WaitForSingleObject</code> Or <code>WaitForMultipleObject</code> • Event has to be deactivated by is_DisableEvent() • Event has to be logged off by is_ExitEvent()
Linux	<ul style="list-style-type: none"> • Event has to be provided by the uc480 API • Event has to be activated by <code>is_EnableEvent()</code> • You have to wait for the event by is_WaitEvent() • Event has to be deactivated by is_DisableEvent()

Input parameters

hCam	Camera handle
which: ID of the event to be released.	
IS_SET_EVENT_AUTOBRIGHTNESS_FINISHED	The automatic brightness control in the run-once mode is completed.
IS_SET_EVENT_CAMERA_MEMORY	In the camera memory mode an image acquisition iteration is finished.
IS_SET_EVENT_CAPTURE_STATUS	There is an information about image capturing available. This information can be requested by is_CaptureStatus() . Note that this event replaces the former <code>IS_SET_EVENT_TRANSFER_FAILED</code> from previous versions.
IS_SET_EVENT_CONNECTIONSPEED_CHANGED	The connection speed of a USB 3 DCx camera changed from USB 2.0 to USB 3.0 or from USB 3.0 to USB 2.0.
IS_SET_EVENT_DEVICE_RECONNECTED	A camera initialized with is_InitCamera() and disconnected afterwards was reconnected.
IS_SET_EVENT_EXTTRIG	An image which was captured following the arrival of a trigger has been transferred completely.

hCam	Camera handle
	This is the earliest possible moment for a new capturing process. The image must then be post-processed by the driver and will be available after the <code>IS_FRAME</code> processing event.
<code>IS_SET_EVENT_FRAME</code>	A new image is available.
<code>IS_SET_EVENT_NEW_DEVICE</code>	A new camera was connected. This is independent of the device handle (hCam is ignored).
<code>IS_SET_EVENT_OVERLAY_DATA_LOST</code>	Direct3D/OpenGL mode: Because of a re-programming the parameters of the overlay are invalid. The overlay must be draw new.
<code>IS_SET_EVENT_REMOVAL</code>	A camera was removed. This is independent of the device handle (hCam is ignored).
<code>IS_SET_EVENT_REMOVE</code>	A camera initialized with is_InitCamera() was disconnected.
<code>IS_SET_EVENT_SEQ</code>	The sequence is completed.
<code>IS_SET_EVENT_STATUS_CHANGED</code>	Linux only: The availability of a camera has changed, e.g. an available camera was opened.
<code>IS_SET_EVENT_STEAL</code>	An image extracted from the overlay is available.
<code>IS_SET_EVENT_WB_FINISHED</code>	The automatic white balance control is completed.

Return values

<code>IS_INVALID_CAMERA_HANDLE</code>	Invalid camera handle
<code>IS_INVALID_PARAMETER</code>	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
<code>IS_NO_SUCCESS</code>	General error message
<code>IS_SUCCESS</code>	Function executed successfully

Related functions

- [is_DisableEvent\(\)](#)
- Windows only: [is_InitEvent\(\)](#)
- Windows only: [is_ExitEvent\(\)](#)
- Linux only: [is_WaitEvent\(\)](#)

Example Windows

```
HANDLE hEvent = CreateEvent(NULL, FALSE, FALSE, NULL);
is_InitEvent(hCam, hEvent, IS_SET_EVENT_FRAME);
is_EnableEvent(hCam, IS_SET_EVENT_FRAME);
is_FreezeVideo(hCam, IS_DONT_WAIT);
DWORD dwRet = WaitForSingleObject(hEvent, 1000);
if (dwRet == WAIT_TIMEOUT)
{
    /* wait timed out */
}
```

```
}  
else if (dwRet == WAIT_OBJECT_0)  
{  
    /* event signalled */  
}  
is_DisableEvent(hCam, IS_SET_EVENT_FRAME);  
is_ExitEvent(hCam, IS_SET_EVENT_FRAME);  
CloseHandle(hEvent);
```



Example Linux

```
is_EnableEvent(hCam, IS_SET_EVENT_FRAME);  
is_FreezeVideo(hCam, IS_DONT_WAIT);  
INT nRet = is_WaitEvent(hCam, IS_SET_EVENT_FRAME, 1000);  
if (nRet == IS_TIMED_OUT)  
{  
    /* wait timed out */  
}  
else if (nRet == IS_SUCCESS)  
{  
    /* event signalled */  
}  
is_DisableEvent(hCam, IS_SET_EVENT_FRAME);
```

Sample programs

- SimpleLive (C++)
- uc480Event (C++)

4.3.22 is_EnableMessage

	
USB 2.0 USB 3.0	-

Syntax

```
INT is_EnableMessage (HCAM hCam, INT which, HWND hWnd)
```

Description

Using `is_EnableMessage()`, you can enable Windows messages. If a particular event occurs, the messages are sent to the application.

Each message is structured as follows:

- **Message:** `IS_uc480_MESSAGE`
- **wParam:** Event (see table)
- **lParam:** DCx camera handle associated with the message

Attention

You have to deactivate Windows messages with `hWnd == NULL` before you free the uc480 API library. Otherwise the application may not close properly.

Input parameters

hCam	Camera handle
<input type="checkbox"/> which: ID of the message to be enabled/disabled	
IS_FRAME	A new image is available.
IS_SEQUENCE	The sequence is completed.
IS_CAPTURE_STATUS	An error occurred during the data transfer, see is_CaptureStatus() . The parameter <code>IS_CAPTURE_STATUS</code> replaces the previous parameter <code>IS_TRANSFER_FAILED</code> . The parameter <code>IS_TRANSFER_FAILED</code> was moved into the new header file <code>uc480_deprecated.h</code> , which contains all obsolete function definitions and constants. If necessary the header file <code>uc480_deprecated.h</code> can be included in addition to the header file <code>uc480.h</code> .
IS_TRIGGER	An image which was captured following the arrival of a trigger has been transferred completely. This is the earliest possible moment for a new capturing process. The image must then be post-processed by the driver and is available after the <code>IS_FRAME</code> message has occurred.
IS_DEVICE_REMOVED	A camera initialized with is_InitCamera() was disconnected.
IS_DEVICE_RECONNECTED	A camera initialized with is_InitCamera() and disconnected afterwards was reconnected.

hCam	Camera handle
IS_NEW_DEVICE	A new camera was connected.
IS_DEVICE_REMOVAL	A camera was removed.
IS_WB_FINISHED	Automatic white balance control is completed (only if this control was started using the <code>IS_SET_AUTO_WB_ONCE</code> function).
IS_AUTOBRIGHTNESS_FINISHED	Automatic brightness control is completed (only if this control was started using the <code>IS_SET_AUTO_BRIGHTNESS_ONCE</code> function).
IS_CAMERA_MEMORY	In the camera memory mode an image acquisition iteration is finished.
IS_CONNECTIONSPEED_CHANGED	The connection speed of a USB 3 DCx camera changed from USB 2.0 to USB 3.0 or from USB 3.0 to USB 2.0.
hWnd	Application window for receiving the message. <code>NULL</code> disables the message designated by the which parameter.



Return values

IS_INVALID_CAMERA_HANDLE	Invalid camera handle
IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
IS_NO_SUCCESS	General error message
IS_SUCCESS	Function executed successfully

Related functions

- [is_InitEvent\(\)](#)

4.3.23 is_ExitCamera

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_ExitCamera (HCAM hCam)
```

Description

`is_ExitCamera()` disables the `hCam` camera handle and releases the data structures and memory areas taken up by the DCx camera. Image memory allocated using the [is_AllocImageMem\(\)](#) function which has not been released yet is automatically released.

Note

We recommend that you call the following functions only from a single thread in order to avoid unpredictable behaviour of the application.

- [is_InitCamera\(\)](#)
- [is_SetDisplayMode\(\)](#)
- [is_ExitCamera\(\)](#)

See also [Programming: Thread programming](#)

Input parameters

hCam	Camera handle
------	---------------



Return values

IS_CANT_OPEN_REGISTRY	Error opening a Windows registry key
IS_CANT_READ_REGISTRY	Error reading settings from the Windows registry
IS_ERROR_CPU_IDLE_STATES_CONFIGURATION	The configuration of the CPU idle has failed.
IS_INVALID_CAMERA_HANDLE	Invalid camera handle
IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
IS_IO_REQUEST_FAILED	An IO request from the uc480 driver failed. Possibly the versions of the <code>uc480.dll</code> (API) and the driver file (<code>uc480_usb.sys</code>) do not match.
IS_NO_IMAGE_MEM_ALLOCATED	The driver could not allocate memory.
IS_NO_SUCCESS	General error message
IS_SUCCESS	Function executed successfully

Related functions

- [is_InitCamera\(\)](#)
- [is_EnableAutoExit\(\)](#)

4.3.24 is_ExitEvent

	
USB 2.0 USB 3.0	-

Syntax

```
INT is_ExitEvent (HCAM hCam, INT which)
```

Description

`is_ExitEvent()` deletes an existing event object. After an event has been deleted, you can no longer enable it by calling the [is_EnableEvent\(\)](#) function.

Input parameters

hCam	Camera handle
which	ID of the event to be deleted. See also is_EnbaleEvent() .

Return values

IS_INVALID_CAMERA_HANDLE	Invalid camera handle
IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
IS_NO_SUCCESS	General error message
IS_SUCCESS	Function executed successfully



Related functions

- [is_EnableEvent\(\)](#)
- [is_InitEvent\(\)](#)

Example

See also [is_ForceTrigger\(\)](#)

4.3.25 is_ExitImageQueue

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_ExitImageQueue (HCAM hCam)
```

Description

`is_ExitImageQueue()` deletes a queue which has been initialized with [is_InitImageQueue\(\)](#) and discards all information about the order of queued images. The image memories will be unlocked. The memory sequence itself persists and can be deleted with [is_ClearSequence\(\)](#).

Input parameters

hCam	Camera handle
------	---------------



Return values

IS_INVALID_CAMERA_HANDLE	Invalid camera handle
IS_NO_SUCCESS	General error message
IS_SUCCESS	Function executed successfully

Related functions

- [is_InitImageQueue\(\)](#)
- [is_WaitForNextImage\(\)](#)

4.3.26 is_Exposure

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_Exposure (HCAM hCam, UINT nCommand, void* pParam, UINT cbSizeOfParam)
```

Description

Using `is_Exposure()` you can query the exposure time ranges available in your camera, and set new exposure times:

- Setting the exposure time
- Exposure time with fine increments
- Setting the long exposure
- Setting the dual exposure

Note on dependencies on other settings

The use of the following functions will affect the exposure time:

- [is_PixelClock\(\)](#)
- [is_SetOptimalCameraTiming\(\)](#)
- [is_SetFrameRate\(\)](#)
- [is_AOI\(\)](#) (if the image size is changed)
- [is_SetSubSampling\(\)](#)
- [is_SetBinning\(\)](#)

Changes made to the image size, the frame rate or the pixel clock frequency also affect the exposure time. For this reason, you need to call `is_Exposure()` again after such changes.

Note on new driver versions

Newer driver versions sometimes allow an extended value range for the exposure time setting. We recommend querying the value range every time and set the exposure time explicitly.

The `nCommand` input parameter is used to select the function mode. The `pParam` input parameter depends on the selected function mode. If you select functions for setting or returning a value, `pParam` contains a pointer to a variable of the `UINT` type. The size of the memory area to which `pParam` refers is specified in the `cbSizeOfParam` input parameter.


Note on older uc480 exposure time functions

The following uc480 API commands are obsolete by the `is_Exposure()` function:

- `is_GetExposureRange()`
- `is_SetExposureTime()`

See also [Obsolete functions](#)

Input parameters

hCam	Camera handle
<div style="border: 1px solid gray; padding: 5px;">  nCommand </div>	

IS_EXPOSURE_CMD_GET_CAPS	Returns the supported function modes. <ul style="list-style-type: none"> ▣ More details • pParam: Pointer to bit mask of type <code>UINT</code> In the bit mask, the status flags from EXPOSURE_CAPS are returned. • nSizeOfParam: 4 Example
pParam	Pointer to a function parameter, whose function depends on nCommand.
cbSizeOfParam	Size (in bytes) of the memory area to which pParam refers.

Contents of the EXPOSURE_CAPS structure

INT	IS_EXPOSURE_CAP_EXPOSURE	The exposure time setting is supported
INT	IS_EXPOSURE_CAP_FINE_INCREMENT	Fine exposure time increments are supported
INT	IS_EXPOSURE_CAP_LONG_EXPOSURE	Long time exposure is supported. Depending on the camera model long time exposure is only supported in trigger mode but not in freerun mode.
INT	IS_EXPOSURE_CAP_DUAL_EXPOSURE	/ and / only: The sensor supports dual exposure. Odd and even lines can be exposed with different exposure times.

Return values

IS_CANT_OPEN_DEVICE	An attempt to initialize or select the camera failed (no camera connected or initialization error).
IS_CANT_COMMUNICATE_WITH_DRIVER	Communication with the driver failed because no driver has been loaded.
IS_INVALID_CAMERA_TYPE	The camera type defined in the .ini file does not match the current camera model.
IS_INVALID_CAMERA_HANDLE	Invalid camera handle
IS_INVALID_MODE	Camera is in standby mode, function not allowed
IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
IS_IO_REQUEST_FAILED	An IO request from the uc480 driver failed. Possibly the versions of the ueye_api.dll (API) and the driver file (ueye_usb.sys or ueye_eth.sys) do not match.
IS_NO_SUCCESS	General error message
IS_NOT_CALIBRATED	The camera does not contain any calibration data.
IS_NOT_SUPPORTED	The camera model used here does not support this function or setting.
IS_SUCCESS	Function executed successfully
IS_TIMED_OUT	A timeout occurred. An image capturing process could not be terminated within the allowable period.

Related functions

- [is_SetFrameRate\(\)](#)
- [is_PixelClock\(\)](#)
- [is_SetOptimalCameraTiming\(\)](#)
- [is_SetAutoParameter\(\)](#)
- [is_AutoParameter\(\)](#)
- [is_SetHardwareGain\(\)](#)

Example



```

UINT nCaps = 0;
INT nRet = is_Exposure(m_hCam, IS_EXPOSURE_CMD_GET_CAPS, (void*)&nCaps, sizeof(nCaps));

if (nRet == IS_SUCCESS)
{
    if (nCaps & IS_EXPOSURE_CAP_FINE_INCREMENT)
    {
        // Fine increment supported
    }
}

```

4.3.26.1 Setting the Exposure Time

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_Exposure (HCAM hCam, UINT nCommand, void* pParam, UINT cbSizeOfParam)
```

Description

Using `is_Exposure()` you can query the exposure time ranges available in your camera, and set new exposure times.

Note on dependencies on other settings

The use of the following functions will affect the exposure time:

- [is_PixelClock\(\)](#)
- [is_SetOptimalCameraTiming\(\)](#)
- [is_SetFrameRate\(\)](#)
- [is_AOI\(\)](#) (if the image size is changed)
- [is_SetSubSampling\(\)](#)
- [is_SetBinning\(\)](#)

Changes made to the image size, the frame rate or the pixel clock frequency also affect the exposure time. For this reason, you need to call `is_Exposure()` again after such changes.

Note on new driver versions

Newer driver versions sometimes allow an extended value range for the exposure time setting. We recommend querying the value range every time and set the exposure time explicitly.

Applying new settings

In freerun mode ([is_CaptureVideo\(\)](#)), any modification of the exposure time will only become effective when image after next is captured. In trigger mode ([is_SetExternalTrigger\(\)](#)), the

modification will be applied to the next image. See also the [Applying new parameters](#) chapter.



In general, the pixel clock is set once when opening the camera and will not be changed. Note that, if you change the pixel clock, the setting ranges for frame rate and exposure time also changes. If you change a parameter, the following order is recommended:

1. Change pixel clock.
2. Query frame rate range and, if applicable, set new value.
3. Query exposure time range and, if applicable, set new value.

If one parameter is changed, the following parameters have to be adjusted due to the dependencies.

Accuracy of the exposure time setting

The increments for setting the exposure time (`IS_EXPOSURE_CMD_GET_EXPOSURE_RANGE_INC`) depend on the sensor's current timing settings (pixel clock, frame rate). The smallest increment usually corresponds to the duration of one pixel row, which is the time it takes the sensor to read out one pixel row.

You can query the actual exposure time setting with the `IS_EXPOSURE_CMD_GET_EXPOSURE` parameter.

Some sensors allow setting the exposure time in smaller increments. Using the `IS_EXPOSURE_CMD_GET_CAPS` parameter, you can check whether your sensor supports this function.

For minimum and maximum exposure times as well as other sensor-based dependencies, please refer to the [Camera and sensor data](#) chapter.

Rounding errors from increments

When calculating a new exposure time based on the returned increment, note that calculations with floating point values in the PC will always be subject to rounding errors. Therefore, an addition or subtraction of the `n*INCREMENT` value might not always produce the exact desired result. In this case, the uc480 API rounds down the floating point value and sets the exposure time to the next lower value.

You can avoid this behavior by additionally adding the value `INCREMENT/2.f` (half increment) when calculating with `n*INCREMENT`. This ensures that the desired value will be set even after rounding.

The `nCommand` input parameter is used to select the function mode. The `pParam` input parameter depends on the selected function mode. If you select functions for setting or returning a value, `pParam` contains a pointer to a variable of the `UINT` type. The size of the memory area to which `pParam` refers is specified in the `cbSizeOfParam` input parameter.



Input parameters

hCam	Camera handle
<div style="border: 1px solid gray; padding: 5px;"> nCommand </div>	

IS_EXPOSURE_CMD_GET_EXPOSURE_DEFAULT	<p>Returns the default setting for the exposure time.</p> <ul style="list-style-type: none"> ▣ More details • pParam: Pointer to variable of type <code>double</code> returning the default value. • nSizeOfParam: 8 <p>If the value 0.0 is returned this means that the exposure time is set to maximum value of 1/frame rate.</p>
IS_EXPOSURE_CMD_GET_EXPOSURE	<p>Returns the currently set exposure time (in ms).</p> <ul style="list-style-type: none"> ▣ More details • pParam: Pointer to variable of type <code>double</code> returning the current value. • nSizeOfParam: 8
IS_EXPOSURE_CMD_GET_EXPOSURE_RANGE_MIN	<p>Returns the minimum exposure time.</p> <ul style="list-style-type: none"> ▣ More details • pParam: Pointer to variable of type <code>double</code> returning the minimum value. • nSizeOfParam: 8
IS_EXPOSURE_CMD_GET_EXPOSURE_RANGE_MAX	<p>Returns the maximum exposure time.</p> <ul style="list-style-type: none"> ▣ More details • pParam: Pointer to variable of type <code>double</code> returning the maximum value. • nSizeOfParam: 8
IS_EXPOSURE_CMD_GET_EXPOSURE_RANGE_INC	<p>Returns the exposure time increment.</p> <ul style="list-style-type: none"> ▣ More details • pParam: Pointer to variable of type <code>double</code> returning the increment. • nSizeOfParam: 8
IS_EXPOSURE_CMD_GET_EXPOSURE_RANGE	<p>Returns the exposure time range.</p> <ul style="list-style-type: none"> ▣ More details • pParam: Pointer to array of type <code>double</code> returning the minimum and maximum values and the increment (in exactly this order). • nSizeOfParam: 24
IS_EXPOSURE_CMD_SET_EXPOSURE	<p>Sets the exposure time (in ms).</p> <ul style="list-style-type: none"> ▣ More details • pParam: Pointer to variable of type <code>double</code> that passes the value to be set. After setting the exposure time this value contains the actually set exposure time. Depending on the sensor the set exposure time may vary slightly from the desired exposure time. • nSizeOfParam: 8 <p>If 0 is passed, the exposure time is set to the maximum value of 1/frame rate.</p>

pParam	Pointer to a function parameter, whose function depends on nCommand.
cbSizeOfParam	Size (in bytes) of the memory area to which pParam refers.

4.3.26.2 Exposure Time with Fine Increments

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_Exposure (HCAM hCam, UINT nCommand, void* pParam, UINT cbSizeOfParam)
```

Description

Using `is_Exposure()` you can query the exposure time ranges with fine increment available in your camera, and set new exposure times.

Note on dependencies on other settings

The use of the following functions will affect the exposure time:

- [is_PixelClock\(\)](#)
- [is_SetOptimalCameraTiming\(\)](#)
- [is_SetFrameRate\(\)](#)
- [is_AOI\(\)](#) (if the image size is changed)
- [is_SetSubSampling\(\)](#)
- [is_SetBinning\(\)](#)

Changes made to the image size, the frame rate or the pixel clock frequency also affect the exposure time. For this reason, you need to call `is_Exposure()` again after such changes.

Note on new driver versions

Newer driver versions sometimes allow an extended value range for the exposure time setting. We recommend querying the value range every time and set the exposure time explicitly.

The `nCommand` input parameter is used to select the function mode. The `pParam` input parameter depends on the selected function mode. If you select functions for setting or returning a value, `pParam` contains a pointer to a variable of the `UINT` type. The size of the memory area to which `pParam` refers is specified in the `cbSizeOfParam` input parameter.

Input parameters

hCam	Camera handle
<div style="border: 1px solid gray; padding: 5px;"> <input type="checkbox"/> nCommand </div>	

IS_EXPOSURE_CMD_GET_FINE_INCREMENT_RANGE_MIN	Returns the minimum exposure time in fine increments for some sensors. <input type="checkbox"/> More details <ul style="list-style-type: none"> • pParam: Pointer to variable of type <code>double</code> returning the minimum value. • nSizeOfParam: 8 Example 1
IS_EXPOSURE_CMD_GET_FINE_INCREMENT_RANGE_MAX	Returns the maximum exposure time in fine increments for some sensors. <input type="checkbox"/> More details <ul style="list-style-type: none"> • pParam: Pointer to variable of type <code>double</code> returning the maximum value. • nSizeOfParam: 8
IS_EXPOSURE_CMD_GET_FINE_INCREMENT_RANGE_INC	Returns the exposure time increment in fine increments for some sensors. <input type="checkbox"/> More details <ul style="list-style-type: none"> • pParam: Pointer to variable of type <code>double</code> returning the increment. • nSizeOfParam: 8
IS_EXPOSURE_CMD_GET_FINE_INCREMENT_RANGE	Returns the exposure time range in fine increments for some sensors. <input type="checkbox"/> More details <ul style="list-style-type: none"> • pParam: Pointer to array of type <code>double</code> returning the minimum and maximum values and the increment (in exactly this order). • nSizeOfParam: 24 Example 2
pParam	Pointer to a function parameter, whose function depends on <code>nCommand</code> .
cbSizeOfParam	Size (in bytes) of the memory area to which <code>pParam</code> refers.

Example 1

```
double dblMin, dblMax, dblInc;

INT nRet = is_Exposure(m_hCam, IS_EXPOSURE_CMD_GET_FINE_INCREMENT_RANGE_MIN,
    (void*)&dblMin, sizeof(dblMin));

INT nRet = is_Exposure(m_hCam, IS_EXPOSURE_CMD_GET_FINE_INCREMENT_RANGE_MAX,
    (void*)&dblMax, sizeof(dblMax));

INT nRet = is_Exposure(m_hCam, IS_EXPOSURE_CMD_GET_FINE_INCREMENT_RANGE_INC,
    (void*)&dblInc, sizeof(dblInc));
```

Example 2

```
double dblRange[3];
double dblMin, dblMax, dblInc;

INT nRet = is_Exposure(m_hCam, IS_EXPOSURE_CMD_GET_FINE_INCREMENT_RANGE,
    (void*)&dblRange, sizeof(dblRange));

if (nRet == IS_SUCCESS)
{
```





```

    dblMin = dblRange[0];
    dblMax = dblRange[1];
    dblInc = dblRange[2];
}

```

4.3.26.3 Setting the Long Exposure

	
USB 2.0	USB 2.0
USB 3.0	USB 3.0

Syntax

```
INT is_Exposure (HCAM hCam, UINT nCommand, void* pParam, UINT cbSizeOfParam)
```

Description

Using `is_Exposure()` you can query the long exposure ranges available in your camera, and set new exposure times.

Note on dependencies on other settings

The use of the following functions will affect the exposure time:

- [is_PixelClock\(\)](#)
- [is_SetOptimalCameraTiming\(\)](#)
- [is_SetFrameRate\(\)](#)
- [is_AOI\(\)](#) (if the image size is changed)
- [is_SetSubSampling\(\)](#)
- [is_SetBinning\(\)](#)

Changes made to the image size, the frame rate or the pixel clock frequency also affect the exposure time. For this reason, you need to call `is_Exposure()` again after such changes.

Note on new driver versions

Newer driver versions sometimes allow an extended value range for the exposure time setting. We recommend querying the value range every time and set the exposure time explicitly.



The `nCommand` input parameter is used to select the function mode. The `pParam` input parameter depends on the selected function mode. If you select functions for setting or returning a value, `pParam` contains a pointer to a variable of the `UINT` type. The size of the memory area to which `pParam` refers is specified in the `cbSizeOfParam` input parameter.

Input parameters

hCam	Camera handle
<div style="border: 1px solid gray; padding: 10px;"> <div style="border-bottom: 1px solid gray; margin-bottom: 10px;"> nCommand </div> </div>	

IS_EXPOSURE_CMD_GET_LONG_EXPOSURE_RANGE_MIN	Returns the minimum long exposure time. <input type="checkbox"/> More information <ul style="list-style-type: none"> • pParam: Pointer to variable of type <code>double</code> returning the minimum value. • nSizeOfParam: 8
IS_EXPOSURE_CMD_GET_LONG_EXPOSURE_RANGE_MAX	Returns the maximum long exposure time. <input type="checkbox"/> More information <ul style="list-style-type: none"> • pParam: Pointer to variable of type <code>double</code> returning the maximum value. • nSizeOfParam: 8
IS_EXPOSURE_CMD_GET_LONG_EXPOSURE_RANGE_INC	Returns the increments for long exposure. <input type="checkbox"/> More information <ul style="list-style-type: none"> • pParam: Pointer to variable of the type <code>double</code> returning the increment. • nSizeOfParam: 8
IS_EXPOSURE_CMD_GET_LONG_EXPOSURE_RANGE	Returns the value range for long exposure. <input type="checkbox"/> More information <ul style="list-style-type: none"> • pParam: Pointer to an array of the type <code>double</code> returning the minimum and maximum values and the increment. • nSizeOfParam: 24
IS_EXPOSURE_CMD_GET_LONG_EXPOSURE_ENABLE	Returns the current settings for long exposure. <input type="checkbox"/> More information pval1: Pointer to a variable of the type <code>uint</code> returning the current setting
IS_EXPOSURE_CMD_SET_LONG_EXPOSURE_ENABLE	Enables/Disables long exposure. <input type="checkbox"/> More information pval1: = Pointer to a variable of the type <code>uint</code> 1 enables control, 0 disables control
pParam	Pointer to a function parameter, whose function depends on nCommand.
cbSizeOfParam	Size (in bytes) of the memory area to which pParam refers.

4.3.26.4 Setting the Dual Exposure

	
USB 3.0	USB 3.0

Syntax

```
INT is_Exposure (HCCAM hCam, UINT nCommand, void* pParam, UINT cbSizeOfParam)
```

Description

Using `is_Exposure()` you can set for the models / and / the dual exposure time.



Dual exposure and AOI

If you want to use dual exposure and AOI simultaneously, you must take care that both

the Y position and the AOI height are a multiple of 4. If this is not the case, the dual exposure cannot be activated. If the dual exposure is activated, you must observe the changed step width.

Note on dependencies on other settings

The use of the following functions will affect the exposure time:

- [is_PixelClock\(\)](#)
- [is_SetOptimalCameraTiming\(\)](#)
- [is_SetFrameRate\(\)](#)
- [is_AOI\(\)](#) (if the image size is changed)
- [is_SetSubSampling\(\)](#)
- [is_SetBinning\(\)](#)

Changes made to the image size, the frame rate or the pixel clock frequency also affect the exposure time. For this reason, you need to call `is_Exposure()` again after such changes.

Note on new driver versions

Newer driver versions sometimes allow an extended value range for the exposure time setting. We recommend querying the value range every time and set the exposure time explicitly.



The `nCommand` input parameter is used to select the function mode. The `pParam` input parameter depends on the selected function mode. If you select functions for setting or returning a value, `pParam` contains a pointer to a variable of the `UINT` type. The size of the memory area to which `pParam` refers is specified in the `cbSizeOfParam` input parameter.



Input parameters

<code>hCam</code>	Camera handle
<input type="checkbox"/> <code>nCommand</code>	
<code>IS_EXPOSURE_CMD_GET_DUAL_EXPOSURE_RATIO_DEFAULT</code>	Returns the default ratio for dual exposure.
<code>IS_EXPOSURE_CMD_GET_DUAL_EXPOSURE_RATIO_RANGE</code>	Returns the range for dual exposure (<code>IS_RANGE_S32</code>).
<code>IS_EXPOSURE_CMD_GET_DUAL_EXPOSURE_RATIO</code>	Returns the ratio between exposure times for dual exposure.
<code>IS_EXPOSURE_CMD_SET_DUAL_EXPOSURE_RATIO</code>	Sets the ratio between exposure times for dual exposure. For model / or / select a percental value between 10 and 100. E.g. 50 means that odd lines are exposed at the selected exposure time and even lines are exposed with 50% of the selected exposure time (Example). Note: The dual exposure cannot be used when the camera is operated with minimum exposure time.
<code>pParam</code>	Pointer to a function parameter, whose function depends on <code>nCommand</code> .
<code>cbSizeOfParam</code>	Size (in bytes) of the memory area to which <code>pParam</code> refers.

Example

```
UINT nRatio = 50;
INT nRet = is_Exposure(m_hCam, IS_EXPOSURE_CMD_SET_DUAL_EXPOSURE_RATIO,
                      (void*)&nRatio, sizeof(nRatio));
```

4.3.27 is_ForceTrigger

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_ForceTrigger (HCAM hCam)
```

Description

You can use `is_ForceTrigger()` to force a software-controlled capture of an image while a capturing process triggered by hardware is in progress. This function can only be used if the triggered capturing process was started using the `IS_DONT_WAIT` parameter.

Input parameters

hCam	Camera handle
------	---------------

Return values

IS_CANT_COMMUNICATE_WITH_DRIVER	Communication with the driver failed because no driver has been loaded.
IS_CANT_OPEN_DEVICE	An attempt to initialize or select the camera failed (no camera connected or initialization error).
IS_INVALID_CAMERA_HANDLE	Invalid camera handle
IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
IS_IO_REQUEST_FAILED	An IO request from the uc480 driver failed. Possibly the versions of the <code>uc480.dll</code> (API) and the driver file (<code>uc480_usb.sys</code>) do not match.
IS_NO_SUCCESS	General error message
IS_SUCCESS	Function executed successfully

Related functions

- [is_FreezeVideo\(\)](#)
- [is_CaptureVideo\(\)](#)
- [is_SetExternalTrigger\(\)](#)

Example

Enable trigger and wait 1 second for the external trigger. If no trigger signal has arrived, force an exception using `is_ForceTrigger()`.



```
HANDLE hEvent = CreateEvent(NULL, TRUE, FALSE, "");
if ( hEvent != NULL )
{
    is_InitEvent(hCam, m_hEvent, IS_SET_EVENT_FRAME);
    is_EnableEvent(hCam, IS_SET_EVENT_FRAME);
```

```
is_SetExternalTrigger(hCam, IS_SET_TRIGGER_HI_LO);
is_FreezeVideo(hCam, IS_DONT_WAIT);

if (WaitForSingleObject(m_hEvent, 1000) != WAIT_OBJECT_0)
{
    // No trigger has been received, so force image capture
    is_ForceTrigger(hCam);
}

is_DisableEvent(hCam, IS_SET_EVENT_FRAME);
is_ExitEvent(hCam, IS_SET_EVENT_FRAME);
}
```

4.3.28 is_FreeImageMem

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_FreeImageMem (HCAM hCam, char* pcImgMem, INT id)
```

Description

`is_FreeImageMem()` releases an image memory that was allocated using [is_AllocImageMem\(\)](#) and removes it from the driver management.

Note

If the memory was not allocated using an SDK function, you need to call `is_FreeImageMem()` as well. Otherwise, there may be errors when the driver keeps trying to access this memory.

This does however not release the memory. So you need to make sure that the memory will be released again.

Input parameters

hCam	Camera handle
pcImgMem	Points to the starting address of the memory (e.g. set in the is_AllocImageMem() function)
id	ID of this memory



Return values

IS_CANT_ADD_TO_SEQUENCE	The image memory is already included in the sequence and cannot be added again.
IS_CANT_CLEANUP_MEMORY	The driver could not release the allocated memory.
IS_CANT_COMMUNICATE_WITH_DRIVER	Communication with the driver failed because no driver has been loaded.
IS_CANT_OPEN_DEVICE	An attempt to initialize or select the camera failed (no camera connected or initialization error).
IS_INVALID_CAMERA_HANDLE	Invalid camera handle
IS_INVALID_MEMORY_POINTER	Invalid pointer or invalid memory ID
IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
IS_IO_REQUEST_FAILED	An IO request from the uc480 driver failed. Possibly the versions of the uc480.dll (API) and the driver file (uc480_usb.sys) do not match.
IS_NO_SUCCESS	General error message
IS_SUCCESS	Function executed successfully

Related functions

- [is_AllocImageMem\(\)](#)

4.3.29 is_FreezeVideo

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

INT is_FreezeVideo (HCAM hCam, INT Wait)

Description

is_FreezeVideo() acquires a single image from the camera. In DIB mode, the image is stored in the active image memory. If ring buffering is used in DIB mode, the captured image is transferred to the next available image memory of the sequence. Once the last available sequence memory has been filled, the sequence event or message will be triggered.

In Direct3D or OpenGL mode, the image is directly copied to the graphics card buffer and then displayed. Image capture will be started by a trigger if you previously enabled the trigger mode using [is_SetExternalTrigger\(\)](#). A hardware triggered image acquisition can be cancelled using [is_StopLiveVideo\(\)](#) if exposure has not started yet.

For further information on the image capture modes of the DCx camera, see the [How to proceed: Image capture](#) section.

Input parameters

hCam	Camera handle
<input type="checkbox"/> Wait	
IS_DONT_WAIT	Timeout value for image capture (see also the How to proceed: Timeout values for image capture section)
IS_WAIT	
Time t	

Return values

IS_BAD_STRUCTURE_SIZE	An internal structure has an incorrect size.
IS_CANT_COMMUNICATE_WITH_DRIVER	Communication with the driver failed because no driver has been loaded.
IS_CANT_OPEN_DEVICE	An attempt to initialize or select the camera failed (no camera connected or initialization error).
IS_CAPTURE_RUNNING	A capturing operation is in progress and must be terminated before you can start another one.
IS_INVALID_BUFFER_SIZE	The image memory has an inappropriate size to store the image in the desired format.
IS_INVALID_CAMERA_TYPE	The camera type defined in the .ini file does not match the current camera model.
IS_INVALID_EXPOSURE_TIME	This setting is not available for the currently set exposure time.
IS_INVALID_CAMERA_HANDLE	Invalid camera handle
IS_INVALID_MEMORY_POINTER	Invalid pointer or invalid memory ID
IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available

	in this mode.
IS_IO_REQUEST_FAILED	An IO request from the uc480 driver failed. Possibly the versions of the uc480.dll (API) and the driver file (uc480_usb.sys) do not match.
IS_NO_ACTIVE_IMG_MEM	No active image memory available. You must set the memory to active using the <code>is_SetImageMem()</code> function or create a sequence using the <code>is_AddToSequence()</code> function.
IS_NO_USB20	The camera is connected to a port which does not support the USB 2.0 high-speed standard.
IS_NO_SUCCESS	General error message
IS_NOT_CALIBRATED	The camera does not contain any calibration data.
IS_NOT_SUPPORTED	The camera model used here does not support this function or setting.
IS_OUT_OF_MEMORY	No memory could be allocated.
IS_TIMED_OUT	A timeout occurred. An image capturing process could not be terminated within the allowable period.
IS_SUCCESS	Function executed successfully
IS_TRANSFER_ERROR	Transfer error. Frequent transfer errors can mostly be avoided by reducing the pixel rate.

Related functions

- [is_HasVideoStarted\(\)](#)
- [is_IsVideoFinish\(\)](#)
- [is_SetExternalTrigger\(\)](#)
- [is_ForceTrigger\(\)](#)
- [is_CaptureVideo\(\)](#)
- [is_SetTimeout\(\)](#)
- [is_CaptureStatus\(\)](#)

Example

Enable trigger mode, set high-active flash mode and capture an image:

```
is_SetExternalTrigger(hCam, IS_SET_TRIGGER_SOFTWARE);



// Set the flash to a high active pulse for each image in the trigger mode
UINT nMode = IO_FLASH_MODE_TRIGGER_HI_ACTIVE;
is_IO(m_hCam, IS_IO_CMD_FLASH_SET_MODE, (void*)&nMode, sizeof(nMode));

is_FreezeVideo(hCam, IS_WAIT);
```

Sample programs

- SimpleAcquire (C++)
- uc480C# Demo (C#)

4.3.30 is_Gamma

	
USB 2.0 USB 3.0 GigE	USB 2.0 USB 3.0 GigE

Syntax

```
INT is_Gamma(HCAM hCam, UINT nCommand, void* pParam, UINT cbSizeOfParam)
```

Description

`is_Gamma()` enables digital gamma correction which applies a gamma characteristic to the image. When hardware color conversion is used on GigE uEye HE cameras the gamma correction is performed in the camera hardware as well. When the color conversion is performed in the PC (software conversion) the gamma correction is performed in software.



When the color format is set to Raw Bayer the gamma correction can not be used.



Typical values for gamma range between 1.6 and 2.2.

The `nCommand` input parameter is used to select the function mode. The `pParam` input parameter depends on the selected function mode. If you select functions for setting or returning a value, `pParam` contains a pointer to a variable of the `UINT` type. The size of the memory area to which `pParam` refers is specified in the `cbSizeOfParam` input parameter.



Note: The following function is obsolete by the `is_Gamma()` function:

- `is_SetGamma()`

See also: [Obsolete functions](#)

Input parameters

<code>hCam</code>	Camera handle
<div style="border: 1px solid gray; padding: 2px;"> <input type="checkbox"/> <code>nCommand</code> </div>	
<code>IS_GAMMA_CMD_SET</code>	Gamma value to be set, multiplied by 100 (Range: 1...1000) Example 1
<code>IS_GAMMA_CMD_GET_DEFAULT</code>	Returns the standard gamma value (Default = 100, corresponds to a gamma value of 1.0) Example 2
<code>IS_GAMMA_CMD_GET</code>	Returns the current set gamma value Example 3
<code>pParam</code>	Pointer to a function parameter, whose function depends on <code>nCommand</code> .
<code>cbSizeOfParam</code>	Size (in bytes) of the memory area to which <code>pParam</code> refers.

Return values

IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
IS_NO_SUCCESS	General error message
IS_NOT_SUPPORTED	The camera model used here does not support this function or setting.
IS_SUCCESS	Function executed successfully

Related functions

- `is_SetHardwareGamma()`

Example 1

```
/* set gamma value to 1.6 */
INT nGamma = 160;
INT nRet = is_Gamma(hCam, IS_GAMMA_CMD_SET, (void*) &nGamma, sizeof(nGamma));
```

Example 2

```
/* set gamma to default value */
INT nGamma;

INT nRet = is_Gamma(hCam, IS_GAMMA_CMD_GET_DEFAULT, (void*)&nGamma, sizeof(nGamma));
if (nRet == IS_SUCCESS)
{
    is_Gamma(hCam, IS_GAMMA_CMD_SET, (void*)&nGamma, sizeof(nGamma));
}
```



Beispiel 3

```
/* read gamma value to nGamma */
INT nGamma;
INT nRet = is_Gamma(hCam, IS_GAMMA_CMD_GET, (void*) &nGamma, sizeof(nGamma));
```

See also:

- Basics: [Characteristics and LUT](#)
- Basics: [Color filter \(Bayer filter\)](#)
- Programming: [is_SetColorConverter\(\)](#)

4.3.31 is_GetActiveImageMem

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_GetActiveImageMem (HCAM hCam, char** ppcMem, INT* pnID)
```

Description

`is_GetActiveImageMem()` returns the pointer to the starting address and the ID number of the active image memory.

If a Direct3D mode is active and image memory was nevertheless allocated, the pointer to the image memory and its ID will be returned. However, in Direct3D mode, the image will not be copied automatically to this image memory.

Input parameters

hCam	Camera handle
ppcMem	Returns the pointer to the starting address of the active image memory.
pnID	Returns the ID of the active image memory.



Return values

IS_CANT_COMMUNICATE_WITH_DRIVER	Communication with the driver failed because no driver has been loaded.
IS_CANT_OPEN_DEVICE	An attempt to initialize or select the camera failed (no camera connected or initialization error).
IS_INVALID_CAMERA_HANDLE	Invalid camera handle
IS_INVALID_MEMORY_POINTER	Invalid pointer or invalid memory ID
IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
IS_IO_REQUEST_FAILED	An IO request from the uc480 driver failed. Possibly the versions of the uc480.dll (API) and the driver file (uc480_usb.sys) do not match.
IS_NO_SUCCESS	General error message
IS_SUCCESS	Function executed successfully

Related functions

- [is_AllocImageMem\(\)](#)
- [is_GetImageMem\(\)](#)
- [is_SetImageMem\(\)](#)
- [is_SetAllocatedImageMem\(\)](#)

4.3.32 is_GetActSeqBuf

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_GetActSeqBuf (HCAM hCam, INT* pnNum,
                    char** ppcMem, char** ppcMemLast);
```

Description

Using `is_GetActSeqBuf()`, you can determine the image memory which is currently used for capturing an image (`ppcMem`) or the image memory that was last used for capturing an image (`ppcMemLast`). This function is only available if you have enabled ring buffering.

Attention

All input parameters of a function have to be initialized with valid values before the function is called; this also applies to parameters that are not used. Variables can be preset with '0', for example. For unused parameters, the `NULL` pointer has to be passed.

Note

This number is not the ID of the image memory that was allocated using the [is_AllocImageMem\(\)](#) function, but the running number from the order in which memory was allocated by the [is_AddToSequence\(\)](#) function.

Input parameters

<code>hCam</code>	Camera handle
<code>pnNum</code>	Contains the number of the image memory currently used for image capturing. If image capturing is already in progress when <code>is_GetActSeqBuf()</code> is called, <code>pnNum</code> will return the value 0 until the sequence arrives at the first image memory again.
<code>ppcMem</code>	Contains the starting address of the image memory currently used for image capturing.
<code>ppcMemLast</code>	Contains the starting address of the image memory last used for image capturing.

Return values



<code>IS_CANT_COMMUNICATE_WITH_DRIVER</code>	Communication with the driver failed because no driver has been loaded.
<code>IS_CANT_OPEN_DEVICE</code>	An attempt to initialize or select the camera failed (no camera connected or initialization error).
<code>IS_INVALID_CAMERA_HANDLE</code>	Invalid camera handle
<code>IS_INVALID_MEMORY_POINTER</code>	Invalid pointer or invalid memory ID
<code>IS_INVALID_PARAMETER</code>	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
<code>IS_IO_REQUEST_FAILED</code>	An IO request from the uc480 driver failed.

	Possibly the versions of the <code>uc480.dll</code> (API) and the driver file (<code>uc480_usb.sys</code>) do not match.
<code>IS_NO_SUCCESS</code>	General error message
<code>IS_SEQUENCE_LIST_EMPTY</code>	The sequence list is empty and cannot be deleted.
<code>IS_SUCCESS</code>	Function executed successfully

Related functions

- [is_AddToSequence\(\)](#)
- [is_GetImageMem\(\)](#)

4.3.33 is_GetAutoInfo

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_GetAutoInfo (HCAM hCam, UC480_AUTO_INFO* pInfo)
```

Description

Using the `is_GetAutoInfo()` function, you can query status information on the automatic image control features. This information is written to the `UC480_AUTO_INFO` structure.

For further information on automatic control, please refer to the [Automatic image control](#) chapter.

Attention

The status information returned in the `UC480_AUTO_INFO` structure is only valid if at least one of the auto control feature has been enabled using [is_SetAutoParameter\(\)](#).

Input parameters

<code>hCam</code>	Camera handle
<code>pinfo</code>	<code>UC480_AUTO_INFO</code> structure (see below)

Contents of the UC480_AUTO_INFO Structure

INT	<code>AutoAbility</code>	Supported auto control features
	<code>AC_SHUTTER</code>	Auto exposure shutter is supported
	<code>AC_SENSOR_SHUTTER</code>	The sensor's internal auto exposure shutter is supported
	<code>AC_FRAMERATE</code>	Auto frame rate is supported
	<code>AC_SENSOR_FRAMERATE</code>	The sensor's internal auto frame rate is supported
	<code>AC_GAIN</code>	Auto gain control is supported
	<code>AC_SENSOR_GAIN</code>	The sensor's internal auto gain control is supported
	<code>AC_SENSOR_AUTO_CONTRAST_CORRECTION</code>	Auto contrast correction for automatic brightness control is supported
	<code>AC_SENSOR_AUTO_CONTRAST_FDT_AOI</code>	Use of face detection as field of view for automatic brightness control is supported
	<code>AC_SENSOR_AUTO_BACKLIGHT_COMP</code>	Backlight compensation for automatic brightness control is supported
	<code>AC_WHITEBAL</code>	Auto white balance is supported
<code>AC_SENSOR_WB</code>	The sensor's internal auto white balance is supported	
<code>AUTO_BRIGHT_STATUS</code>	<code>sBrightCtrlStatus</code>	Status of automatic brightness control, see below
<code>AUTO_WB_STATUS</code>	<code>sWBCtrlStatus</code>	Status of auto white balance, see below
DWORD	<code>AShutterPhotomCaps</code>	Returns a bit mask containing all supported photometry settings (fields of view) for auto

		exposure shutter, see below.
DWORD	AGainPhotomCaps	Returns a bit mask containing all supported photometry settings (fields of view) for auto gain control, see below.
DWORD	AAntiFlickerCaps	Returns a bit mask containing all supported anti flicker settings for automatic control, see below.
DWORD	SensorWBModeCaps	Returns a bit mask containing all supported settings for the sensor's auto white balance, see below.
DWORD	reserved[8]	Reserved space for extensions

Contents of the UC480_AUTO_INFO::AUTO_BRIGHT_STATUS Structure

INT	curValue	Current average brightness of the image (actual value); the following rule applies independently of the image bit depth: 0 = black 255 = white
INT	curError	Current control deviation (error)
INT	curController	Current parameter value
	AC_SHUTTER	Exposure time (shutter)
	AC_GAIN	Gain
INT	curCtrlStatus	Current control status
	ACS_ADJUSTING	Control is active.
	ACS_FINISHED	Control is completed.
	ACS_DISABLED	Control is disabled.

Contents of the UC480_AUTO_INFO::AUTO_WB_STATUS Structure

INT	curController	Current white balance control
	AC_WB_RED_CHANNEL	Value of the red channel
	AC_WB_GREEN_CHANNEL	Value of the green channel
	AC_WB_BLUE_CHANNEL	Value of the blue channel
AUTO_WB_CHANNEL_STATUS	RedChannel	See AUTO_WB_CHANNEL_STATUS
AUTO_WB_CHANNEL_STATUS	GreenChannel	See AUTO_WB_CHANNEL_STATUS
AUTO_WB_CHANNEL_STATUS	BlueChannel	See AUTO_WB_CHANNEL_STATUS

Contents of the UC480_AUTO_INFO::AUTO_WB_STATUS::AUTO_WB_CHANNEL_STATUS Structure

INT	curValue	Current average grayscale value (actual value)
INT	curError	Current control deviation (error)
INT	curCtrlStatus	Current control status
	ACS_ADJUSTING	Control is active.
	ACS_FINISHED	Control is completed.
	ACS_DISABLED	Control is disabled.

Status Flags in UC480_AUTO_INFO::AShutterPhotomCaps

AS_PM_NONE	The entire field of view is used for metering.
AS_PM_SENS_CENTER_WEIGHTED	Metering is based on the entire field of view, but gives greater emphasis to the center area of the image.
AS_PM_SENS_CENTER_SPOT	Only a small area in the image center is used for metering.
AS_PM_SENS_PORTRAIT	Metering is based on that part of the field of view that corresponds to the portrait format.
AS_PM_SENS_LANDSCAPE	Metering is based on that part of the field of view that corresponds to the landscape format.

Status Flags in UC480_AUTO_INFO::AGainPhotomCaps

AG_PM_NONE	The entire field of view is used for metering.
AG_PM_SENS_CENTER_WEIGHTED	Metering is based on the entire field of view, but gives greater emphasis to the center area of the image.
AG_PM_SENS_CENTER_SPOT	Only a small area in the image center is used for metering.
AG_PM_SENS_PORTRAIT	Metering is based on that part of the field of view that corresponds to the portrait format.
AG_PM_SENS_LANDSCAPE	Metering is based on that part of the field of view that corresponds to the landscape format.

Status Flags in UC480_AUTO_INFO::AAntiFlickerCaps

ANTIFLCK_MODE_OFF	Anti flicker function disabled.
ANTIFLCK_MODE_SENS_AUTO	The anti flicker mode is selected automatically (50 or 60 Hz).
ANTIFLCK_MODE_SENS_50_FIXED	The anti flicker mode is set to a fixed value of 50 Hz.
ANTIFLCK_MODE_SENS_60_FIXED	The anti flicker mode is set to a fixed value of 60 Hz.

Status Flags in UC480_AUTO_INFO::SensorWBModeCaps

WB_MODE_DISABLE	Disables the sensor's auto white balance
WB_MODE_AUTO	Sensor automatically determines auto white balance
WB_MODE_ALL_PULLIN	Sensor automatically determines auto white balance using the Gray World algorithm. This algorithm assumes that the average color value in the scene is gray.
WB_MODE_INCANDESCENT_LAMP	Sensor sets auto white balance to incandescent light
WB_MODE_FLUORESCENT_DL	Sensor sets auto white balance to fluorescent light (daylight type)

WB_MODE_OUTDOOR_CLEAR_SKY	Sensor sets auto white balance to direct daylight
WB_MODE_OUTDOOR_CLOUDY	Sensor sets auto white balance to cloudy sky



Return values

IS_CANT_COMMUNICATE_WITH_DRIVER	Communication with the driver failed because no driver has been loaded.
IS_CANT_OPEN_DEVICE	An attempt to initialize or select the camera failed (no camera connected or initialization error).
IS_INVALID_CAMERA_HANDLE	Invalid camera handle
IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
IS_IO_REQUEST_FAILED	An IO request from the uc480 driver failed. Possibly the versions of the uc480.dll (API) and the driver file (uc480_usb.sys) do not match.
IS_NO_SUCCESS	General error message
IS_NOT_SUPPORTED	The camera model used here does not support this function or setting.
IS_SUCCESS	Function executed successfully

Related functions

- [is_SetAutoParameter\(\)](#)

4.3.34 is_GetBusSpeed

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_GetBusSpeed (HCAM hCam)
```

Description

Using `is_GetBusSpeed()`, you can query whether a camera is connected to a USB 2.0 or USB 3.0 host controller. You can see in the uc480 Camera Manager below ["General Information"](#) which kind of USB host controller are available on your PC.



Input parameters

hCam	Camera handle
------	---------------

Return values

IS_INVALID_CAMERA_HANDLE	Invalid camera handle
IS_NO_SUCCESS	General error message
IS_SUCCESS	Function executed successfully
IS_USB_10	The controller to which the camera is connected does not support USB 2.0.
IS_USB_20	The camera is connected to a USB 2.0 controller.
IS_USB_30	The camera is connected to a USB 3.0 controller.

4.3.35 is_GetCameraInfo

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_GetCameraInfo (HCAM hCam, CAMINFO* pInfo)
```

Description

is_GetCameraInfo() reads out the data hard-coded in the EEPROM and writes it to the data structure that pInfo points to.

Note

The serial number or model name should not be used to find a specific camera (e.g. in order to control this specific camera). If you use the serial number, the software may not find the serial number after exchanging the camera. The model name can be changed when updating the camera driver.

Instead, we recommend identifying a camera by a fixed camera ID, the camera type or by the sensor ID (see [is_GetCameraList\(\)](#)). The advantage of the camera ID is that you can set it manually. That means if you exchange a camera, you can set the same camera ID for the new camera.

Attention

For technical reasons, the following values for CAMINFO::Type are internally redirected to the same value:

IS_CAMERA_TYPE_UC480_USB_SE and IS_CAMERA_TYPE_UC480_USB_RE

You can use the parameter strSensorName of the [is_GetSensorInfo\(\)](#) function to discern the camera models DCU223x, DCU224x and DCC1240x.

Input parameters

hCam	Camera handle
pInfo	Pointer to a CAMINFO data structure

Contents of the CAMINFO Structure

char	SerNo[12]	Serial number of the camera
char	ID[20]	Manufacturer of the camera
char	Version[10]	For USB cameras, this value indicates the USB board hardware version (e.g. v2.10)
char	Date[12]	System date of the final quality check (e.g. 01.08.2011 (DD.MM.YYYY))
unsigned char	Select	Camera ID
unsigned char	Type	Camera type: <ul style="list-style-type: none"> IS_CAMERA_TYPE_UEYE_USB_SE: DCU223x, DCU224x and DCC1240x IS_CAMERA_TYPE_UEYE_USB_LE: DCC1545M / DCC1645C IS_CAMERA_TYPE_UEYE_USB3_CP: DCC3240x

char	Reserved[8]	Reserved
------	-------------	----------



Return values

IS_INVALID_CAMERA_HANDLE	Invalid camera handle
IS_NO_SUCCESS	General error message
IS_SUCCESS	Function executed successfully

Related functions

- [is_CameraStatus\(\)](#)
- [is_GetSensorInfo\(\)](#)

4.3.36 is_GetCameraList

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_GetCameraList (UC480_CAMERA_LIST* pucl)
```

Description

Using `is_GetCameraList()`, you can query information about the connected cameras. To get all information that is available, you need to adjust the field size to the number of connected cameras. The following tables explain the structures used for that purpose.

Input parameters

<code>pucl</code>	Handle to the <code>UC480_CAMERA_LIST</code> structure
-------------------	--

Contents of the UC480_CAMERA_LIST Structure

ULONG	<code>dwCount</code>	Has to be initialized with the number of cameras connected to the system. This value can be read out with is_GetNumberOfCameras() .
UC480_CAMERA_INFO	<code>uci[1]</code>	Placeholder for 1 ... n UC480_CAMERA_INFO structures

Contents of the UC480_CAMERA_LIST::UC480_CAMERA_INFO Structure

DWORD	<code>dwCameraID</code>	Customizable camera ID. This ID is stored in the camera and is persistent.
DWORD	<code>dwDeviceID</code>	Internal device ID. This ID is generated by the driver depending on order of connection and camera type. The device ID is not persistent.
DWORD	<code>dwSensorID</code>	Sensor ID
DWORD	<code>dwInUse</code>	1 = camera is being used. 0 = camera is not being used.
Char	<code>SerNo[16]</code>	Serial number of the camera
Char	<code>Model[16]</code>	Camera model
DWORD	<code>dwStatus</code>	Information for the status of the camera
DWORD	<code>dwReserved[15]</code>	Reserved for later use

Note

The serial number or model name should not be used to find a specific camera (e.g. in order to control this specific camera). If you use the serial number, the software may not find the serial number after exchanging the camera. The model name can be changed when updating the camera driver.

Instead, we recommend identifying a camera by a fixed camera ID, the camera type or by the sensor ID. The advantage of the camera ID is that you can set it manually. That means if you exchange a camera, you can set the same camera ID for the new camera.

Return values

IS_ACCESS_VIOLATION	Not enough memory allocated for the UC480_CAMERA_LIST structure
IS_CANT_OPEN_DEVICE	An attempt to initialize or select the camera failed (no camera connected or initialization error).
IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
IS_IO_REQUEST_FAILED	An IO request from the uc480 driver failed. Possibly the versions of the uc480.dll (API) and the driver file (uc480_usb.sys) do not match.
IS_NO_SUCCESS	General error message
IS_SUCCESS	Function executed successfully

Related functions



- [is_GetNumberOfCameras\(\)](#)

Example

```
// At least one camera must be available
INT nNumCam;
if( is_GetNumberOfCameras( &nNumCam ) == IS_SUCCESS) {
    if( nNumCam >= 1 ) {
        // Create new list with suitable size
        UC480_CAMERA_LIST* pucl;
        pucl = (UC480_CAMERA_LIST*) new BYTE [sizeof (DWORD) + nNumCam * sizeof (UC480_CAMERA_INFO)];
        pucl->dwCount = nNumCam;

        //Retrieve camera info
        if (is_GetCameraList(pucl) == IS_SUCCESS) {
            int iCamera;
            for (iCamera = 0; iCamera < (int)pucl->dwCount; iCamera++) {
                //Test output of camera info on the screen
                printf("Camera %i Id: %d", iCamera,
                    pucl->uci[iCamera].dwCameraID);
            }
        }
    }
    delete [] pucl;
}
```

4.3.37 is_GetCameraLUT

	
USB 3.0	USB 3.0

Syntax

```
INT is_GetCameraLUT (HCAM hCam,
                    UINT Mode, UINT NumberOfEntries,
                    double* pRed_Grey,
                    double* pGreen,
                    double* pBlue)
```

Description

`is_GetCameraLUT()` returns the current LUT values. Using the `is_SetCameraLUT()` function, you can select a different LUT for the camera.

Note

The `is_SetCameraLUT()` function is only supported by DCC3240 cameras.



Input parameters

<code>hCam</code>	Camera handle
Mode	
<code>IS_GET_CAMERA_LUT_USER</code>	Returns the LUT values set by the user without modifications.
<code>IS_GET_CAMERA_LUT_COMPLETE</code>	Returns the LUT values set by the user after the gamma, contrast and brightness values have been taken into account.
NumberOfEntries	
<code>IS_CAMERA_LUT_64</code>	LUT with 64 values
<code>pRed_Grey</code>	Pointer to the array to which the red channel values or the gray scale value of the LUT are written.
<code>pGreen</code>	Pointer to the array to which the green channel values of the LUT are written.
<code>pBlue</code>	Pointer to the array to which the blue channel values of the LUT are written.

Return values

<code>IS_INVALID_CAMERA_HANDLE</code>	Invalid camera handle
<code>IS_NO_SUCCESS</code>	General error message
<code>IS_NOT_SUPPORTED</code>	The camera model used here does not support this function or setting.
<code>IS_SUCCESS</code>	Function executed successfully

4.3.38 is_GetColorConverter

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_GetColorConverter (HCAM hCam,
                        INT ColorMode,
                        INT* pCurrentConvertMode,
                        INT* pDefaultConvertMode,
                        INT* pSupportedConvertModes)
```

Description

For color cameras, `is_GetColorConverter()` returns the set mode or all available Bayer conversion modes for the specified color mode. The return value depends on the selected color mode. For further information, please refer to the [Appendix: color and memory formats](#) section.

Input parameters

<code>hCam</code>	Camera handle
<code>ColorMode</code>	Color mode for which the converter is to be returned For a list of all available color formats and the associated input parameters, see the Appendix: Color and memory formats section.
<code>pCurrentConvertMode</code>	Currently selected converter for this color mode
<code>pDefaultConvertMode</code>	Default converter for this color mode
<code>pSupportedConvertModes</code>	All converters supported for this color mode. Possible converters are: <code>IS_CONV_MODE_NONE</code> <code>IS_CONV_MODE_SOFTWARE_3X3</code> <code>IS_CONV_MODE_SOFTWARE_5X5</code> <code>IS_CONV_MODE_HARDWARE_3X3</code>



Return values

<code>IS_INVALID_COLOR_FORMAT</code>	Invalid color format
<code>IS_INVALID_CAMERA_HANDLE</code>	Invalid camera handle
<code>IS_NO_SUCCESS</code>	General error message
<code>IS_SUCCESS</code>	Function executed successfully

Related functions

- [is_SetColorConverter\(\)](#)
- [is_SetColorMode\(\)](#)

4.3.39 is_GetColorDepth

	
USB 2.0 USB 3.0	-

Syntax

```
INT is_GetColorDepth(HCAM hCam, INT* pnCol, INT* pnColMode)
```

Description

`is_GetColorDepth()` retrieves the current Windows Desktop color setting and returns the bit depth per pixel and the matching uc480 color mode. The color mode can be passed directly to the [is_SetColorMode\(\)](#) function. You need to pass the bit depth when allocating an image memory.

Input parameters

hCam	Camera handle
pnCol	Returns the bit depth of the color setting.
pnColMode	Returns the uc480 color mode that corresponds to pnCol. For a list of all available color formats and the associated input parameters, see the Appendix: color and Memory Formats section.



Return values

IS_INVALID_CAMERA_HANDLE	Invalid camera handle
IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
IS_NO_SUCCESS	General error message
IS_SUCCESS	Function executed successfully

Related functions

- [is_SetColorMode\(\)](#)
- [is_AllocImageMem\(\)](#)

4.3.40 is_GetDLLVersion

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_GetDLLVersion()
```

Description

Returns the version of the `uc480.dll`.

Input parameters

<none>

Return values

The return value contains the version number which is coded as follows:

Bits 31-24:	Major version
Bits 23-16:	Minor version
Bits 15-0:	Build version



Related functions

- [is_GetOsVersion\(\)](#)

Example

```
int version = is_GetDLLVersion();
int build = version & 0xFFFF;
version = version >> 16;
int minor = version & 0xFF;
version = version >> 8;
int major = version & 0xFF;
printf("API version %d.%d.%d \n\n", major, minor, build);
```

4.3.41 is_GetError

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_GetError (HCAM hCam, INT* pErr, IS_CHAR** ppcErr)
```

Description

`is_GetError()` queries the last error that occurred and returns the associated error code and message. We recommend to use this function after an error has occurred that returned `IS_NO_SUCCESS`. Each error message will be overwritten when a new error occurs.

Input parameters

<code>hCam</code>	Camera handle
<code>pErr</code>	Pointer to the variable containing the error code
<code>ppcErr</code>	Pointer to the string containing the error text



Return values

<code>IS_INVALID_CAMERA_HANDLE</code>	Invalid camera handle
<code>IS_INVALID_PARAMETER</code>	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
<code>IS_NO_SUCCESS</code>	General error message
<code>IS_SUCCESS</code>	Function executed successfully

Related functions

- [is_CaptureStatus\(\)](#)
- [is_SetErrorReport\(\)](#)
- [is_CameraStatus\(\)](#)

4.3.42 is_GetFramesPerSecond

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_GetFramesPerSecond (HCAM hCam, double* dblFPS)
```

Description

In live capture mode started by [is_CaptureVideo\(\)](#), the `is_GetFramesPerSecond()` function returns the number of frames actually captured per second.

Input parameters

hCam	Camera handle
dblFPS	Returns the current frame rate.



Return values

IS_INVALID_CAMERA_HANDLE	Invalid camera handle
IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
IS_NO_SUCCESS	General error message
IS_SUCCESS	Function executed successfully

Related functions

- [is_GetFrameTimeRange\(\)](#)
- [is_SetFrameRate\(\)](#)
- [is_PixelClock\(\)](#)
- [is_Exposure\(\)](#)

4.3.43 is_GetFrameTimeRange

	
USB 2.0	USB 2.0

Syntax

```
INT is_GetFrameTimeRange (HCAM hCam,
                          double* min, double* max, double* intervall)
```

Description

Using `is_GetFrameTimeRange()`, you can read out the frame rate settings which are available for the current pixel clock setting. The returned values indicate the minimum and maximum frame duration in seconds. You can set the frame duration between `min` and `max` in increments defined by the `intervall` parameter.

The following applies:

$$fps_{min} = \frac{1}{max}$$

$$fps_{max} = \frac{1}{min}$$

$$fps_n = \frac{1}{(min + n * intervall)}$$

The call of this function makes only sense in the [freerun mode](#).

Note

The use of the following functions will affect the frame duration:

- [is_PixelClock\(\)](#)
- [is_SetOptimalCameraTiming\(\)](#)
- [is_AOI\(\)](#) (if the image size is changed)
- [is_SetSubSampling\(\)](#)
- [is_SetBinning\(\)](#)

Changes made to the window size, the frame rate or the read-out timing (pixel clock frequency) also affect the defined frame duration. For this reason, you need to call `is_GetFrameTimeRange()` again after such changes.

Input parameters

hCam	Camera handle
min	Returns the minimum available frame duration.
max	Returns the maximum available frame duration.
intervall	Returns the increment you can use to change the frame duration.



Return values

IS_INVALID_CAMERA_HANDLE	Invalid camera handle
IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
IS_NO_SUCCESS	General error message
IS_SUCCESS	Function executed successfully

Related functions

- [is_GetFramesPerSecond\(\)](#)
- [is_SetFrameRate\(\)](#)
- [is_PixelClock\(\)](#)
- [is_Exposure\(\)](#)

4.3.44 is_GetImageHistogram

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_GetImageHistogram (HCAM hCam,
                        int nID, INT ColorMode, DWORD* pHistoMem)
```

Description

`is_GetImageHistogram()` computes the histogram of the submitted image. The histogram always contains 256 values per channel. For color modes with a bit depth of more than 8 bits, the system evaluates the 8 most significant bits (MSBs).

Input parameters

<code>hCam</code>	Camera handle
<code>nID</code>	Memory ID
<code>ColorMode</code>	Color mode of the image with the <code>nID</code> memory ID For a list of all available color formats and the associated input parameters, see the Appendix: Color and memory formats section.
<code>pHistoMem</code>	Pointer to a <code>DWORD</code> array The array must be allocated in such a way that it can accommodate 3*256 values for color formats and in raw Bayer mode. In monochrome mode, the array must be able to accommodate 1*256 values.

Return values

<code>IS_CANT_COMMUNICATE_WITH_DRIVER</code>	Communication with the driver failed because no driver has been loaded.
<code>IS_CANT_OPEN_DEVICE</code>	An attempt to initialize or select the camera failed (no camera connected or initialization error).
<code>IS_INVALID_COLOR_FORMAT</code>	Invalid color format
<code>IS_INVALID_CAMERA_HANDLE</code>	Invalid camera handle
<code>IS_INVALID_MEMORY_POINTER</code>	Invalid pointer or invalid memory ID
<code>IS_INVALID_PARAMETER</code>	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
<code>IS_IO_REQUEST_FAILED</code>	An IO request from the <code>uc480</code> driver failed. Possibly the versions of the <code>uc480.dll</code> (API) and the driver file (<code>uc480_usb.sys</code>) do not match.
<code>IS_NO_SUCCESS</code>	General error message
<code>IS_NULL_POINTER</code>	Invalid array
<code>IS_SUCCESS</code>	Function executed successfully

Related functions

- [is_SetColorMode\(\)](#)

Example

```
char * pcSource;
INT nIDSource;
is_AllocImageMem (hCam, 256, 256, 24, &pcSource, &nIDSource);

int nX, nY, nBits, nPitch;
is_InquireImageMem (hCam, pcSource, nIDSource, &nX, &nY, &nBits, &nPitch);



//Create RGB test image
for (int j = 0; j < nY; j++)
{
    for (int i = 0; i < nX*3; i += 3)
    {
        pcSource[i + j*nPitch] = 0; // Blue pixels
        pcSource[i + j*nPitch + 1] = i/3; // Green pixels
        pcSource[i + j*nPitch + 2] = 255; // Red pixels
    }
}

// Create memory for RGB histogram
DWORD bgrBuffer [256*3];

//Create pointer for each histogram color
DWORD * pBlueHisto = bgrBuffer;
DWORD * pGreenHisto = bgrBuffer + 256;
DWORD * pRedHisto = bgrBuffer + 512;

//Retrieve histogram and release memory
is_GetImageHistogram (hCam, nIDSource, IS_CM_RGB8_PACKED, bgrBuffer);
is_FreeImageMem (hCam, pcSource, nIDSource);
```

4.3.45 is_GetImageInfo

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_GetImageInfo (HCAM hCam, INT nImageBufferID, UEYEIMAGEINFO* pImageInfo, INT nImageInfoSize)
```

- [Description](#)
- [Input parameters](#)
- [Contents of the UEYEIMAGEINFO structure](#)
- [Status flags in UEYEIMAGEINFO::dwIoStatus](#)
- [Contents of the UEYETIME structure](#)
- [Return values](#)
- [Related functions](#)
- [Code sample](#)

Description

`is_GetImageInfo()` provides additional information on the images you take. The function returns a timestamp indicating the time of image capture, and the states of the camera I/Os at that point in time. To get information on the last image that was taken, call `is_GetImageInfo()` directly after receiving the `IS_FRAME` event.

Using the function with USB DCx Cameras

The `u64TimestampDevice` timestamp returns the time when image data transfer to the PC was completed.

The `UEYETIME` structure returns the timestamp (with a resolution of 1 ms) synchronized with the PC system time.

Attention

Image buffers that are part of a sequence need to be locked using [is_LockSeqBuf\(\)](#). This is important to ensure correct assignment between image data and image information. Otherwise, it may happen that an image buffer is filled with new image data. In this case, the image information will not match the image data any more.

Input parameters

<code>hCam</code>	Camera handle
<code>nImageBufferID</code>	ID of the image buffer for which information is requested
<code>pImageInfo</code>	Pointer to a <code>UC480IMAGEINFO</code> type structure to which the information will be written
<code>nImageInfoSize</code>	Size of the structure

Contents of the UC480IMAGEINFO structure

DWORD	<code>dwFlags</code>	Internal status flags (currently not used)
-------	----------------------	--

BYTE	byReserved1[4]	Reserved
unsigned long long	u64TimestampDevice	Internal timestamp of image capture (tick count of the camera in 0.1 μ s steps)
UEYETIME	TimestampSystem	Structure with timestamp information in PC system time format, see UEYETIME below
DWORD	dwIoStatus	With DCC3240x Cameras: Returns the states of the digital I/Os at the time of image capture: <ul style="list-style-type: none"> • Digital input (trigger): Pending signal • GPIO as input: Pending signal • GPIO as output: Set level With all other cameras, dwIoStatus is empty. See dwIoStatus below.
WORD	wAOIIndex	AOI index (only AOI sequence mode of DCC1240x)
WORD	wAOICycle	Readout cycles (only AOI sequence mode of DCC1240x)
unsigned long long	u64FrameNumber	Internal image number
DWORD	dwImageBuffers	Number of image buffers existing in the camera
DWORD	dwImageBuffersInUse	Number of image buffers in use in the camera
DWORD	dwReserved3	Reserved
DWORD	dwImageHeight	Image height
DWORD	dwImageWidth	Image width

Status flags in UC480IMAGEINFO::dwIoStatus

Bit combination	State of digital input	State of GPIO 1	State of GPIO 2
000	0	0	0
001	0	0	1
010	0	1	0
011	0	1	1
100	1	0	0
101	1	0	1
110	1	1	0
111	1	1	1

Contents of the UC480IMAGEINFO::UEYETIME structure

WORD	wYear	Timestamp year
WORD	wMonth	Timestamp month
WORD	wDay	Timestamp day
WORD	wHour	Timestamp hour
WORD	wMinute	Timestamp minute
WORD	wSecond	Timestamp second
WORD	wMilliseconds	Timestamp millisecond
WORD	wReserved[2]	Reserved

Return values

IS_INVALID_CAMERA_HANDLE	Invalid camera handle
IS_INVALID_PARAMETER	<p>One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.</p> <p>This may happen when e.g.:</p> <ul style="list-style-type: none"> • more memory is allocated than the <code>UC480IMAGEINFO</code> structure needs • <code>nImageBufferID <= 0</code> • <code>pImageInfo == NULL</code> • <code>nImageInfoSize <= 0</code>
IS_NO_SUCCESS	General error message
IS_SUCCESS	Function executed successfully

Related functions

- [is_CaptureStatus\(\)](#)
- [is_LockSeqBuf\(\)](#)
- [is_UnlockSeqBuf\(\)](#)
- [is_SetImageMem\(\)](#)
- [is_IO\(\)](#)

Example

```
UEYEIMAGEINFO ImageInfo;
// Read out camera timestamp
INT nRet = is_GetImageInfo( m_hCam,
                           m_lMemoryId, &ImageInfo, sizeof(ImageInfo));
if (nRet == IS_SUCCESS)
{
    unsigned long long u64TimestampDevice;
    u64TimestampDevice = ImageInfo.u64TimestampDevice;



    CString Str; // Read out timestamp in system time
    Str.Format( "%02d.%02d.%04d, %02d:%02d:%02d:%03d",
               ImageInfo.TimestampSystem.wDay,
               ImageInfo.TimestampSystem.wMonth,
               ImageInfo.TimestampSystem.wYear,
               ImageInfo.TimestampSystem.wHour,
               ImageInfo.TimestampSystem.wMinute,
               ImageInfo.TimestampSystem.wSecond,
               ImageInfo.TimestampSystem.wMilliseconds);

    DWORD dwTotalBuffers = ImageInfo.dwImageBuffers;
    DWORD dwUsedBuffers = ImageInfo.dwImageBuffersInUse;
}
```

Sample Program

- uc480Timestamp (C++)

4.3.46 is_GetImageMem

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_GetImageMem (HCAM hCam, VOID** pMem)
```

Description

`is_GetImageMem()` returns the pointer to the starting address of the active image memory. If you use ring buffering, `is_GetImageMem()` returns the starting address of the image memory last used for image capturing.

Input parameters

<code>hCam</code>	Camera handle
<code>pMem</code>	Pointer to the starting address of the image memory

Return values

<code>IS_CANT_COMMUNICATE_WITH_DRIVER</code>	Communication with the driver failed because no driver has been loaded.
<code>IS_CANT_OPEN_DEVICE</code>	An attempt to initialize or select the camera failed (no camera connected or initialization error).
<code>IS_INVALID_CAMERA_HANDLE</code>	Invalid camera handle
<code>IS_INVALID_MEMORY_POINTER</code>	Invalid pointer or invalid memory ID
<code>IS_INVALID_PARAMETER</code>	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
<code>IS_IO_REQUEST_FAILED</code>	An IO request from the <code>uc480</code> driver failed. Possibly the versions of the <code>uc480.dll</code> (API) and the driver file (<code>uc480_usb.sys</code>) do not match.
<code>IS_NO_SUCCESS</code>	General error message
<code>IS_SUCCESS</code>	Function executed successfully



Related functions

- [is_GetImageMemPitch\(\)](#)
- [is_AllocImageMem\(\)](#)
- [is_AddToSequence\(\)](#)
- [is_SetImageMem\(\)](#)
- [is_SetAllocatedImageMem\(\)](#)

Sample programs

- `uc480PixelPeek` (C++)

4.3.47 is_GetImageMemPitch

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_GetImageMemPitch (HCAM hCam, INT* pPitch)
```

Description

`is_GetImageMemPitch()` returns the line increment (in bytes). The line increment is defined as the number of bytes from the beginning of a line to the beginning of the next line. It may be greater than suggested by the parameters passed when calling [is_AllocImageMem\(\)](#). The line increment is always a number that can be divided by 4.

The line increment is calculated as:

```
line = width * int[(bitapixel + 7) / 8]
```

```
lineinc = line + adjust
```

```
adjust = 0 - if line can be divided by 4 without remainder
```

```
adjust = 4 - rest(line / 4) if line cannot be divided by 4 without remainder
```

Input parameters

hCam	Camera handle
pPitch	Pointer to the variable containing the line increment

Return values



IS_CANT_COMMUNICATE_WITH_DRIVER	Communication with the driver failed because no driver has been loaded.
IS_CANT_OPEN_DEVICE	An attempt to initialize or select the camera failed (no camera connected or initialization error).
IS_INVALID_CAMERA_HANDLE	Invalid camera handle
IS_INVALID_MEMORY_POINTER	Invalid pointer or invalid memory ID
IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
IS_IO_REQUEST_FAILED	An IO request from the uc480 driver failed. Possibly the versions of the <code>uc480.dll</code> (API) and the driver file (<code>uc480_usb.sys</code>) do not match.
IS_NO_SUCCESS	General error message
IS_SUCCESS	Function executed successfully

Related functions

- [is_GetImageMem\(\)](#)
- [is_AllocImageMem\(\)](#)
- [is_AddToSequence\(\)](#)

- [is_SetImageMem\(\)](#)
- [is_SetAllocatedImageMem\(\)](#)

4.3.48 is_GetNumberOfCameras

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_GetNumberOfCameras (INT* pNumCams)
```

Description

`is_GetNumberOfCameras()` returns the number of DCx Cameras connected to the PC.

Input parameters

pNumCams	Returns the number of connected cameras.
----------	--



Return values

IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
IS_SUCCESS	Function executed successfully

Related functions

- [is_GetCameraList\(\)](#)
- [is_DeviceInfo\(\)](#)

4.3.49 is_GetOsVersion

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_GetOsVersion ()
```

Description

is_GetOsVersion() returns the operating system type used at runtime.

Input parameters

<none>



Return values

IS_OS_WIN2000	Windows 2000 operating system
IS_OS_WINXP	Windows XP operating system
IS_OS_WINSERVER2003	Windows Server 2003 operating system
IS_OS_WINVISTA	Windows Vista operating system
IS_OS_WIN7	Windows 7 operating system
IS_OS_LINUX26	Linux 2.6 operating system
IS_OS_UNDETERMINED	Unknown operating system

Related functions

- [is_GetDLLVersion\(\)](#)

4.3.50 is_GetSensorInfo

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_GetSensorInfo (HCAM hCam, SENSORINFO* pInfo)
```

Description

Using `is_GetSensorInfo()`, you can query information about the sensor type used in the camera. The information contained in the `SENSORINFO` structure is listed in the table below.

The [ueye.h](#) file provides a complete up-to-date list of all supported sensor types. To quickly locate the list, search the file for the keyword "Sensor types".

Input parameters

hCam	Camera handle
pInfo	Pointer to the <code>SENSORINFO</code> Structure

Contents of the `SENSORINFO` structure

WORD	SensorID	Returns the sensor type (e.g.: <code>IS_SENSOR_XXXXX_X</code>).
Char	strSensorName[32]	Returns the camera model (e.g.: <code>XXXXXXXX-X</code>).
Char	nColorMode	Returns the sensor color mode. <code>IS_COLORMODE_BAYER</code> <code>IS_COLORMODE_MONOCHROME</code>
DWORD	nMaxWidth	Returns the maximum image width
DWORD	nMaxHeight	Returns the maximum image height
BOOL	bMasterGain	Indicates whether the sensor provides analog master gain
BOOL	bRGain	Indicates whether the sensor provides analog red channel gain
BOOL	bGGain	Indicates whether the sensor provides analog green channel gain
BOOL	bBGain	Indicates whether the sensor provides analog blue channel gain
BOOL	bGlobShutter	Indicates whether the sensor has a global shutter. <code>TRUE</code> = global shutter <code>FALSE</code> = rolling shutter
WORD	wPixelSize	Returns the pixel size in μm (e.g. 465 is equivalent to 4.65 μm)
Char	Reserved[14]	Reserved

Return values



<code>IS_INVALID_CAMERA_HANDLE</code>	Invalid camera handle
<code>IS_INVALID_PARAMETER</code>	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.

IS_NO_SUCCESS	General error message
IS_NOT_SUPPORTED	The camera model used here does not support this function or setting.
IS_SUCCESS	Function executed successfully

Related functions

- [is_GetCameraInfo\(\)](#)
- [is_CameraStatus\(\)](#)

4.3.51 is_GetSensorScalerInfo

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_GetSensorScalerInfo (HCAM hCam,
                           SENSORSCALERINFO* pSensorScalerInfo,
                           INT nSensorScalerInfoSize)
```

Description

Using `is_GetSensorScalerInfo()` you can query information on the image scaling which is supported by some sensors.

Note

Internal image scaling is only supported by DCC1240x and DCC3240x series sensors.

Input parameters

<code>hCam</code>	Camera handle
<code>pSensorScalerInfo</code>	Pointer to a <code>SENSORSCALERINFO</code> type structure to which the information will be written
<code>nSensorScalerInfoSize</code>	Size of the structure

Contents of the SENSORSCALERINFO structure

INT	<code>nCurrMode</code>	Returns the current mode
INT	<code>nNumberOfSteps</code>	Returns the number of steps for the scaling factor
double	<code>dblFactorIncrement</code>	Returns the increment for the scaling factor
double	<code>dblMinFactor</code>	Returns the minimum scaling factor
double	<code>dblMaxFactor</code>	Returns the maximum scaling factor
double	<code>dblCurrFactor</code>	Returns the current scaling factor
INT	<code>nSupportedModes</code>	Returns the supported function modes, see is_SetSensorScaler()
BYTE	<code>bReserved[84]</code>	Reserved



Return values

<code>IS_INVALID_CAMERA_HANDLE</code>	Invalid camera handle
<code>IS_INVALID_PARAMETER</code>	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
<code>IS_NO_SUCCESS</code>	General error message
<code>IS_NOT_SUPPORTED</code>	The camera model used here does not support this function or setting.
<code>IS_SUCCESS</code>	Function executed successfully

Related functions

- [is_SetSensorScaler\(\)](#)

4.3.52 is_GetSupportedTestImages

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_GetSupportedTestImages (HCAM hCam, INT* SupportedTestImages)
```

Description

`is_GetSupportedTestImages()` returns all test images supported by the camera. You can enable the sensor test image feature using [is_SetSensorTestImage\(\)](#).

Input parameters

hCam	Camera handle
SupportedTestImages	Returns a bit mask of all test images supported by the camera.
IS_TEST_IMAGE_NONE	No test image
IS_TEST_IMAGE_WHITE	White image
IS_TEST_IMAGE_BLACK	Black image
IS_TEST_IMAGE_HORIZONTAL_GREYSCALE	Horizontal grayscale
IS_TEST_IMAGE_VERTICAL_GREYSCALE	Vertical grayscale
IS_TEST_IMAGE_DIAGONAL_GREYSCALE	Diagonal grayscale
IS_TEST_IMAGE_WEDGE_GRAY_SENSOR	Gray wedges, generated by the sensor
IS_TEST_IMAGE_WEDGE_COLOR	Color wedges
IS_TEST_IMAGE_ANIMATED_WEDGE_GRAY_SENSOR	Gray wedges, animated, generated by the sensor
IS_TEST_IMAGE_ANIMATED_WEDGE_COLOR	Color wedges, animated
IS_TEST_IMAGE_COLOR_BARS1	Color bars
IS_TEST_IMAGE_GREY_AND_COLOR_BARS	Gray and color bars
IS_TEST_IMAGE_MOVING_GREY_AND_COLOR_BARS	Gray and color bars, animated
IS_TEST_IMAGE_ANIMATED_LINE	Line, animated
IS_TEST_IMAGE_ALTERNATE_PATTERN	Alternating pattern (raw Bayer mode only)
IS_TEST_IMAGE_RAMPING_PATTERN	Diagonal color pattern
IS_TEST_IMAGE_MONOCHROME_HORIZONTAL_BARS	Monochrome bars, horizontal
IS_TEST_IMAGE_MONOCHROME_VERTICAL_BARS	Monochrome bars, vertical
IS_TEST_IMAGE_COLDPIXEL_GRID	Camera image overlaid with a grid of blue dots
IS_TEST_IMAGE_HOTPIXEL_GRID	Camera image overlaid with a grid of red dots
IS_TEST_IMAGE_VARIABLE_GREY	Adjustable grayscale image
IS_TEST_IMAGE_VARIABLE_RED_PART	Image with adjustable red content
IS_TEST_IMAGE_VARIABLE_GREEN_PART	Image with adjustable green content
IS_TEST_IMAGE_VARIABLE_BLUE_PART	Image with adjustable blue content



Return values

IS_INVALID_CAMERA_HANDLE	Invalid camera handle
IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
IS_NO_SUCCESS	General error message
IS_NOT_SUPPORTED	The camera model used here does not support this function or setting.
IS_SUCCESS	Function executed successfully

Related functions

- [is_SetSensorTestImage\(\)](#)
- [is_GetTestImageValueRange\(\)](#)

4.3.53 is_GetTestImageValueRange

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_GetTestImageValueRange (HCAM hCam,
                               INT TestImage,
                               INT* TestImageValueMin,
                               INT* TestImageValueMax)
```

Description

Using `is_GetTestImageValueRange()`, you can query the value range of the additional parameter required for some camera test images. You can enable the sensor test image feature using [is_SetSensorTestImage\(\)](#).

Input parameters

hCam	Camera handle
TestImage	Test image for which the value range is queried
TestImageValueMin	Minimum value
TestImageValueMax	Maximum value



Return values

IS_INVALID_CAMERA_HANDLE	Invalid camera handle
IS_NO_SUCCESS	General error message
IS_NOT_SUPPORTED	The camera model used here does not support this function or setting. In this case, the <code>TestImageValueMin</code> and <code>TestImageValueMax</code> parameters are equal to 0.
IS_SUCCESS	Function executed successfully

Related functions

- [is_GetSupportedTestImages\(\)](#)
- [is_SetSensorTestImage\(\)](#)

4.3.54 is_GetTimeout

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_GetTimeout (HCAM hCam, UINT nMode, UINT* pTimeout)
```

Description

Using `is_GetTimeout()`, you can read out user-defined timeout values from the uc480 API.

For further information, please refer to the [How to proceed: Timeout values for image capture](#) section.

Input parameters

hCam	Camera handle
nMode	Selects the timeout value to be returned
IS_TRIGGER_TIMEOUT	Returns the timeout value in steps of 10 ms for triggered image capture
pTimeout	Pointer to the variable that holds the timeout value. Returns 0 if the default value of the uc480 API is used.

Return values

IS_INVALID_CAMERA_HANDLE	Invalid camera handle
IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
IS_NO_SUCCESS	General error message
IS_NOT_SUPPORTED	The camera model used here does not support this function or setting.
IS_SUCCESS	Function executed successfully



Related functions

- [is_SetTimeout\(\)](#)
- [is_CaptureVideo\(\)](#)
- [is_FreezeVideo\(\)](#)
- [is_SetExternalTrigger\(\)](#)

Example

```
// Return user-defined timeout
UINT nTimeout;
INT ret = is_GetTimeout(hCam, IS_TRIGGER_TIMEOUT, &nTimeout);
```

4.3.55 is_GetUsedBandwidth

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_GetUsedBandwidth (HCAM hCam)
```

Description

`is_GetUsedBandwidth()` returns the bus bandwidth (in MByte/s) currently used by all initialized or selected cameras. This is an approximate value which is calculated based on the pixel clock that has been set and the data format (bits per pixel). The actual data load on the bus can slightly deviate from this value.

Input parameters

hCam	Camera handle
------	---------------



Return values

INT value	The total current bus bandwidth (in MByte/s)
-----------	--

Related functions

- [is_PixelClock\(\)](#)

4.3.56 is_GetVsyncCount

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_GetVsyncCount (HCAM hCam, long* pIntr, long* pActIntr)
```

Description

`is_GetVsyncCount()` reads out the VSYNC counter. It will be incremented by 1 each time the sensor starts capturing an image.

Input parameters

<code>hCam</code>	Camera handle
<code>pIntr</code>	Current VSYNC count
<code>pActIntr</code>	Current Frame SYNC count



Return values

<code>IS_INVALID_CAMERA_HANDLE</code>	Invalid camera handle
<code>IS_NO_SUCCESS</code>	General error message
<code>IS_SUCCESS</code>	Function executed successfully

Related functions

- [is_GetFramesPerSecond\(\)](#)

4.3.57 is_HasVideoStarted

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_HasVideoStarted (HCAM hCam, BOOL* pbo)
```

Description

Using `is_HasVideoStarted()`, you can check whether the image digitizing process has started. This function is helpful when the [is_FreezeVideo\(\)](#) function was called with the `IS_DONT_WAIT` parameter.

Input parameters

hCam	Camera handle
pbo	Returns the digitizing status: 0 = Image capturing has not started yet. 1 = Image capturing has started.



Return values

IS_INVALID_CAMERA_HANDLE	Invalid camera handle
IS_NO_SUCCESS	General error message
IS_SUCCESS	Function executed successfully

Related functions

- [is_FreezeVideo\(\)](#)
- [is_IsVideoFinish\(\)](#)

4.3.58 is_HotPixel

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_HotPixel (HCAM hCam, UINT nCommand, void* pParam, UINT nSizeOfParam)
```

Description

`is_HotPixel()` configures the correction of sensor hot pixels. The correction is performed by the software. The hot pixel list is stored in the camera's non-volatile EEPROM. Some sensor models can also correct hot pixels directly in the sensor.

For further information on hot pixel correction, please refer to [Basics: Hot pixels](#).

The `nCommand` input parameter is used to select the function mode. The `pParam` input parameter depends on the selected function mode. If you select functions for setting or returning a value, `pParam` contains a pointer to a variable of the `UINT` type. The size of the memory area to which `pParam` refers is specified in the `nSizeOfParam` input parameter.

Attention

This correction will not work with subsampling or with binning factors greater than 2.

Note

Previous hot pixel functions

The `is_HotPixel()` function comprises all hot pixel correction functions. The following uc480 API commands are therefore obsolete:

- `is_SetBadPixelCorrection()`
- `is_SetBadPixelCorrectionTable()`
- `is_LoadBadPixelCorrectionTable()`
- `is_SaveBadPixelCorrectionTable()`

See also [Obsolete functions](#)

Input parameters

hCam	Camera handle
<input type="checkbox"/> nCommand	
IS_HOTPIXEL_DISABLE_CORRECTION	Disables hot pixel correction (Example 1)
IS_HOTPIXEL_ENABLE_CAMERA_CORRECTION	Enables hot pixel correction using the hot pixel list(s) stored in the camera EEPROM.
IS_HOTPIXEL_ENABLE_SOFTWARE_USER_CORRECTION	Enables hot pixel correction using the user's hot pixel list stored in the computer. This requires the user's hot pixel list to be set (<code>IS_HOTPIXEL_SET_SOFTWARE_USER_LIST</code>)
IS_HOTPIXEL_ENABLE_SENSOR_CORRECTION	Enables sensor's own hot pixel correction function (if available).
IS_HOTPIXEL_DISABLE_SENSOR_CORRECTION	Disables the sensor's own hot pixel correction function.

hCam	Camera handle
IS_HOTPIXEL_GET_CORRECTION_MODE	Returns the currently set hot pixel correction mode (Example 2)
IS_HOTPIXEL_GET_SUPPORTED_CORRECTION_MODES	Returns the supported hot pixel correction modes. The return value is a bitmask with the following constants (combined by OR): <ul style="list-style-type: none"> IS_HOTPIXEL_ENABLE_CAMERA_CORRECTION: Hot pixel correction is possible via the hot pixel list in the camera EEPROM. IS_HOTPIXEL_ENABLE_SOFTWARE_USER_CORRECTION: Hot pixel correction is possible via the user-defined hot pixel list. IS_HOTPIXEL_ENABLE_SENSOR_CORRECTION: Hot pixel correction is possible via the sensor-internal hot pixel correction.
IS_HOTPIXEL_GET_SOFTWARE_USER_LIST_EXISTS	Indicates whether the user-defined hot pixel list exists in the computer (Example 3)
IS_HOTPIXEL_GET_SOFTWARE_USER_LIST_NUMBER	Returns the number of hot pixels in the user-defined hot pixel list stored in the computer.
IS_HOTPIXEL_GET_SOFTWARE_USER_LIST	Returns the user-defined hot pixel list stored in the computer.
IS_HOTPIXEL_SET_SOFTWARE_USER_LIST	Sets the user-defined hot pixel list that is stored in the computer.
IS_HOTPIXEL_SAVE_SOFTWARE_USER_LIST IS_HOTPIXEL_SAVE_SOFTWARE_USER_LIST_UNICODE	Saves the user-defined hot pixel list to a file. The function can also be used with Unicode file names. (Example 4)
IS_HOTPIXEL_LOAD_SOFTWARE_USER_LIST IS_HOTPIXEL_LOAD_SOFTWARE_USER_LIST_UNICODE	Loads the user-defined hot pixel list from a file. The function can also be used with Unicode file names.
IS_HOTPIXEL_GET_CAMERA_FACTORY_LIST_EXISTS	Indicates whether the factory-set hot pixel list exists.
IS_HOTPIXEL_GET_CAMERA_FACTORY_LIST_NUMBER	Returns the number of hot pixels in the factory-set hot pixel list.
IS_HOTPIXEL_GET_CAMERA_FACTORY_LIST	Returns the factory-set hot pixel list.
IS_HOTPIXEL_GET_CAMERA_USER_LIST_EXISTS	Indicates whether the user-defined hot pixel list exists in the camera EEPROM.
IS_HOTPIXEL_GET_CAMERA_USER_LIST_NUMBER	Returns the number of hot pixels in the user-defined hot pixel list stored in the camera EEPROM.
IS_HOTPIXEL_GET_CAMERA_USER_LIST	Returns the user-defined hot pixel list stored in the camera EEPROM.
IS_HOTPIXEL_SET_CAMERA_USER_LIST	Sets the user-defined hot pixel list stored in the camera EEPROM (Example 5)
IS_HOTPIXEL_DELETE_CAMERA_USER_LIST	Deletes the user-defined hot pixel list from the camera EEPROM.
IS_HOTPIXEL_GET_CAMERA_USER_LIST_MAX_NUMBER	Returns the maximum number of hot pixels that the user can store in the camera EEPROM.

hCam	Camera handle
IS_HOTPIXEL_GET_MERGED_CAMERA_LIST_NUMB ER	Returns the number of hot pixels in a merged list that combines the entries from the factory-set hot pixel list with those of the user-defined hot pixels list stored in the camera EEPROM.
IS_HOTPIXEL_GET_MERGED_CAMERA_LIST	Returns the merged list (Example 6)
pParam	Pointer to a function parameter, whose function depends on nCommand.
nSizeOfParam	Size (in bytes) of the memory area to which pParam refers.

Return values

IS_CANT_COMMUNICATE_WITH_DRIVER	Communication with the driver failed because no driver has been loaded.
IS_CANT_OPEN_DEVICE	An attempt to initialize or select the camera failed (no camera connected or initialization error).
IS_INVALID_CAMERA_TYPE	The camera type defined in the .ini file does not match the current camera model.
IS_INVALID_CAMERA_HANDLE	Invalid camera handle
IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
IS_IO_REQUEST_FAILED	An IO request from the uc480 driver failed. Possibly the versions of the uc480.dll (API) and the driver file (uc480_usb.sys) do not match.
IS_NO_SUCCESS	General error message
IS_NOT_CALIBRATED	The camera does not contain any calibration data.
IS_NOT_SUPPORTED	The camera model used here does not support this function or setting.
IS_NULL_POINTER	Invalid array
IS_OUT_OF_MEMORY	No memory could be allocated.
IS_SUCCESS	Function executed successfully
IS_TIMED_OUT	A timeout occurred. An image capturing process could not be terminated within the allowable period.

Example 1

```
// Enable/disable correction
INT nRet = is_HotPixel(hCam, IS_HOTPIXEL_DISABLE_CORRECTION, NULL, NULL);
nRet = is_HotPixel(hCam, IS_HOTPIXEL_ENABLE_CAMERA_CORRECTION, NULL, NULL);
nRet = is_HotPixel(hCam, IS_HOTPIXEL_ENABLE_SOFTWARE_USER_CORRECTION, NULL, NULL);
nRet = is_HotPixel(hCam, IS_HOTPIXEL_ENABLE_SENSOR_CORRECTION, NULL, NULL);
nRet = is_HotPixel(hCam, IS_HOTPIXEL_DISABLE_SENSOR_CORRECTION, NULL, NULL);
```

Example 2

```
// Read out current mode
INT nMode = 0;
```



```

INT nRet = is_HotPixel(hCam, IS_HOTPIXEL_GET_CORRECTION_MODE,
                      (void*)&nMode, sizeof(nMode));
// Query supported modes
INT nRet = is_HotPixel(hCam, IS_HOTPIXEL_GET_SUPPORTED_CORRECTION_MODES,
                      (void*)&nMode, sizeof(nMode));

```

Example 3

```

// Query user-defined hot pixel list
INT nRet = is_HotPixel(hCam, IS_HOTPIXEL_GET_SOFTWARE_USER_LIST_EXISTS, NULL, NULL);
if (nRet == IS_SUCCESS)
{
    // Query the number of hot pixels in the user-defined list
    INT nNumber = 0;
    nRet = is_HotPixel(hCam, IS_HOTPIXEL_GET_SOFTWARE_USER_LIST_NUMBER,
                      (void*)&nNumber, sizeof(nNumber));
    if (nRet == IS_SUCCESS)
    {
        // Allocate sufficient memory. Each hot pixel needs two WORDS
        // memory space.
        // Additional memory space of one WORD per hot pixel is required for numbering.
        WORD *pList = new WORD[1 + 2 * nNumber];
        nRet = is_HotPixel(hCam, IS_HOTPIXEL_GET_SOFTWARE_USER_LIST,
                          (void*)pList, (1 + 2 * nNumber) * sizeof(WORD));

        // Change a value and save the list.
        // The number of the hot pixel has to be specified in pList[0]
        pList[1] = 100;
        nRet = is_HotPixel(hCam, IS_HOTPIXEL_SET_SOFTWARE_USER_LIST,
                          (void*)pList, (1 + 2 * nNumber) * sizeof(WORD));

        // Delete unneeded list
        delete [] pList;
    }
}

```

Example 4

```
// Save user-defined list to file
char File1[100];
ZeroMemory(File1, sizeof(File1));
strcpy(File1, "c:\\test.txt");

nRet = is_HotPixel(hCam, IS_HOTPIXEL_LOAD_SOFTWARE_USER_LIST, (void*)File1, 0);

nRet = is_HotPixel(hCam, IS_HOTPIXEL_SAVE_SOFTWARE_USER_LIST, (void*)File1, 0);

// Unicode
wchar_t File2[100];
ZeroMemory(File2, sizeof(File2));
wcsncpy(File2, L"c:\\test.txt");

nRet = is_HotPixel(hCam, IS_HOTPIXEL_LOAD_SOFTWARE_USER_LIST_UNICODE, (void*)File2, 0);

nRet = is_HotPixel(hCam, IS_HOTPIXEL_SAVE_SOFTWARE_USER_LIST_UNICODE, (void*)File2, 0);
```

Example 5

```
// Save user-defined list to the camera EEPROM
INT nNumber = 0;
INT nRet = is_HotPixel(hCam, IS_HOTPIXEL_GET_CAMERA_USER_LIST_MAX_NUMBER,
                      (void*)&nNumber, sizeof(nNumber));
if (nRet == IS_SUCCESS)
{
    // Write the maximum number of hot pixels to EEPROM
    WORD *pList = new WORD[1 + 2 * nNumber];
    pList[0] = nNumber;
    for (int i = 0; i < nNumber; i++)
    {
        pList[1 + 2 * i] = x_value;
        pList[2 + 2 * i] = y_value;
    }

    nRet = is_HotPixel(hCam, IS_HOTPIXEL_SET_CAMERA_USER_LIST,
                      (void*)pList, (1 + 2 * nNumber) * sizeof(WORD));

    delete [] pList;



    // Delete user-defined EEPROM list
    nRet = is_HotPixel(hCam, IS_HOTPIXEL_DELETE_CAMERA_USER_LIST, NULL, NULL);
}
}
```

Example 6

```
// Return combined list
INT nNumber = 0;
INT nRet = is_HotPixel(hCam, IS_HOTPIXEL_GET_MERGED_CAMERA_LIST_NUMBER,
                      (void*)&nNumber, sizeof(nNumber));
if (nRet == IS_SUCCESS)
{
    // Allocate sufficient memory. Each hot pixel needs two WORDS
    // memory space.
    // Additional memory space of one WORD per hot pixel is required for numbering.
    WORD *pList = new WORD[1 + 2 * nNumber];
    nRet = is_HotPixel(hCam, IS_HOTPIXEL_GET_MERGED_CAMERA_LIST,
                      (void*)pList, (1 + 2 * nNumber) * sizeof(WORD));

    // Delete unneeded list
    delete [] pList;
}
}
```

4.3.59 is_ImageFile

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_ImageFile (HCAM hCam, UINT nCommand, void* pParam, UINT cbSizeOfParam)
```

Description

`is_ImageFile()` loads and save an image from or to a file. The image must be BMP, JPEG or PNG format. The image is loaded into the active image memory or read-out from the active image memory.

Note

When saving an image [is_FreezeVideo\(\)](#) should not be called with the `IS_DONT_WAIT` parameter, because the image acquisition might not be completed at the time of saving.

The bitmap is stored with the color depth that was used when allocating the image memory (in DIB mode) or that was set for the current color mode (in Direct3D mode). You can save images with a bit depth of more than 8 bit in the PNG format. 12 bit formats are converted into 16 bit. JPEG files are always saved with a color depth of 8 or 24 bits.

Note

In Direct3D or OpenGL mode, overlay data is not saved.

The `nCommand` input parameter is used to select the function mode. The `pParam` input parameter depends on the selected function mode. If you select functions for setting or returning a value, `pParam` contains a pointer to a variable of the `UINT` type. The size of the memory area to which `pParam` refers is specified in the `cbSizeOfParam` input parameter.

Note

The following functions are obsolete by the `is_ImageFile()` function:

- `is_LoadImage()`
- `is_LoadImageMem()`
- `is_SaveImage()`
- `is_SaveImageMem()`
- `is_SaveImageEx()`
- `is_SaveImageMemEx()`

See also: [Obsolete functions](#)

Input parameters

hCam	Camera handle
<input type="checkbox"/> nCommand	

IS_IMAGE_FILE_CMD_LOAD	Loads an image file (bmp, jpg, png) (Example 1) The function can be used with UNICODE file names.
IS_IMAGE_FILE_CMD_SAVE	Saves an image file (bmp, jpg, png) (Example 2) The function can be used with UNICODE file names.
pParam	Pointer to a function parameter, whose function depends on nCommand.
cbSizeOfParam	Size (in bytes) of the memory area to which pParam refers.

Contents of the IMAGE_FILE_PARAMS structure

wchar_t	pwchFileName	Name of the file to be loaded/saved (Unicode). If <code>NULL</code> is passed, the "Open file"/"Save as" dialog opens.
UINT	nFileType	File type to be saved: <ul style="list-style-type: none"> • <code>IS_IMG_BMP</code> • <code>IS_IMG_JPG</code> • <code>IS_IMG_PNG</code>
UINT	nQuality	Sets the image quality for JPEG and PNG (and therefore the compression). The higher the value, the better the quality is: <ul style="list-style-type: none"> • 100 = maximum quality with minimum compression • If the parameter is set to 0, the the default value of 75 is used. For BMP the parameter is ignored.
char	ppcImageMem	When loading: Pointer to an image memory and pointer to the corresponding ID. If both pointers are <code>NULL</code> the image is loaded into the active image memory. If both pointers are valid a new memory is allocated. This memory must be released with is_FreeImageMem() .
UINT	pnImageID	When saving: Pointer to an image memory and pointer to the corresponding ID. If both pointers are <code>NULL</code> the image is saved from the active image memory. If both pointers are valid the corresponding memory is used.
BYTE	reserved[32]	reserved

Return values

<code>IS_FILE_READ_INVALID_BMP_ID</code>	The specified file is not a valid bitmap file.
<code>IS_FILE_READ_OPEN_ERROR</code>	The file cannot be opened.
<code>IS_INVALID_PARAMETER</code>	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
<code>IS_NO_SUCCESS</code>	General error message
<code>IS_NOT_SUPPORTED</code>	The camera model used here does not support this function or setting.
<code>IS_SUCCESS</code>	Function executed successfully

Related functions

- [is_GetImageMem\(\)](#)
- [is_SetImageMem\(\)](#)

Example 1

```
IMAGE_FILE_PARAMS ImageFileParams;

ImageFileParams.pwchFileName = NULL;
ImageFileParams.pnImageID = NULL;
ImageFileParams.ppcImageMem = NULL;
ImageFileParams.nQuality = 0;

// Load bitmap into active memory (with file open dialog)
ImageFileParams.nFileType = IS_IMG_BMP;
INT nRet = is_ImageFile(m_hCam, IS_IMAGE_FILE_CMD_LOAD, (void*)&ImageFileParams,
                       sizeof(ImageFileParams));

// Load jpeg into active memory (with file open dialog)
ImageFileParams.nFileType = IS_IMG_JPG;
nRet = is_ImageFile(m_hCam, IS_IMAGE_FILE_CMD_LOAD, (void*)&ImageFileParams,
                   sizeof(ImageFileParams));

// Alloc image memory and load bitmap into it (without file open dialog)
char* pcMemory = NULL;
UINT nID = 0;
ImageFileParams.pwchFileName = L"c:\\test.bmp";
ImageFileParams.pnImageID = &nID;
ImageFileParams.ppcImageMem = &pcMemory;
ImageFileParams.nFileType = IS_IMG_BMP;
nRet = is_ImageFile(m_hCam, IS_IMAGE_FILE_CMD_LOAD, (void*)&ImageFileParams,
                   sizeof(ImageFileParams));
```

Example 2

```
IMAGE_FILE_PARAMS ImageFileParams;



ImageFileParams.pwchFileName = NULL;
ImageFileParams.pnImageID = NULL;
ImageFileParams.ppcImageMem = NULL;
ImageFileParams.nQuality = 0;

// Save bitmap from active memory to file (with file open dialog)
ImageFileParams.nFileType = IS_IMG_BMP;
INT nRet = is_ImageFile(m_hCam, IS_IMAGE_FILE_CMD_SAVE, (void*)&ImageFileParams,
                       sizeof(ImageFileParams));

// Save jpeg from active memory with quality 80 (without file open dialog)
ImageFileParams.pwchFileName = L"c:\\test.jpg";
ImageFileParams.nFileType = IS_IMG_JPG;
ImageFileParams.nQuality = 80;
nRet = is_ImageFile(m_hCam, IS_IMAGE_FILE_CMD_SAVE, (void*)&ImageFileParams,
                   sizeof(ImageFileParams));

// Save png from special memory with quality 50 (with file open dialog)
ImageFileParams.pwchFileName = NULL; ImageFileParams.pnImageID = &nID; // valid ID
ImageFileParams.ppcImageMem = &pcMemory; // valid buffer
ImageFileParams.nFileType = IS_IMG_PNG;
ImageFileParams.nQuality = 50;
nRet = is_ImageFile(m_hCam, IS_IMAGE_FILE_CMD_SAVE, (void*)&ImageFileParams,
                   sizeof(ImageFileParams));
```

4.3.60 is_ImageFormat

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_ImageFormat (HCAM hCam,
                  UINT nCommand,
                  void *pParam,
                  UINT nSizeOfParam)
```

Description

Using `is_ImageFormat()`, you can query a list of possible image sizes and set a new image format if supported by your DCx Camera model. This is useful for sensors that do not support a free selection of the area of interest or image format. Using the AOI, binning/subsampling or scaling functions, the driver sets the selected image format to achieve the best possible image quality. For a complete list of available image formats see [table "Image formats"](#) below.

The `nCommand` input parameter is used to select the function mode. The `pParam` input parameter depends on the selected function mode. If you select functions for setting or returning a value, `pParam` contains a pointer to a variable of the `UINT` type. The size of the memory area to which `pParam` refers is specified in the `nSizeOfParam` input parameter.

Input parameters

hCam	Camera handle
<input type="checkbox"/> nCommand	
IMGFRMT_CMD_GET_LIST	Returns a list of all image formats supported by the sensor. <input type="checkbox"/> More details You can query the number of entries in the list with <code>IMGFRMT_CMD_GET_NUM_ENTRIES</code> . <ul style="list-style-type: none"> • <code>pParam</code>: Pointer to list of type IMAGE_FORMAT_LIST. The list must be preallocated as specified below. • <code>nSizeOfParam</code>: Size of the list Size of <code>(IMAGE_FORMAT_LIST + (number of list entries - 1) * Size of (IMAGE_FORMAT_INFO))</code>
IMGFRMT_CMD_GET_NUM_ENTRIES	Returns the number of entries in the list. <input type="checkbox"/> More details <ul style="list-style-type: none"> • <code>pParam</code>: Pointer to variable of type <code>UINT</code> returning the number of list entries. • <code>nSizeOfParam</code>: 4
IMGFRMT_CMD_SET_FORMAT	Sets the desired image format. <input type="checkbox"/> More details <ul style="list-style-type: none"> • <code>pParam</code>: Pointer to variable of type <code>UINT</code> passing the format ID of the desired image format. • <code>nSizeOfParam</code>: 4

hCam	Camera handle
IMGFRMT_CMD_GET_ARBITRARY_AOI_SUPPORTED	Returns if the sensor supports a free selection of the area of interest (AOI). <ul style="list-style-type: none"> ▣ More details • pParam: Pointer to variable of type <code>UINT</code> indicating if free AOI selection is supported: 0 = No free AOI supported 1 = Free AOI supported • nSizeOfParam: 4
pParam	Pointer to a function parameter, whose function depends on <code>nCommand</code> .
nSizeOfParam	Size (in bytes) of the memory area to which <code>pParam</code> refers.

Contents of the `IMAGE_FORMAT_LIST` list of image formats

<code>UINT</code>	<code>nSizeOfListEntry</code>	Must be preset with the size of a list entry in bytes
<code>UINT</code>	<code>nNumListElements</code>	Must be preset with the number of list entries (from <code>IMGFRMT_CMD_GET_NUM_ENTRIES</code>)
<code>UINT</code>	<code>nReserved[4]</code>	Reserved
<code>IMAGE_FORMAT_INFO</code>	<code>FormatInfo[0]</code>	First entry in the list. After having been filled by <code>IMGFRMT_CMD_GET_LIST</code> , the list contains additional entries <code>FormatInfo[1]...</code> <code>FormatInfo[nNumListElements-1]</code> .

Contents of the list entry `IMAGE_FORMAT_INFO`

<code>INT</code>	<code>nFormatID</code>	Format ID of the specified image format (see table "Image formats" below)
<code>UINT</code>	<code>nWidth</code>	Width of the area of interest
<code>UINT</code>	<code>nHeight</code>	Height of the area of interest
<code>UINT</code>	<code>nX0</code>	Start point of the area of interest (X)
<code>UINT</code>	<code>nY0</code>	Start point of the area of interest (Y)
<code>UINT</code>	<code>nSupportedCaptureModes</code>	Image capture modes supported for this format (see table below)
<code>UINT</code>	<code>nBinningMode</code>	Binning mode used
<code>UINT</code>	<code>nSubsamplingMode</code>	Subsampling mode used
<code>IS_CHAR</code>	<code>strFormatName[64]</code>	Description of the format
<code>double</code>	<code>dSensorScalerFactor</code>	Scaling factor used (only sensors that support scaling).
<code>UINT</code>	<code>nReserved[24]</code>	Reserved

Possible values for `CAPTUREMODE`

<code>CAPTMODE_SINGLE</code>	Freerun mode, single frame (freerun snap)
<code>CAPTMODE_FREERUN</code>	Freerun mode, continuous (freerun live)
<code>CAPTMODE_TRIGGER_SOFT_SINGLE</code>	Software triggered mode, single frame

CAPTMODE_TRIGGER_SOFT_CONTINUOUS	Software triggered mode, continuous
CAPTMODE_TRIGGER_HW_SINGLE	Hardware triggered mode, single frame
CAPTMODE_TRIGGER_HW_CONTINUOUS	Hardware triggered mode, continuous

For further information on the image capture modes, see also in the Basics: [Operating modes](#) chapter.

Image formats of CMOS sensors

Format ID	Resolution	Name	Camera model		
			DCC1645C	DCC1545M	DCC1240x
1	3264x2448	(8M)			
2	3264x2176	(8M 3:2)			
3	3264x1836	(8M 16:9)			
4	2592x1944	(5M)			
5	2048x1536	(3M)			
6	1920x1080	(Full HD 16:9)			
7	1632x1224	(2M)			
8	1280x960	(1.2M 4:3)	X	X	X
9	1280x720	(HD 16:9)	X	X	X
11	960x480	(WVGA 2:1)	X	X	X
12	800x480	(WVGA)	X	X	X
13	640x480	(VGA)	X	X	X
14	640x360	(VGA 16:9)	X	X	X
15	400x240	(WQVGA)	X	X	X
16	352x288	(CIF)	X	X	X
17	288x352	(CIF Portrait)	X	X	X
18	320x240	(QVGA)	X	X	X
19	240x320	(QVGA Portrait)	X	X	X
20	1600x1200	(UXGA)			
21	3840x2748	(10M)			
22	1920x1080	(Full HD 16:9, HQ)			
23	2560x1920	(5M)			
24	768x576	(CCIR)	X	X	X
25	1280x1024	(1.3M SXGA)	X	X	X
26	2448x2048	(5M)			
27	1024x768	(XGA)	X	X	X
28	1024x1024	(1M)	X	X	
29	800x600	(SVGA)	X	X	X
30	1360x1024	(1.4M 4:3)			

Image formats of CCD sensors

Format ID	Resolution	Name	Camera model	
			DCU-223x	DCU-224x
1	3264x2448	(8M)		
2	3264x2176	(8M 3:2)		
3	3264x1836	(8M 16:9)		
4	2592x1944	(5M)		
5	2048x1536	(3M)		
6	1920x1080	(Full HD 16:9)		
7	1632x1224	(2M)		
8	1280x960	(1.2M 4:3)		X
9	1280x720	(HD 16:9)		X
11	960x480	(WVGA 2:1)	X	X
12	800x480	(WVGA)	X	X
13	640x480	(VGA)	X	X
14	640x360	(VGA 16:9)	X	
15	400x240	(WQVGA)	X	
16	352x288	(CIF)	X	
17	288x352	(CIF Portrait)	X	
18	320x240	(QVGA)	X	
19	240x320	(QVGA Portrait)	X	
20	1600x1200	(UXGA)		
21	3840x2748	(10M)		
22	1920x1080	(Full HD 16:9, HQ)		
23	2560x1920	(5M)		
24	768x576	(CCIR)	X	X
25	1280x1024	(1.3M SXGA)		X
26	2448x2048	(5M)		
27	1024x768	(XGA)	X	X
28	1024x1024	(1M)		X
29	800x600	(SVGA)	X	X
30	1360x1024	(1.4M 4:3)		

Return values

IS_CANT_ADD_TO_SEQUENCE	The image memory is already included in the sequence and cannot be added again.
IS_BAD_STRUCTURE_SIZE	An internal structure has an incorrect size.
IS_CANT_COMMUNICATE_WITH_DRIVER	Communication with the driver failed because no driver has been loaded.

IS_CANT_OPEN_DEVICE	An attempt to initialize or select the camera failed (no camera connected or initialization error).
IS_CAPTURE_RUNNING	A capturing operation is in progress and must be terminated before you can start another one.
IS_DR_CANNOT_CREATE_SURFACE	The image surface or overlay surface could not be created.
IS_DR_CANNOT_CREATE_TEXTURE	The texture could not be created.
IS_DR_CANNOT_CREATE_VERTEX_BUFFER	The vertex buffer could not be created.
IS_DR_DEVICE_OUT_OF_MEMORY	Not enough graphics memory available.
IS_DR_LIBRARY_NOT_FOUND	The DirectRenderer library could not be found.
IS_INVALID_BUFFER_SIZE	The image memory has an inappropriate size to store the image in the desired format.
IS_INVALID_CAMERA_TYPE	The camera type defined in the .ini file does not match the current camera model.
IS_INVALID_CAPTURE_MODE	The function can not be executed in the current camera operating mode (free run, trigger or standby).
IS_INVALID_CAMERA_HANDLE	Invalid camera handle
IS_INVALID_MEMORY_POINTER	Invalid pointer or invalid memory ID
IS_INVALID_MODE	Camera is in standby mode, function not allowed
IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
IS_INVALID_PIXEL_CLOCK	This setting is not available for the currently set pixel clock frequency.
IS_IO_REQUEST_FAILED	An IO request from the uc480 driver failed. Possibly the versions of the uc480.dll (API) and the driver file (uc480_usb.sys) do not match.
IS_NO_ACTIVE_IMG_MEM	No active image memory available. You must set the memory to active using the is_SetImageMem() function or create a sequence using the is_AddToSequence() function.
IS_NO_SUCCESS	General error message
IS_NOT_CALIBRATED	The camera does not contain any calibration data.
IS_NOT_SUPPORTED	The camera model used here does not support this function or setting.
IS_NULL_POINTER	Invalid array
IS_OUT_OF_MEMORY	No memory could be allocated.
IS_SEQUENCE_BUF_ALREADY_LOCKED	The memory could not be locked. The pointer to the buffer is invalid.
IS_SUCCESS	Function executed successfully
IS_TIMED_OUT	A timeout occurred. An image capturing process could not be terminated within the allowable period.
IS_TRIGGER_ACTIVATED	The function cannot be used because the camera is waiting for a trigger signal.

Related functions

- [is_AOI\(\)](#)
- [is_SetBinning\(\)](#)
- [is_SetSubSampling\(\)](#)

Example

```
HIDS hCam;
char strCamFileName[256];
int nRet;

// Get number of available formats and size of list
UINT count;
UINT bytesNeeded = sizeof(IMAGE_FORMAT_LIST);
nRet = is_ImageFormat(hCam, IMGFRMT_CMD_GET_NUM_ENTRIES, &count, 4);
bytesNeeded += (count - 1) * sizeof(IMAGE_FORMAT_INFO);
void* ptr = malloc(bytesNeeded);

// Create and fill list
IMAGE_FORMAT_LIST* pformatList = (IMAGE_FORMAT_LIST*) ptr;
pformatList->nSizeOfListEntry = sizeof(IMAGE_FORMAT_INFO);
pformatList->nNumListElements = count;
nRet = is_ImageFormat(hCam, IMGFRMT_CMD_GET_LIST, pformatList, bytesNeeded);

// Activate trigger mode for capturing high resolution images (USB uEye XS)
nRet = is_StopLiveVideo(hCam, IS_WAIT);
nRet = is_SetExternalTrigger(hCam, IS_SET_TRIGGER_SOFTWARE);



// Prepare for creating image buffers
char* pMem = NULL;
int memID = 0;

// Set each format and then capture an image
IMAGE_FORMAT_INFO formatInfo;
for (int i = 0; i < count; i++)
{
    formatInfo = pformatList->FormatInfo[i];
    int width = formatInfo.nWidth;
    int height = formatInfo.nHeight;

    // Allocate image mem for current format, set format
    nRet = is_AllocImageMem(hCam, width, height, 24, &pMem, &memID);
    nRet = is_SetImageMem(hCam, pMem, memID);
    nRet = is_ImageFormat(hCam, IMGFRMT_CMD_SET_FORMAT, &formatInfo.nFormatID, 4);

    // Capture image
    nRet = is_FreezeVideo(hCam, IS_WAIT);
}
}
```

4.3.61 is_InitCamera

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_InitCamera (HCAM* phf, HWND hWnd)
```

Description

is_InitCamera() starts the driver and establishes the connection to the camera. After successful initialization, this function assigns the camera handle. All subsequent functions require this handle as the first parameter.

When using Direct3D or OpenGL for image display, you can pass a handle to the output window.

Notes

1. Multi-camera environments

When using multiple cameras in parallel operation on a single system, you should assign a unique camera ID to each camera. To initialize or select a camera with is_InitCamera(), the phCam handle must previously have been set to the desired camera ID.

To initialize or select the next available camera without specifying a camera ID, phCam has to be preset with 0.

2. Thread safety

We recommend that you call the following functions exclusively from a single thread in order to avoid unpredictable behavior of the application.

- [is_InitCamera\(\)](#)
- [is_SetDisplayMode\(\)](#)
- [is_ExitCamera\(\)](#)

See also [General: Thread programming](#)

Input parameters

phCam	<p>Pointer to the camera handle</p> <p>When you call this function, the pointer value has the following meaning:</p> <p>0: The first available camera will be initialized or selected.</p> <p>1-254: The camera with the specified camera ID will be initialized or selected.</p>
*phCam IS_USE_DEVICE_ID	<p>The camera is opened using the device ID instead of the camera ID. For details on device ID please refer to the is_GetCameraList() chapter.</p>
*phCam IS_ALLOW_STARTER_FW_UPLOAD	<p>During initialization of the camera, this parameter checks whether a new version of the starter firmware is required. If it is, the new starter firmware is updated automatically (only GigE uEye SE/RE/CP cameras).</p> <p>To ensure backward compatibility of applications, always call is_InitCamera() without the IS_ALLOW_STARTER_FW_UPLOAD parameter first. Only if an error occurs, call the function with this parameter set (see</p>

	Example below).
hWnd	Pointer to the window where the Direct3D image will be displayed If hWnd = NULL, DIB mode will be used for image display.

Return values

IS_ALL_DEVICES_BUSY	All cameras are in use
IS_BAD_STRUCTURE_SIZE	An internal structure has an incorrect size.
IS_CANT_ADD_TO_SEQUENCE	The image memory is already included in the sequence and cannot be added again.
IS_CANT_COMMUNICATE_WITH_DRIVER	Communication with the driver failed because no driver has been loaded.
IS_CANT_OPEN_DEVICE	An attempt to initialize or select the camera failed (no camera connected or initialization error).
IS_CANT_OPEN_REGISTRY	Error opening a Windows registry key
IS_CANT_READ_REGISTRY	Error reading settings from the Windows registry
IS_CAPTURE_RUNNING	A capturing operation is in progress and must be terminated before you can start another one.
IS_CRC_ERROR	A CRC error-correction problem occurred while reading the settings.
IS_DEVICE_ALREADY_PAIRIED	The device is already paired.
IS_DEVICE_NOT_COMPATIBLE	The device is not compatible to the drivers.
IS_DR_CANNOT_CREATE_SURFACE	The image surface or overlay surface could not be created.
IS_DR_CANNOT_CREATE_TEXTURE	The texture could not be created.
IS_DR_CANNOT_CREATE_VERTEX_BUFFER	The vertex buffer could not be created.
IS_DR_DEVICE_OUT_OF_MEMORY	Not enough graphics memory available.
IS_DR_LIBRARY_NOT_FOUND	The DirectRenderer library could not be found.
IS_ERROR_CPU_IDLE_STATES_CONFIGURATION	The configuration of the CPU idle has failed.
IS_FILE_WRITE_OPEN_ERROR	File cannot be opened for writing or reading.
IS_INCOMPATIBLE_SETTING	Because of other incompatible settings the function is not possible.
IS_INVALID_BUFFER_SIZE	The image memory has an inappropriate size to store the image in the desired format.
IS_INVALID_CAMERA_TYPE	The camera type defined in the .ini file does not match the current camera model.
IS_INVALID_CAPTURE_MODE	The function can not be executed in the current camera operating mode (free run, trigger or standby).
IS_INVALID_DEVICE_ID	The device ID is invalid. Valid IDs start from 1 for USB cameras.
IS_INVALID_EXPOSURE_TIME	This setting is not available for the currently set exposure time.
IS_INVALID_CAMERA_HANDLE	Invalid camera handle

IS_INVALID_MEMORY_POINTER	Invalid pointer or invalid memory ID
IS_INVALID_MODE	Camera is in standby mode, function not allowed
IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
IS_INVALID_PIXEL_CLOCK	This setting is not available for the currently set pixel clock frequency.
IS_IO_REQUEST_FAILED	An IO request from the uc480 driver failed. Possibly the versions of the uc480.dll (API) and the driver file (uc480_usb.sys) do not match.
IS_NO_ACTIVE_IMG_MEM	No active image memory available. You must set the memory to active using the <code>is_SetImageMem()</code> function or create a sequence using the <code>is_AddToSequence()</code> function.
IS_NO_IMAGE_MEM_ALLOCATED	The driver could not allocate memory.
IS_NO_IR_FILTER	No IR filter available
IS_NO_SUCCESS	General error message
IS_NOT_CALIBRATED	The camera does not contain any calibration data.
IS_NOT_SUPPORTED	The camera model used here does not support this function or setting.
IS_NULL_POINTER	Invalid array
IS_OUT_OF_MEMORY	No memory could be allocated.
IS_SEQUENCE_BUF_ALREADY_LOCKED	The memory could not be locked. The pointer to the buffer is invalid.
IS_STARTER_FW_UPLOAD_NEEDED	The camera's starter firmware is not compatible with the driver and needs to be updated.
IS_SUCCESS	Function executed successfully
IS_TIMED_OUT	A timeout occurred. An image capturing process could not be terminated within the allowable period.
IS_TRIGGER_ACTIVATED	The function cannot be used because the camera is waiting for a trigger signal.

Related functions

- [is_ExitCamera\(\)](#)
- [is_EnableAutoExit\(\)](#)
- [is_GetCameraList\(\)](#)
- [is_SetCameraID\(\)](#)
- [is_GetCameraInfo\(\)](#)

Example

```
//Open camera with ID 1
HIDS hCam = 1;
INT nRet = is_InitCamera (&hCam, NULL);
```



```
if (nRet != IS_SUCCESS)
{
    //Check if GigE uEye SE needs a new starter firmware
    if (nRet == IS_STARTER_FW_UPLOAD_NEEDED)
    {
        //Calculate time needed for updating the starter firmware
        INT nTime;
        is_GetDuration (hCam, IS_SE_STARTER_FW_UPLOAD, &nTime);
        /*
        e.g. have progress bar displayed in separate thread
        */



        //Upload new starter firmware during initialization
        hCam = hCam | IS_ALLOW_STARTER_FW_UPLOAD;
        nRet = is_InitCamera (&hCam, NULL);

        /*
        end progress bar
        */
    }
}
```

Sample programs

- uc480MultipleCameraScan (C++)
- uc480Console (C++)
- uc480C# Demo (C#)

4.3.62 is_InitEvent

	
USB 2.0 USB 3.0	-

Syntax

```
INT is_InitEvent (HCAM hCam, HANDLE hEv, INT which)
```

Description

is_InitEvent() initializes the event handle for the specified event object. This registers the event object in the uc480 kernel driver.

Attention

Using USB cameras under Windows

The following events require a Windows message loop. This message loop has to be executed by the thread that loads the uc480 API. The message loop is usually provided by the application window. In some cases, the message loop might not be created automatically (e.g. in console applications). In this case you will need to implement the message loop yourself.

This applies to the following uc480 events:

- IS_SET_EVENT_REMOVE
- IS_SET_EVENT_REMOVAL
- IS_SET_EVENT_DEVICE_RECONNECTED
- IS_SET_EVENT_NEW_DEVICE

If no message loop exists, a USB camera will not be automatically detected after reconnecting.

Input parameters

hCam	Camera handle
hEv	Event handle created by the <code>CreateEvent()</code> Windows API function.
which	ID of the event to be initialized (see is_EnableEvent())

Return values

IS_INVALID_CAMERA_HANDLE	Invalid camera handle
IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
IS_NO_SUCCESS	General error message
IS_SUCCESS	Function executed successfully

Related functions

- [is_EnableEvent\(\)](#)
- [is_DisableEvent\(\)](#)
- [is_ExitEvent\(\)](#)

Example



```
HANDLE hEvent = CreateEvent(NULL, FALSE, FALSE, NULL);

//Enable frame event, start image capture and wait for event
is_InitEvent(hCam, hEvent, IS_SET_EVENT_FRAME);
is_EnableEvent(hCam, IS_SET_EVENT_FRAME);
is_FreezeVideo(hCam, IS_DONT_WAIT);
DWORD dwRet = WaitForSingleObject(hEvent, 1000);
if (dwRet == WAIT_TIMEOUT)
{
    /* wait timed out */
}
else if (dwRet == WAIT_OBJECT_0)
{
    /* event signalled */
}
is_DisableEvent(hCam, IS_SET_EVENT_FRAME);
is_ExitEvent(hCam, IS_SET_EVENT_FRAME);
CloseHandle(hEvent);
```

Sample programs

- SimpleLive (C++)
- uc480Event (C++)

4.3.63 is_InitImageQueue

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_InitImageQueue (HCAM hCam, INT nMode)
```

Description

`is_InitImageQueue()` enables the queue mode for existing image memory sequences. New images will be added to the end of the queue on arrival (FIFO principle). The image memory sequence has to be created with [is_AddToSequence\(\)](#) prior to calling `is_InitImageQueue()`. With [is_WaitForNextImage\(\)](#) you can query the pointer and sequence ID of the first (i.e. oldest) image in the sequence.

Note

Image memory sequences can also be used without queue mode. In this case the current image memory has to be queried with [is_GetActSeqBuf\(\)](#) on every frame event. Disadvantage of this proceeding is that at very high frame rates it may happen that additional images arrive between the frame event and accessing/locking the memory. The images arriving in this period will be skipped when you query the current image.

When the queue mode is used (`is_InitImageQueue()`), however, you can be sure to always receive the oldest image which has not yet been queried. In addition, image memories are automatically locked immediately after receiving the image. This prevents images from being overwritten when very high frame rates and few image memories are used.

Input parameters

hCam	Camera handle
nMode	Queue mode. Currently only <code>nMode = 0</code> is supported.

Return values

IS_INVALID_CAMERA_HANDLE	Invalid camera handle
IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
IS_NO_SUCCESS	General error message
IS_SUCCESS	Function executed successfully

Related functions

- [is_ExitImageQueue\(\)](#)
- [is_WaitForNextImage\(\)](#)
- [is_AddToSequence\(\)](#)



Example

```
// A previously initialized camera continuously captures images
// until a timeout or transfer error occurs.
// Note: image memories have to be allocated before this
```

```
is_InitImageQueue (m_hCam, 0);
INT nMemID = 0;
char *pBuffer = NULL;

while (IS_SUCCESS == is_WaitForNextImage(m_hCam, 1000, &pBuffer, &nMemID))
{
    is_SaveImageMem (m_hCam, "image.bmp", pBuffer, nMemID);
    is_UnlockSeqBuf (m_hCam, nMemID, pBuffer);
}
is_ExitImageQueue (m_hCam);
```

4.3.64 is_InquireImageMem

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_InquireImageMem (HCAM hCam, char* pcMem, int nID,
                       int* pnX, int* pnY,
                       int* pnBits, int* pnPitch);
```

Description

`is_InquireImageMem()` reads out the properties of an allocated image memory.

Input parameters

<code>hCam</code>	Camera handle
<code>pMem</code>	Pointer to the starting address of the image memory as allocated by is_AllocImageMem()
<code>NID</code>	ID of the image memory as allocated by is_AllocImageMem()
<code>pnX</code>	Returns the width used to define the image memory. You can also pass <code>NULL</code> instead.
<code>pnY</code>	Returns the height used to define the image memory. You can also pass <code>NULL</code> instead.
<code>pnBits</code>	Returns the bit width used to define the image memory. You can also pass <code>NULL</code> instead.
<code>pnPitch</code>	Returns the line increment of the image memory. You can also pass <code>NULL</code> instead.

Return values



<code>IS_CANT_COMMUNICATE_WITH_DRIVER</code>	Communication with the driver failed because no driver has been loaded.
<code>IS_CANT_OPEN_DEVICE</code>	An attempt to initialize or select the camera failed (no camera connected or initialization error).
<code>IS_INVALID_CAMERA_HANDLE</code>	Invalid camera handle
<code>IS_INVALID_MEMORY_POINTER</code>	Invalid pointer or invalid memory ID
<code>IS_INVALID_PARAMETER</code>	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
<code>IS_IO_REQUEST_FAILED</code>	An IO request from the <code>uc480</code> driver failed. Possibly the versions of the <code>uc480.dll</code> (API) and the driver file (<code>uc480_usb.sys</code>) do not match.
<code>IS_NO_SUCCESS</code>	General error message
<code>IS_SUCCESS</code>	Function executed successfully

Related functions

- [is_AllocImageMem\(\)](#)

- [is_SetImageMem\(\)](#)
- [is_SetAllocatedImageMem\(\)](#)
- [is_GetColorDepth\(\)](#)

4.3.65 is_IO

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_IO(HCAM hCam, UINT nCommand, void* pParam, UINT cbSizeOfParam)
```

Description

With the `is_IO()` function you control all [flash and trigger functions](#) and the additional digital outputs (GPIOs) of some DCx Camera models. For information on GPIO wiring, please refer to the [Electrical Specifications](#) chapter.

Additionally you can toggle the color of the LED on the back of the DCU22x and DCC1240x camera housing.

Note

GPIOs are available only for DCC3240x cameras. The GPIOs are not provided with optocouplers and use TTL/LVCMOS voltages. For information on GPIO wiring, please refer to the [Electrical specifications](#) chapter.

- Rolling shutter cameras:

Using `is_IO()`, you can determine the times required to implement a global flash function for rolling shutter cameras. This way, a rolling shutter camera can also be used as a global shutter camera provided that no ambient light falls on the sensor outside the flash period.

If the exposure time is set too short so that no global flash operation is possible, the function returns `IS_NO_SUCCESS`.

Note

To use a rolling shutter camera with the global start function, first call the `is_SetGlobalShutter()` function. Otherwise, incorrect values will be returned for `Delay` and `Duration`.

- Global shutter cameras:

In freerun mode, the exposure of global shutter cameras is delayed if the exposure time is not set to the maximum value. `is_IO()` determines the required delay in order to synchronize exposure and flash operation. In triggered mode, the return values for delay and flash duration are 0, since no delay is necessary before exposure starts.

For further information, please refer to the chapters Camera basics: [Shutter methods](#), [Digital input/output \(trigger/flash\)](#) and [Operating modes](#).

Attention

Accuracy of flash synchronization

The following parameters have an influence on the camera's internal timing:

- [Image geometry](#) (CMOS and CCD sensors)
- [Pixel clock](#) (CMOS and CCD sensors)
- [Exposure time](#) (CCD sensors)

If you change any of these parameters, you will have to set the flash duration and flash delay parameters once again.

The `nCommand` input parameter is used to select the function mode. The `pParam` input parameter depends on the selected function mode. If you select functions for setting or returning a value, `pParam` contains a pointer to a variable of the `UINT` type. The size of the memory area to which

pParam refers is specified in the cbSizeOfParam input parameter.

Note

The following functions are obsolete by the is_IO() function:

- is_GetGlobalFlashDelays()
- is_SetFlashDelay()
- is_SetFlashStrobe()
- is_SetIO()
- is_SetIOMask()
- is_SetLED()

See also: [Obsolete functions](#)

Input parameters

hCam	Camera handle
<input type="checkbox"/> nCommand	
GPIO	
IS_IO_CMD_GPIOS_GET_SUPPORTED	Returns the supported GPIO ports (Example 1) <ul style="list-style-type: none"> • IO_FLASH_GPIO_PORT_MASK (IO_FLASH_MODE_GPIO_1 IO_FLASH_MODE_GPIO_2)
IS_IO_CMD_GPIOS_GET_SUPPORTED_INPUTS	Returns the supported GPIO inputs (Example 1)
IS_IO_CMD_GPIOS_GET_SUPPORTED_OUTPUTS	Returns the supported GPIO outputs (Example 1)
IS_IO_CMD_GPIOS_GET_DIRECTION	Returns the input/output mask of the GPIOs (Example 2)
IS_IO_CMD_GPIOS_SET_DIRECTION	Set the GPIO on input/output (Example 2) <ul style="list-style-type: none"> • IO_FLASH_MODE_GPIO_1: Sets GPIO 1 as output. • IO_FLASH_MODE_GPIO_2: Sets GPIO 2 as output.
IS_IO_CMD_GPIOS_GET_STATE	Returns the state of the GPIO (High, Low) (Example 2)
IS_IO_CMD_GPIOS_SET_STATE	Sets the state of the GPIOs if they are defined as output (High, Low) (Example 2)
IS_IO_CMD_GPIOS_GET_CONFIGURATION	Returns the configuration of a GPIO port (Example 10)
IS_IO_CMD_GPIOS_SET_CONFIGURATION	Sets the configuration of a GPIO port (Example 11 , Example 12)
Flash	
IS_IO_CMD_FLASH_GET_SUPPORTED_GPIOS	Returns the GPIOs which can be used for flash output (Example 6)
IS_IO_CMD_FLASH_GET_MODE	Returns the current flash mode, see below (Example 7)
IS_IO_CMD_FLASH_SET_MODE	Sets the flash mode (Example 7) <ul style="list-style-type: none"> • IO_FLASH_MODE_OFF: Disables the digital output. • IO_FLASH_MODE_TRIGGER_LO_ACTIVE: Enables

hCam	Camera handle
	<p>the flash strobe in trigger mode. The digital output is set to low level for the flash duration.</p> <ul style="list-style-type: none"> • <code>IO_FLASH_MODE_TRIGGER_HI_ACTIVE</code>: Enables the flash strobe in trigger mode. The digital output is set to high level for the flash duration. • <code>IO_FLASH_MODE_CONSTANT_HIGH</code>: Statically sets the digital output to high level (HIGH). • <code>IO_FLASH_MODE_CONSTANT_LOW</code>: Statically sets the digital output to low level (LOW). • <code>IO_FLASH_MODE_FREERUN_LO_ACTIVE</code>: Enables the flash strobe in freerun mode. The digital output is set to low level for the flash duration. • <code>IO_FLASH_MODE_FREERUN_HI_ACTIVE</code>: Enables the flash strobe in freerun mode. The digital output is set to high level for the flash duration.
<code>IS_IO_CMD_FLASH_GET_GLOBAL_PARAMS</code>	Returns the parameters for the global exposure window (Example 4)
<code>IS_IO_CMD_FLASH_APPLY_GLOBAL_PARAMS</code>	Returns the parameters for the global exposure window and sets them as flash parameters (Example 4)
<code>IS_IO_CMD_FLASH_GET_PARAMS</code>	Returns the current values for flash delay and duration (Example 5)
<code>IS_IO_CMD_FLASH_SET_PARAMS</code>	Sets the current values for flash delay and duration (Example 5)
<code>IS_IO_CMD_FLASH_GET_PARAMS_MIN</code>	Returns the minimum possible values for flash delay and duration (Example 5)
<code>IS_IO_CMD_FLASH_GET_PARAMS_MAX</code>	Returns the maximum possible values for flash delay and duration (Example 5)
<code>IS_IO_CMD_FLASH_GET_PARAMS_INC</code>	Returns the increments for flash delay and duration (Example 5)
<code>IS_IO_CMD_FLASH_GET_GPIO_PARAMS_MIN</code>	Returns the minimum possible parameters for the GPIOs as shorter flash delay and flash duration are possible when using the GPIOs for flash (Example 13).
<code>IS_IO_CMD_FLASH_SET_GPIO_PARAMS</code>	<p>Sets the flash delay and flash duration and allows the minimum values for GPIOs.</p> <p>Attention: For values below 20 μs an unpredictable behavior can occur when flashing is done via the normal flash pin (Example 13).</p>
Pulse-width modulation	
<code>IS_IO_CMD_PWM_GET_SUPPORTED_GPIOs</code>	Returns the GPIOs which can be used for pulse-width modulation (PWM) (Example 6)
<code>IS_IO_CMD_PWM_GET_MODE</code>	Returns the current PWM mode (Example 9)
<code>IS_IO_CMD_PWM_SET_MODE</code>	<p>Sets the current PWM mode (Example 9)</p> <ul style="list-style-type: none"> • <code>IS_FLASH_MODE_PWM</code>: Sets the flash output as output for PWM mode.

hCam	Camera handle
	<ul style="list-style-type: none"> • IO_GPIO_1: Sets GPIO 1 as output. • IO_GPIO_2: Sets GPIO 2 as output.
IS_IO_CMD_PWM_GET_PARAMS	Returns the current values of the PWM parameters (Example 8)
IS_IO_CMD_PWM_SET_PARAMS	Sets the current values of the PWM parameters (Example 8)
IS_IO_CMD_PWM_GET_PARAMS_MIN	Returns the minimum possible values for PWM parameters (Example 8)
IS_IO_CMD_PWM_GET_PARAMS_MAX	Returns the maximum possible values for PWM parameters (Example 8)
IS_IO_CMD_PWM_GET_PARAMS_INC	Returns the increments of the PWM parameters (Example 8)
LED	
IS_IO_CMD_LED_GET_STATE	Returns the state of the LED (Example 3)
IS_IO_CMD_LED_SET_STATE	Sets the state of the LED (Example 3) <ul style="list-style-type: none"> • IO_LED_STATE_1: Sets LED to orange. • IO_LED_STATE_2: Sets LED to green.
IS_IO_CMD_LED_TOGGLE_STATE	Toggles between the LED states (Example 3)
pParam	Pointer to a function parameter, whose function depends on nCommand.
cbSizeOfParam	Size (in bytes) of the memory area to which pParam refers.

Contents of the IO_FLASH_PARAMS structure

INT	s32Delay	Flash delay (in μ s)
UINT	u32Duration	Flash duration (in μ s) If 0 is passed, the flash output will be active until the end of the exposure time. For sensors with Global Start Shutter this is the time until the end of exposure of the first sensor row.

Contents of the IO_PWM_PARAMS structure

double	dblFrequency_Hz	Frequency of the pulse-width modulation (PWM) 1.0...10 000 Hz
double	dbl_DutyCycle	Duty cycle of the pulse-width modulation 0.0...1.0 (1.0 corresponds to 100 %)

Contents of the IO_GPIO_CONFIGURATION structure

UINT	u32Gpio	Sets the GPIO whose configuration is to be read or set (<code>IO_GPIO_1</code> , <code>IO_GPIO_2</code>). So this value must be initialized before the GPIO configuration is read or set.
UINT	u32Caps	When reading the configuration: ORed bitmask of the supported GPIO modes (<code>IS_GPIO_INPUT</code> <code>IS_GPIO_OUTPUT...</code>).
UINT	u32Configuration	<ul style="list-style-type: none"> When reading the configuration: returns the current set configuration When setting the configuration: sets the configuration
UINT	u32State	<ul style="list-style-type: none"> When reading the configuration: returns the current state of the GPIO (0 = Low, 1 = High). When setting the configuration: sets the state of the GPIO (0 = Low, 1 = High).
UINT	u32Reserved[12]	Reserved

Return values

<code>IS_CANT_COMMUNICATE_WITH_DRIVER</code>	Communication with the driver failed because no driver has been loaded.
<code>IS_CANT_OPEN_DEVICE</code>	An attempt to initialize or select the camera failed (no camera connected or initialization error).
<code>IS_INVALID_CAMERA_HANDLE</code>	Invalid camera handle
<code>IS_INVALID_PARAMETER</code>	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
<code>IS_IO_REQUEST_FAILED</code>	An IO request from the uc480 driver failed. Possibly the versions of the <code>uc480.dll</code> (API) and the driver file (<code>uc480_usb.sys</code>) do not match.
<code>IS_NO_SUCCESS</code>	General error message
<code>IS_NOT_SUPPORTED</code>	The camera model used here does not support this

	function or setting.
IS_SUCCESS	Function executed successfully
IS_TRIGGER_ACTIVATED	The function cannot be used because the camera is waiting for a trigger signal.
IS_TRIGGER_NOT_ACTIVATED	The function is not possible as trigger is disabled.

Example 1

```
INT nRet = IS_SUCCESS;

UINT nSupportedIOs = 0;
nRet = is_IO(m_hCam, IS_IO_CMD_GPIOS_GET_SUPPORTED,
             (void*)&nSupportedIOs, sizeof(nSupportedIOs));

UINT nSupportedInputs = 0;
nRet = is_IO(m_hCam, IS_IO_CMD_GPIOS_GET_SUPPORTED_INPUTS,
             (void*)&nSupportedInputs, sizeof(nSupportedInputs));

UINT nSupportedOutputs = 0;
nRet = is_IO(m_hCam, IS_IO_CMD_GPIOS_GET_SUPPORTED_OUTPUTS,
             (void*)&nSupportedOutputs, sizeof(nSupportedOutputs));
```

Example 2

```
INT nRet = IS_SUCCESS;

UINT nDirection = 0;

// Get direction
nRet = is_IO(m_hCam, IS_IO_CMD_GPIOS_GET_DIRECTION,
             (void*)&nDirection, sizeof(nDirection));

// Set GPIO1 and GPIO2 to output
nDirection = IO_GPIO_1 | IO_GPIO_2;
nRet = is_IO(m_hCam, IS_IO_CMD_GPIOS_SET_DIRECTION,
             (void*)&nDirection, sizeof(nDirection));

// Get the current state of the GPIOs
UINT nCurrentState = 0;
nRet = is_IO(m_hCam, IS_IO_CMD_GPIOS_GET_STATE,
             (void*)&nCurrentState, sizeof(nCurrentState));

// Set GPIO1 to high, GPIO2 to low
nCurrentState = IO_GPIO_1;
nRet = is_IO(m_hCam, IS_IO_CMD_GPIOS_SET_STATE,
             (void*)&nCurrentState, sizeof(nCurrentState));
```

Example 3

```
INT nRet = IS_SUCCESS;

UINT nLED = 0;

// Get the current state of the LED
nCurrentState = 0;
nRet = is_IO(m_hCam, IS_IO_CMD_LED_GET_STATE,
             (void*)&nCurrentState, sizeof(nCurrentState));

// Set LED to state 1 (red)
nCurrentState = IO_LED_STATE_1;
nRet = is_IO(m_hCam, IS_IO_CMD_LED_SET_STATE,
             (void*)&nCurrentState, sizeof(nCurrentState));

// Toggle LED state to green
nRet = is_IO(m_hCam, IS_IO_CMD_LED_TOGGLE_STATE, NULL, 0);
```

Example 4

```
INT nRet = IS_SUCCESS;

// Read the global flash params
IO_FLASH_PARAMS flashParams;
```

```
INT nRet = is_IO(m_hCam, IS_IO_CMD_FLASH_GET_GLOBAL_PARAMS,
                (void*)&flashParams, sizeof(flashParams));
if (nRet == IS_SUCCESS)
{
    INT nDelay      = flashParams.s32Delay;
    UINT nDuration = flashParams.u32Duration;
}

// Read the global flash params and set the flash params to these values
INT nRet = is_IO(m_hCam, IS_IO_CMD_FLASH_APPLY_GLOBAL_PARAMS, NULL, 0);
```

Example 5

```
INT nRet = IS_SUCCESS;

IO_FLASH_PARAMS flashParams;

// Get the minimum values for flash delay and flash duration
nRet = is_IO(m_hCam, IS_IO_CMD_FLASH_GET_PARAMS_MIN,
             (void*)&flashParams, sizeof(flashParams));
if (nRet == IS_SUCCESS)
{
    INT nFlashDelayMin    = flashParams.s32Delay;
    UINT nFlashDurationMin = flashParams.u32Duration;
}

// Get the maximum values for flash delay and flash duration
nRet = is_IO(m_hCam, IS_IO_CMD_FLASH_GET_PARAMS_MAX,
             (void*)&flashParams, sizeof(flashParams));
if (nRet == IS_SUCCESS)
{
    INT nFlashDelayMax    = flashParams.s32Delay;
    UINT nFlashDurationMax = flashParams.u32Duration;
}

// Get the increment for flash delay and flash duration
nRet = is_IO(m_hCam, IS_IO_CMD_FLASH_GET_PARAMS_INC,
             (void*)&flashParams, sizeof(flashParams));
if (nRet == IS_SUCCESS)
{
    UINT nFlashDelayInc    = flashParams.s32Delay;
    UINT nFlashDurationInc = flashParams.u32Duration;
}

// Get the current values for flash delay and flash duration
nRet = is_IO(m_hCam, IS_IO_CMD_FLASH_GET_PARAMS,
             (void*)&flashParams, sizeof(flashParams));
if (nRet == IS_SUCCESS)
{
    INT nCurrentFlashDelay    = flashParams.s32Delay;
    UINT nCurrentFlashDuration = flashParams.u32Duration;
}

// Set the current values for flash delay and flash duration
nRet = is_IO(m_hCam, IS_IO_CMD_FLASH_SET_PARAMS,
             (void*)&flashParams, sizeof(flashParams));
```

Example 6

```
INT nRet = IS_SUCCESS;

// Get all GPIOs that can be used as flash output
UINT nGPIOs_Flash = 0;
INT nRet = is_IO(m_hCam, IS_IO_CMD_FLASH_GET_SUPPORTED_GPIO,
                (void*)&nGPIOs_Flash, sizeof(nGPIOs_Flash));

// Get all GPIOs that can be used for the PWM
UINT nGPIOs_PWM = 0;
INT nRet = is_IO(m_hCam, IS_IO_CMD_PWM_GET_SUPPORTED_GPIO,
                (void*)&nGPIOs_PWM, sizeof(nGPIOs_PWM));
```

Example 7

```
INT nRet = IS_SUCCESS;

// Disable flash
UINT nMode = IO_FLASH_MODE_OFF;
nRet = is_IO(m_hCam, IS_IO_CMD_FLASH_SET_MODE, (void*)&nMode, sizeof(nMode));
```



```
// Set the flash to a constant low output
nMode = IO_FLASH_MODE_CONSTANT_LOW;
nRet = is_IO(m_hCam, IS_IO_CMD_FLASH_SET_MODE, (void*)&nMode, sizeof(nMode));

// Set the flash to a high active pulse for each image in the trigger mode
nMode = IO_FLASH_MODE_TRIGGER_HI_ACTIVE;
nRet = is_IO(m_hCam, IS_IO_CMD_FLASH_SET_MODE, (void*)&nMode, sizeof(nMode));

// Get the current flash mode
nRet = is_IO(m_hCam, IS_IO_CMD_FLASH_GET_MODE, (void*)&nMode, sizeof(nMode));
```

Example 8

```
INT nRet = IS_SUCCESS;

IO_PWM_PARAMS m_pwmParams;

// Get the minimum values of the PWM parameters
nRet = is_IO(m_hCam, IS_IO_CMD_PWM_GET_PARAMS_MIN,
             (void*)&m_pwmParams, sizeof(m_pwmParams));
if (nRet == IS_SUCCESS)
{
    double dblFrequencyMin = m_pwmParams.dblFrequency_Hz;
    double dblDutyCycleMin = m_pwmParams.dblDutyCycle;
}

// Get the maximum values of the PWM parameters
nRet = is_IO(m_hCam, IS_IO_CMD_PWM_GET_PARAMS_MAX,
             (void*)&m_pwmParams, sizeof(m_pwmParams));
if (nRet == IS_SUCCESS)
{
    double dblFrequencyMax = m_pwmParams.dblFrequency_Hz;
    double dblDutyCycleMax = m_pwmParams.dblDutyCycle;
}

// Get the increment of the PWM parameters
nRet = is_IO(m_hCam, IS_IO_CMD_PWM_GET_PARAMS_INC,
             (void*)&m_pwmParams, sizeof(m_pwmParams));
if (nRet == IS_SUCCESS)
{
    double dblFrequencyInc = m_pwmParams.dblFrequency_Hz;
    double dblDutyCycleInc = m_pwmParams.dblDutyCycle;
}

// Get the current values of the PWM parameters
nRet = is_IO(m_hCam, IS_IO_CMD_PWM_GET_PARAMS,
             (void*)&m_pwmParams, sizeof(m_pwmParams));
if (nRet == IS_SUCCESS)
{
    double dblFrequency = m_pwmParams.dblFrequency_Hz;
    double dblDutyCycle = m_pwmParams.dblDutyCycle;
}

// Set the current values of the PWM parameters (1 KHz, 50% duty cycle)
m_pwmParams.dblFrequency_Hz = 1000;
m_pwmParams.dblDutyCycle = 0.5;
nRet = is_IO(m_hCam, IS_IO_CMD_PWM_SET_PARAMS,
             (void*)&m_pwmParams, sizeof(m_pwmParams));
```

Example 9

```
INT nRet = IS_SUCCESS;

// Set GPIO1 as PWM output
UINT nMode = IO_GPIO_1;
nRet = is_IO(m_hCam, IS_IO_CMD_PWM_SET_MODE,
             (void*)&nMode, sizeof(nMode));

// Set GPIO1, GPIO2 and the flash pin as PWM output
nMode = IO_GPIO_1 | IO_GPIO_2 | IS_FLASH_MODE_PWM;
nRet = is_IO(m_hCam, IS_IO_CMD_PWM_SET_MODE, (void*)&nMode, sizeof(nMode));

// Get the current PWM mode
nRet = is_IO(m_hCam, IS_IO_CMD_PWM_GET_MODE, (void*)&nMode, sizeof(nMode));
```

Example 10

```

INT nRet = IS_SUCCESS;

IO_GPIO_CONFIGURATION gpioConfiguration;

// Read information about GPIO1
gpioConfiguration.u32Gpio = IO_GPIO_1;

nRet = is_IO(hCam, IS_IO_CMD_GPIOS_GET_CONFIGURATION, (void*)&gpioConfiguration,
            sizeof(gpioConfiguration) );

if (nRet == IS_SUCCESS)
{
    if ((gpioConfiguration.u32Caps & IS_GPIO_PWM) != 0)
    {
        // GPIO1 supports PWM
    }
    if ((gpioConfiguration.u32Caps & IS_GPIO_FLASH) != 0)
    {
        // GPIO1 supports Flash
    }
    if (gpioConfiguration.u32Configuration == IS_GPIO_OUTPUT)
    {
        // GPIO1 is currently configured as output
        if (gpioConfiguration.u32State == 1)
        {
            // GPIO1 is currently output HIGH
        }
    }
}
}

```

Example 11

```

INT nRet = IS_SUCCESS;

IO_GPIO_CONFIGURATION gpioConfiguration;

// Set configuration of GPIO1 (OUTPUT LOW)
gpioConfiguration.u32Gpio = IO_GPIO_1;
gpioConfiguration.u32Configuration = IS_GPIO_OUTPUT;
gpioConfiguration.u32State = 0;

nRet = is_IO(hCam, IS_IO_CMD_GPIOS_SET_CONFIGURATION, (void*)&gpioConfiguration,
            sizeof(gpioConfiguration));

```

Example 12

```

INT nRet = IS_SUCCESS;

IO_GPIO_CONFIGURATION gpioConfiguration;

// Set configuration of GPIO1 (COM-port TX)
// GPIO1 configured as RX is not supported!
gpioConfiguration.u32Gpio = IO_GPIO_1;
gpioConfiguration.u32Configuration = IS_GPIO_COMPORT_TX;

// GPIO2 will be configured as IS_GPIO_COMPORT_RX automatically!
nRet = is_IO(hCam, IS_IO_CMD_GPIOS_SET_CONFIGURATION, (void*)&gpioConfiguration,
            sizeof(gpioConfiguration));

// The following code leads to the same setting
// Set configuration of GPIO2 (COM-port RX)
gpioConfiguration.u32Gpio = IO_GPIO_2;
gpioConfiguration.u32Configuration = IS_GPIO_COMPORT_RX;

// GPIO1 will be configured as IS_GPIO_COMPORT_TX automatically!
nRet = is_IO(hCam, IS_IO_CMD_GPIOS_SET_CONFIGURATION, (void*)&gpioConfiguration,
            sizeof(gpioConfiguration));

```



Example 13

```
INT nRet = IS_SUCCESS;

IO_FLASH_PARAMS flashParams;
// Get the minimum values for the GPIO flash delay and flash duration
nRet = is_IO(m_hCam, IS_IO_CMD_FLASH_GET_GPIO_PARAMS_MIN, (void*)&flashParams,
             sizeof(flashParams));

// Set the minimum values for flash delay and flash duration. Be careful: The normal flash does not work with
nRet = is_IO(m_hCam, IS_IO_CMD_FLASH_SET_GPIO_PARAMS, (void*)&flashParams,
             sizeof(flashParams));
```

4.3.66 is_IsVideoFinish

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_IsVideoFinish (HCAM hCam, INT* pbo)
```

Description

Using `is_IsVideoFinish()`, you can check whether an image has been captured and stored completely in the image memory. This function is helpful if the [is_FreezeVideo\(\)](#) function was called with the `IS_DONT_WAIT` parameter.

By setting the `*pbo==IS_CAPTURE_STATUS` parameter before calling `is_IsVideoFinish()`, you can also check whether a transfer or post-processing error occurred.

Input parameters

hCam	Camera handle
pbo	<p>By setting <code>*pbo != IS_CAPTURE_STATUS</code> before calling the function, <code>pbo</code> contains the following digitizing status:</p> <ul style="list-style-type: none"> • <code>IS_VIDEO_NOT_FINISH</code> = Digitizing of the image is not completed yet. • <code>IS_VIDEO_FINISH</code> = Digitizing of the image is completed. <p>By setting <code>*pbo == IS_CAPTURE_STATUS</code> before calling the function, <code>pbo</code> contains the following digitizing status:</p> <ul style="list-style-type: none"> • <code>IS_VIDEO_NOT_FINISH</code> = Digitizing of the image is not completed yet. • <code>IS_VIDEO_FINISH</code> = Digitizing of the image is completed. • <code>IS_CAPTURE_STATUS</code> = Transfer error or conversion problem (e.g. destination memory is invalid) <p>The parameter <code>IS_CAPTURE_STATUS</code> replaces the previous parameter <code>IS_TRANSFER_FAILED</code>.</p> <p>The parameter <code>IS_TRANSFER_FAILED</code> was moved into the new header file <code>uc480_deprecated.h</code>, which contains all obsolete function definitions and constants. If necessary the header file <code>uc480_deprecated.h</code> can be included in addition to the header file <code>uc480.h</code>.</p>

Return values



<code>IS_INVALID_CAMERA_HANDLE</code>	Invalid camera handle
<code>IS_INVALID_PARAMETER</code>	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
<code>IS_NO_SUCCESS</code>	General error message
<code>IS_SUCCESS</code>	Function executed successfully

Related functions

- [is_FreezeVideo\(\)](#)

- [is_HasVideoStarted\(\)](#)

4.3.67 is_LockSeqBuf

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_LockSeqBuf (HCAM hCam, INT nNum, char* pcMem)
```

Description

`is_LockSeqBuf()` locks write access to an image memory within a sequence. In the capturing process, locked image memories will be skipped in the sequence list of image memories to be used. This way, you can avoid that image data which are required for further processing will be overwritten by newly captured data. Full access to the image memory is still guaranteed. You can lock any number of image memories at the same time.

Using the [is_UnlockSeqBuf\(\)](#) function, you can re-enable write access to the image memory.

Input parameters

hCam	Camera handle
nNum	Number of the image memory to be locked (1...max) or IS_IGNORE_PARAMETER: The image memory will be identified by its starting address only.
pcMem	Starting address of the image memory to be locked

Attention

`nNum` indicates the location in the sequence list, not the memory ID assigned using [is_AllocImageMem\(\)](#).



Return values

IS_BAD_STRUCTURE_SIZE	An internal structure has an incorrect size.
IS_INVALID_CAMERA_HANDLE	Invalid camera handle
IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
IS_NO_SUCCESS	General error message
IS_SUCCESS	Function executed successfully

Related functions

- [is_UnlockSeqBuf\(\)](#)
- [is_AddToSequence\(\)](#)
- [is_SetImageMem\(\)](#)
- [is_SetAllocatedImageMem\(\)](#)

4.3.68 is_LUT

	
USB 2.0 USB 3.0 GigE	USB 2.0 USB 3.0 GigE

Syntax

```
INT is_LUT(HCAM hCam, UINT nCommand, void* pParam, UINT cbSizeOfParams);
```

Description

Using `is_LUT()`, you can enable a hardware or software LUT for uEye cameras. This LUT which will be applied to the image in the camera. A number of predefined LUTs are available. Alternatively, you define your own LUT. It is possible to define a LUT without enabling it at the same time. You can also query the current LUT used by the camera.

Each lookup table (LUT) for the uEye contains modification values for the image brightness and contrast parameters. When a LUT is used, each brightness value in the image will be replaced by a value from the table. LUTs are typically used to enhance the image contrast or the gamma curve. The values must be in the range between 0.0 and 1.0. A linear LUT containing 64 equidistant values between 0.0 and 1.0 has no effect on the image.

For further information on LUTs, please refer to the LUT properties section.



Note: The following functions are obsolete by the `is_LUT()` function:

- `is_GetCameraLUT()`
- `is_SetCameraLUT()`

See also: [Obsolete functions](#)

Input parameters

hCam	Camera handle
<div style="border: 1px solid #ccc; padding: 5px;"> <div style="display: flex; align-items: flex-start;"> <div style="width: 20px; text-align: center; border-right: 1px solid #ccc; padding-right: 5px;"> <input type="checkbox"/> </div> <div style="padding-left: 5px;">nCommand</div> </div> </div>	

IS_LUT_CMD_SET_ENABLED	<p>Enables/disables the LUT (Example 1):</p> <ul style="list-style-type: none"> IS_LUT_DISABLED: LUT is disabled. IS_LUT_ENABLED: LUT is enabled. If no other LUT is defined the linear LUT IS_CAMERA_LUT_IDENTITY is set.
IS_LUT_CMD_SET_MODE	<p>Sets the calculation mode (Example 2):</p> <ul style="list-style-type: none"> IS_LUT_MODE_ID_DEFAULT: Default mode: If supported the LUT/gamma is done on hardware otherwise in software IS_LUT_MODE_ID_FORCE_HARDWARE: LUT/gamma in hardware (if not possible: IS_LUT_STATE_ID_NOT_SUPPORTED) IS_LUT_MODE_ID_FORCE_SOFTWARE: LUT/gamma in software (if not possible: IS_LUT_STATE_ID_NOT_SUPPORTED)
IS_LUT_CMD_GET_STATE	<p>Returns the current status information on the LUT (Example 3):</p> <ul style="list-style-type: none"> IS_LUT_STATE_ID_INACTIVE: LUT is inactive (No gamma and no LUT is set) IS_LUT_STATE_ID_NOT_SUPPORTED: LUT is not supported with the current settings IS_LUT_STATE_ID_HARDWARE_LUT: LUT is calculated in hardware IS_LUT_STATE_ID_HARDWARE_GAMMA: Gamma is calculated in hardware IS_LUT_STATE_ID_HARDWARE_LUTANDGAMMA: LUT and gamma are calculated in hardware IS_LUT_STATE_ID_SOFTWARE_LUT: LUT is calculated in software IS_LUT_STATE_ID_SOFTWARE_GAMMA: Gamma is calculated in software IS_LUT_STATE_ID_SOFTWARE_LUTANDGAMMA: LUT and gamma are calculated in software
IS_LUT_CMD_GET_SUPPORT_INFO	<p>Returns the current support information of the LUT (Example 3)</p>
IS_LUT_CMD_SET_USER_LUT	<p>Sets the user-defined values of the LUT (Example 4)</p>
IS_LUT_CMD_GET_USER_LUT	<p>Returns the user-defined values of the LUT (unchanged) (Example 5)</p>
IS_LUT_CMD_GET_COMPLETE_LUT	<p>Returns the user-defined values of the LUT after the values for gamma, contrast and brightness are included (Example 5)</p>
IS_LUT_CMD_GET_PRESET_LUT	<p>Returns the predefined LUT (Example 6):</p> <ul style="list-style-type: none"> IS_LUT_PRESET_ID_IDENTITY: linear LUT, no image modifications IS_LUT_PRESET_ID_NEGATIVE: inverts the image IS_LUT_PRESET_ID_GLOW1: false-color display of the image IS_LUT_PRESET_ID_GLOW2: false-color display of the image

pParam	Pointer to a function parameter, whose function depends on nCommand.
cbSizeOfParam	Size (in bytes) of the memory area to which pParam refers.

Contents of the IS_LUT_CONFIGURATION_PRESET_64 structure

IS_LUT_PRESET	predefinedLutID	ID of the predefined LUT
IS_LUT_CONFIGURATION_64	lutConfiguration	Predefined LUT

Contents of the IS_LUT_CONFIGURATION_64 structure

DOUBLE	dblValues[3][IS_LUT_64]	Defines a LUT with 64 knee points. This results in 32 sections with a start and end point each.
BOOL	bAllChannelsAreEqual	If TRUE, the same LUT is applied to all three channels.

Contents of the IS_LUT_STATE structure

BOOL	bLUTEnabled	LUT is enabled
INT	nLUTStateID	ID of the LUT status information
INT	nLUTModeID	ID of the LUT mode
INT	nLUTBits	Used bits of the LUT

Contents of the IS_LUT_SUPPORT_INFO structure

BOOL	bSupportLUTHardware	Hardware LUT is supported
BOOL	bSupportLUTSoftware	Software LUT is supported
INT	nBitsHardware	Used bits of the hardware LUT
INT	nBitsSoftware	Used bits of the software LUT
INT	nChannelsHardware	Supported channels for hardware LUT
INT	nChannelsSoftware	Supported channels for software LUT

Return values

IS_FILE_READ_OPEN_ERROR	The file cannot be opened.
IS_FILE_WRITE_OPEN_ERROR	File cannot be opened for writing or reading.
IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
IS_NO_SUCCESS	General error message
IS_NOT_SUPPORTED	The camera model used here does not support this function or setting.
IS_SUCCESS	Function executed successfully

Example 1

```
/* Enable the last set LUT */
IS_LUT_ENABLED_STATE nLutEnabled = IS_LUT_ENABLED;
INT nRet = is_LUT(hCam, IS_LUT_CMD_SET_ENABLED, (void*) &nLutEnabled, sizeof(nLutEnabled));
```

Example 2

```
/* Force LUT and gamma to be included in software (former gamma behavior for cameras with USB3) */
IS_LUT_MODE nLutMode = IS_LUT_MODE_ID_FORCE_SOFTWARE;
INT nRet = is_LUT(hCam, IS_LUT_CMD_SET_MODE, &nLutMode , sizeof(nLutMode));
```

Example 3

```
/* Readout the current LUT state */
IS_LUT_STATE lutState;
INT nRet = is_LUT(hCam, IS_LUT_CMD_GET_STATE, (void*) &lutState, sizeof(lutState));

/* Readout the current LUT support information */
IS_LUT_SUPPORT_INFO lutSupportInfo;
nRet = is_LUT(hCam, IS_LUT_CMD_GET_SUPPORT_INFO, (void*) &lutSupportInfo, sizeof(lutSupportInfo));
```

Example 4

```
/* Set user-defined LUT */
IS_LUT_CONFIGURATION_64 userdefinedLUT;
...
... setze Werte in userdefinedLUT.dblValues[][] ...
...
INT nRet = is_LUT(hCam, IS_LUT_CMD_SET_USER_LUT, (void*) &userdefinedLUT, sizeof(userdefinedLUT));
```

Example 5

```
/* Readout of the set LUT */
/* User-defined LUT (without included gamma) */
IS_LUT_CONFIGURATION_64 userLUT;
INT nRet = is_LUT(hCam, IS_LUT_CMD_GET_USER_LUT, (void*) &userLUT, sizeof(userLUT));

/* Complet LUT (with included gamma) */
IS_LUT_CONFIGURATION_64 completeLUT;
nRet = is_LUT(hCam, IS_LUT_CMD_GET_COMPLETE_LUT, (void*) &completeLUT, sizeof(completeLUT));
```

Example 6

```
/* Set predefined LUT "Glow1" */
/* Firstly load GLOW1 */
IS_LUT_CONFIGURATION_PRESET_64 presetLUT;
presetLUT.predefinedLutID = IS_LUT_PRESET_ID_GLOW1;
INT nRet = is_LUT(hCam, IS_LUT_CMD_GET_PRESET_LUT, (void*) &presetLUT, sizeof(presetLUT));



/* Secondly set GLOW1 */
if (IS_SUCCESS == nRet)
{
    nRet = is_LUT(hCam, IS_LUT_CMD_SET_USER_LUT, (void*) &presetLUT, sizeof(presetLUT));
}
```

Example 7

```
/* Load the "lutFile.xml" file and set LUT */
wchar_t* pFilename = L"lutFile.xml";
INT nRet = is_LUT(hCam, IS_LUT_CMD_LOAD_FILE, (void*) pFilename , NULL);

/* Save the current set LUT into the "lutFile2.xml" file */
wchar_t* pFilename2 = L"lutFile2.xml";
nRet = is_LUT(hCam, IS_LUT_CMD_SAVE_FILE, (void*) pFilename2 , NULL);
```

4.3.69 is_Measure

	
USB 2.0	USB 2.0
USB 3.0	USB 3.0

Syntax

```
INT is_Measure(HCAM hCam, UINT nCommand, void* pParam, UINT cbSizeOfParam)
```

Description

The function allows the measurement of the sharpness in a defined AOI of the current image. To get a sharpness value the edges in the image are evaluated. The sharpness can only be indicated as a relative value as it depends on the edges in the current image. An image with less edges will reach the sharpness value of an image with a lot of edges.

The higher the value, the better the sharpness. The value can be used in comparative measurements to detect changes in the image acquisition of the same object, e.g. caused by readjusted lenses.

The `nCommand` input parameter is used to select the function mode. The `pParam` input parameter depends on the selected function mode. If you select functions for setting or returning a value, `pParam` contains a pointer to a variable of the `UINT` type. The size of the memory area to which `pParam` refers is specified in the `cbSizeOfParam` input parameter.

Input parameters

<code>hCam</code>	Camera handle
<div style="border: 1px solid gray; padding: 2px;"> <input type="checkbox"/> <code>nCommand</code> </div>	
<code>IS_MEASURE_CMD_SHARPNESS_AOI_SET</code>	Sets an AOI in which the sharpness is measured. In the image are up to 5 AOIs possible. These AOIs can also overlap. (Example 1)
<code>IS_MEASURE_CMD_SHARPNESS_AOI_INQUIRE</code>	Returns information of the AOI, e.g. the sharpness(Example 2)
<code>IS_MEASURE_CMD_SHARPNESS_AOI_SET_PRESET</code>	Sets different predefined AOIs in the image (Example 3)
<code>pParam</code>	Pointer to a function parameter, whose function depends on <code>nCommand</code> .
<code>cbSizeOfParam</code>	Size (in bytes) of the memory area to which <code>pParam</code> refers.

Content of the MEASURE_SHARPNESS_AOI_INFO structure

<code>UINT</code>	<code>u32NumberAOI</code>	ID of the AOI
<code>UINT</code>	<code>u32SharpnessValue</code>	Relative sharpness value in the defined AOI
<code>IS_RECT</code>	<code>rcAOI</code>	Position and size of the AOI: <ul style="list-style-type: none"> • <code>s32X</code>: X position • <code>s32Y</code>: Y position • <code>s32Width</code>: AOI width • <code>s32Height</code>: AOI height

Content of the MEASURE_SHARPNESS_AOI_PRESETS enumeration

<code>IS_MEASURE_SHARPNESS_AOI_PRESET_1</code>	Predefined AOI for the sharpness measurement (in each of the four image corner and in the center)
--	---

Return values

IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
IS_NO_ACTIVE_IMG_MEM	No active image memory available. You must set the memory to active using the <code>is_SetImageMem()</code> function or create a sequence using the <code>is_AddToSequence()</code> function.
IS_NO_SUCCESS	General error message
IS_NOT_SUPPORTED	The camera model used here does not support this function or setting.
IS_SUCCESS	Function executed successfully

Example 1

```

INT nRet = IS_SUCCESS;

/* Create info object */
MEASURE_SHARPNESS_AOI_INFO measureSharpnessInfo;

/* Set values of AOI_0: The position and size of the AOI equals the whole image */
measureSharpnessInfo.u32NumberAOI = 0;
measureSharpnessInfo.rcAOI.s32X = 0;
measureSharpnessInfo.rcAOI.s32Y = 0;
measureSharpnessInfo.rcAOI.s32Width = m_s32MaxImageWidth;
measureSharpnessInfo.rcAOI.s32Height = m_s32MaxImageHeight;

/* Set AOI_0 */
nRet = is_Measure(m_hCam, IS_MEASURE_CMD_SHARPNESS_AOI_SET, (void*)&measureSharpnessInfo,
                 sizeof(measureSharpnessInfo));

/* Set values of AOI_1 */
measureSharpnessInfo.u32NumberAOI = 1;
measureSharpnessInfo.rcAOI.s32X = 50;
measureSharpnessInfo.rcAOI.s32Y = 50;
measureSharpnessInfo.rcAOI.s32Width = 200;
measureSharpnessInfo.rcAOI.s32Height = 200;

/* Set AOI_1 */
nRet = is_Measure(m_hCam, IS_MEASURE_CMD_SHARPNESS_AOI_SET, (void*)&measureSharpnessInfo,
                 sizeof(measureSharpnessInfo));

```

Example 2

```

INT nRet = IS_SUCCESS;

/* Create info object */
MEASURE_SHARPNESS_AOI_INFO measureSharpnessInfo;

/* Get values of AOI_0 */
measureSharpnessInfo.u32NumberAOI = 0;



nRet = is_Measure(m_hCam, IS_MEASURE_CMD_SHARPNESS_AOI_INQUIRE, (void*)&measureSharpnessInfo,
                 sizeof(measureSharpnessInfo));
if (nRet == IS_SUCCESS)
{
    UINT s32Sharpness = measureSharpnessInfo.u32SharpnessValue;
}

```

Example 3

```
INT nRet = IS_SUCCESS;  
  
/* Set preset */  
UINT nPreset = IS_MEASURE_SHARPNESS_AOI_PRESET_1;  
  
nRet = is_Measure(m_hCam, IS_MEASURE_CMD_SHARPNESS_AOI_SET_PRESET, (void*)&nPreset,  
                sizeof(nPreset));
```

4.3.70 is_ParameterSet

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_ParameterSet(HCAM hCam, UINT nCommand, void* pParam, UINT cbSizeOfParam)
```

Description

The `is_ParameterSet()` function saves the current camera parameters to a file or to the EEPROM of the camera and loads the parameter set from a file or the EEPROM.

The `nCommand` input parameter is used to select the function mode. The `pParam` input parameter depends on the selected function mode. If you select functions for setting or returning a value, `pParam` contains a pointer to a variable of the `UINT` type. The size of the memory area to which `pParam` refers is specified in the `cbSizeOfParam` input parameter.

Only camera-specific ini files can be loaded. The [uc480_parameter file](#) section in the appendix describes the structure of a uc480 ini file.

Attention

When loading an ini file, make sure that the image size (AOI) and color depth parameters in the ini file match those in the allocated memory. Otherwise, display errors may occur.

Note

The following functions are obsolete by the `is_ParameterSet()` function:

- `is_SaveParameters()`
- `is_LoadParameters()`

See also: [Obsolete functions](#)

Input parameter

hCam	Camera handle
<input type="checkbox"/> nCommand	

IS_PARAMETERSET_CMD_LOAD_EEPROM	<p>Loads a camera parameter set from the EEPROM (Example 1)</p> <p>The parameter sets in the EEPROM of the camera can be loaded via special file names:</p> <ul style="list-style-type: none"> • \\cam\\set1 oder /cam/set1
IS_PARAMETERSET_CMD_LOAD_FILE	<p>Loads a camera parameter set from a file (Example 2)</p> <p>You must pass the path to the ini file as Unicode string. You can pass either a relative or an absolute path. If you pass <code>NULL</code> the "Open file" dialog opens.</p>
IS_PARAMETERSET_CMD_SAVE_EEPROM	<p>Saves a camera parameter set in the EEPROM (Example 3)</p> <p>The parameter sets in the EEPROM of the camera can be saved with special file names:</p> <ul style="list-style-type: none"> • \\cam\\set1 oder /cam/set1
IS_PARAMETERSET_CMD_SAVE_FILE	<p>Saves a camera parameter set in a file (Example 4)</p> <p>You must pass the path to the ini file as Unicode string. You can pass either a relative or an absolute path. If you pass <code>NULL</code> the "Save as" dialog opens.</p>
IS_PARAMETERSET_CMD_GET_NUMBER_SUPPORTED	<p>Returns the number of supported parameter sets in the camera EEPROM (Example 5). At the moment this is "1" for all cameras.</p>
pParam	<p>Pointer to a function parameter, whose function depends on <code>nCommand</code>.</p>
cbSizeOfParam	<p>Size (in bytes) of the memory area to which <code>pParam</code> refers.</p>

Return values

IS_INVALID_CAMERA_TYPE	The camera type defined in the .ini file does not match the current camera model.
IS_NO_SUCCESS	General error message
IS_SUCCESS	Function executed successfully

Related functions

- [is_CameraStatus\(\)](#)

Example 1

```
INT nRet = is_ParameterSet(m_hCam, IS_PARAMETERSET_CMD_LOAD_EEPROM, NULL, NULL);
```

Example 2

```
// Load parameters from file (open filebox)
INT nRet = is_ParameterSet(m_hCam, IS_PARAMETERSET_CMD_LOAD_FILE, NULL, NULL);

// Load parameters from specified file
nRet = is_ParameterSet(m_hCam, IS_PARAMETERSET_CMD_LOAD_FILE, L"file.ini", NULL);
```


Example 3

```
INT nRet = is_ParameterSet(m_hCam, IS_PARAMETERSET_CMD_SAVE_EEPROM, NULL, NULL);
```

Example 4



```
// Save parameters to file (open filebox)
INT nRet = is_ParameterSet(m_hCam, IS_PARAMETERSET_CMD_SAVE_FILE, NULL, NULL);

// Save parameters to specified file
nRet = is_ParameterSet(m_hCam, IS_PARAMETERSET_CMD_SAVE_FILE, L"file.ini", NULL);
```

Example 5

```
// Get the number of supported parameter sets in the camera EEPROM
UINT nNumber;
INT nRet = is_ParameterSet(m_hCam, IS_PARAMETERSET_CMD_GET_NUMBER_SUPPORTED, (void*)
    &nNumber, sizeof(nNumber));
```

4.3.71 is_PixelClock

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_PixelClock(HCAM hCam, UINT nCommand, void* pParam, UINT cbSizeOfParam)
```

Description

The function returns the adjustable pixel clock range sets the pixel clock. Due to an excessive pixel clock for USB cameras, images may get lost during the transfer. If you change the pixel clock on-the-fly, the current image capturing process will be aborted.

The pixel clock limit values can vary, depending on the camera model and operating mode. For detailed information on the pixel clock range of a specific camera model, please refer to the [Camera and sensor data](#) chapter.

The `nCommand` input parameter is used to select the function mode. The `pParam` input parameter depends on the selected function mode. If you select functions for setting or returning a value, `pParam` contains a pointer to a variable of the `UINT` type. The size of the memory area to which `pParam` refers is specified in the `cbSizeOfParam` input parameter.

Note

Note: The following functions are obsolete by the `is_PixelClock()` function:

- `is_SetPixelClock()`
- `is_GetPixelClockRange()`

See also: [Obsolete functions](#)

Input parameter

<code>hCam</code>	Camera handle
<input type="checkbox"/> <code>nCommand</code>	
<code>IS_PIXELCLOCK_CMD_GET_NUMBER</code>	Returns the number of discrete pixel clock which are supported by the camera (Example 1).
<code>IS_PIXELCLOCK_CMD_GET_LIST</code>	Returns the list with discrete pixel clocks.
<code>IS_PIXELCLOCK_CMD_GET_RANGE</code>	Returns the range for the pixel clock (Example 2)
<code>IS_PIXELCLOCK_CMD_GET_DEFAULT</code>	Returns the default pixel clock (Example 4)
<code>IS_PIXELCLOCK_CMD_GET</code>	Returns the current set pixel clock in MHz (Example 3)
<code>IS_PIXELCLOCK_CMD_SET</code>	Sets the pixel clock in MHz (Example 4)
<code>pParam</code>	Pointer to a function parameter, whose function depends on <code>nCommand</code> .
<code>cbSizeOfParam</code>	Size (in bytes) of the memory area to which <code>pParam</code> refers.

Return values

<code>IS_INVALID_MODE</code>	Camera is in standby mode, function not allowed
------------------------------	---

IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
IS_NO_SUCCESS	General error message
IS_NOT_SUPPORTED	The camera model used here does not support this function or setting.
IS_SUCCESS	Function executed successfully

Related functions

- [is_GetFramesPerSecond\(\)](#)
- [is_GetFrameTimeRange\(\)](#)
- [is_Exposure\(\)](#)
- [is_SetOptimalCameraTiming\(\)](#)
- [is_SetFrameRate\(\)](#)
- [is_SetAutoParameter\(\)](#)
- [is_SetBinning\(\)](#)
- [is_SetSubSampling\(\)](#)
- [is_AOI\(\)](#)

Example 1

```

UINT nNumberOfSupportedPixelClocks = 0;
INT nRet = is_PixelClock(hCam, IS_PIXELCLOCK_CMD_GET_NUMBER,
                        (void*)&nNumberOfSupportedPixelClocks,
                        sizeof(nNumberOfSupportedPixelClocks));
if ((nRet == IS_SUCCESS) && (m_nNumberOfSupportedPixelClocks > 0))
{
    // No camera has more than 150 different pixel clocks.
    // Of course, the list can be allocated dynamically
    UINT nPixelClockList[150];
    ZeroMemory(&nPixelClockList, sizeof(nPixelClockList));

    nRet = is_PixelClock(hCam, IS_PIXELCLOCK_CMD_GET_LIST,
                        (void*)nPixelClockList,
                        nNumberOfSupportedPixelClocks * sizeof(UINT));
}

```

Example 2

```

UINT nRange[3];
ZeroMemory(nRange, sizeof(nRange));

// Get pixel clock range
INT nRet = is_PixelClock(hCam, IS_PIXELCLOCK_CMD_GET_RANGE, (void*)nRange, sizeof(nRange));
if (nRet == IS_SUCCESS)
{
    UINT nMin = nRange[0];
    UINT nMax = nRange[1];
    UINT nInc = nRange[2];
}

```

Example 3

```

UINT nPixelClock;

// Get current pixel clock
nRet = is_PixelClock(hCam, IS_PIXELCLOCK_CMD_GET, (void*)&nPixelClock, sizeof(nPixelClock));

```



Example 4

```
UINT nPixelClockDefault;

// Get default pixel clock
INT nRet = is_PixelClock(hCam, IS_PIXELCLOCK_CMD_GET_DEFAULT,
                        (void*)&nPixelClockDefault, sizeof(nPixelClockDefault));

if (nRet == IS_SUCCESS)
{
    // Set this pixel clock
    nRet = is_PixelClock(hCam, IS_PIXELCLOCK_CMD_SET,
                        (void*)&nPixelClockDefault, sizeof(nPixelClockDefault));
}
```

4.3.72 is_ReadEEPROM

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_ReadEEPROM (HCAM hCam, INT Adr, char* pcString, INT Count)
```

Description

Using `is_ReadEEPROM()`, you can read the contents of the camera EEPROM. Besides the hard-coded factory information, the EEPROM of the DCx Camera can hold 64 bytes of user data.

Input parameters

<code>hCam</code>	Camera handle
<code>Adr</code>	Starting address for data reads Value range: 0...63
<code>pcString</code>	Pointer to the buffer for the data to read (min. size = <code>Count</code>)
<code>Count</code>	Number of characters to read

Return values

<code>IS_BAD_STRUCTURE_SIZE</code>	An internal structure has an incorrect size.
<code>IS_CANT_COMMUNICATE_WITH_DRIVER</code>	Communication with the driver failed because no driver has been loaded.
<code>IS_CANT_OPEN_DEVICE</code>	An attempt to initialize or select the camera failed (no camera connected or initialization error).
<code>IS_CRC_ERROR</code>	A CRC error-correction problem occurred while reading the settings.
<code>IS_INVALID_CAMERA_HANDLE</code>	Invalid camera handle
<code>IS_INVALID_PARAMETER</code>	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
<code>IS_IO_REQUEST_FAILED</code>	An IO request from the uc480 driver failed. Possibly the versions of the <code>uc480.dll</code> (API) and the driver file (<code>uc480_usb.sys</code>) do not match.
<code>IS_NO_SUCCESS</code>	General error message
<code>IS_NOT_CALIBRATED</code>	The camera does not contain any calibration data.
<code>IS_NOT_SUPPORTED</code>	The camera model used here does not support this function or setting.
<code>IS_SUCCESS</code>	Function executed successfully



Related functions

- [is_WriteEEPROM\(\)](#)

Example

```
char buffer[64];  
is_ReadEEPROM( hCam, 0x00, buffer, 64 );
```

4.3.73 is_RenderBitmap

	
USB 2.0 USB 3.0	-

Syntax

```
INT is_RenderBitmap (HCAM hCam, INT nMemID, HWND hwnd, INT nMode)
```

Description

Using `is_RenderBitmap()`, you can output an image from an image memory in the specified window. For the display, Windows bitmap functionality is used. The image is displayed in the format you specified when allocating the image memory.

The `bitspixel` parameter of the [is_AllocImageMem\(\)](#) function defines the color depth and display type. RGB16 and RGB15 require the same amount of memory but can be distinguished by the `bitspixel` parameter.

Attention

`is_RenderBitmap()` can render Y8 and RGB formats. For displaying YUV/YCbCr formats please use the [is_DirectRenderer\(\)](#) function (see also [Color and memory formats](#)).

Input parameters

hCam	Camera handle
nMemID	ID of the image memory whose contents is to be displayed
hwnd	Output window handle
<input type="checkbox"/> nMode	
IS_RENDER_NORMAL	The image is rendered normally. It will be displayed in 1:1 scale as stored in the image memory.
IS_RENDER_FIT_TO_WINDOW	The image size is adjusted to fit the output window.
IS_RENDER_DOWNSCALE_1_2	Displays the image at 50 % of its original size.
IS_RENDER_PLANAR_COLOR_RED	Renders the red color component of the planar format in red.
IS_RENDER_PLANAR_COLOR_GREEN	Renders the green color component of the planar format in green.
IS_RENDER_PLANAR_COLOR_BLUE	Renders the blue color component of the planar format in blue.
IS_RENDER_PLANAR_MONO_RED	Renders the red color component of the planar format in gray shades.
IS_RENDER_PLANAR_MONO_GREEN	Renders the green color component of the planar format in gray shades.
IS_RENDER_PLANAR_MONO_BLUE	Renders the blue color component of the planar format in gray shades.
<input type="checkbox"/> The following options can be linked by a logical OR using the nMode parameter:	

IS_RENDER_MIRROR_UPDOWN	Mirrors the displayed image along the horizontal axis.
-------------------------	--

Return values

IS_CANT_COMMUNICATE_WITH_DRIVER	Communication with the driver failed because no driver has been loaded.
IS_CANT_OPEN_DEVICE	An attempt to initialize or select the camera failed (no camera connected or initialization error).
IS_INVALID_CAMERA_HANDLE	Invalid camera handle
IS_INVALID_MEMORY_POINTER	Invalid pointer or invalid memory ID
IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
IS_IO_REQUEST_FAILED	An IO request from the uc480 driver failed. Possibly the versions of the uc480.dll (API) and the driver file (uc480_usb.sys) do not match.
IS_NO_SUCCESS	General error message
IS_SUCCESS	Function executed successfully

Related functions

- [is_AllocImageMem\(\)](#)
- [is_SetColorMode\(\)](#)
- [is_SetDisplayMode\(\)](#)
- [is_DirectRender\(\)](#)

Example



Fit image to window and display it upside down:

```
is_RenderBitmap (hCam, nMemID, hwnd, IS_RENDER_FIT_TO_WINDOW | IS_RENDER_MIRROR_UPDOWN);
```

Sample programs

- SimpleAcquire (C++)
- SimpleLive (C++)

4.3.74 is_ResetToDefault

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_ResetToDefault (HCAM hCam)
```

Description

`is_ResetToDefault()` resets all parameters to the camera-specific defaults as specified by the driver. By default, the camera uses full resolution, a medium speed and color level gain values adapted to daylight exposure.

Input parameters

hCam	Camera handle
------	---------------

Return values



IS_CANT_ADD_TO_SEQUENCE	The image memory is already included in the sequence and cannot be added again.
IS_CANT_COMMUNICATE_WITH_DRIVER	Communication with the driver failed because no driver has been loaded.
IS_CANT_OPEN_DEVICE	An attempt to initialize or select the camera failed (no camera connected or initialization error).
IS_DR_CANNOT_CREATE_SURFACE	The image surface or overlay surface could not be created.
IS_DR_CANNOT_CREATE_TEXTURE	The texture could not be created.
IS_DR_CANNOT_CREATE_VERTEX_BUFFER	The vertex buffer could not be created.
IS_DR_DEVICE_OUT_OF_MEMORY	Not enough graphics memory available.
IS_DR_LIBRARY_NOT_FOUND	The DirectRenderer library could not be found.
IS_INVALID_CAMERA_TYPE	The camera type defined in the .ini file does not match the current camera model.
IS_INVALID_CAPTURE_MODE	The function can not be executed in the current camera operating mode (free run, trigger or standby).
IS_INVALID_CAMERA_HANDLE	Invalid camera handle
IS_INVALID_MEMORY_POINTER	Invalid pointer or invalid memory ID
IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
IS_IO_REQUEST_FAILED	An IO request from the uc480 driver failed. Possibly the versions of the <code>uc480.dll</code> (API) and the driver file (<code>uc480_usb.sys</code>) do not match.
IS_NO_SUCCESS	General error message
IS_NOT_CALIBRATED	The camera does not contain any calibration data.

IS_NOT_SUPPORTED	The camera model used here does not support this function or setting.
IS_NULL_POINTER	Invalid array
IS_OUT_OF_MEMORY	No memory could be allocated.
IS_SEQUENCE_BUF_ALREADY_LOCKED	The memory could not be locked. The pointer to the buffer is invalid.
IS_SUCCESS	Function executed successfully
IS_TIMED_OUT	A timeout occurred. An image capturing process could not be terminated within the allowable period.

Related functions

- [is_ParameterSet\(\)](#)

4.3.75 is_SetAllocatedImageMem

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_SetAllocatedImageMem (HCAM hCam, INT width, INT height, INT bitspixel,
                             char* pcImgMem, int* pid)
```

Description

Using `is_SetAllocatedImageMem()`, you can make a memory allocated by a user the active memory for storing digitized images in it. The allocated memory must be large enough and must always be locked globally.

Depending on the selected image format you need more than one byte per pixel for image memory:

```
unsigned int uBytesPerPixel = bitspixel/8;
if (uImageSize % bitspixel != 0)
{
    uBytesPerPixel++;
}
unsigned int uImageSize = width * height * uBytesPerPixel;
```

You can call the [is_AddToSequence\(\)](#) function to add a memory which was set using `is_SetAllocatedImageMem()` to a sequence.

The address of this memory will be passed to the uc480 driver. For this, you can use the `is_SetAllocatedImageMem()` function. In addition, you need to specify the image size, just as you do when calling [is_AllocImageMem\(\)](#). The returned memory ID is required by other functions for memory access.

The memory area must be removed from the driver management again using the [is_FreeImageMem\(\)](#) function. Please note that this does not release the memory. You then need to make sure that the memory will be released again.

After `is_SetAllocatedImageMem` you must call `is_SetImageMem` or `is_AddToSequence` in order that the image caption can be carried out in the image memory.

Input parameters

hCam	Camera handle
width	Image width
height	Image height
bitspixel	Image color depth (bits per pixel)
pcImgMem	Pointer to the starting address of the allocated memory
pid	Returns the ID of this memory.

Return values

IS_CANT_ADD_TO_SEQUENCE	The image memory is already included in the sequence and cannot be added again.
IS_INVALID_CAMERA_HANDLE	Invalid camera handle
IS_INVALID_MEMORY_POINTER	Invalid pointer or invalid memory ID
IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid

	range or is not supported for this sensor or is not available in this mode.
IS_NO_SUCCESS	General error message
IS_OUT_OF_MEMORY	No memory could be allocated.
IS_SEQUENCE_BUF_ALREADY_LOCKED	The memory could not be locked. The pointer to the buffer is invalid.
IS_SUCCESS	Function executed successfully

Related functions

- [is_AllocImageMem\(\)](#)
- [is_FreeImageMem\(\)](#)
- [is_AddToSequence\(\)](#)
- [is_SetImageMem\(\)](#)
- [is_GetColorDepth\(\)](#)
- [is_GetImgMemPitch\(\)](#)



Example Windows

```
HANDLE hMem = GlobalAlloc(0, uImageSize);
char* pcMem = (char*)GlobalLock(hMem);
INT nRet = is_SetAllocatedImageMem(hCam, uWidth, uHeight, uBitspixel, pcMem, &iMemID);
[...]
nRet = is_FreeImageMem(hCam, pcMem, iMemID);
GlobalUnlock(hMem);
GlobalFree(hMem);
```

Example Linux

```
char* pcMem = (char*)malloc(uImageSize);
int iRet = mlock(pcMem, uImageSize);
INT nRet = is_SetAllocatedImageMem(hCam, uWidth, uHeight, uBitspixel, pcMem, &iMemID);
[...]
nRet = is_FreeImageMem(hCam, pcMem, iMemID);
iRet = munlock(pcMem, uImageSize);
free(pcMem);
```

4.3.76 is_SetAutoParameter

	
USB 2.0 USB 3.0 GigE	USB 2.0 USB 3.0 GigE

Syntax

```
INT is_SetAutoParameter (HCAM hCam, INT param, double* pval1, double* pval2)
```

Description

Using `is_SetAutoParameter()`, you can control the automatic gain, exposure shutter, frame rate and white balance control values.

For further information on automatic control, please refer to the [Automatic image control](#) chapter.

- Control is only active as long as the camera is capturing images.
- A manual change of the exposure time and gain settings disables the auto functions.
- When the **auto exposure shutter** function is enabled, you cannot modify the pixel clock frequency.
- The **auto frame rate** function is only available when the auto shutter control is on. Auto frame rate and auto gain control cannot be used simultaneously.
- The **auto gain** function can only be used for cameras with master gain control. Auto white balance is only available for cameras with hardware RGB gain control.
- The **sensor's internal auto features** are not supported by the sensors of DCx camera models.

Attention

Automatic controls when using very high frame rates

Using very high frame rates can cause that too many control commands are sent to the camera. When using frame rates higher than 100 fps you should increase the value for `IS_SET_AUTO_SKIPFRAMES`. Thus, less image will be used for the automatic controls which takes load off the camera.

Input parameters

hCam	Camera handle
param	Configure auto function
<input type="checkbox"/> Enabling auto functions and querying the status	
IS_SET_ENABLE_AUTO_GAIN	Enables/disables the auto gain control function <input type="checkbox"/> Control parameters pval1 = 1 enables control, 0 disables control
IS_GET_ENABLE_AUTO_GAIN	Returns the current auto gain setting or white level adjustment <input type="checkbox"/> Control parameters pval1: returns the current setting
IS_SET_ENABLE_AUTO_SENSOR_GAIN	Enables/disables the internal auto gain control function of the sensor ¹ <input type="checkbox"/> Control parameters

hCam	Camera handle
	pval1 = 1 enables control, 0 disables control
IS_GET_ENABLE_AUTO_SENSOR_GAIN	Returns the current auto gain setting of the sensor* ¹ <input type="checkbox"/> Control parameters pval1: returns the current setting
IS_SET_ENABLE_AUTO_SHUTTER	Enables/disables the auto exposure shutter function. <input type="checkbox"/> Control parameters pval1 = 1 enables control, 0 disables control
IS_GET_ENABLE_AUTO_SHUTTER	Returns the current auto exposure shutter setting. <input type="checkbox"/> Control parameters pval1: returns the current setting
IS_SET_ENABLE_AUTO_SENSOR_SHUTTER	Enables/disables the sensor's internal auto exposure shutter function* ¹ <input type="checkbox"/> Control parameters pval1 = 1 enables control, 0 disables control
IS_GET_ENABLE_AUTO_SENSOR_SHUTTER	Returns the sensor's current auto exposure shutter setting* ¹ <input type="checkbox"/> Control parameters pval1: returns the current setting
IS_SET_ENABLE_AUTO_WHITEBALANCE	Enables/disables the auto white balance function. <input type="checkbox"/> Control parameters pval1 = 1 enables control, 0 disables control
IS_GET_ENABLE_AUTO_WHITEBALANCE	Returns the current auto white balance setting. <input type="checkbox"/> Control parameters pval1: returns the current setting
IS_SET_ENABLE_AUTO_SENSOR_WHITEBALANCE	Enables/disables the sensor's internal auto white balance function* ¹ <input type="checkbox"/> Control parameters pval1: white balance mode (see is_GetAutoInfo()): <ul style="list-style-type: none"> • WB_MODE_DISABLE • WB_MODE_AUTO • WB_MODE_ALL_PULLIN • WB_MODE_INCANDESCENT_LAMP • WB_MODE_FLUORESCENT_DL • WB_MODE_OUTDOOR_CLEAR_SKY • WB_MODE_OUTDOOR_CLOUDY
IS_GET_ENABLE_AUTO_SENSOR_WHITEBALANCE	Returns the sensor's current auto white balance setting* ¹ <input type="checkbox"/> Control parameters pval1: returns the current setting
IS_SET_ENABLE_AUTO_FRAMERATE	Enables/disables the auto frame rate function.

hCam	Camera handle
	<p>☐ Control parameters pval1 = 1 enables control, 0 disables control</p>
IS_GET_ENABLE_AUTO_FRAMERATE	<p>Returns the current auto frame rate setting.</p> <p>☐ Control parameters pval1: returns the current setting</p>
IS_SET_ENABLE_AUTO_SENSOR_FRAMERATE	<p>Enables/disables the sensor's internal auto frame rate function*¹</p> <p>☐ Control parameters pval1 = 1 enables control, 0 disables control</p>
IS_GET_ENABLE_AUTO_SENSOR_FRAMERATE	<p>Returns the sensor's current auto frame rate setting*¹</p> <p>☐ Control parameters pval1: returns the current setting</p>
<p>☐ Adjusting auto gain control/auto exposure shutter</p>	
IS_SET_AUTO_REFERENCE	<p>Sets the setpoint for auto gain control/auto exposure shutter.</p> <p>☐ Control parameters pval1: defines the setpoint (average brightness of the image); the following rule applies independently of the image bit depth:</p> <ul style="list-style-type: none"> • 0 = black • 128 = 50% gray (default) • 255 = white <p>☐ Note on the sensor's internal control functionality</p> <p>When using the sensor's internal control functionality, you can only used values in a range between [44...235]. The increment in this range is 4. Smaller values are automatically set to 44, larger values to 235.</p>
IS_GET_AUTO_REFERENCE	<p>Returns the set point for auto gain control/auto exposure shutter.</p> <p>☐ Control parameters pval1: returns the current setting</p>
IS_SET_AUTO_GAIN_MAX	<p>Sets the upper limit for auto gain control.</p> <p>☐ Control parameters pval1: valid value for gain (0...100)</p>
IS_GET_AUTO_GAIN_MAX	<p>Returns the upper limit for auto gain control.</p> <p>☐ Control parameters pval1: returns the current setting</p>
IS_SET_AUTO_SHUTTER_MAX	<p>Sets the upper limit for auto exposure shutter.</p> <p>☐ Control parameters pval1: valid exposure value (0 sets the value continuously to max. exposure)</p>

hCam	Camera handle
IS_GET_AUTO_SHUTTER_MAX	Returns the upper limit for auto exposure shutter. <ul style="list-style-type: none"> ☐ Control parameters pval1: returns the current setting
IS_SET_AUTO_BRIGHTNESS_ONCE	Enables/disables the automatic disable for automatic brightness control (gain and exposure time) ² <ul style="list-style-type: none"> ☐ Control parameters pval1 = 1 enables control, 0 disables control
IS_GET_AUTO_BRIGHTNESS_ONCE	Returns the automatic disable status ² <ul style="list-style-type: none"> ☐ Control parameters pval1: returns the current setting
Speed and hysteresis	
IS_SET_AUTO_SPEED	Sets the speed value for the auto function. <ul style="list-style-type: none"> ☐ Control parameters pval1: defines the control speed (0...100)
IS_GET_AUTO_SPEED	Returns the speed value for the auto function. <ul style="list-style-type: none"> ☐ Control parameters pval1: returns the current setting
IS_SET_AUTO_SKIPFRAMES	Sets the number of frames to be skipped during automatic control. <ul style="list-style-type: none"> ☐ Control parameters pval1: defines the number of frames to be skipped during automatic control (default: 4)
IS_GET_AUTO_SKIPFRAMES	Returns the number of frames to be skipped during automatic control. <ul style="list-style-type: none"> ☐ Control parameters pval1: returns the current setting
IS_GET_AUTO_SKIPFRAMES_RANGE	Returns the permissible range for the number of frames to be skipped. <ul style="list-style-type: none"> ☐ Control parameters <ul style="list-style-type: none"> • pval1: returns the minimum permitted value • pval2: returns the maximum permitted value
IS_SET_AUTO_HYSTERESIS	Sets the hysteresis value for auto exposure shutter and auto gain control. <ul style="list-style-type: none"> ☐ Control parameters pval1: defines the hysteresis value (default: 2)
IS_GET_AUTO_HYSTERESIS	Returns the hysteresis value for auto exposure shutter and auto gain control. <ul style="list-style-type: none"> ☐ Control parameters pval1: returns the current setting
IS_GET_AUTO_HYSTERESIS_RANGE	Returns the permissible range for the hysteresis value. <ul style="list-style-type: none"> ☐ Control parameters <ul style="list-style-type: none"> • pval1: returns the minimum permitted value

hCam	Camera handle
	<ul style="list-style-type: none"> • pval2: returns the maximum permitted value
<p>☐ Photometric settings for auto gain control/auto exposure shutter</p>	
IS_SET_SENS_AUTO_SHUTTER_PHOTOM	<p>Sets the photometry mode for auto exposure shutter.</p> <p>☐ Control parameters</p> <p>pval1: defines which fields of view are used for auto exposure shutter (see is_GetAutoInfo()):</p> <ul style="list-style-type: none"> • AS_PM_NONE • AS_PM_SENS_CENTER_WEIGHTED • AS_PM_SENS_CENTER_SPOT • AS_PM_SENS PORTRAIT • AS_PM_SENS LANDSCAPE
IS_GET_SENS_AUTO_SHUTTER_PHOTOM	<p>Returns the photometry mode for auto exposure shutter.</p> <p>☐ Control parameters</p> <p>pval1: returns the current setting</p>
IS_GET_SENS_AUTO_SHUTTER_PHOTOM_DEF	<p>Returns the default photometry mode for auto exposure shutter.</p> <p>☐ Control parameters</p> <p>pval1: returns the default</p>
IS_SET_SENS_AUTO_GAIN_PHOTOM	<p>Sets the photometry mode for auto gain control.</p> <p>☐ Control parameters</p> <p>pval1: defines which fields of view are used for auto gain control (see is_GetAutoInfo()):</p> <ul style="list-style-type: none"> • AG_PM_NONE • AG_PM_SENS_CENTER_WEIGHTED • AG_PM_SENS_CENTER_SPOT • AG_PM_SENS PORTRAIT • AG_PM_SENS LANDSCAPE
IS_GET_SENS_AUTO_GAIN_PHOTOM	<p>Returns the photometry mode for auto gain control.</p> <p>☐ Control parameters</p> <p>pval1: returns the current setting</p>
IS_GET_SENS_AUTO_GAIN_PHOTOM_DEF	<p>Returns the default photometry mode for auto gain control.</p> <p>☐ Control parameters</p> <p>pval1: returns the default</p>
<p>☐ Adjusting auto white balance</p>	
IS_SET_AUTO_WB_OFFSET	<p>Sets the offset values for the red and blue channels.</p> <p>☐ Control parameters</p> <ul style="list-style-type: none"> • pval1: defines the red level offset (-50...50) • pval2: defines the blue level offset (-50...50)
IS_GET_AUTO_WB_OFFSET	<p>Returns the offset values for the red and blue</p>

hCam	Camera handle
	channels. <ul style="list-style-type: none"> ☐ Control parameters <ul style="list-style-type: none"> • pval1: returns the red level offset (-50...50) • pval2: returns the blue level offset (-50...50)
IS_SET_AUTO_WB_GAIN_RANGE	Sets the color gain limits for auto white balance. <ul style="list-style-type: none"> ☐ Control parameters <ul style="list-style-type: none"> • pval1: sets the lowest gain value • pval2: sets the highest gain value
IS_GET_AUTO_WB_GAIN_RANGE	Returns the color gain limits for auto white balance. <ul style="list-style-type: none"> ☐ Control parameters <ul style="list-style-type: none"> • pval1: returns the minimum permitted gain value • pval2: returns the maximum permitted gain value
IS_SET_AUTO_WB_ONCE	Sets automatic disabling of auto white balance ^{*2} <ul style="list-style-type: none"> ☐ Control parameters <ul style="list-style-type: none"> pval1 = 1 enables control, 0 disables control
IS_GET_AUTO_WB_ONCE	Returns the automatic disable status ^{*2} <ul style="list-style-type: none"> ☐ Control parameters <ul style="list-style-type: none"> pval1: returns the current setting
Speed and hysteresis	
IS_SET_AUTO_WB_SPEED	Sets the speed for auto white balance. <ul style="list-style-type: none"> ☐ Control parameters <ul style="list-style-type: none"> pval1: defines the control speed (0...100)
IS_GET_AUTO_WB_SPEED	Returns the speed for auto white balance. <ul style="list-style-type: none"> ☐ Control parameters <ul style="list-style-type: none"> pval1: returns the current setting
IS_SET_AUTO_WB_HYSTERESIS	Sets the hysteresis value for auto white balance. <ul style="list-style-type: none"> ☐ Control parameters <ul style="list-style-type: none"> pval1: defines the hysteresis value (default: 2)
IS_GET_AUTO_WB_HYSTERESIS	Returns the hysteresis value for auto white balance. <ul style="list-style-type: none"> ☐ Control parameters <ul style="list-style-type: none"> pval1: returns the current setting
IS_GET_AUTO_WB_HYSTERESIS_RANGE	Returns the permissible range for the hysteresis value. <ul style="list-style-type: none"> ☐ Control parameters <ul style="list-style-type: none"> • pval1: returns the minimum permitted value • pval2: returns the maximum permitted value
IS_SET_AUTO_WB_SKIPFRAMES	Sets the number of frames to be skipped during automatic control. <ul style="list-style-type: none"> ☐ Control parameters

hCam	Camera handle
	pval1: defines the number (default: 4)
IS_GET_AUTO_WB_SKIPFRAMES	Returns the number of frames to be skipped during automatic control. <input type="checkbox"/> Control parameters pval1: returns the current setting
IS_GET_AUTO_WB_SKIPFRAMES_RANGE	Returns the permissible range for the number of frames to be skipped. <input type="checkbox"/> Control parameters <ul style="list-style-type: none"> • pval1: returns the minimum permitted value • pval2: returns the maximum permitted value
<input type="checkbox"/> Default values for auto gain control/auto exposure shutter	
NULL is passed for the pval1, pval2 parameters.	
IS_DEFAULT_AUTO_BRIGHT_REFERENCE	Default setpoint for auto gain control and auto exposure shutter.
IS_MIN_AUTO_BRIGHT_REFERENCE	Minimum setpoint for auto gain control and auto exposure shutter.
IS_MAX_AUTO_BRIGHT_REFERENCE	Maximum setpoint for auto gain control and auto exposure shutter.
IS_DEFAULT_AUTO_SPEED	Default value for auto speed.
IS_MAX_AUTO_SPEED	Maximum value for auto speed
<input type="checkbox"/> Default values for Auto White Balance	
NULL is passed for the pval1, pval2 parameters.	
IS_MIN_WB_OFFSET	Minimum value for auto white balance offset.
IS_MAX_WB_OFFSET	Maximum value for auto white balance offset.
IS_DEFAULT_AUTO_WB_SPEED	Default value for auto white balance speed.
IS_MIN_AUTO_WB_SPEED	Minimum value for auto white balance speed.
IS_MAX_AUTO_WB_SPEED	Maximum value for auto white balance speed.
pval1	Control parameter, can have a variable value depending on the corresponding auto function
pval2	Control parameter, can have a variable value depending on the corresponding auto function

*1 Not all sensors support this feature (see information box)

*2 Not with use of sensor's internal control functionality

Return values

IS_CANT_COMMUNICATE_WITH_DRIVER	Communication with the driver failed because no driver has been loaded.
IS_CANT_OPEN_DEVICE	An attempt to initialize or select the camera failed (no camera connected or initialization error).
IS_INCOMPATIBLE_SETTING	Because of other incompatible settings the function is not possible.

IS_INVALID_CAMERA_TYPE	The camera type defined in the .ini file does not match the current camera model.
IS_INVALID_CAMERA_HANDLE	Invalid camera handle
IS_INVALID_MODE	Camera is in standby mode, function not allowed
IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
IS_INVALID_WB_BINNING_MODE	Mono binning/mono sub-sampling do not support automatic white balance.
IS_IO_REQUEST_FAILED	An IO request from the uc480 driver failed. Possibly the versions of the uc480.dll (API) and the driver file (uc480_usb.sys) do not match.
IS_NO_SUCCESS	General error message
IS_NOT_CALIBRATED	The camera does not contain any calibration data.
IS_NOT_SUPPORTED	The camera model used here does not support this function or setting.
IS_NULL_POINTER	Invalid array
IS_SUCCESS	Function executed successfully
IS_TIMED_OUT	A timeout occurred. An image capturing process could not be terminated within the allowable period.

Related functions

- [is_ParameterSet\(\)](#)
- [is_GetAutoInfo\(\)](#)
- [is_SetHardwareGain\(\)](#)
- [is_SetHWGainFactor\(\)](#)
- [is_Exposure\(\)](#)
- [is_SetFrameRate\(\)](#)
- [is_AOI\(\)](#)



Examples

```
//Enable auto gain control:
double dEnable = 1;
int ret = is_SetAutoParameter (hCam, IS_SET_ENABLE_AUTO_GAIN, &dEnable, 0);

//Set brightness setpoint to 128:
double nominal = 128;
int ret = is_SetAutoParameter (hCam, IS_SET_AUTO_REFERENCE, &nominal, 0);

//Return shutter control limit:
double maxShutter;
int ret = is_SetAutoParameter (hCam, IS_GET_AUTO_SHUTTER_MAX, &maxShutter, 0);
```

4.3.77 is_SetBinning

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

INT is_SetBinning (HCAM hCam, INT mode)

Description

Using `is_SetBinning()`, you can enable the binning mode both in horizontal and in vertical direction. This way, the image size in the binning direction can be reduced without scaling down the area of interest. Depending on the sensor used, the sensitivity or the frame rate can be increased while binning is enabled.

To enable horizontal and vertical binning at the same time, you can link the horizontal and vertical binning parameters by a logical OR.

The adjustable binning factors of each sensor are listed in the [Camera and sensor data](#) chapter.

Note

Some sensors allow a higher pixel clock setting if binning or subsampling has been activated. If you set a higher pixel clock and then reduce the binning/subsampling factors again, the driver will automatically select the highest possible pixel clock for the new settings.

Attention

Changes to the image geometry or pixel clock affect the value ranges of the frame rate and exposure time. After executing `is_SetBinning()`, calling the following functions is recommended in order to keep the defined camera settings:

- [is_SetFrameRate\(\)](#)
- [is_Exposure\(\)](#)
- If you are using the DCx Camera's flash function: [is_IO\(\)](#)

Note

For the models DCC1240x Binning can be used only combined for the horizontal and the vertical direction. Please see also the information in section [DCC1240x / DCC3240x Application Notes](#).

Input parameters

hCam	Camera handle
mode	
IS_BINNING_DISABLE	Disables binning.
IS_BINNING_2X_VERTICAL	Enables vertical binning with factor 2.
IS_BINNING_3X_VERTICAL	Enables vertical binning with factor 3.
IS_BINNING_4X_VERTICAL	Enables vertical binning with factor 4.
IS_BINNING_5X_VERTICAL	Enables vertical binning with factor 5.
IS_BINNING_6X_VERTICAL	Enables vertical binning with factor 6.
IS_BINNING_8X_VERTICAL	Enables vertical binning with factor 8.
IS_BINNING_16X_VERTICAL	Enables vertical binning with factor 16.

hCam	Camera handle
IS_BINNING_2X_HORIZONTAL	Enables horizontal binning with factor 2.
IS_BINNING_3X_HORIZONTAL	Enables horizontal binning with factor 3.
IS_BINNING_4X_HORIZONTAL	Enables horizontal binning with factor 4.
IS_BINNING_5X_HORIZONTAL	Enables horizontal binning with factor 5.
IS_BINNING_6X_HORIZONTAL	Enables horizontal binning with factor 6.
IS_BINNING_8X_HORIZONTAL	Enables horizontal binning with factor 8.
IS_BINNING_16X_HORIZONTAL	Enables horizontal binning with factor 16.
IS_GET_BINNING	Returns the current setting.
IS_GET_BINNING_FACTOR_VERTICAL	Returns the vertical binning factor.
IS_GET_BINNING_FACTOR_HORIZONTAL	Returns the horizontal binning factor.
IS_GET_SUPPORTED_BINNING	Returns the supported binning modes.
IS_GET_BINNING_TYPE	Indicates whether the camera uses color-proof binning (IS_BINNING_COLOR) or not (IS_BINNING_MONO)

Return values



When used with	Current setting
IS_GET_BINNING	
IS_GET_BINNING_FACTOR_VERTICAL IS_GET_BINNING_FACTOR_HORIZONTAL	Current setting: Returns the current factor as integer value (2, 3, 4, 5, 6, 8, 16)
IS_GET_BINNING_TYPE	Returns IS_BINNING_COLOR if the camera uses color-proof binning; otherwise, IS_BINNING_MONO is returned.
IS_GET_SUPPORTED_BINNING	Returns the supported binning modes linked by logical ORs.
IS_BAD_STRUCTURE_SIZE	An internal structure has an incorrect size.
IS_CANT_ADD_TO_SEQUENCE	The image memory is already included in the sequence and cannot be added again.
IS_CANT_COMMUNICATE_WITH_DRIVER	Communication with the driver failed because no driver has been loaded.
IS_CANT_OPEN_DEVICE	An attempt to initialize or select the camera failed (no camera connected or initialization error).
IS_CAPTURE_RUNNING	A capturing operation is in progress and must be terminated before you can start another one.
IS_DR_CANNOT_CREATE_SURFACE	The image surface or overlay surface could not be created.
IS_DR_CANNOT_CREATE_TEXTURE	The texture could not be created.
IS_DR_CANNOT_CREATE_VERTEX_BUFFER	The vertex buffer could not be created.
IS_DR_DEVICE_OUT_OF_MEMORY	Not enough graphics memory available.
IS_DR_LIBRARY_NOT_FOUND	The DirectRenderer library could not be found.
IS_INVALID_BUFFER_SIZE	The image memory has an inappropriate size to

	store the image in the desired format.
IS_INVALID_CAMERA_TYPE	The camera type defined in the .ini file does not match the current camera model.
IS_INVALID_CAPTURE_MODE	The function can not be executed in the current camera operating mode (free run, trigger or standby).
IS_INVALID_CAMERA_HANDLE	Invalid camera handle
IS_INVALID_MEMORY_POINTER	Invalid pointer or invalid memory ID
IS_INVALID_MODE	Camera is in standby mode, function not allowed
IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
IS_INVALID_PIXEL_CLOCK	This setting is not available for the currently set pixel clock frequency.
IS_IO_REQUEST_FAILED	An IO request from the uc480 driver failed. Possibly the versions of the uc480.dll (API) and the driver file (uc480_usb.sys) do not match.
IS_NO_ACTIVE_IMG_MEM	No active image memory available. You must set the memory to active using the <code>is_SetImageMem()</code> function or create a sequence using the <code>is_AddToSequence()</code> function.
IS_NO_SUCCESS	General error message
IS_NOT_CALIBRATED	The camera does not contain any calibration data.
IS_NOT_SUPPORTED	The camera model used here does not support this function or setting.
IS_NULL_POINTER	Invalid array
IS_OUT_OF_MEMORY	No memory could be allocated.
IS_SEQUENCE_BUF_ALREADY_LOCKED	The memory could not be locked. The pointer to the buffer is invalid.
IS_SUCCESS	Function executed successfully
IS_TIMED_OUT	A timeout occurred. An image capturing process could not be terminated within the allowable period.
IS_TRIGGER_ACTIVATED	The function cannot be used because the camera is waiting for a trigger signal.

Related functions

- [is_SetSubSampling\(\)](#)
- [is_AOI\(\)](#)
- [is_SetImagePos\(\)](#)
- [is_PixelClock\(\)](#)

4.3.78 is_SetCameraID

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_SetCameraID (HCAM hCam, INT nID)
```

Description

Using `is_SetCameraID()`, you can assign a unique camera ID to a camera. Thus, it is possible to access the camera directly with the [is_InitCamera\(\)](#) function.

The camera ID is stored in the non-volatile memory of the camera. The factory default camera ID is 1. The camera ID can also be changed in the uc480 Camera Manager.

Input parameters

hCam	Camera handle
<input type="checkbox"/> nID	
1...254	New camera ID
IS_GET_CAMERA_ID	Returns the current ID.

Return values



When used together with	Current ID
IS_GET_CAMERA_ID	
IS_ACCESS_VIOLATION	An internal error has occurred. The camera ID cannot be changed because the camera is running in the boot-boost mode.
IS_CANT_COMMUNICATE_WITH_DRIVER	Communication with the driver failed because no driver has been loaded.
IS_CANT_OPEN_DEVICE	An attempt to initialize or select the camera failed (no camera connected or initialization error).
IS_INVALID_CAMERA_TYPE	The camera type defined in the .ini file does not match the current camera model.
IS_INVALID_CAMERA_HANDLE	Invalid camera handle
IS_INVALID_MODE	Camera is in standby mode, function not allowed
IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
IS_IO_REQUEST_FAILED	An IO request from the uc480 driver failed. Possibly the versions of the <code>uc480.dll</code> (API) and the driver file (<code>uc480_usb.sys</code>) do not match.
IS_NO_SUCCESS	General error message
IS_NOT_CALIBRATED	The camera does not contain any calibration data.
IS_NOT_SUPPORTED	The camera model used here does not support

	this function or setting.
IS_NULL_POINTER	Invalid array
IS_SUCCESS	Function executed successfully
IS_TIMED_OUT	A timeout occurred. An image capturing process could not be terminated within the allowable period.

Related functions

- [is_InitCamera\(\)](#)
- [is_GetCameraInfo\(\)](#)
- [is_CameraStatus\(\)](#)

4.3.79 is_SetColorConverter

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_SetColorConverter (HCAM hCam, INT ColorMode, INT ConvertMode)
```

Description

Using `is_SetColorConverter()`, you can select the type of Bayer conversion for color cameras. Software conversion is done on the PC. The use of a larger filter mask results in a higher image quality, but increases the computational load. For further information, please refer to the [Camera basics: Color filters](#) chapter.

Note

Software conversion with the large filter mask should only be used for sensors whose green pixels have the same sensitivity. This applies to all DCU22xX CCD cameras.

For all other sensors, we recommend using the standard filter mask.

Attention

While free run mode is active, you cannot change the color conversion type. To do so, you must first stop the capturing process using [is_StopLiveVideo\(\)](#) or set the camera to trigger mode (see [is_SetExternalTrigger\(\)](#)).

Input parameters

hCam	Camera handle
ColorMode	Color mode for which the converter is to be set. For a list of all available color formats and the associated input parameters, see the Appendix: Color and memory formats section.
ConvertMode	Conversion mode selection
IS_CONV_MODE_NONE	No conversion
IS_CONV_MODE_SOFTWARE	Only for monochrome cameras, if you want to add a gamma
IS_CONV_MODE_SOFTWARE_3X3	Software conversion using the standard filter mask (default)
IS_CONV_MODE_SOFTWARE_5X5	Software conversion using a large filter mask
IS_CONV_MODE_HARDWARE_3X3	(Not applicable to DCx Cameras)
IS_CONV_MODE_OPENCL_3X3	Software conversion using the standard filter mask, but conversion is done on the graphic board
IS_CONV_MODE_OPENCL_5X5	(Not applicable to DCx Cameras)

Return values



IS_INVALID_COLOR_FORMAT	Parameter <code>ColorMode</code> invalid or not supported
IS_CANT_COMMUNICATE_WITH_DRIVER	Communication with the driver failed because

	no driver has been loaded.
IS_CANT_OPEN_DEVICE	An attempt to initialize or select the camera failed (no camera connected or initialization error).
IS_CAPTURE_RUNNING	A capturing operation is in progress and must be terminated before you can start another one.
IS_INVALID_COLOR_FORMAT	Invalid color format
IS_INVALID_CAMERA_HANDLE	Invalid camera handle
IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
IS_IO_REQUEST_FAILED	An IO request from the uc480 driver failed. Possibly the versions of the uc480.dll (API) and the driver file (uc480_usb.sys) do not match.
IS_NO_IR_FILTER	No IR filter available
IS_NO_SUCCESS	General error message
IS_NOT_SUPPORTED	The camera model used here does not support this function or setting.
IS_OUT_OF_MEMORY	No memory could be allocated.
IS_SUCCESS	Function executed successfully

Related functions

- [is_GetColorConverter\(\)](#)
- [is_SetColorMode\(\)](#)
- [is_Convert\(\)](#)

4.3.80 is_SetColorCorrection

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_SetColorCorrection (HCAM hCam, INT nEnable, double* factors)
```

Description

For color cameras, `is_SetColorCorrection()` enables color correction in the uc480 driver. This enhances the rendering of colors for cameras with color sensors. Color correction is a digital correction based on a color matrix which is adjusted individually for each sensor.

Note

After changing this parameter, perform manual or automatic white balancing in order to obtain correct color rendering (see also [is_SetAutoParameter\(\)](#)).

Input parameters

<code>hCam</code>	Camera handle
<input type="checkbox"/> <code>nEnable</code>	
<code>IS_CCOR_ENABLE_NORMAL</code>	Enables simple color correction. This parameter replaces <code>IS_CCOR_ENABLE</code> .
<code>IS_CCOR_ENABLE_BG40_ENHANCED</code>	Enables color correction for cameras with optical IR filter glasses of the BG40 type.
<code>IS_CCOR_ENABLE_HQ_ENHANCED</code>	Enables color correction for cameras with optical IR filter glasses of the HQ type.
<code>IS_CCOR_SET_IR_AUTOMATIC</code>	Enables color correction for cameras with optical IR filter glasses. The glass type is set automatically as specified in the camera EEPROM.
<code>IS_CCOR_DISABLE</code>	Disables color correction.
<code>IS_GET_CCOR_MODE</code>	Returns the current setting.
<code>IS_GET_SUPPORTED_CCOR_MODE</code>	Returns all supported color correction modes. See the Return values section.
<code>IS_GET_DEFAULT_CCOR_MODE</code>	Returns the default color correction mode.
<code>factors</code>	Sets the strength of the color correction between 0.0 (no correction) and 1.0 (strong correction).

Return values



When used together with <code>IS_GET_CCOR_MODE</code>	Current setting
When used together with <code>IS_GET_SUPPORTED_CCOR_MODE</code>	When used for color cameras and together with <code>IS_GET_SUPPORTED_CCOR_MODE</code> , this parameter returns the supported values linked by a logical OR: <ul style="list-style-type: none"> <code>IS_CCOR_ENABLE_NORMAL</code> <code>IS_CCOR_ENABLE_BG40_ENHANCED</code>

	<ul style="list-style-type: none"> IS_CCOR_ENABLE_HQ_ENHANCED <p>When used for monochrome cameras, the system returns 0.</p>
When used together with IS_GET_DEFAULT_CCOR_MODE	<p>When used for color cameras and together with IS_GET_DEFAULT_CCOR_MODE, this parameter returns the default color correction mode:</p> <ul style="list-style-type: none"> IS_CCOR_ENABLE_NORMAL IS_CCOR_ENABLE_HQ_ENHANCED <p>When used for monochrome cameras, the system returns 0.</p>
IS_CANT_COMMUNICATE_WITH_DRIVER	Communication with the driver failed because no driver has been loaded.
IS_CANT_OPEN_DEVICE	An attempt to initialize or select the camera failed (no camera connected or initialization error).
IS_INVALID_CAMERA_HANDLE	Invalid camera handle
IS_INVALID_MODE	Camera is in standby mode, function not allowed
IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
IS_IO_REQUEST_FAILED	An IO request from the uc480 driver failed. Possibly the versions of the uc480.dll (API) and the driver file (uc480_usb.sys) do not match.
IS_NO_IR_FILTER	No IR filter available
IS_NO_SUCCESS	General error message
IS_NOT_SUPPORTED	The camera model used here does not support this function or setting.
IS_SUCCESS	Function executed successfully

Related functions

- [is_SetColorConverter\(\)](#)
- [is_SetColorMode\(\)](#)
- [is_SetAutoParameter\(\)](#)

4.3.81 is_SetColorMode

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_SetColorMode (HCAM hCam, INT Mode)
```

Description

`is_SetColorMode()` sets the color mode to be used when image data are saved or displayed by the graphics card. For this purpose, the allocated image memory must be large enough to accommodate the data with the selected color mode. When images are transferred directly to the graphics card memory, make sure that the display settings match the color mode settings. Otherwise, the images will be displayed with altered colors or are not clearly visible.

Notes

Display Modes

This function is only supported in the bitmap (DIB) display mode. Use the [is_SetDisplayMode\(\)](#) function to display other color formats in Direct3D or OpenGL mode.

Bit Depth

Color formats with a bit depth of more than 8 bits per channel are not supported by DCx camera models.

RGB15/16

For the RGB16 and RGB15 data formats, the MSBs of the internal 8-bit R, G and B colors are used.

Input parameters

hCam	Camera handle
Mode: Color mode to be set. For a list of all available color formats and the associated input parameters, see the Appendix: Color and memory formats section.	
IS_CM_MONO16	Grayscale (16), for monochrome and color cameras, LUT/gamma active
IS_CM_MONO12	Grayscale (12), for monochrome and color cameras, LUT/gamma active
IS_CM_MONO8	Grayscale (8), for monochrome and color cameras, LUT/gamma active
IS_CM_SENSOR_RAW16	Raw sensor data (16), for monochrome and color cameras, LUT/gamma active
IS_CM_SENSOR_RAW12	Raw sensor data (12), for monochrome and color cameras, LUT/gamma active
IS_CM_SENSOR_RAW8	Raw sensor data (8), for monochrome and color cameras, LUT/gamma active
IS_CM_RGB12_PACKED	RGB36 (12 12 12), for monochrome and color cameras, LUT/gamma active.
IS_CM_RGB10_PACKED	RGB30 (10 10 10), for monochrome and color cameras, LUT/gamma active

hCam	Camera handle
IS_CM_RGB8_PACKED	RGB24 (8 8 8), for monochrome and color cameras, LUT/gamma active
IS_CM_RGBA12_PACKED	RGB48 (12 12 12), for monochrome and color cameras, LUT/gamma active
IS_CM_RGBA8_PACKED	RGB32 (8 8 8), for monochrome and color cameras, LUT/gamma active
IS_CM_RGBY8_PACKED	RGBY (8 8 8 8), for monochrome and color cameras, LUT/gamma active
IS_CM_BGR12_PACKED	BGR36 (12 12 12), for monochrome and color cameras, LUT/gamma active
IS_CM_BGR10_PACKED	BGR30 (10 10 10), for monochrome and color cameras, LUT/gamma active
IS_CM_BGR8_PACKED	BGR24 (8 8 8), for monochrome and color cameras, LUT/gamma active
IS_CM_BGRA12_PACKED	BGR48 (12 12 12), for monochrome and color cameras, LUT/gamma active
IS_CM_BGRA8_PACKED	BGR32 (8 8 8), for monochrome and color cameras, LUT/gamma active
IS_CM_BGRY8_PACKED	BGRY (8 8 8), for monochrome and color cameras, LUT/gamma active
IS_CM_RGB8_PLANAR	Planar RGB (8) for monochrome and color cameras, LUT/gamma active
IS_CM_BGR565_PACKED	BGR16 (5 6 5), for monochrome and color cameras, LUT/gamma active
IS_CM_BGR5_PACKED	BGR15 (5 5 5), for monochrome and color cameras, LUT/gamma active
IS_CM_UYVY_PACKED	YUV 4:2:2 (8 8), for monochrome and color cameras, LUT/gamma active
IS_CM_CBYCRY_PACKED	YCbCr, 4:2:2 (8 8), for monochrome and color cameras, LUT/gamma active
IS_GET_COLOR_MODE	Returns the current setting.

Return values

When used together with	Current setting
IS_GET_COLOR_MODE	
IS_CANT_COMMUNICATE_WITH_DRIVER	Communication with the driver failed because no driver has been loaded.
IS_CANT_OPEN_DEVICE	An attempt to initialize or select the camera failed (no camera connected or initialization error).
IS_CAPTURE_RUNNING	A capturing operation is in progress and must be terminated before you can start another one.
IS_INVALID_CAMERA_TYPE	The camera type defined in the .ini file does not match the current camera model.
IS_INVALID_COLOR_FORMAT	Invalid color format

IS_INVALID_CAMERA_HANDLE	Invalid camera handle
IS_INVALID_MODE	Camera is in standby mode, function not allowed
IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
IS_IO_REQUEST_FAILED	An IO request from the uc480 driver failed. Possibly the versions of the uc480.dll (API) and the driver file (uc480_usb.sys) do not match.
IS_NO_IR_FILTER	No IR filter available
IS_NO_SUCCESS	General error message
IS_NOT_CALIBRATED	The camera does not contain any calibration data.
IS_NOT_SUPPORTED	The camera model used here does not support this function or setting.
IS_NULL_POINTER	Invalid array
IS_OUT_OF_MEMORY	No memory could be allocated.
IS_SUCCESS	Function executed successfully
IS_TIMED_OUT	A timeout occurred. An image capturing process could not be terminated within the allowable period.

Related functions

- [is_SetDisplayMode\(\)](#)
- [is_SetColorConverter\(\)](#)
- [is_SetColorCorrection\(\)](#)
- [is_GetColorDepth\(\)](#)
- [is_AllocImageMem\(\)](#)
- [is_RenderBitmap\(\)](#)



Note on obsolete parameters

The following parameters for color formats are obsolete. Only the new parameters should be used:

Old parameter	New parameter
IS_SET_CM_RGB32	IS_CM_BGRA8_PACKED
IS_SET_CM_RGB24	IS_CM_BGR8_PACKED
IS_SET_CM_RGB16	IS_CM_BGR565_PACKED
IS_SET_CM_RGB15	IS_CM_BGR555_PACKED
IS_SET_CM_Y8	IS_CM_MONO8
IS_SET_CM_BAYER	IS_CM_BAYER_RG8
IS_SET_CM_UYVY	IS_CM_UYVY_PACKED
IS_SET_CM_UYVY_MONO	IS_CM_UYVY_MONO_PACKED
IS_SET_CM_UYVY_BAYER	IS_CM_UYVY_BAYER_PACKED
IS_SET_CM_CBYCRY	IS_CM_CBYCRY_PACKED
IS_SET_CM_RGBY	IS_CM_BGRY8_PACKED
IS_SET_CM_RGB30	IS_CM_BGR10V2_PACKED

Old parameter	New parameter
IS_SET_CM_Y12	IS_CM_MONO12
IS_SET_CM_BAYER12	IS_CM_BAYER_RG12
IS_SET_CM_Y16	IS_CM_MONO16
IS_SET_CM_BAYER16	IS_CM_BAYER_RG16
IS_CM_BGR10V2_PACKED	IS_CM_BGR10_PACKED
IS_CM_RGB10V2_PACKED	IS_CM_RGB10_PACKED
IS_CM_BGR555_PACKED	IS_CM_BGR5_PACKED
IS_CM_BAYER_RG8	IS_CM_SENSOR_RAW8
IS_CM_BAYER_RG12	IS_CM_SENSOR_RAW12
IS_CM_BAYER_RG16	IS_CM_SENSOR_RAW16

4.3.82 is_SetDisplayMode

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_SetDisplayMode (HCAM hCam, INT Mode)
```

Description

Using `is_SetDisplayMode()`, you can set the way in which images will be displayed on the screen.

For live videos including overlays, you can use the Direct3D or OpenGL mode. These modes are not supported by all graphics cards. The graphics card must have sufficient extended memory because the overlay mode requires additional memory up to the size needed for the current screen resolution.

For further information on the display modes of the DCx camera, see the [How to proceed: Image display](#) section.

Attention

The Direct3D display mode is not available on Linux operating systems.

Note

We recommend that you call the following functions exclusively from a single thread in order to avoid unpredictable behaviour of the application.

- [is_InitCamera\(\)](#)
- [is_SetDisplayMode\(\)](#)
- [is_ExitCamera\(\)](#)

See also [Programming: Thread programming](#)

Input parameters

hCam	Camera handle
Mode	
IS_SET_DM_DIB	Captures an image in system memory (RAM). Using is_RenderBitmap() , you can define the image display (default).
IS_SET_DM_DIRECT3D	Image display in Direct3D mode
IS_SET_DM_DIRECT3D IS_SET_DM_MONO	Monochrome image display in Direct3D mode
IS_SET_DM_DIRECT3D IS_SET_DM_BAYER	Raw Bayer format image display in Direct3D mode
IS_SET_DM_OPENGL	Image display in OpenGL mode
IS_SET_DM_OPENGL IS_SET_DM_MONO	Monochrome image display in OpenGL mode
IS_SET_DM_OPENGL IS_SET_DM_BAYER	Raw Bayer format image display in OpenGL mode
IS_GET_DISPLAY_MODE	Returns the current setting.

Attention

The new Direct3D mode completely replaces the "BackBuffer" and "Overlay Surface" display modes from DirectDraw. It is advisable not to use these modes any longer (see also [Obsolete functions](#)). To activate the obsolete modes, do the following:

IS_SET_DM_DIRECTDRAW IS_SET_DM_BACKBUFFER	Image display in DirectDraw BackBuffer mode
IS_SET_DM_DIRECTDRAW IS_SET_DM_ALLOW_OVERLAY	Image display in DirectDraw Overlay Surface mode
IS_SET_DM_ALLOW_SCALING	Real-time scaling in Overlay Surface mode

Return values

When used with	Current setting
IS_GET_DISPLAY_MODE	
IS_CANT_ADD_TO_SEQUENCE	The image memory is already included in the sequence and cannot be added again.
IS_CANT_COMMUNICATE_WITH_DRIVER	Communication with the driver failed because no driver has been loaded.
IS_CANT_OPEN_DEVICE	An attempt to initialize or select the camera failed (no camera connected or initialization error).
IS_CAPTURE_RUNNING	A capturing operation is in progress and must be terminated before you can start another one.
IS_DR_CANNOT_CREATE_SURFACE	The image surface or overlay surface could not be created.
IS_DR_CANNOT_CREATE_TEXTURE	The texture could not be created.
IS_DR_CANNOT_CREATE_VERTEX_BUFFER	The vertex buffer could not be created.
IS_DR_DEVICE_OUT_OF_MEMORY	Not enough graphics memory available.
IS_DR_LIBRARY_NOT_FOUND	The DirectRenderer library could not be found.
IS_INVALID_CAMERA_TYPE	The camera type defined in the .ini file does not match the current camera model.
IS_INVALID_COLOR_FORMAT	Invalid color format
IS_INVALID_CAMERA_HANDLE	Invalid camera handle
IS_INVALID_MEMORY_POINTER	Invalid pointer or invalid memory ID
IS_INVALID_MODE	Camera is in standby mode, function not allowed
IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
IS_IO_REQUEST_FAILED	An IO request from the uc480 driver failed. Possibly the versions of the uc480.dll (API) and the driver file (uc480_usb.sys) do not match.
IS_NO_IR_FILTER	No IR filter available
IS_NO_SUCCESS	General error message
IS_NOT_CALIBRATED	The camera does not contain any calibration data.

IS_NOT_SUPPORTED	The camera model used here does not support this function or setting.
IS_NULL_POINTER	Invalid array
IS_OUT_OF_MEMORY	No memory could be allocated.
IS_SEQUENCE_BUF_ALREADY_LOCKED	The memory could not be locked. The pointer to the buffer is invalid.
IS_SUCCESS	Function executed successfully
IS_TIMED_OUT	A timeout occurred. An image capturing process could not be terminated within the allowable period.

Related functions

- [is_RenderBitmap\(\)](#)
- [is_SetColorMode\(\)](#)
- [is_DirectRenderer\(\)](#)

Example

```
is_SetDisplayMode (hCam, Mode);
```



```
//Bitmap mode (images are digitized and stored in system memory):
```

```
Mode = IS_SET_DM_DIB
```

```
//Direct3D mode
```

```
Mode = IS_SET_DM_DIRECT3D
```

4.3.83 is_SetDisplayPos

	
USB 2.0 USB 3.0	-

Syntax

```
INT is_SetDisplayPos (HCAM hCam, INT x, INT y)
```

Description

`is_SetDisplayPos()` allows you to move an area of interest when rendering images using [is_RenderBitmap\(\)](#). The function moves the camera image by the selected offset within the output window. The image memory remains unchanged.

Note

To set the size and position of an area of interest in memory, use the [is_AOI\(\)](#) functions.

Input parameters

hCam	Camera handle
x	Offset in x direction, measured from the top left corner of the output window
y	Offset in y direction, measured from the top left corner of the output window



Return values

IS_INVALID_CAMERA_HANDLE	Invalid camera handle
IS_NO_SUCCESS	General error message
IS_SUCCESS	Function executed successfully

Related functions

- [is_AOI\(\)](#)
- [is_RenderBitmap\(\)](#)
- [is_SetDisplayMode\(\)](#)

4.3.84 is_SetErrorReport

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_SetErrorReport (HCAM hCam, INT Mode)
```

Description

Using `is_SetErrorReport()`, you can enable/disable error event logging. If error reporting is enabled, errors will automatically be displayed in a dialog box. Cancelling the dialog box disables the error report. Even with disabled error reporting, you can still query errors using the [is_GetError\(\)](#) function.

Note

`is_SetErrorReport()` can be called before calling [is_InitCamera\(\)](#).

You only need to enable the `is_SetErrorReport()` function once for all cameras in the application.

Input parameters

hCam	Camera handle Or 0 if no camera has been initialized yet
Mode	
IS_DISABLE_ERR_REP	Disables error reporting.
IS_ENABLE_ERR_REP	Enables error reporting.
IS_GET_ERR_REP_MODE	Current status of error reporting.



Return values

When used with IS_GET_ERR_REP_MODE	Current setting
IS_SUCCESS	Function executed successfully

Related functions

- [is_GetError\(\)](#)
- [is_CaptureStatus\(\)](#)
- [is_CameraStatus\(\)](#)

4.3.85 is_SetExternalTrigger

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_SetExternalTrigger (HCAM hCam, INT nTriggerMode)
```

Description

Using `is_SetExternalTrigger()`, you can activate the trigger mode. If the camera is in standby mode, it quits this mode and activates trigger mode.

In hardware trigger mode, image capture is delayed for each function call until the selected trigger event has occurred.

In software trigger mode, an image is captured immediately when [is_FreezeVideo\(\)](#) is called, or a continuous triggered capture is started when [is_CaptureVideo\(\)](#) is called. In hardware trigger mode, you can use the [is_ForceTrigger\(\)](#) command to trigger an image capture even if no electric signal is present.

When you disable the trigger functionality, you can query the signal level at the trigger input. This option causes the camera to change to freerun mode.

For further information on the image capture modes of the DCx camera, see [How to proceed: Image capture](#).

Input parameters

hCam	Camera handle	
<input type="checkbox"/> nTriggerMode		
	Trigger mode	Trigger event
IS_SET_TRIGGER_OFF	Off	-
IS_SET_TRIGGER_HI_LO	Hardware trigger	Falling signal edge
IS_SET_TRIGGER_LO_HI	Hardware trigger	Rising signal edge
IS_SET_TRIGGER_PRE_HI_LO	(Not supported by DCx Cameras)	
IS_SET_TRIGGER_PRE_LO_HI	(Not supported by DCx Cameras)	
IS_SET_TRIGGER_HI_LO_SYNC	(Not supported by DCx Cameras)	
IS_SET_TRIGGER_LO_HI_SYNC	(Not supported by DCx Cameras)	
IS_SET_TRIGGER_SOFTWARE	Software trigger	Call of is_FreezeVideo() (single frame mode) Call of is_CaptureVideo() (continuous mode)
IS_GET_EXTERNALTRIGGER	Returns the trigger mode setting	
IS_GET_TRIGGER_STATUS	Returns the current signal level at the trigger input	
IS_GET_SUPPORTED_TRIGGER_MODE	Returns the supported trigger modes	

Return values

When used with IS_GET_EXTERNALTRIGGER	Returns the current setting
When used with IS_GET_TRIGGER_STATUS	Returns the current signal level at the trigger input
When used with IS_GET_SUPPORTED_TRIGGER_MODE	Returns the supported modes linked by logical ORs
IS_INVALID_CAPTURE_MODE	The function can not be executed in the current camera operating mode (free run, trigger or standby).
IS_INVALID_CAMERA_HANDLE	Invalid camera handle
IS_NO_SUCCESS	General error message
IS_SUCCESS	Function executed successfully

Related functions

- [is_CaptureVideo\(\)](#)
- [is_FreezeVideo\(\)](#)
- [is_ForceTrigger\(\)](#)
- [is_SetTriggerCounter\(\)](#)
- [is_SetTriggerDelay\(\)](#)
- [is_IO\(\)](#)

Example

```
//Enable trigger mode and set high-active flash mode.
is_SetExternalTrigger(hCam, IS_SET_TRIGGER_SOFTWARE);



// Set the flash to a high active pulse for each image in the trigger mode
UINT nMode = IO_FLASH_MODE_TRIGGER_HI_ACTIVE;
is_IO(m_hCam, IS_IO_CMD_FLASH_SET_MODE, (void*)&nMode, sizeof(nMode));

is_FreezeVideo(hCam, IS_WAIT);
```

Sample programs

- uc480 Simple Trigger (C++)
- uc480 IO (C++)

4.3.86 is_SetFrameRate

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_SetFrameRate (HCAM hCam, double FPS, double* newFPS)
```

Description

Using `is_SetFrameRate()`, you can set the sensor frame rate in freerun mode (live mode). Since this value depends on the sensor timing, the exposure time actually used may slightly deviate from the value set here. After you have called the function, the actual frame rate is returned through the `newFPS` parameter.

If the frame rate is set too high, it might not be possible to transfer every single frame. In this case, the effective frame rate may vary from the set value.

For minimum and maximum frame rates as well as other sensor-based dependencies, please refer to [Camera and sensor data](#) chapter.

Note

Newer driver versions sometimes allow an extended value range for the frame rate setting. We recommend to query the value range every time and set the frame rate explicitly.

Changes to the frame rate affect the value ranges of the exposure time. After executing `is_SetFrameRate()`, calling the function [is_Exposure\(\)](#) is recommended in order to keep the defined camera settings.

Attention

The use of the following functions will affect the frame rate:

- [is_PixelClock\(\)](#)
- [is_SetOptimalCameraTiming\(\)](#)
- [is_AOI\(\)](#) (if the image size is changed)
- [is_SetSubSampling\(\)](#)
- [is_SetBinning\(\)](#)

Changes made to the window size or the read-out timing (pixel clock frequency) also affect the defined frame rate. For this reason, you need to call `is_SetFrameRate()` again after such changes.

Attention

To be able to set the default frame rate, you have to set a pixel clock equal to or higher than the default pixel clock.

Input parameters

hCam	Camera handle
FPS	Desired frame rate in frames per second (fps)
IS_GET_FRAMERATE	Returns the set frame rate in the <code>newFPS</code> parameter. To query the frame rate actually reached by the camera, use is_GetFramesPerSecond() .
IS_GET_DEFAULT_FRAMERATE	Returns the default frame rate.
newFPS	Returns the frame rate actually set.



Return values

IS_INVALID_CAMERA_HANDLE	Invalid camera handle
IS_INVALID_MODE	Camera is in standby mode, function not allowed
IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
IS_NO_SUCCESS	General error message
IS_SUCCESS	Function executed successfully

Related functions

- [is_GetFramesPerSecond\(\)](#)
- [is_GetFrameTimeRange\(\)](#)
- [is_PixelClock\(\)](#)
- [is_SetOptimalCameraTiming\(\)](#)
- [is_Exposure\(\)](#)
- [is_SetAutoParameter\(\)](#)
- [is_AOI\(\)](#)
- [is_SetSubSampling\(\)](#)
- [is_SetBinning\(\)](#)
- [is_CaptureVideo\(\)](#)

4.3.87 is_SetGainBoost

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_SetGainBoost (HCAM hCam, INT mode)
```

Description

In some cameras, `is_SetGainBoost()` enables an additional analog hardware gain boost feature on the sensor.

Input parameters

hCam	Camera handle
<input type="checkbox"/> mode	
IS_GET_GAINBOOST	Returns the current state of the gain boost function.
IS_SET_GAINBOOST_ON	Enables the gain boost function.
IS_SET_GAINBOOST_OFF	Disables the gain boost function.
IS_GET_SUPPORTED_GAINBOOST	Indicates whether the camera supports a gain boost feature or not.



Return values

Current setting when used together with IS_GET_GAINBOOST	Returns IS_SET_GAINBOOST_ON if the function is enabled, otherwise it returns IS_SET_GAINBOOST_OFF.
Current setting when used together with IS_GET_SUPPORTED_GAINBOOST	Returns IS_SET_GAINBOOST_ON if the function is supported, otherwise it returns IS_SET_GAINBOOST_OFF.
IS_INVALID_CAMERA_HANDLE	Invalid camera handle
IS_INVALID_MODE	Camera is in standby mode, function not allowed
IS_NO_SUCCESS	General error message
IS_SUCCESS	Function executed successfully

Related functions

- [is_SetHardwareGain\(\)](#)
- [is_SetHWGainFactor\(\)](#)
- [is_SetAutoParameter\(\)](#)

4.3.88 is_SetGamma

	
USB 2.0 USB 3.0 GigE	USB 2.0 USB 3.0 GigE

Syntax

INT is_SetGamma (HCAM hCam, INT nGamma)

Description

is_SetGamma() enables digital gamma correction which applies a gamma characteristic to the image. When hardware color conversion is used on GigE uEye HE cameras the gamma correction is performed in the camera hardware as well. When the color conversion is performed in the PC (software conversion) the gamma correction is performed in software.

Notes

1. When the color format is set to Raw Bayer the gamma correction can not be used.
2. Typical values for gamma range between 1.6 and 2.2.

Input parameters

hCam	Camera handle
nGamma	Gamma value to be set, multiplied by 100 (Range: 1...1000. Default = 100, corresponds to a gamma value of 1.0)
IS_GET_GAMMA	Returns the current setting.



Return values

When used with IS_GET_GAMMA	Current setting
IS_CANT_COMMUNICATE_WITH_DRIVER	Communication with the driver failed because no driver has been loaded.
IS_CANT_OPEN_DEVICE	An attempt to initialize or select the camera failed (no camera connected or initialization error).
IS_INVALID_CAMERA_HANDLE	Invalid camera handle
IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
IS_IO_REQUEST_FAILED	An IO request from the uc480 driver failed. Possibly the versions of the uc480.dll (API) and the driver file (uc480_usb.sys) do not match.
IS_NO_SUCCESS	General error message
IS_NOT_SUPPORTED	The camera model used here does not support this function or setting.
IS_SUCCESS	Function executed successfully

See also:

- Basics: [Characteristics and LUT](#)
- Basics: [Color filter \(Bayer filter\)](#)
- Programming: [is_SetColorConverter\(\)](#)

4.3.89 is_SetHardwareGain

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_SetHardwareGain (HCAM hCam, INT nMaster, INT nRed, INT nGreen, INT nBlue)
```

Description

`is_SetHardwareGain()` controls the sensor gain channels. These can be set between 0 % and 100 % independently of each other. The actual gain factor obtained for the value 100 % depends on the sensor and is specified in [Camera and sensor data](#) chapter.

You can use the [is_GetSensorInfo\(\)](#) function to query the available gain controls.

Notes

1. Sensor Gain

A signal gain will also result in a noise gain. High gain settings are therefore not recommended.

We suggest the following gain settings:

1. Enable the gain boost function [is_SetGainBoost\(\)](#).
2. If required, adjust the gain setting with `is_SetHardwareGain()`

New gain settings might only become effective when the next image is captured. This depends on the time when the gain settings are changed.

2. Linearity of Sensor Gain

On DCx Cameras, you can set the gain factor in increments from 0 to 100. These increments are not graduated linearly throughout the range due to the sensor. The increments will typically be greater in the upper range than in the lower range.

The maximum gain factor settings also vary from sensor to sensor.

3. Default Settings for RGB gains

The default setting values for the red, green and blue channel gain factors depend on the color correction matrix that has been set. If you select a different color correction matrix, the returned default values might change (see also [is_SetColorCorrection\(\)](#)).

Input parameters

<code>hCam</code>	Camera handle
<code>nMaster</code>	Sets the overall gain factor (0...100).
<code>IS_IGNORE_PARAMETER</code>	The master gain factor will not be changed.
<code>IS_GET_MASTER_GAIN</code>	Returns the master gain factor.
<code>IS_GET_RED_GAIN</code>	Returns the red channel gain factor.
<code>IS_GET_GREEN_GAIN</code>	Returns the green channel gain factor.
<code>IS_GET_BLUE_GAIN</code>	Returns the blue channel gain factor.
<code>IS_GET_DEFAULT_MASTER</code>	Returns the default master gain factor.
<code>IS_GET_DEFAULT_RED</code>	Returns the default red channel gain factor.
<code>IS_GET_DEFAULT_GREEN</code>	Returns the default green channel gain factor.
<code>IS_GET_DEFAULT_BLUE</code>	Returns the default blue channel gain factor.

hCam	Camera handle
IS_SET_ENABLE_AUTO_GAIN	Enables the auto gain functionality (see also is_SetAutoParameter()). You can disable the auto gain functionality by setting a value for nMaster.
nRed	Sets the red channel gain factor (0...100).
IS_IGNORE_PARAMETER	The red channel gain factor will not be changed.
nGreen	Sets the green channel gain factor (0...100).
IS_IGNORE_PARAMETER	The green channel gain factor will not be changed.
nBlue	Sets the blue channel gain factor (0...100).
IS_IGNORE_PARAMETER	The blue channel gain factor will not be changed.

Return values



When used with	Current setting
IS_GET_MASTER_GAIN IS_GET_RED_GAIN IS_GET_GREEN_GAIN IS_GET_BLUE_GAIN	
IS_CANT_COMMUNICATE_WITH_DRIVER	Communication with the driver failed because no driver has been loaded.
IS_CANT_OPEN_DEVICE	An attempt to initialize or select the camera failed (no camera connected or initialization error).
IS_INVALID_CAMERA_TYPE	The camera type defined in the .ini file does not match the current camera model.
IS_INVALID_CAMERA_HANDLE	Invalid camera handle
IS_INVALID_MODE	Camera is in standby mode, function not allowed
IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
IS_IO_REQUEST_FAILED	An IO request from the uc480 driver failed. Possibly the versions of the uc480.dll (API) and the driver file (uc480_usb.sys) do not match.
IS_NO_SUCCESS	General error message
IS_NOT_CALIBRATED	The camera does not contain any calibration data.
IS_NOT_SUPPORTED	The camera model used here does not support this function or setting.
IS_NULL_POINTER	Invalid array
IS_SUCCESS	Function executed successfully
IS_TIMED_OUT	A timeout occurred. An image capturing process could not be terminated within the allowable period.

Related functions

- [is_SetHWGainFactor\(\)](#)

- [is_GetSensorInfo\(\)](#)
- [is_SetGainBoost\(\)](#)
- [is_SetAutoParameter\(\)](#)

4.3.90 is_SetHWGainFactor

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_SetHWGainFactor (HCAM hCam, INT nMode, INT nFactor)
```

Description

`is_SetHWGainFactor()` uses gain factors to control sensor gain channels. These channels can be set independently of each other. The [is_SetHardwareGain\(\)](#) does not use factors for setting the gain channels, but standardized values between 0 and 100. The actual gain factor is sensor-dependent and can be found in [Camera and sensor data](#) chapter.

You can use the [is_GetSensorInfo\(\)](#) function to query the available gain controls.

Depending on the time when the gain settings are changed, these changes might only become effective when the next image is captured.

Input parameters

<code>hCam</code>	Camera handle
<input type="checkbox"/> <code>nMode</code>	
<code>IS_GET_MASTER_GAIN_FACTOR</code>	Returns the master gain factor.
<code>IS_GET_RED_GAIN_FACTOR</code>	Returns the red channel gain factor.
<code>IS_GET_GREEN_GAIN_FACTOR</code>	Returns the green channel gain factor.
<code>IS_GET_BLUE_GAIN_FACTOR</code>	Returns the blue channel gain factor.
<code>IS_SET_MASTER_GAIN_FACTOR</code>	Sets the master gain factor.
<code>IS_SET_RED_GAIN_FACTOR</code>	Sets the red channel gain factor.
<code>IS_SET_GREEN_GAIN_FACTOR</code>	Sets the green channel gain factor.
<code>IS_SET_BLUE_GAIN_FACTOR</code>	Sets the blue channel gain factor.
<code>IS_GET_DEFAULT_MASTER_GAIN_FACTOR</code>	Returns the default master gain factor.
<code>IS_GET_DEFAULT_RED_GAIN_FACTOR</code>	Returns the default red channel gain factor.
<code>IS_GET_DEFAULT_GREEN_GAIN_FACTOR</code>	Returns the default green channel gain factor.
<code>IS_GET_DEFAULT_BLUE_GAIN_FACTOR</code>	Returns the default blue channel gain factor.
<code>IS_INQUIRE_MASTER_GAIN_FACTOR</code>	Converts the index value for the master gain factor.
<code>IS_INQUIRE_RED_GAIN_FACTOR</code>	Converts the index value for the red channel gain factor.
<code>IS_INQUIRE_GREEN_GAIN_FACTOR</code>	Converts the index value for the green channel gain factor.
<code>IS_INQUIRE_BLUE_GAIN_FACTOR</code>	Converts the index value for the blue channel gain factor.
<code>nFactor</code>	Gain value (100 = gain factor 1, i. e. no effect)

For converting a gain value from the [is_SetHardwareGain\(\)](#) function, you can set the `nMode` parameter to one of the `IS_INQUIRE_x_FACTOR` values. In this case, the value range for `nFactor` is

between 0 and 100.

To set the gain using `IS_SET...GAIN_FACTOR`, you must set the `nFactor` parameter to an integer value in the range from 100 to the maximum value. By calling `IS_INQUIRE_x_FACTOR` and specifying the value 100 for `nFactor`, you can query the maximum value. A gain value of 100 means no gain, a gain value of 200 means gain to the double level (factor 2), etc.

Return values

When used with <code>IS_GET_MASTER_GAIN_FACTOR</code> <code>IS_GET_RED_GAIN_FACTOR</code> <code>IS_GET_GREEN_GAIN_FACTOR</code> <code>IS_GET_BLUE_GAIN_FACTOR</code>	Current setting
When used with <code>IS_SET_MASTER_GAIN_FACTOR</code> <code>IS_SET_RED_GAIN_FACTOR</code> <code>IS_SET_GREEN_GAIN_FACTOR</code> <code>IS_SET_BLUE_GAIN_FACTOR</code>	Defined setting
When used with <code>IS_GET_DEFAULT_MASTER_GAIN_FACTOR</code> <code>IS_GET_DEFAULT_RED_GAIN_FACTOR</code> <code>IS_GET_DEFAULT_GREEN_GAIN_FACTOR</code> <code>IS_GET_DEFAULT_BLUE_GAIN_FACTOR</code>	Default setting
When used with <code>IS_INQUIRE_MASTER_GAIN_FACTOR</code> <code>IS_INQUIRE_RED_GAIN_FACTOR</code> <code>IS_INQUIRE_GREEN_GAIN_FACTOR</code> <code>IS_INQUIRE_BLUE_GAIN_FACTOR</code>	Converted gain index
<code>IS_INVALID_CAMERA_HANDLE</code>	Invalid camera handle
<code>IS_INVALID_MODE</code>	Camera is in standby mode, function not allowed
<code>IS_NO_SUCCESS</code>	General error message
<code>IS_SUCCESS</code>	Function executed successfully

Related functions



- [is_SetHardwareGain\(\)](#)
- [is_SetGainBoost\(\)](#)
- [is_SetAutoParameter\(\)](#)
- [is_GetSensorInfo\(\)](#)

Example

```
//Set master gain factor to 3.57:
INT ret = is_SetHWGainFactor (hCam, IS_SET_MASTER_GAIN_FACTOR, 357);
//ret has the value 363 for the UI-1460-C

//Query the maximum gain factor for the red channel:
ret = is_SetHWGainFactor (hCam, IS_INQUIRE_RED_GAIN_FACTOR, 100);
//ret has the value 725 for the UI-1460-C
```

4.3.91 is_SetImageMem

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_SetImageMem (HCAM hCam, char* pcImgMem, INT id)
```

Description

`is_SetImageMem()` makes the specified image memory the active memory. Only an active image memory can receive image data. When you call [is_FreezeVideo\(\)](#), the captured image is stored in the image buffer designated by `pcImgMem` and `id`. For `pcImgMem`, you must pass a pointer which was created by [is_AllocImageMem\(\)](#), passing any other pointer will result in an error message. You may pass the same pointer multiple times.

Notes

- In the Direct3D or OpenGL modes, there is no need to set an image memory.
- If you want the application to be compatible with the FALCON SDK, make sure to call `is_SetImageSize()` after `is_SetImageMem()`.

Input parameters

<code>hCam</code>	Camera handle
<code>pcImgMem</code>	Pointer to the starting position in the memory.
<code>id</code>	ID of this memory.



Return values

<code>IS_INVALID_CAMERA_HANDLE</code>	Invalid camera handle
<code>IS_NO_SUCCESS</code>	General error message
<code>IS_SUCCESS</code>	Function executed successfully

Related functions

- [is_AllocImageMem\(\)](#)
- [is_FreeImageMem\(\)](#)
- [is_AddToSequence\(\)](#)
- [is_SetAllocatedImageMem\(\)](#)
- [is_GetColorDepth\(\)](#)
- [is_GetImageMem\(\)](#)
- [is_GetImageMemPitch\(\)](#)

4.3.92 is_SetOptimalCameraTiming

	
USB 2.0 USB 3.0	—

Syntax

```
INT is_SetOptimalCameraTiming (HCAM hCam, INT Mode, INT Timeout,
                              INT* pMaxPxlClk, double* pMaxFrameRate)
```

Description

Using `is_SetOptimalCameraTiming()`, you can determine the highest possible pixel clock frequency for the current configuration. This function sets the pixel clock for which no transfer errors will occur during the Timeout period. Moreover, it returns the highest frame rate available for this pixel clock frequency.

`is_SetOptimalCameraTiming()` can only be executed in free-run mode ([is_CaptureVideo\(\)](#)). If the return value is \neq `IS_SUCCESS`, no clock setting will be made.

Attention

The function should be executed in a separate thread and run in the background to allow for the computational load caused by additional color conversions, etc. Otherwise, it will not be able to return the optimum values.

Changes to the image geometry or pixel clock affect the value ranges of the frame rate and exposure time. After executing `is_SetOptimalCameraTiming()`, calling the following functions is recommended in order to keep the defined camera settings:

- [is_SetFrameRate\(\)](#)
- [is_Exposure\(\)](#)
- If you are using the DCx Camera's flash function: [is_IO\(\)](#)

Input parameters

hCam	Camera handle
Mode	
IS_BEST_PCLK_RUN_ONCE	The function makes one attempt to determine the optimum pixel clock and returns immediately.
Timeout [4000...20000]	Sets the period (in milliseconds) during which no transfer error may occur. The adjustable range is between 4 and 20 seconds. The higher the value you set for this parameter, the more stable the determined pixel clock value will be. This, in turn, increases the runtime of the function correspondingly.
pMaxPxlClk	Returns the maximum pixel clock frequency (in MHz).
pMaxFrameRate	Returns the maximum frame rate (in fps).

Return values



IS_AUTO_EXPOSURE_RUNNING	This setting cannot be changed while automatic exposure time control is enabled.
--------------------------	--

IS_INVALID_IMAGE_SIZE	Invalid image size This value is returned if e.g. the function is called with active binning or subsampling
IS_INVALID_MODE	Camera is in standby mode, function not allowed
IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
IS_NOT_SUPPORTED	The camera model used here does not support this function or setting.
IS_TRIGGER_ACTIVATED	The function cannot be used because the camera is waiting for a trigger signal.

Related functions

- [is_PixelClock\(\)](#)
- [is_SetFrameRate\(\)](#)
- [is_SetAutoParameter\(\)](#)
- [is_CaptureVideo\(\)](#)

4.3.93 is_SetRopEffect

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_SetRopEffect (HCAM hCam, INT effect, INT param, INT reserved)
```

Description

`is_SetRopEffect()` enables functions for real-time image geometry modification (Rop = raster operation).

Input parameters

hCam	Camera handle
<input type="checkbox"/> effect	
IS_SET_ROP_MIRROR_UPDOWN	Mirrors the image along the horizontal axis.
IS_SET_ROP_MIRROR_LEFTRIGHT	Mirrors the image along the vertical axis. Depending on the sensor, this operation is performed in the camera or in the PC software.
IS_GET_ROP_EFFECT	Returns the current settings.
param	Turns the Rop effect on/off. 0 = Turn off 1 = Turn on
reserved	Reserved. 0 must be passed.

Return values



When used with IS_GET_ROP_EFFECT	Current setting
IS_CANT_COMMUNICATE_WITH_DRIVER	Communication with the driver failed because no driver has been loaded.
IS_CANT_OPEN_DEVICE	An attempt to initialize or select the camera failed (no camera connected or initialization error).
IS_INVALID_CAMERA_TYPE	The camera type defined in the .ini file does not match the current camera model.
IS_INVALID_CAPTURE_MODE	The function can not be executed in the current camera operating mode (free run, trigger or standby).
IS_INVALID_CAMERA_HANDLE	Invalid camera handle
IS_INVALID_MODE	Camera is in standby mode, function not allowed
IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
IS_IO_REQUEST_FAILED	An IO request from the uc480 driver failed. Possibly the versions of the uc480.dll (API) and the driver file (uc480_usb.sys) do not match.

IS_NO_SUCCESS	General error message
IS_NOT_CALIBRATED	The camera does not contain any calibration data.
IS_NOT_SUPPORTED	The camera model used here does not support this function or setting.
IS_NULL_POINTER	Invalid array
IS_SUCCESS	Function executed successfully
IS_TIMED_OUT	A timeout occurred. An image capturing process could not be terminated within the allowable period.

Related functions

- [is_SetBinning\(\)](#)
- [is_SetSubSampling\(\)](#)
- [is_AOI\(\)](#)

4.3.94 is_SetSaturation

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_SetSaturation (HCAM hCam, INT ChromU, INT ChromV)
```

Description

Using `is_SetSaturation()`, you can set the software color saturation.

Note

In the YUV format, color information (i.e. the color difference signals) is provided by the U and V channels. In the U channel, this information results from the difference between the blue level and Y (luminance), in the V channel from the difference between the red level and Y.

For use in other color formats than YUV, U and V are converted using a driver matrix.

Input parameters

hCam	Camera handle
ChromU	U saturation: value multiplied by 100. Range: [IS_MIN_SATURATION ... IS_MAX_SATURATION]
IS_GET_SATURATION_U	Returns the current value for the U saturation.
ChromV	V saturation: value multiplied by 100. Range: [IS_MIN_SATURATION ... IS_MAX_SATURATION]
IS_GET_SATURATION_V	Returns the current value for the V saturation.



Return values

When used with IS_GET_SATURATION_U IS_GET_SATURATION_V	Current setting
IS_INVALID_CAMERA_HANDLE	Invalid camera handle
IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
IS_NO_SUCCESS	General error message
IS_SUCCESS	Function executed successfully

Related functions

- [is_SetColorMode\(\)](#)
- [is_SetColorCorrection\(\)](#)
- [is_SetColorConverter\(\)](#)

4.3.95 is_SetSensorScaler

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_SetSensorScaler (HCAM hCam, UINT nMode, double dblFactor)
```

Description

`is_SetSensorScaler()` enables internal image scaling for some sensors. This allows to reduce the image resolution by adjustable factors. Thus, the amount of data from high resolution sensors can be reduced.

Note

- Internal image scaling is only supported [DCC1240x and DCC3240x](#) series cameras.
- The use of the internal scaler has no effect on the attainable frame rate.

Input parameters

<code>hCam</code>	Camera handle
<input type="checkbox"/> <code>nMode</code> : Function mode	
<code>IS_ENABLE_SENSOR_SCALER</code>	Enable image scaling
<code>IS_ENABLE_SENSOR_SCALER IS_ENABLE_ANTI_ALIASING</code>	Enable image scaling with smoothed edges (anti-aliasing effect)
<code>dblFactor</code>	Scaling factor

Return values

<code>IS_CANT_COMMUNICATE_WITH_DRIVER</code>	Communication with the driver failed because no driver has been loaded.
<code>IS_CANT_OPEN_DEVICE</code>	An attempt to initialize or select the camera failed (no camera connected or initialization error).
<code>IS_INVALID_CAMERA_TYPE</code>	The camera type defined in the .ini file does not match the current camera model.
<code>IS_INVALID_CAMERA_HANDLE</code>	Invalid camera handle
<code>IS_INVALID_PARAMETER</code>	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
<code>IS_IO_REQUEST_FAILED</code>	An IO request from the uc480 driver failed. Possibly the versions of the <code>uc480.dll</code> (API) and the driver file (<code>uc480_usb.sys</code>) do not match.
<code>IS_NO_SUCCESS</code>	General error message
<code>IS_NOT_CALIBRATED</code>	The camera does not contain any calibration data.
<code>IS_NOT_SUPPORTED</code>	The camera model used here does not support this function or setting.
<code>IS_NULL_POINTER</code>	Invalid array
<code>IS_SUCCESS</code>	Function executed successfully

IS_TIMED_OUT	A timeout occurred. An image capturing process could not be terminated within the allowable period.
--------------	---

Related functions

- [is_GetSensorScalerInfo\(\)](#)



Example

```
SENSORSCALERINFO Info;
INT nRet;
double dblNewFactor;

// Query information on image scaling
nRet = is_GetSensorScalerInfo (hCam, &Info,
                              sizeof(Info));

// Enable scaling with anti aliasing
dblNewFactor = Info.dblMinFactor + Info.dblFactorIncrement;
nRet = is_SetSensorScaler (hCam, IS_ENABLE_SENSOR_SCALER |
                          IS_ENABLE_ANTI_ALIASING, dblNewFactor);
```

4.3.96 is_SetSensorTestImage

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_SetSensorTestImage (HCAM hCam, INT TestImage, INT Param)
```

Description

`is_SetSensorTestImage()` enables a test image function in the sensor. You can select different test images. The test images supported by a particular camera can be queried using the [is_GetSupportedTestImages\(\)](#) function. For some test images, the `Param` parameter provides additional options. If the test image does not support additional parameters, `Param` will be ignored.

Attention

Manually changing the pixel clock will disable the test image mode.

Input parameters

<code>hCam</code>	Camera handle
<code>TestImage</code>	The test image to be set. See also is_GetSupportedTestImages() .
<code>Param</code>	Additional parameter for used to modify the test image. Not available for all test images.

Return values



<code>IS_BAD_STRUCTURE_SIZE</code>	An internal structure has an incorrect size.
<code>IS_CANT_COMMUNICATE_WITH_DRIVER</code>	Communication with the driver failed because no driver has been loaded.
<code>IS_CANT_OPEN_DEVICE</code>	An attempt to initialize or select the camera failed (no camera connected or initialization error).
<code>IS_CAPTURE_RUNNING</code>	A capturing operation is in progress and must be terminated before you can start another one.
<code>IS_INVALID_BUFFER_SIZE</code>	The image memory has an inappropriate size to store the image in the desired format.
<code>IS_INVALID_CAMERA_TYPE</code>	The camera type defined in the .ini file does not match the current camera model.
<code>IS_INVALID_CAMERA_HANDLE</code>	Invalid camera handle
<code>IS_INVALID_MEMORY_POINTER</code>	Invalid pointer or invalid memory ID
<code>IS_INVALID_PARAMETER</code>	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
<code>IS_IO_REQUEST_FAILED</code>	An IO request from the uc480 driver failed. Possibly the versions of the <code>uc480.dll</code> (API) and the driver file (<code>uc480_usb.sys</code>) do not match.
<code>IS_NO_ACTIVE_IMG_MEM</code>	No active image memory available. You must set the memory to active using the <code>is_SetImageMem()</code> function or create a

	sequence using the <code>is_AddToSequence()</code> function.
<code>IS_NO_SUCCESS</code>	General error message
<code>IS_NOT_CALIBRATED</code>	The camera does not contain any calibration data.
<code>IS_NOT_SUPPORTED</code>	The camera model used here does not support this function or setting.
<code>IS_NULL_POINTER</code>	Invalid array
<code>IS_SUCCESS</code>	Function executed successfully
<code>IS_TIMED_OUT</code>	A timeout occurred. An image capturing process could not be terminated within the allowable period.
<code>IS_TRIGGER_ACTIVATED</code>	The function cannot be used because the camera is waiting for a trigger signal.

Related functions

- [is_GetSupportedTestImages\(\)](#)
- [is_GetTestImageValueRange\(\)](#)

4.3.97 is_SetSubSampling

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_SetSubSampling (HCAM hCam, INT mode)
```

Description

Using `is_SetSubSampling()`, you can enable sub-sampling mode both in horizontal and in vertical directions. This allows you to reduce the image size in the sub-sampling direction without scaling down the area of interest. In order to simultaneously enable horizontal and vertical sub-sampling, the horizontal and vertical sub-sampling parameters can be linked by a logical OR.

Some monochrome sensors are limited by their design to mere color sub-sampling. In case of fine image structures, this can result in slight artifacts.

The adjustable sub-sampling factors of each sensor are listed in [Camera and sensor data](#) chapter.

Note

Some sensors allow a higher pixel clock setting when binning or subsampling is activated. If you set a higher pixel clock and then reduce the binning/subsampling factors again, the driver will automatically select the highest possible pixel clock for the new settings.

Attention

Changes to the image geometry or pixel clock affect the value ranges of the frame rate and exposure time. After executing `is_SetSubSampling()`, calling the following functions is recommended in order to keep the defined camera settings:

- [is_SetFrameRate\(\)](#)
- [is_Exposure\(\)](#)
- If you are using the DCx Camera's flash function: [is_IO\(\)](#)

Input parameters

hCam	Camera handle
<input type="checkbox"/> mode	
IS_SUBSAMPLING_DISABLE	Disables sub-sampling.
IS_SUBSAMPLING_2X_VERTICAL	Enables vertical sub-sampling with factor 2.
IS_SUBSAMPLING_3X_VERTICAL	Enables vertical sub-sampling with factor 3.
IS_SUBSAMPLING_4X_VERTICAL	Enables vertical sub-sampling with factor 4.
IS_SUBSAMPLING_5X_VERTICAL	Enables vertical sub-sampling with factor 5.
IS_SUBSAMPLING_6X_VERTICAL	Enables vertical sub-sampling with factor 6.
IS_SUBSAMPLING_8X_VERTICAL	Enables vertical sub-sampling with factor 8.
IS_SUBSAMPLING_16X_VERTICAL	Enables vertical sub-sampling with factor 16.
IS_SUBSAMPLING_2X_HORIZONTAL	Enables horizontal sub-sampling with factor 2.
IS_SUBSAMPLING_3X_HORIZONTAL	Enables horizontal sub-sampling with factor 3.
IS_SUBSAMPLING_4X_HORIZONTAL	Enables horizontal sub-sampling with factor 4.

hCam	Camera handle
IS_SUBSAMPLING_5X_HORIZONTAL	Enables horizontal sub-sampling with factor 5.
IS_SUBSAMPLING_6X_HORIZONTAL	Enables horizontal sub-sampling with factor 6.
IS_SUBSAMPLING_8X_HORIZONTAL	Enables horizontal sub-sampling with factor 8.
IS_SUBSAMPLING_16X_HORIZONTAL	Enables horizontal sub-sampling with factor 16.
IS_GET_SUBSAMPLING	Returns the current setting.
IS_GET_SUBSAMPLING_FACTOR_VERTICAL	Returns the vertical sub-sampling factor
IS_GET_SUBSAMPLING_FACTOR_HORIZONTAL	Returns the horizontal sub-sampling factor
IS_GET_SUBSAMPLING_TYPE	Indicates whether the camera uses color-proof sub-sampling.
IS_GET_SUPPORTED_SUBSAMPLING	Returns the supported sub-sampling modes.

Return values

When used with IS_GET_SUBSAMPLING	Current setting: Returns an ORing of the defined constants from <code>ueye.h</code> , e.g. <code>IS_SUBSAMPLING_2X_HORIZONTAL</code>
When used with IS_GET_SUBSAMPLING_FACTOR_VERTICAL IS_GET_SUBSAMPLING_FACTOR_HORIZONTAL	Current setting: Returns the current factor as integer value (2, 3, 4, 5, 6, 8, 16)
When used with IS_GET_SUBSAMPLING_TYPE	Current setting: Returns <code>IS_SUBSAMPLING_COLOR</code> if the camera uses color-proof sub-sampling, else <code>IS_SUBSAMPLING_MONO</code>
In Verbindung mit IS_GET_SUPPORTED_SUBSAMPLING	Returns the supported sub-sampling modes linked by logical ORs
IS_BAD_STRUCTURE_SIZE	An internal structure has an incorrect size.
IS_CANT_ADD_TO_SEQUENCE	The image memory is already included in the sequence and cannot be added again.
IS_CANT_COMMUNICATE_WITH_DRIVER	Communication with the driver failed because no driver has been loaded.
IS_CANT_OPEN_DEVICE	An attempt to initialize or select the camera failed (no camera connected or initialization error).
IS_CAPTURE_RUNNING	A capturing operation is in progress and must be terminated before you can start another one.
IS_DR_CANNOT_CREATE_SURFACE	The image surface or overlay surface could not be created.
IS_DR_CANNOT_CREATE_TEXTURE	The texture could not be created.
IS_DR_CANNOT_CREATE_VERTEX_BUFFER	The vertex buffer could not be created.
IS_DR_DEVICE_OUT_OF_MEMORY	Not enough graphics memory available.
IS_DR_LIBRARY_NOT_FOUND	The DirectRenderer library could not be found.



IS_INVALID_BUFFER_SIZE	The image memory has an inappropriate size to store the image in the desired format.
IS_INVALID_CAMERA_TYPE	The camera type defined in the .ini file does not match the current camera model.
IS_INVALID_CAPTURE_MODE	The function can not be executed in the current camera operating mode (free run, trigger or standby).
IS_INVALID_CAMERA_HANDLE	Invalid camera handle
IS_INVALID_MEMORY_POINTER	Invalid pointer or invalid memory ID
IS_INVALID_MODE	Camera is in standby mode, function not allowed
IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
IS_INVALID_PIXEL_CLOCK	This setting is not available for the currently set pixel clock frequency.
IS_IO_REQUEST_FAILED	An IO request from the uc480 driver failed. Possibly the versions of the uc480.dll (API) and the driver file (uc480_usb.sys) do not match.
IS_NO_ACTIVE_IMG_MEM	No active image memory available. You must set the memory to active using the <code>is_SetImageMem()</code> function or create a sequence using the <code>is_AddToSequence()</code> function.
IS_NO_SUCCESS	General error message
IS_NOT_CALIBRATED	The camera does not contain any calibration data.
IS_NOT_SUPPORTED	The camera model used here does not support this function or setting.
IS_NULL_POINTER	Invalid array
IS_OUT_OF_MEMORY	No memory could be allocated.
IS_SEQUENCE_BUF_ALREADY_LOCKED	The memory could not be locked. The pointer to the buffer is invalid.
IS_SUCCESS	Function executed successfully
IS_TIMED_OUT	A timeout occurred. An image capturing process could not be terminated within the allowable period.
IS_TRIGGER_ACTIVATED	The function cannot be used because the camera is waiting for a trigger signal.

Related functions

- [is_SetBinning\(\)](#)

- [is_AOI\(\)](#)
- [is_ImageFormat\(\)](#)
- [is_PixelClock\(\)](#)

4.3.98 is_SetTimeout

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_SetTimeout (HCAM hCam, UINT nMode, UINT Timeout)
```

Description

Using `is_SetTimeout()`, you can change user-defined timeout values of the uc480 API. If no user-defined timeout is set, the default value of the uc480 API is used for the relevant timeout.

For further information, please refer to the [How to proceed: Timeout values for image capture](#) section.

Note

The user-defined timeout only applies to the specified camera at runtime of the program.

Input parameters

hCam	Camera handle
nMode	Selects the timeout value to be set
IS_TRIGGER_TIMEOUT	Sets the timeout value for triggered image capture
Timeout	Timeout value in 10 ms. Value range [0; 4...429496729] (corresponds to 40 ms to approx. 1193 hours) 0 = use default value of the uc480 API For 1...3, the value 4 is used.

Return values

IS_CANT_COMMUNICATE_WITH_DRIVER	Communication with the driver failed because no driver has been loaded.
IS_CANT_OPEN_DEVICE	An attempt to initialize or select the camera failed (no camera connected or initialization error).
IS_INVALID_CAMERA_HANDLE	Invalid camera handle
IS_IO_REQUEST_FAILED	An IO request from the uc480 driver failed. Possibly the versions of the <code>uc480.dll</code> (API) and the driver file (<code>uc480_usb.sys</code>) do not match.
IS_NO_SUCCESS	General error message
IS_NOT_SUPPORTED	The camera model used here does not support this function or setting.
IS_SUCCESS	Function executed successfully



Related functions

- [is_GetTimeout\(\)](#)
- [is_CaptureVideo\(\)](#)
- [is_FreezeVideo\(\)](#)
- [is_SetExternalTrigger\(\)](#)

Example

```
// Set user-defined timeout to 120 seconds  
is_SetTimeout(hCam, IS_TRIGGER_TIMEOUT, 12000);
```

4.3.99 is_SetTriggerCounter

	
USB 2.0	USB 2.0

Syntax

```
INT is_SetTriggerCounter (HCAM hCam, INT nValue)
```

Description

`is_SetTriggerCounter()` returns the number of images captured in hardware or software trigger mode.

Note

In freerun mode, the counter always returns 0 even when images were captured.

Input parameters

hCam	Camera handle
<input type="checkbox"/> nValue	
IS_GET_TRIGGER_COUNTER	Returns the current count for triggered image captures
Other values	Resets the counter for triggered image captures



Return values

IS_CANT_COMMUNICATE_WITH_DRIVER	Communication with the driver failed because no driver has been loaded.
IS_CANT_OPEN_DEVICE	An attempt to initialize or select the camera failed (no camera connected or initialization error).
IS_INVALID_CAMERA_HANDLE	Invalid camera handle
IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
IS_IO_REQUEST_FAILED	An IO request from the uc480 driver failed. Possibly the versions of the <code>uc480.dll</code> (API) and the driver file (<code>uc480_usb.sys</code>) do not match.
IS_NO_SUCCESS	General error message
IS_SUCCESS	Function executed successfully

Related functions

- [is_SetExternalTrigger\(\)](#)
- [is_CameraStatus\(\)](#)

4.3.100 is_SetTriggerDelay

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_SetTriggerDelay (HCAM hCam, INT nTriggerDelay)
```

Description

Using `is_SetTriggerDelay()`, you can set the delay time between the arrival of a trigger signal and the start of exposure. The trigger signal can be initiated by hardware or by software.

The delay time set here adds to the delay caused by the sensor. The delay times of each sensor are listed in [Camera and sensor data](#) chapter.

Input parameters

hCam	Camera handle
nTriggerDelay	Time by which the image capture is delayed (in μ s) 0 = deactivate trigger delay
IS_GET_TRIGGER_DELAY	Returns the currently set delay time.
IS_GET_MIN_TRIGGER_DELAY	Returns the minimum adjustable value.
IS_GET_MAX_TRIGGER_DELAY	Returns the maximum adjustable value.
IS_GET_TRIGGER_DELAY_GRANULARITY	Returns the resolution of the adjustable delay time.



Return values

When used with	Current setting
IS_GET_TRIGGER_DELAY	
IS_CANT_COMMUNICATE_WITH_DRIVER	Communication with the driver failed because no driver has been loaded.
IS_CANT_OPEN_DEVICE	An attempt to initialize or select the camera failed (no camera connected or initialization error).
IS_INVALID_CAMERA_HANDLE	Invalid camera handle
IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
IS_IO_REQUEST_FAILED	An IO request from the uc480 driver failed. Possibly the versions of the <code>uc480.dll</code> (API) and the driver file (<code>uc480_usb.sys</code>) do not match.
IS_NO_SUCCESS	General error message
IS_SUCCESS	Function executed successfully

Related functions

- [is_IO\(\)](#)
- [is_SetExternalTrigger\(\)](#)

4.3.101 is_StopLiveVideo

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_StopLiveVideo (HCAM hCam, INT Wait)
```

Description

`is_StopLiveVideo()` stops live mode or cancels a hardware triggered image capture in case the exposure has not yet started.

Input parameters

hCam	Camera handle
<input type="checkbox"/> Wait	
IS_WAIT	The function waits until the image save is complete.
IS_DONT_WAIT	The function returns immediately. Digitizing the image is completed in the background.
IS_FORCE_VIDEO_STOP	Digitizing is stopped immediately.



Return values

IS_INVALID_CAMERA_HANDLE	Invalid camera handle
IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
IS_NO_SUCCESS	General error message
IS_SUCCESS	Function executed successfully
IS_TIMED_OUT	A timeout occurred. An image capturing process could not be terminated within the allowable period.

Related functions

- [is_FreezeVideo\(\)](#)
- [is_CaptureVideo\(\)](#)
- [is_SetDisplayMode\(\)](#)

4.3.102 is_Trigger

	
GigE USB 3.0	GigE USB 3.0

Syntax

```
INT is_Trigger (HCAM hCam, UINT nCommand, void* pParam, UINT cbSizeOfParam)
```

Description

`is_Trigger()` activates the burst trigger mode in GigE and USB 3 uEye cameras. In burst trigger mode, the camera captures a series of images in rapid succession on receipt of a single trigger signal. The trigger signal can be generated by the software ([is_FreezeVideo\(\)](#)) or transmitted via the digital input of the camera. The burst images are captured and transferred at maximum speed. The maximum speed depends on the pixel clock parameter (see [is_PixelClock\(\)](#)) and the exposure time parameter (see [is_Exposure\(\)](#)). `is_Transfer()` allows adjusting the latency of image data transfer.

The `nCommand` input parameter is used to select the function mode. The `pParam` input parameter depends on the selected function mode. If you select functions for setting or returning a value, `pParam` contains a pointer to a variable of the `UINT` type. The size of the memory area to which `pParam` refers is specified in the `cbSizeOfParam` input parameter.



This function is currently only supported by the GigE and USB 3 uEye camera series.



Note on the uEye trigger functions

In a future version of the uEye API, all trigger functions will be combined in a single `is_Trigger()` function. The following functions will then be obsolete:

- [is_SetExternalTrigger\(\)](#)
- [is_ForceTrigger\(\)](#)
- [is_SetTriggerCounter\(\)](#)
- [is_SetTriggerDelay\(\)](#)
- [is_TriggerDebounce\(\)](#)



Note on trigger delay in burst trigger mode

If you set a trigger delay with the [is_SetTriggerDelay\(\)](#) function, the delay will only apply to the first image after each trigger signal.

Input parameters

hCam	Camera handle
<div style="border: 1px solid gray; padding: 5px;"> <input type="checkbox"/> nCommand </div>	

IS_TRIGGER_CMD_GET_BURST_SIZE_SUPPORTED	<p>Returns if the camera supports the burst trigger mode.</p> <ul style="list-style-type: none"> ▣ More details • pParam: Pointer to variable of type <code>UINT</code>. • nSizeOfParam: 4 <p>Example 1</p>
IS_TRIGGER_CMD_GET_BURST_SIZE_RANGE	<p>Returns the value range, the default value and the increment for the number of images in a burst.</p> <ul style="list-style-type: none"> ▣ More details • pParam: Pointer to variable of type RANGE_OF_VALUES_U32 returning the values. • nSizeOfParam: <code>sizeof (RANGE_OF_VALUES_U32)</code>
IS_TRIGGER_CMD_GET_BURST_SIZE	<p>Returns the currently set number of images in a burst.</p> <ul style="list-style-type: none"> ▣ More details • pParam: Pointer to variable of type <code>UINT</code> returning the current value. • nSizeOfParam: 4 <p>Example 2</p>
IS_TRIGGER_CMD_SET_BURST_SIZE	<p>Sets the number of images in a burst.</p> <ul style="list-style-type: none"> ▣ More details • pParam: Pointer to variable of type <code>UINT</code> that passes the value to be set. • nSizeOfParam: 4
pParam	Pointer to a function parameter, whose function depends on <code>nCommand</code> .
cbSizeOfParam	Size (in bytes) of the memory area to which <code>pParam</code> refers.

Contents of the RANGE_OF_VALUES_U32 structure

UINT	u32Minimum	Minimum value
UINT	u32Maximum	Maximum value
UINT	u32Increment	Increment
UINT	u32Default	Default value

Return values

IS_INVALID_CAMERA_HANDLE	Invalid camera handle
IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
IS_NO_SUCCESS	General error message
IS_NOT_SUPPORTED	The camera model used here does not support this function or setting.
IS_SUCCESS	Function executed successfully

Related functions

- [is_SetExternalTrigger\(\)](#)
- [is_ForceTrigger\(\)](#)
- [is_SetTriggerCounter\(\)](#)
- [is_SetTriggerDelay\(\)](#)
- [is_TriggerDebounce\(\)](#)

Example 1

```

UINT nTriggerBurstSizeSupported = 0;
INT nRet = is_Trigger(m_hCam, IS_TRIGGER_CMD_GET_BURST_SIZE_SUPPORTED,
                    (void*)&nTriggerBurstSizeSupported, sizeof(nTriggerBurstSizeSupported));

if (nRet == IS_SUCCESS) {
    // Burst size supported
    if (nTriggerBurstSizeSupported == 1) {
        RANGE_OF_VALUES_U32 rangeBurstSize;
        nRet = is_Trigger(m_hCam, IS_TRIGGER_CMD_GET_BURST_SIZE_RANGE,
                        (void*)&rangeBurstSize, sizeof(rangeBurstSize));

        if (nRet == IS_SUCCESS) {
            UINT nMin, nMax;
            nMin = rangeBurstSize.u32Minimum;
            nMax = rangeBurstSize.u32Maximum;
        }
    }
}

```

Example 2



```

UINT nTriggerBurstSize = 0;
INT nRet = is_Trigger(m_hCam,
                    IS_TRIGGER_CMD_GET_BURST_SIZE,
                    (void*)&nTriggerBurstSize,
                    sizeof(nTriggerBurstSize)
                    );

nRet = is_Trigger(m_hCam,
                IS_TRIGGER_CMD_SET_BURST_SIZE,
                (void*)&nTriggerBurstSize,
                sizeof(nTriggerBurstSize)
                );

```

4.3.103 is_TriggerDebounce

	
GigE	GigE

Syntax

```

INT is_TriggerDebounce (HCAM hCam, UINT nCommand,
                      void *pParam, UINT nSizeOfParam)

```

Description

Using `is_TriggerDebounce()`, you can suppress disturbances at the digital input when you are running a GigE uEye camera in trigger mode. The signal at the digital input is only recognized as a trigger if the signal level remains constant at the target level for a user-selectable time. The signal edge and a delay (`DELAY_TIME`) can be set as parameters. It is recommend to use automatic signal edge selection.

Example: Mode set to "rising edge" (`TRIGGER_DEBOUNCE_MODE_RISING_EDGE`) and delay set to 50 μ s. The camera will not trigger the image capture on the rising edge until the digital signal has remained at the high level for longer than 50 μ s without interruption. If this is not the case, the signal

is regarded as a disturbance and ignored.

The `nCommand` input parameter is used to select the function mode. The `pParam` input parameter depends on the selected function mode. If you select functions for setting or returning a value, `pParam` contains a pointer to a variable of the `UINT` type. The size of the memory area to which `pParam` refers is specified in the `nSizeOfParam` input parameter.



`is_TriggerDebounce()` delays the start of a triggered image capture by the selected time (`DELAY_TIME`).



This function is currently only supported by the GigE uEye camera series.

Input parameters

hCam	Camera handle
<p data-bbox="82 766 248 792">[-] nCommand</p>	

Setting the function mode	
TRIGGER_DEBOUNCE_CMD_GET_SUPPORTED_MODES	<p>Returns the function modes supported by the camera.</p> <ul style="list-style-type: none"> ▣ More details • pParam: Pointer to a bit mask of type <code>UINT</code>. The bit mask returns the supported modes, linked by logical ORs (see Function modes table). • nSizeOfParam: 4
TRIGGER_DEBOUNCE_CMD_SET_MODE	<p>Sets a function mode.</p> <ul style="list-style-type: none"> ▣ More details • pParam: Pointer to variable of type <code>UINT</code> that passes the value to be set. • nSizeOfParam: 4
TRIGGER_DEBOUNCE_CMD_GET_MODE	<p>Returns the set function mode.</p> <ul style="list-style-type: none"> ▣ More details • pParam: Pointer to variable of type <code>UINT</code> returning the current value. • nSizeOfParam: 4
TRIGGER_DEBOUNCE_CMD_GET_MODE_DEFAULT	<p>Returns the default mode.</p> <ul style="list-style-type: none"> ▣ More details • pParam: Pointer to variable of type <code>UINT</code> returning the default value. • nSizeOfParam: 4
Setting the delay time	
TRIGGER_DEBOUNCE_CMD_SET_DELAY_TIME	<p>Sets a delay time (in μs).</p> <ul style="list-style-type: none"> ▣ More details • pParam: Pointer to variable of type <code>UINT</code> that passes the value to be set. • nSizeOfParam: 4
TRIGGER_DEBOUNCE_CMD_GET_DELAY_TIME	<p>Returns the set delay time (in μs).</p> <ul style="list-style-type: none"> ▣ More details • pParam: Pointer to variable of type <code>UINT</code> returning the current value. • nSizeOfParam: 4
TRIGGER_DEBOUNCE_CMD_GET_DELAY_TIME_MIN	<p>Returns the minimum value for the delay (in μs).</p> <ul style="list-style-type: none"> ▣ More details • pParam: Pointer to variable of type <code>UINT</code> returning the minimum value. • nSizeOfParam: 4
TRIGGER_DEBOUNCE_CMD_GET_DELAY_TIME_MAX	<p>Returns the maximum value for the delay (in μs).</p> <ul style="list-style-type: none"> ▣ More details • pParam: Pointer to variable of type <code>UINT</code> returning the maximum value. • nSizeOfParam: 4
TRIGGER_DEBOUNCE_CMD_GET_DELAY_TIME_INC	<p>Returns the increment for setting the delay.</p> <ul style="list-style-type: none"> ▣ More details • pParam: Pointer to variable of type <code>UINT</code> returning the increment.

pParam	Pointer to a function parameter, whose function depends on nCommand.
nSizeOfParam	Size (in bytes) of the memory area to which pParam refers.

Function modes

TRIGGER_DEBOUNCE_MODE_NONE	Disables debouncing the digital input.
TRIGGER_DEBOUNCE_MODE_FALLING_EDGE	Debounces the digital input for falling edge signals.
TRIGGER_DEBOUNCE_MODE_RISING_EDGE	Debounces the digital input for rising edge signals.
TRIGGER_DEBOUNCE_MODE_BOTH_EDGES	Debounces the digital input for rising or falling edge signals.
TRIGGER_DEBOUNCE_MODE_AUTOMATIC	Debounces the digital input with automatic edge selection (recommended). The edge is selected based on the set trigger edge (see is_SetExternalTrigger()).

Return values

IS_CANT_COMMUNICATE_WITH_DRIVER	Communication with the driver failed because no driver has been loaded.
IS_CANT_OPEN_DEVICE	An attempt to initialize or select the camera failed (no camera connected or initialization error).
IS_INVALID_CAMERA_HANDLE	Invalid camera handle
IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
IS_IO_REQUEST_FAILED	An IO request from the uc480 driver failed. Possibly the versions of the uc480.dll (API) and the driver file (uc480_usb.sys) do not match.
IS_NO_SUCCESS	General error message
IS_NOT_SUPPORTED	The camera model used here does not support this function or setting.
IS_SUCCESS	Function executed successfully



Related functions

- [is_SetExternalTrigger\(\)](#)
- [is_SetTriggerDelay\(\)](#)

Example

```
INT value = TRIGGER_DEBOUNCE_MODE_AUTOMATIC;
is_TriggerDebounce(hCam, TRIGGER_DEBOUNCE_CMD_SET_MODE, (void*)&value, sizeof(value));
```

4.3.104 is_UnlockSeqBuf

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_UnlockSeqBuf (HCAM hCam, INT nNum, char* pcMem)
```

Description

Using `is_UnlockSeqBuf()`, you unlock a previously locked image memory in order to make it available again for storing captured images. The image memory is re-inserted at its previous position in the sequence list.

Input parameters

hCam	Camera handle
nNum	Number of the image memory to unlock. When you pass <code>IS_IGNORE_PARAMETER</code> , the image memory is only identified by its starting address. nNum identifies the position in the sequence list, not the memory ID assigned with is_AllocImageMem() .
pcMem	Starting address of the image memory



Return values

IS_BAD_STRUCTURE_SIZE	An internal structure has an incorrect size.
IS_INVALID_CAMERA_HANDLE	Invalid camera handle
IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
IS_NO_SUCCESS	General error message
IS_SUCCESS	Function executed successfully

Related functions

- [is_LockSeqBuf\(\)](#)
- [is_WaitForNextImage\(\)](#)

4.3.105 is_WaitEvent

	
-	USB 2.0 USB 3.0

Syntax

```
INT is_WaitEvent (HCAM hCam, INT which, INT nTimeout)
```

Description

`is_WaitEvent()` allows waiting for uc480 events. The function indicates successful execution when the event has occurred within the specified timeout.

Input parameters

hCam	Camera handle
which	ID of the event (see is_EnableEvent())
nTimeout	Time (in ms) that the function will wait for an event to occur. Using the constant <code>INFINITE</code> you can set the time for the timeout to infinity.

Return values

IS_NO_SUCCESS	General error message
IS_SUCCESS	Function executed successfully
IS_TIMED_OUT	A timeout occurred. An image capturing process could not be terminated within the allowable period. That means the specified timeout expired without the event having occurred.

Related functions



- [is_EnableEvent\(\)](#)
- [is_DisableEvent\(\)](#)

Example

```
//Activate and initialize FRAME event
is_EnableEvent(hCam, IS_SET_EVENT_FRAME);

//Start image capture and wait 1000 ms for event to occur
is_FreezeVideo(hCam, IS_DONT_WAIT);
INT nRet = is_WaitEvent(hCam, IS_SET_EVENT_FRAME, 1000);
if (nRet == IS_TIMED_OUT)
{
    /* wait timed out */
}
else if (nRet == IS_SUCCESS)
{
    /* event signalled */
}
is_DisableEvent(hCam, IS_SET_EVENT_FRAME);
```


4.3.106 is_WaitForNextImage

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_WaitForNextImage(HCAM hCam, UINT timeout, char** ppcMem, INT* imageID)
```

Description

`is_WaitForNextImage()` returns the pointer and sequence ID of the first (i.e. oldest) image in a memory sequence. The queue mode has to be enabled for the memory sequence (see [is_InitImageQueue\(\)](#)). If the sequence does not contain images, `is_WaitForNextImage()` waits until a new image arrives or until the specified time has elapsed.

Note

Note that also image capture errors are added to the `ImageQueue` like images. If a call of `is_WaitForNextImage()` returns the `IS_CAPTURE_STATUS` return value then you can check by a new call of the function, if any further images were enqueued into the `ImageQueue` after the error.

Attention

Image memories in a sequence with queue mode are automatically locked. The image memories will have to be unlocked with [is_UnlockSeqBuf\(\)](#) in order to be re-used in the sequence.

Input parameters

<code>hCam</code>	Camera handle
<code>timeout</code>	Timeout in ms. Range $0 \dots 2^{32}-1$ If no images are in the sequence and no image arrives during the timeout, the function returns <code>IS_TIMED_OUT</code> .
<code>ppcMem</code>	Pointer to a variable which will receive the address of the last image in the sequence.
<code>imageID</code>	Pointer to a variable which will receive the sequence ID of the oldest image in the sequence.

Return values

<code>IS_CANT_COMMUNICATE_WITH_DRIVER</code>	Communication with the driver failed because no driver has been loaded.
<code>IS_CANT_OPEN_DEVICE</code>	An attempt to initialize or select the camera failed (no camera connected or initialization error).
<code>IS_INVALID_CAMERA_HANDLE</code>	Invalid camera handle
<code>IS_INVALID_MEMORY_POINTER</code>	Invalid pointer or invalid memory ID
<code>IS_INVALID_PARAMETER</code>	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
<code>IS_IO_REQUEST_FAILED</code>	An IO request from the <code>uc480</code> driver failed. Possibly the versions of the <code>uc480.dll</code> (API) and the driver file (<code>uc480_usb.sys</code>) do not match.

IS_NO_SUCCESS	General error message
IS_SUCCESS	Function executed successfully
IS_TIMED_OUT	A timeout occurred. An image capturing process could not be terminated within the allowable period.
IS_CAPTURE_STATUS	A transfer error occurred or no image memory was available for saving. The parameter <code>IS_CAPTURE_STATUS</code> replaces the previous parameter <code>IS_TRANSFER_FAILED</code> . The parameter <code>IS_TRANSFER_FAILED</code> was moved into the new header file <code>uc480_deprecated.h</code> , which contains all obsolete function definitions and constants. If necessary the header file <code>uc480_deprecated.h</code> can be included in addition to the header file <code>uc480.h</code> .



Related functions

- [is_InitImageQueue\(\)](#)
- [is_ExitImageQueue\(\)](#)
- [is_UnlockSeqBuf\(\)](#)

See also:

- How to proceed: [Image memory sequences](#)
- How to proceed: [Allocating image memory](#)
- Transfer error: uc480 Viewer Image infos
- Transfer error: [is_GetImageInfo\(\)](#)

4.3.107 is_WriteEEPROM

	
USB 2.0 USB 3.0	USB 2.0 USB 3.0

Syntax

```
INT is_WriteEEPROM (HCAM hCam, INT Adr, char* pcString, INT Count)
```

Description

Using `is_WriteEEPROM()`, you can write data to the EEPROM of the camera. Besides the hard-coded factory information, the EEPROM of the DCx Camera can hold 64 bytes of user data.

Input parameters

hCam	Camera handle
Adr	Starting address for data writes (0...63)
pcString	Pointer to the string containing the data to be written
Count	Number of characters to be written (1...64)

Return values

IS_CANT_COMMUNICATE_WITH_DRIVER	Communication with the driver failed because no driver has been loaded.
IS_CANT_OPEN_DEVICE	An attempt to initialize or select the camera failed (no camera connected or initialization error).
IS_INVALID_CAMERA_TYPE	The camera type defined in the .ini file does not match the current camera model.
IS_INVALID_CAMERA_HANDLE	Invalid camera handle
IS_INVALID_MODE	Camera is in standby mode, function not allowed
IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
IS_IO_REQUEST_FAILED	An IO request from the uc480 driver failed. Possibly the versions of the <code>uc480.dll</code> (API) and the driver file (<code>uc480_usb.sys</code>) do not match.
IS_NO_SUCCESS	General error message
IS_NOT_CALIBRATED	The camera does not contain any calibration data.
IS_NOT_SUPPORTED	The camera model used here does not support this function or setting.
IS_NULL_POINTER	Invalid array
IS_SUCCESS	Function executed successfully
IS_TIMED_OUT	A timeout occurred. An image capturing process could not be terminated within the allowable period.

Related functions

- [is_ReadEEPROM\(\)](#)

4.4 AVI Function Descriptions

The functions of the `uc480_tools.dll` enable you to save images captured with the DCx Camera as sequences to an AVI file. The [How to proceed: Capturing AVIs](#) chapter shows the command sequence for capturing an AVI video.

Notes



- **Video Resolution**

If the width or height of the video resolution is greater than 2048 pixels, some media players may not be able to play the video.

- **Video Compression**

To reduce the file size, the single frames are stored in the AVI container using an adjustable JPEG compression. It is possible to extract single frames from the AVI file.

4.4.1 isavi_AddFrame

	
USB 2.0 USB 3.0	-

Syntax

```
INT isavi_AddFrame (INT nAviID, char* pcImageMem)
```

Description

`isavi_AddFrame()` adds a new frame to an AVI sequence.

Input parameters

<code>nAviID</code>	Instance ID set by the isavi_InitAVI() function
<code>pcImageMem</code>	Pointer to the memory containing the image.



Return values

<code>IS_AVI_NO_ERR</code>	Function executed successfully.
<code>IS_AVI_ERR_INVALID_ID</code>	The indicated AVI instance could not be found. Either the AVI ID is invalid or the instance has already been deleted using isavi_ExitAVI() .
<code>IS_AVI_ERR_COMPRESS</code>	The last image compression failed.
<code>IS_AVI_ERR_COMPRESSION_RUN</code>	The current image could not be processed since compression is still in progress.
<code>IS_AVI_ERR_INVALID_FILE</code>	The file has no valid AVI format.

Related functions

- [isavi_InitAVI\(\)](#)

4.4.2 isavi_CloseAVI

	
USB 2.0 USB 3.0	-

Syntax

```
INT isavi_CloseAVI (INT nAviID)
```

Description

isavi_CloseAVI() closes an AVI file which was opened using [isavi_OpenAVI\(\)](#).

Input parameters

nAviID	Instance ID set by the isavi_InitAVI() function.
--------	--



Return values

IS_AVI_NO_ERR	Function executed successfully.
IS_AVI_ERR_INVALID_ID	The indicated AVI instance could not be found. Either the AVI ID is invalid or the instance has already been deleted using isavi_ExitAVI() .

Related functions

- [isavi_OpenAVI\(\)](#)
- [isavi_InitAVI\(\)](#)
- [isavi_ExitAVI\(\)](#)

4.4.3 isavi_DisableEvent

	
USB 2.0 USB 3.0	-

Syntax

```
INT isavi_DisableEvent (INT nAviID, INT which)
```

Description

isavi_DisableEvent() disables the specified event. The disabled event no longer triggers an event notification when the associated event occurs.

Input parameters

nAviID	Instance ID set by the isavi_InitAVI() function.
which	Name of the event to be disabled.
IS_AVI_SET_EVENT_FRAME_SAVED	A new frame was saved to the AVI file.

Return values



IS_AVI_NO_ERR	Function executed successfully.
IS_AVI_ERR_INVALID_ID	The indicated AVI instance could not be found.

	Either the AVI ID is invalid or the instance has already been deleted using isavi_ExitAVI() .
IS_AVI_ERR_PARAMETER	One of the submitted parameters is outside the valid range.

Related functions

- [isavi_EnableEvent\(\)](#)

4.4.4 isavi_EnableEvent

	
USB 2.0 USB 3.0	-

Syntax

```
INT isavi_EnableEvent (INT nAviID, INT which)
```

Description

`isavi_EnableEvent()` enables the specified event. The enabled event triggers an event notification when the associated event occurs.

Input parameters

nAviID	Instance ID set by the isavi_InitAVI() function.
which	Name of the event to be enabled.
IS_AVI_SET_EVENT_FRAME_SAVED	A new frame was saved to the AVI file.



Return values

IS_AVI_NO_ERR	Function executed successfully.
IS_AVI_ERR_INVALID_ID	The indicated AVI instance could not be found. Either the AVI ID is invalid or the instance has already been deleted using isavi_ExitAVI() .
IS_AVI_ERR_PARAMETER	One of the submitted parameters is outside the valid range.

Related functions

- [isavi_DisableEvent\(\)](#)

4.4.5 isavi_ExitAVI

	
USB 2.0 USB 3.0	-

Syntax

```
INT isavi_ExitAVI (INT nAviID)
```

Description

`isavi_ExitAVI()` terminates and deletes the instance of the uc480 AVI interface.

Input parameters

nAviID	Instance ID set by the isavi_InitAVI() function.
--------	--



Return values

IS_AVI_NO_ERR	Function executed successfully.
IS_AVI_ERR_INVALID_ID	The indicated AVI instance could not be found. Either the AVI ID is invalid or the instance has already been deleted using isavi_ExitAVI() .
IS_AVI_ERR_INVALID_FILE	The file has no valid AVI format.

Related functions

- [isavi_InitAVI\(\)](#)
- [isavi_OpenAVI\(\)](#)
- [isavi_CloseAVI\(\)](#)

4.4.6 isavi_ExitEvent

	
USB 2.0 USB 3.0	-

Syntax

```
INT isavi_ExitEvent (INT nAviID, INT which)
```

Description

`isavi_ExitEvent()` deletes the specified event. The deleted event no longer occurs and needs to be re-created using [isavi_InitEvent\(\)](#) before it can be enabled and used.

Input parameters

<code>nAviID</code>	Instance ID set by the isavi_InitAVI() function.
<code>which</code>	Name of the event to be deleted.
<code>IS_AVI_SET_EVENT_FRAME_SAVED</code>	A new frame was saved to the AVI file.



Return values

<code>IS_AVI_NO_ERR</code>	Function executed successfully.
<code>IS_AVI_ERR_INVALID_ID</code>	The indicated AVI instance could not be found. Either the AVI ID is invalid or the instance has already been deleted using isavi_ExitAVI() .
<code>IS_AVI_ERR_PARAMETER</code>	One of the submitted parameters is outside the valid range.

Related functions

- [isavi_InitEvent\(\)](#)
- [isavi_EnableEvent\(\)](#)
- [isavi_DisableEvent\(\)](#)

4.4.7 isavi_GetAVIFilename

	
USB 2.0 USB 3.0	-

Syntax

```
INT isavi_GetAVIFilename (INT nAviID, char* strName)
```

Description

Using `isavi_GetAVIFilename()`, you can read out the filename of the current AVI file. This function is helpful if an AVI file was opened with the [isavi_OpenAVI\(\)](#) function and a Null parameter was specified.

Note

The functions [isavi_OpenAVI\(\)](#) and `isavi_GetAVIFilename()` do not support UNICODE strings. For this purpose you use the functions [isavi_OpenAVIW\(\)](#) and [isavi_GetAVIFilenameW\(\)](#).

If the AVI file was created using a UNICODE string, only the `isavi_GetAVIFilenameW()` function can return the right file string.

Input parameters

nAviID	Instance ID set by the isavi_InitAVI() function.
strName	Pointer to the memory location where the filename is written to. The allocated memory must be large enough to accommodate the full file path.
NULL	When NULL is passed the function returns the length of the filename.



Return values

IS_AVI_NO_ERR	Function executed successfully.
IS_AVI_ERR_INVALID_ID	The indicated AVI instance could not be found. Either the AVI ID is invalid or the instance has already been deleted using isavi_ExitAVI() .

Related functions

- [isavi_GetAVISize\(\)](#)
- [isavi_OpenAVI\(\)](#)

4.4.8 isavi_GetAVIFilenameW

	
USB 2.0 USB 3.0	-

Syntax

```
INT isavi_GetAVIFilenameW (INT nAviID, wchar_t* strName)
```

Description

Using `isavi_GetAVIFilenameW()`, you can read out the filename of the current AVI file as UNICODE string. This function is helpful if an AVI file was opened with the [isavi_OpenAVIW\(\)](#) function and a `NULL` parameter was specified.

Input parameters

<code>nAviID</code>	Instance ID set by the isavi_InitAVI() function.
<code>strName</code>	Pointer to the memory location where the filename is written to. The allocated memory must be large enough to accommodate the full file path as UNICODE string.
<code>NULL</code>	When <code>NULL</code> is passed the function returns the length of the filename.



Return values

<code>IS_AVI_NO_ERR</code>	Function executed successfully.
<code>IS_AVI_ERR_INVALID_ID</code>	The indicated AVI instance could not be found. Either the AVI ID is invalid or the instance has already been deleted using isavi_ExitAVI() .

Related functions

- [isavi_GetAVISize\(\)](#)
- [isavi_OpenAVIW\(\)](#)

4.4.9 isavi_GetAVISize

	
USB 2.0 USB 3.0	-

Syntax

```
INT isavi_GetAVISize (INT nAviID, float* size)
```

Description

Use `isavi_GetAVISize()` to retrieve the size of the frame sequence saved to the current AVI file.

Input parameters

nAviID	Instance ID set by the isavi_InitAVI() function
size	The size in kBytes



Return values

IS_AVI_NO_ERR	Function executed successfully.
IS_AVI_ERR_INVALID_ID	The indicated AVI instance could not be found. Either the AVI ID is invalid or the instance has already been deleted using isavi_ExitAVI() .

Related functions

- [isavi_GetAVIFileName\(\)](#)

4.4.10 isavi_GetnCompressedFrames

	
USB 2.0 USB 3.0	-

Syntax

```
INT isavi_GetnCompressedFrames (INT nAviID, unsigned long* nFrames)
```

Description

Using `isavi_GetnCompressedFrames()`, you can read out the number of frames saved to the current AVI file.

Input parameters

nAviID	Instance ID set by the isavi_InitAVI() function.
nFrames	The number of frames



Return values

IS_AVI_NO_ERR	Function executed successfully.
IS_AVI_ERR_INVALID_ID	The indicated AVI instance could not be found. Either the AVI ID is invalid or the instance has already been deleted using isavi_ExitAVI() .

Related functions

- [isavi_GetnLostFrames\(\)](#)
- [isavi_ResetFrameCounters\(\)](#)

4.4.11 isavi_GetnLostFrames

	
USB 2.0 USB 3.0	-

Syntax

```
INT isavi_GetnLostFrames (INT nAviID, unsigned long* nFrames)
```

Description

Using `isavi_GetnLostFrames()`, you can read out the number of frames that have been discarded. A frame will be discarded if it cannot be processed because a compression operation is still in progress.

Input parameters

nAviID	Instance ID set by the isavi_InitAVI() function.
nFrames	The number of frames



Return values

IS_AVI_NO_ERR	Function executed successfully.
IS_AVI_ERR_INVALID_ID	The indicated AVI instance could not be found. Either the AVI ID is invalid or the instance has already been deleted using isavi_ExitAVI() .

Related functions

- [isavi_GetnCompressedFrames\(\)](#)
- [isavi_ResetFrameCounters\(\)](#)

4.4.12 isavi_InitAVI

	
USB 2.0 USB 3.0	-

Syntax

```
INT isavi_InitAVI (INT* pnAviID, HIDS hCam)
```

Description

`isavi_InitAVI()` initializes an instance of the uc480 AVI interface. Multiple instances can be created simultaneously.

Input parameters

<code>pnAviID</code>	Pointer. Returns the instance ID which is needed for calling the other uc480 AVI functions.
<code>hCam</code>	Handle of a selected or initialized DCx camera.



Return values

<code>IS_AVI_NO_ERR</code>	Function executed successfully.
<code>IS_AVI_ERR_PARAMETER</code>	One of the submitted parameters is outside the valid range.
<code>IS_AVI_ERR_NO_CODEEC_AVAIL</code>	The maximum number of instances allowed in this system has been reached. It is not possible to create another instance.
<code>IS_AVI_ERR_INVALID_UEYE</code>	No DCx camera was found.

Related functions

- [isavi_ExitAVI\(\)](#)
- [isavi_OpenAVI\(\)](#)
- [isavi_CloseAVI\(\)](#)

4.4.13 isavi_InitEvent

	
USB 2.0 USB 3.0	-

Syntax

```
INT isavi_InitEvent (INT nAviID, INT which)
```

Description

`isavi_InitEvent()` creates the specified event. This includes registering the event object in the uc480 AVI interface and creating an event handler. Before you can use a new event, you must enable it by calling [isavi_EnableEvent\(\)](#).

Input parameters

nAviID	Instance ID set by the isavi_InitAVI() function.
which	Name of the event to be created.
IS_AVI_SET_EVENT_FRAME_SAVED	A new frame was saved to the AVI file.

Return values

IS_AVI_NO_ERR	Function executed successfully.
IS_AVI_ERR_INVALID_ID	The indicated AVI instance could not be found. Either the AVI ID is invalid or the instance has already been deleted using isavi_ExitAVI() .
IS_AVI_ERR_EVENT_FAILED	The event could not be initialized. The Windows <code>SetEvent</code> function failed.
IS_AVI_ERR_PARAMETER	One of the submitted parameters is outside the valid range.

Related functions

- [isavi_ExitEvent\(\)](#)
- [isavi_EnableEvent\(\)](#)
- [isavi_DisableEvent\(\)](#)



Example

Create and enable an event object for the "Frame saved" event:

```
HANDLE hEvent = CreateEvent( NULL, TRUE, FALSE, "" );
if ( hEvent != NULL )
{
    isavi_InitEvent( AviDest, hEvent, IS_AVI_SET_EVENT_FRAME_SAVED );
    isavi_EnableEvent( AviDest, IS_AVI_SET_EVENT_FRAME_SAVED );

    if ( WaitForSingleObject( hEvent, 1000 ) == WAIT_OBJECT_0 )
    {
        //Frame was captured successfully...
    }
    isavi_DisableEvent( AviDest, IS_AVI_SET_EVENT_FRAME_SAVED );
    isavi_ExitEvent( AviDest, IS_AVI_SET_EVENT_FRAME_SAVED );
}
```


4.4.14 isavi_OpenAVI

	
USB 2.0 USB 3.0	-

Syntax

```
INT isavi_OpenAVI (INT nAviID, const char* strFileName)
```

Description

isavi_OpenAVI() opens a new or existing AVI file.

Note

The functions [isavi_OpenAVI\(\)](#) and [isavi_GetAVIFileName\(\)](#) do not support UNICODE strings. For this purpose you use the functions [isavi_OpenAVIW\(\)](#) and [isavi_GetAVIFileNameW\(\)](#).

Input parameters

nAviID	Instance ID set by the isavi_InitAVI() function.
pFileName	Pointer to the name to be assigned to the AVI file. If <code>NULL</code> is passed, the "Open File" dialog is displayed.



Return values

IS_AVI_NO_ERR	Function executed successfully.
IS_AVI_ERR_INVALID_ID	The indicated AVI instance could not be found. Either the AVI ID is invalid or the instance has already been deleted using isavi_ExitAVI() .
IS_AVI_ERR_CAPTURE_RUNNING	Another capturing operation is in progress or an AVI file is still open.
IS_AVI_ERR_INVALID_FILE	The file has no valid AVI format.
IS_AVI_ERR_NEW_FAILED	No memory could be allocated for the AVI file.
IS_AVI_ERR_CREATESTREAM	No AVI stream could be created.

Related functions

- [isavi_GetAVIFileName\(\)](#)
- [isavi_CloseAVI\(\)](#)
- [isavi_InitAVI\(\)](#)
- [isavi_ExitAVI\(\)](#)

4.4.15 isavi_OpenAVIW

	
USB 2.0 USB 3.0	-

Syntax

```
INT isavi_OpenAVIW (INT nAviID, const wchar_t* strFileName)
```

Description

isavi_OpenAVIW() opens a new or existing AVI file.

Input parameters

nAviID	Instance ID set by the isavi_InitAVI() function.
pFileName	Pointer to the name to be assigned to the AVI file. The file name is passed as UNICODE string. If <code>NULL</code> is passed, the "Open File" dialog is displayed.



Return values

IS_AVI_NO_ERR	Function executed successfully.
IS_AVI_ERR_INVALID_ID	The indicated AVI instance could not be found. Either the AVI ID is invalid or the instance has already been deleted using isavi_ExitAVI() .
IS_AVI_ERR_CAPTURE_RUNNING	Another capturing operation is in progress or an AVI file is still open.
IS_AVI_ERR_INVALID_FILE	The file has no valid AVI format.
IS_AVI_ERR_NEW_FAILED	No memory could be allocated for the AVI file.
IS_AVI_ERR_CREATESTREAM	No AVI stream could be created.

Related functions

- [isavi_GetAVIFileNameW\(\)](#)
- [isavi_CloseAVI\(\)](#)
- [isavi_InitAVI\(\)](#)
- [isavi_ExitAVI\(\)](#)

4.4.16 isavi_ResetFrameCounters

	
USB 2.0 USB 3.0	-

Syntax

```
INT isavi_ResetFrameCounters (INT nAviID)
```

Description

`isavi_ResetFrameCounters()` resets the counters for saved and discarded images.

Input parameters

nAviID	Instance ID set by the isavi_InitAVI() function.
--------	--



Return values

IS_AVI_NO_ERR	Function executed successfully.
IS_AVI_ERR_INVALID_ID	The indicated AVI instance could not be found. Either the AVI ID is invalid or the instance has already been deleted using isavi_ExitAVI() .

Related functions

- [isavi_GetnCompressedFrames\(\)](#)
- [isavi_GetnLostFrames\(\)](#)

4.4.17 isavi_SetFrameRate

	
USB 2.0 USB 3.0	-

Syntax

```
INT isavi_SetFrameRate (INT nAviID, double fr)
```

Description

`isavi_SetFrameRate()` sets the frame rate for AVI capturing. You can set the frame rate after opening the AVI file. This value does not have to be equal to the frame rate set for the DCx camera.

Input parameters

nAviID	Instance ID set by the isavi_InitAVI() function.
fr	The frame rate to be set. Default = 25.0



Return values

IS_AVI_NO_ERR	Function executed successfully.
IS_AVI_ERR_INVALID_ID	The indicated AVI instance could not be found. Either the AVI ID is invalid or the instance has already been deleted using isavi_ExitAVI() .
IS_AVI_ERR_WRITE_INFO	The AVI file could not be modified.
IS_AVI_ERR_INVALID_FILE	The file has no valid AVI format.

Related functions

- [isavi_SetImageQuality\(\)](#)
- [isavi_SetImageSize\(\)](#)

4.4.18 isavi_SetImageQuality

	
USB 2.0 USB 3.0	-

Syntax

```
INT isavi_SetImageQuality (INT nAviID, INT q)
```

Description

`isavi_SetImageQuality()` indicates the quality for the frames to be compressed. You can change the image quality at any time; it then applies to all subsequent frames. For compression, the system uses the JPEG algorithm.

Input parameters

nAviID	Instance ID set by the isavi_InitAVI() function.
q	Image quality [1 = lowest ... 100 = highest]



Return values

IS_AVI_NO_ERR	Function executed successfully.
IS_AVI_ERR_INVALID_ID	The indicated AVI instance could not be found. Either the AVI ID is invalid or the instance has already been deleted using isavi_ExitAVI() .
IS_AVI_ERR_INVALID_VALUE	The <code>q</code> parameter is outside the range of 1...100.
IS_AVI_ERR_INVALID_FILE	The file has no valid AVI format.

Related functions

- [isavi_SetFrameRate\(\)](#)
- [isavi_SetImageSize\(\)](#)

4.4.19 isavi_SetImageSize

	
USB 2.0 USB 3.0	-

Syntax

```
INT isavi_SetImageSize (INT nAviID, INT cMode, INT Width, INT Height,
                       INT PosX, INT PosY, INT LineOffset)
```

Description

`isavi_SetImageSize()` sets the size and position of the area of interest which will be saved to the AVI file. Only the defined area of interest of each frame will be saved. In addition, this function specifies the input color format of the frames. You define these settings only once for the entire video.

Note

The supported input color formats are RGB32, RGB24, Y8 and raw Bayer. The output file will always be in RGB24 format, regardless of the input data format. For further information on the structure of the different color formats, see the [Appendix: Color and memory formats](#) section.

Attention

When an area of interest is used, the width (`width`) and height (`Height`) of the AOI must be at least 16 pixel. The AOI width must be a multiple of 8.

Input parameters

<code>nAviID</code>	Instance ID set by the isavi_InitAVI() function.
<code>cMode</code>	Color format of the input frames captured by the DCx Camera.
<code>Width</code>	Width of the entire frame or of the area of interest.
<code>Height</code>	Height of the entire frame or of the area of interest.
<code>PosX</code>	X position (offset) of the area of interest.
<code>PosY</code>	Y position (offset) of the area of interest.
<code>LineOffset</code>	Line increment. The line increment is the difference between the width of the entire frame (in pixel) and the area of interest (in pixel).

Return values

<code>IS_AVI_NO_ERR</code>	Function executed successfully.
<code>IS_AVI_ERR_INVALID_ID</code>	The indicated AVI instance could not be found. Either the AVI ID is invalid or the instance has already been deleted using isavi_ExitAVI() .
<code>IS_AVI_ERR_INVALID_FILE</code>	The file has no valid AVI format.
<code>IS_AVI_ERR_CAPTURE_RUNNING</code>	Another capturing operation is in progress or an AVI file is still open.
<code>IS_AVI_ERR_ALLOC_MEMORY</code>	No memory could be allocated.
<code>IS_AVI_ERR_INVALID_CM</code>	The submitted color mode is not supported for AVI capturing.

IS_AVI_ERR_INVALID_SIZE	The submitted size is invalid.
IS_AVI_ERR_INVALID_POSITION	The submitted position is invalid.

Related functions

- [isavi_SetFrameRate\(\)](#)
- [isavi_SetImageQuality\(\)](#)

Example



```
// Query image buffer geometry
int nWidth, nHeight, nBits, nPitch;
is_InquireImageMem ( hCam, pLast, nImageID,
                    &nWidth, &nHeight,
                    &nBits, &nPitch);

INT nOffsetX = is_SetImagePos ( hCam, IS_GET_IMAGE_POS_X_ABS, 0);
INT nOffsetY = is_SetImagePos ( hCam, IS_GET_IMAGE_POS_Y_ABS, 0);

// Derive pixel pitch from buffer byte pitch
INT nPitchPx=0;
nPitchPx = (nPitch * 8 ) / nBits;

INT nAwiWidth = nWidth / 8 * 8; // Width must be multiple of 8
INT LineOffsetPx = nPitchPx - nAwiWidth ;
isavi_SetImageSize( nAwiId, m_cMode,
                   nAwiWidth, nHeight,
                   nOffsetX, nOffsetY,
                   LineOffsetPx);
```

4.4.20 isavi_StartAVI

	
USB 2.0 USB 3.0	-

Syntax

```
INT isavi_StartAVI (INT nAwiID)
```

Description

isavi_StartAVI() starts the image capture thread.

Input parameters

nAwiID	Instance ID set by the isavi_InitAVI() function.
--------	--



Return values

IS_AVI_NO_ERR	Function executed successfully.
IS_AVI_ERR_INVALID_ID	The indicated AVI instance could not be found. Either the AVI ID is invalid or the instance has already been deleted using isavi_ExitAVI() .
IS_AVI_ERR_INVALID_FILE	The file has no valid AVI format.
IS_AVI_ERR_PLAY_RUNNING	A playback is already running.

Related functions

- [isavi_StopAVI\(\)](#)
- [isavi_InitEvent\(\)](#)
- [isavi_ExitAVI\(\)](#)

4.4.21 isavi_StopAVI

	
USB 2.0 USB 3.0	-

Syntax

```
INT isavi_StopAVI (INT nAviID)
```

Description

`isavi_StopAVI()` stops the image capture thread. Subsequent calls of [isavi_AddFrame\(\)](#) will be ignored.

Input parameters

nAviID	Instance ID set by the isavi_InitAVI() function.
--------	--

Return values

IS_AVI_NO_ERR	Function executed successfully.
IS_AVI_ERR_CAPTURE_NOT_RUNNING	No capturing operation is running or no AVI file is opened.
IS_AVI_ERR_INVALID_ID	The indicated AVI instance could not be found. Either the AVI ID is invalid or the instance has already been deleted using isavi_ExitAVI() .
IS_AVI_ERR_INVALID_FILE	The file has no valid AVI format.
IS_AVI_ERR_PLAY_NOT_RUNNING	No playback is running.



Related functions

- [isavi_StartAVI\(\)](#)
- [isavi_InitEvent\(\)](#)
- [isavi_ExitAVI\(\)](#)

4.5 RAW function descriptions

The functions of the `uEye_tools.dll` enable you to save images captured with the uEye camera as an image sequences to an RAW file or open it.

4.5.1 `israw_AddFrame`

	
USB 2.0 USB 3.0 GigE	-

Syntax

```
INT israw_AddFrame(UINT unFileID, const char* pcData, UINT64 unTimestampDevice)
```

Description

`israw_AddFrame()` adds a new frame to an image sequence.

Input parameters

<code>unFileID</code>	Instance ID set by the israw_InitFile() function
<code>pcData</code>	Pointer to the image data
<code>unTimestampDevice</code>	Data timestamp



Return values

<code>IS_AVI_NO_ERR</code>	No error, function executed successfully.
<code>IS_AVI_ERR_EXCEPTION</code>	An error has occurred.
<code>IS_AVI_ERR_FILE_NOT_OPEN</code>	The file is not open.
<code>IS_AVI_ERR_GENERIC</code>	Generic error
<code>IS_AVI_ERR_INVALID_ID</code>	The indicated instance could not be found. Either the ID is invalid or the instance has already been deleted using israw_ExitFile() .
<code>IS_AVI_ERR_NOT_SUPPORTED</code>	The file was opened in "read-only" mode (see israw_InitFile()).

Related functions

- [israw_GetFrame\(\)](#)

4.5.2 `israw_CloseFile`

	
USB 2.0 USB 3.0 GigE	-

Syntax

```
INT israw_CloseFile(UINT unFileID)
```

Description

`israw_CloseFile()` closes a RAW file which was opened using [israw_OpenFile\(\)](#).

Input parameters

unFileID	Instance ID set by the israw_InitFile() function
----------	--



Return values

IS_AVI_NO_ERR	No error, function executed successfully.
IS_AVI_ERR_EXCEPTION	An error has occurred.
IS_AVI_ERR_GENERIC	Generic error
IS_AVI_ERR_INVALID_ID	The indicated instance could not be found. Either the ID is invalid or the instance has already been deleted using israw_ExitFile() .

Related functions

- [israw_OpenFile\(\)](#)

4.5.3 israw_ExitFile

	
USB 2.0 USB 3.0 GigE	-

Syntax

```
INT israw_ExitFile(UINT unFileID)
```

Description

`israw_ExitFile()` terminates and deletes the instance of the uEye RAW interface.

Input parameters

unFileID	Instance ID set by the israw_InitFile() function
----------	--



Return values

IS_AVI_NO_ERR	No error, function executed successfully.
IS_AVI_ERR_EXCEPTION	An error has occurred.
IS_AVI_ERR_GENERIC	Generic error

Related functions

- [israw_InitFile\(\)](#)

4.5.4 israw_GetFrame

	
USB 2.0 USB 3.0 GigE	-

Syntax

```
INT israw_GetFrame(UINT unFileID, char* pData, UINT64* punTimestampDevice)
```

Description

`israw_GetFrame()` returns the image data from an image sequence.

Input parameters

<code>unFileID</code>	Instance ID set by the israw_InitFile() function
<code>pData</code>	Pointer to the image data
<code>punTimestampDevice</code>	Data timestamp



Return values

<code>IS_AVI_NO_ERR</code>	No error, function executed successfully.
<code>IS_AVI_ERR_EXCEPTION</code>	An error has occurred.
<code>IS_AVI_ERR_FILE_NOT_OPEN</code>	The file is not open.
<code>IS_AVI_ERR_GENERIC</code>	Generic error
<code>IS_AVI_ERR_INVALID_ID</code>	The indicated instance could not be found. Either the ID is invalid or the instance has already been deleted using israw_ExitFile() .
<code>IS_AVI_ERR_NOT_SUPPORTED</code>	The file was opened in "read-only" mode (see israw_InitFile()).
<code>IS_AVI_ERR_READ_FAILED</code>	The file could not be read.

Related functions

- [israw_SeekFrame\(\)](#)

4.5.5 israw_GetImageInfo

	
USB 2.0 USB 3.0 GigE	-

Syntax

```
INT israw_GetImageInfo(UINT unFileID, UINT* punWidth, UINT* punHeight, UINT* punBitsPerPixel)
```

Description

`israw_GetImageInfo()` returns the width, height, and bit depth of the data which are saved in the RAW file.

Input parameters

<code>unFileID</code>	Instance ID set by the israw_InitFile() function
<code>punWidth</code>	Image width
<code>punHeight</code>	Image height
<code>punBitsPerPixel</code>	Bit depth (bits per pixel)

Return values



<code>IS_AVI_NO_ERR</code>	No error, function executed successfully.
<code>IS_AVI_ERR_EXCEPTION</code>	An error has occurred.

IS_AVI_ERR_FILE_NOT_OPEN	The file is not open.
IS_AVI_ERR_GENERIC	Generic error
IS_AVI_ERR_INVALID_ID	The indicated instance could not be found. Either the ID is invalid or the instance has already been deleted using israw_ExitFile() .
IS_AVI_ERR_PARAMETER	One of the submitted parameters is outside the valid range.

Related functions

- [israw_SetImageInfo\(\)](#)

4.5.6 israw_GetSize

	
USB 2.0 USB 3.0 GigE	-

Syntax

```
INT israw_GetSize(UINT unFileID, float* pfSize)
```

Description

`israw_GetSize()` returns the size of the saved image sequence in the current RAW file.

Input parameters

unFileID	Instance ID set by the israw_InitFile() function
pfSize	Size in kBytes



Return values

IS_AVI_NO_ERR	No error, function executed successfully.
IS_AVI_ERR_EXCEPTION	An error has occurred.
IS_AVI_ERR_FILE_NOT_OPEN	The file is not open.
IS_AVI_ERR_GENERIC	Generic error
IS_AVI_ERR_INVALID_ID	The indicated instance could not be found. Either the ID is invalid or the instance has already been deleted using israw_ExitFile() .

Related functions

- [israw_OpenFile\(\)](#)

4.5.7 israw_InitFile

	
USB 2.0 USB 3.0 GigE	-

Syntax

```
INT israw_InitFile(UINT* punFileID, INT nAccessMode)
```

Description

`israw_InitFile()` initializes an instance of the uEye RAW interface. Multiple instances can be created simultaneously.

Input parameters

<code>punFileID</code>	Pointer in which the instance ID is returned. This ID is needed for calling other functions.
<code>nAccessMode</code>	File access mode: <ul style="list-style-type: none"> • Read: <code>IS_FILE_ACCESS_MODE_READ</code> • Write: <code>IS_FILE_ACCESS_MODE_WRITE</code>



Return values

<code>IS_AVI_NO_ERR</code>	No error, function executed successfully.
<code>IS_AVI_ERR_EXCEPTION</code>	An error has occurred.
<code>IS_AVI_ERR_GENERIC</code>	Generic error
<code>IS_AVI_ERR_PARAMETER</code>	One of the submitted parameters is outside the valid range.
<code>IS_AVI_ERR_NEW_FAILED</code>	The maximum number of instances have been reached.

Related functions

- [israw_ExitFile\(\)](#)

4.5.8 israw_OpenFile

	
USB 2.0 USB 3.0 GigE	-

Syntax

```
INT israw_OpenFile(UINT unFileID, const char* strFileName)
```

Description

`israw_OpenFile()` opens a new or existing RAW file.

Input parameters

<code>unFileID</code>	Instance ID set by the israw_InitFile() function
<code>strFileName</code>	RAW file name

Return values



<code>IS_AVI_NO_ERR</code>	No error, function executed successfully.
<code>IS_AVI_ERR_EXCEPTION</code>	An error has occurred.
<code>IS_AVI_ERR_FILE_NOT_OPEN</code>	The file is not open.
<code>IS_AVI_ERR_GENERIC</code>	Generic error
<code>IS_AVI_ERR_INVALID_ID</code>	The indicated instance could not be found. Either the ID is invalid or the instance has

	already been deleted using israw_ExitFile() .
IS_AVI_ERR_INVALID_VALUE	No valid file name was selected.

Related functions

- [israw_CloseFile\(\)](#)

4.5.9 israw_SeekFrame

	
USB 2.0 USB 3.0 GigE	-

Syntax

```
INT israw_SeekFrame(UINT unFileID, UINT unFrame)
```

Decsription

`israw_SeekFrame()` jumps to a specific image within an image sequence.

Input parameters

unFileID	Instance ID set by the israw_InitFile() function
unFrame	Image to be jumped to (index starts with "0")



Return values

IS_AVI_NO_ERR	No error, function executed successfully.
IS_AVI_ERR_EXCEPTION	An error has occurred.
IS_AVI_ERR_FILE_NOT_OPEN	The file is not open.
IS_AVI_ERR_GENERIC	Generic error
IS_AVI_ERR_INVALID_ID	The indicated instance could not be found. Either the ID is invalid or the instance has already been deleted using israw_ExitFile() .
IS_AVI_ERR_NOT_SUPPORTED	The file was opened in "read-only" mode (see israw_InitFile()).
IS_AVI_ERR_SEEK_FAILED	The image could not be found.

Related functions

- [israw_GetFrame\(\)](#)

4.5.10 israw_SetImageInfo

	
USB 2.0 USB 3.0 GigE	-

Syntax

```
INT israw_SetImageInfo(UINT unFileID, UINT unWidth, UINT unHeight, UINT unBitsPerPixel)
```

Description

`israw_SetImageInfo()` sets the width, height, and bit depth of the data that is saved in the RAW file (see [is_raw_AddFrame\(\)](#)).

Input parameters

<code>unFileID</code>	Instance ID set by the israw_InitFile() function
<code>Width</code>	Image width
<code>Height</code>	Image height
<code>BitsPerPixel</code>	Bit depth (bits per pixel)

Return values

<code>IS_AVI_NO_ERR</code>	No error, function executed successfully.
<code>IS_AVI_ERR_EXCEPTION</code>	An error has occurred.
<code>IS_AVI_ERR_GENERIC</code>	Generic error
<code>IS_AVI_ERR_INVALID_ID</code>	The indicated instance could not be found. Either the ID is invalid or the instance has already been deleted using israw_ExitFile() .
<code>IS_AVI_ERR_INVALID_VALUE</code>	No valid file name was selected.
<code>IS_AVI_ERR_NOT_SUPPORTED</code>	The file was opened in "read-only" mode (see israw_InitFile()).

Related functions

- [israw_AddFrame\(\)](#)
- [israw_GetImageInfo\(\)](#)

4.6 Obsolete functions

We are continuously extending and enhancing the uc480 API. The resulting product upgrades sometimes require replacing obsolete functions with new ones. We recommend against using the obsolete functions. They will continue to be supported for reasons of backward compatibility, but they will not be documented any longer.



Note on older functions

If it is necessary to continue working with the older functions, it is possible to add the `uc480_deprecated.h` header file additionally to the `uc480.h` header file. The `uc480_deprecated.h` header file contains all obsolete function definitions and constants which are no longer part of the `uc480.h` header file.

The following table lists the obsolete functions and indicates the recommended alternatives. See also [History of API functions](#).

Obsolete function	Recommended alternative	No longer documented since version	New function valid from version
<code>is_ConvertImage()</code>	is_Convert()	4.50	4.00
<code>is_DisableDDOverlay()</code>	is_DirectRenderer()	4.50	3.40
<code>is_EnableDDOverlay()</code>		4.50	
<code>is_GetCameraLUT</code>	is_LUT()	4.50	4.40
<code>is_GetCameraType()</code>	is_GetCameraInfo()	4.50	2.00

<code>is_GetCaptureErrorInfo()</code>	is_CaptureStatus()	4.50	3.90
<code>is_GetDC()</code>	is_DirectRenderer()	4.50	3.40
<code>is_GetDDOvlSurface()</code>		4.50	
<code>is_GetEthDeviceInfo</code>	is_DeviceInfo()	4.50	4.00
<code>is_GetExposureRange</code>	is_Exposure()	4.50	3.90
<code>is_GetGlobalFlashDelays()</code>	is_IO()	4.50	3.90
<code>is_GetLastMemorySequence()</code>	The uEye memory board is not supported any longer (see below).	3.30	
<code>is_GetMemorySequenceWindow()</code>			
<code>is_GetNumberOfMemoryImages()</code>			
<code>is_GetOsVersion</code>	-	4.50	
<code>is_GetPixelClockRange()</code>	is_PixelClock()	4.50	4.00
<code>is_GetRevisionInfo()</code>	is_GetCameraInfo()	3.20	2.00
<code>is_GetHWGain()</code>	is_SetAutoParameter()	3.31	2.20
<code>is_HideDDOverlay()</code>	is_DirectRenderer()	4.50	3.40
<code>is_IsMemoryBoardConnected()</code>	The uEye memory board is not supported any longer (see below).	3.30	
<code>is_LoadBadPixelCorrectionTable()</code>	is_HotPixel()	4.50	3.80
<code>is_LoadImage()</code>	is_ImageFile()	4.50	4.00
<code>is_LoadImageMem()</code>		4.50	
<code>is_LoadParameters()</code>	is_ParameterSet()	4.50	4.00
<code>is_LockDDMem()</code>	is_DirectRenderer()	4.50	3.40
<code>is_LockDDOverlayMem()</code>		4.50	
<code>is_MemoryFreezeVideo()</code>	The uEye memory board is not supported any longer (see below).	3.30	
<code>is_ReleaseDC()</code>	is_DirectRenderer()	4.50	3.40
<code>is_ResetCaptureErrorInfo()</code>	is_CaptureStatus()	4.50	3.90
<code>is_ResetMemory()</code>	The uEye memory board is not supported any longer (see below).	3.30	
<code>is_SaveBadPixelCorrectionTable()</code>	is_HotPixel()	4.50	3.80
<code>is_SaveImage()</code>	is_ImageFile()	4.50	4.00
<code>is_SaveImageEx()</code>		4.50	
<code>is_SaveImageMem()</code>		4.50	
<code>is_SaveImageMemEx()</code>		4.50	
<code>is_SaveParameters()</code>	is_ParameterSet()	4.50	4.00
<code>is_SetAOI()</code>	is_AOI()	4.50	3.80
<code>is_SetAutoCfgIpSetup()</code>	<code>is_IpConfig()</code>	4.50	3.82
<code>is_SetBadPixelCorrection()</code>	is_HotPixel()	4.50	3.80
<code>is_SetBadPixelCorrectionTable</code>		4.50	

()				
is_SetBayerConversion()	is_SetColorConverter()	4.50	3.30	
is_SetBlCompensation()	is_Blacklevel()	4.50	4.01	
is_SetBrightness()	is_Exposure() is_SetHardwareGain()	3.40	3.90 2.00	
is_SetCameraLUT	is_LUT()	4.50	4.40	
is_SetContrast()	is_Gamma() is_SetHardwareGamma() is_Blacklevel()	3.40	4.40 2.00 4.01	
is_SetConvertParam()	is_Convert()	4.50	4.00	
is_SetDDUpdateTime()	is_DirectRenderer()	4.50	3.40	
is_SetEdgeEnhancement()	is_EdgeEnhancement()	4.50	4.00	
is_SetExposureTime()	is_Exposure()	4.50	3.90	
is_SetFlashDelay()	is_IO()	4.50	3.90	
is_SetFlashStrobe()		4.50		
is_SetGamma()	is_Gamma()	4.50	4.40	
is_SetHwnd()	is_DirectRenderer()	4.50	3.40	
is_SetImageAOI()		4.50	3.80	
is_SetImageSize()		is_AOI()		4.50
is_SetImagePos()		4.50		
is_SetIO()	is_IO()	4.50	3.90	
is_SetIOMask()		4.50		
is_SetKeyColor()	is_DirectRenderer()	4.50	3.40	
is_SetLED()	is_IO()	4.50	3.90	
is_SetMemoryMode()	The uEye memory board is not supported any longer (see below).	3.30		
is_SetPixelClock()	is_PixelClock()		4.00	
is_SetWhiteBalance()	is_SetAutoParameter()	3.31	2.20	
is_SetWhiteBalanceMultipliers()				
is_SetPersistentIpCfg()	is_IpConfig()	4.50	3.82	
is_ShowDDOverlay()	is_DirectRenderer()	4.50	3.40	
is_StealVideo()		4.50		
is_TransferImage()	The uEye memory board is not supported any longer (see below).	3.30		
is_TransferMemorySequence()				
is_UnlockDDMem()		4.50	3.40	
is_UnlockDDOverlayMem()	is_DirectRenderer()	4.50		
is_UpdateDisplay()		4.50		



The `is_SetWhiteBalance()` and `is_SetWhiteBalanceMultipliers()` functions have been completely replaced by the [is_SetAutoParameter\(\)](#) function and are no longer supported by the uc480 API.

4.7 Programming Notes

Note

Parameter Validity

Functions that refer to an initialized camera have the camera handle `HIDS hCam` as the first parameter. All parameters that are set using these functions remain valid for as long as the handle is valid, that is, until you close the corresponding camera or exit the program. The next time you open the camera, it is initialized with the defaults again.

Attention

All input parameters of a function have to be initialized with valid values before the function is called; this also applies to parameters that are not used. Variables can be preset with '0', for example. For unused parameters, the NULL pointer has to be passed.

The uc480.h header file

The `uc480.h` header file contains all the definitions and constants needed for the uc480 API. After the installation of the uc480 drivers you will find this file in the directory:

- Windows: `C:\Program Files\Thorlabs\DCx Cameras\Develop\include`
- Linux: `/usr/include`

Note

Note on older functions

If it is necessary to continue working with the older functions, it is possible to add the `uc480_deprecated.h` header file additionally to the `uc480.h` header file. The `uc480_deprecated.h` header file contains all obsolete function definitions and constants which are no longer part of the `uc480.h` header file.

See also:

- [Programming in C / C++](#)
- [Programming in C#](#)
- [Programming in VB.NET](#)
- [Programming in Delphi](#)
- [Programming with ActiveX](#)
- [Thread programming](#)

4.7.1 Programming in C/C++

For programming with the uc480 API, we suggest to use the C/C++ programming language. This programming language offers efficient access to all functions of the uc480 API. Enabling access to image memory contents through pointers, C/C++ is especially suitable for image processing applications.

Most of the uc480 sample programs were created in Microsoft Visual Studio using the C++ programming language.

Hint

We suggest that you keep the function libraries (DLL, AX and OCX files) in the default directory. After the installation, these files reside e.g. under Window (32 bit) in `C:\Windows\System32\`.

Copying these files to other locations may result in version conflicts.

Required Files

In order to access the uc480 API, make sure to include the following files in your project:

- Header file: `uc480.h`
- Lib file: `uc480.lib`
- Function library (DLL): `uc480.dll`

In order to access the DCx Camera AVI functions, make sure to include the following files in your project:

- Header file: `uc480_tools.h`
- Lib file: `uc480_tools.lib`
- Function library (DLL): `uc480_tools.dll`

In order to access the uc480 DirectShow functions, make sure to include the following files in your project:

- Header file: `uc480CaptureInterface.h`
- DirectShow interface: `uc480capture.ax`

Note

Programming under Linux

In order to access the uc480 API, make sure to include the following files in your project:

- Header file: `uc480.h`
- Library: `libuc480.so`

4.7.2 Programming in C#

We suggest to use the C# programming language for the creation of visualization applications. While it is possible to access image memory contents, doing so is more tedious than in C/C++ due to the 'managed code'. To access image memory contents in C#, you can use 'unsafe code' or the 'Marshal class'. Some system-level functions, such as Windows event handling, can be integrated using the Windows API.

The uc480 SDK includes sample programs for programming with Microsoft Visual Studio in the C# programming language.

Required Files

In order to access the uc480 API in C#, make sure to include the following files in your project:

- Header file: `uc480.cs`
- Function library (DLL): `uc480.dll`

In order to access the uc480 AVI functions in C#, make sure to include the following files in your project:

- Header file: `uc480_tools.cs`
- Function library (DLL): `uc480_tools.dll`

Hint

We suggest that you keep the function libraries (DLL, AX and OCX files) in the default directory. After the installation, these files reside e.g. under Window (32 bit) in `C:\Windows\System32\`.

Copying these files to other locations may result in version conflicts.

4.7.3 Programming in VB.NET

We suggest to use the Visual Basic programming language for the creation of applications which are exclusively used for visualization purposes. The access to image memory contents is extremely tedious due to the missing pointer arithmetics.

We suggest to use the [uc480 ActiveX component](#) when programming in VB.Net. The uc480 SDK includes a sample program for programming with Microsoft Visual Studio in the VB.NET programming language using the uc480 ActiveX component.

The constants can be looked up in the `uc480.h` file.

4.7.4 Programming in Delphi

The uc480 SDK does not provide direct integration of the uc480 API for the Delphi programming language. In order to use the uc480 API in Delphi, you need to create separate header files. We suggest to use the uc480 ActiveX component (see also [Programming with ActiveX](#)) when programming in Delphi.

Attention

To use the `uc480_api.dll` in Delphi, the `cdecl` calling convention has to be used.

Hint

We suggest that you keep the function libraries (DLL, AX and OCX files) in the default directory. After the installation, these files reside e.g. under Window (32 bit) in `C:\Windows\System32\`.

Copying these files to other locations may result in version conflicts.

4.7.5 Programming with ActiveX

The uc480 SDK comes with an ActiveX component that allows you to use almost all functions of the DCx Camera. Programming the uc480 ActiveX component is described in the uc480 ActiveX Manual. After the installation, you will find this manual in the `C:\Programs\uc480\Help` directory.

Note

ActiveX is only available on Windows systems.

Required Files

In order to access the uc480 ActiveX component, make sure to include the following file in your project:

- ActiveX control: `uc480Cam.ocx`

Hint

We suggest that you keep the function libraries (DLL, AX and OCX files) in the default directory. After the installation, these files reside e.g. under Window (32 bit) in `C:\Windows\System32\`.

Copying these files to other locations may result in version conflicts.

4.7.6 Thread Programming

In general, the uc480 API is thread-safe. This means that the uc480 API can be accessed by multiple threads simultaneously. Simultaneous attempts to call the same function are recognized and prevented by the driver.

Note

Multi-threading

We recommend that you call the following functions from one thread per camera in order to avoid unpredictable behavior of the application.

- [is_InitCamera\(\)](#)
- [is_SetDisplayMode\(\)](#)
- [is_ExitCamera\(\)](#)

Attention

Using USB cameras under Windows

The following events require a Windows message loop. This message loop has to be executed by the thread that loads the uc480 API. The message loop is usually provided by the application window. In some cases, the message loop might not be created automatically (e.g. in console applications). In this case you will need to implement the message loop yourself.

This applies to the following uc480 events:

- IS_SET_EVENT_REMOVED
- IS_SET_EVENT_REMOVAL
- IS_SET_EVENT_DEVICE_RECONNECTED
- IS_SET_EVENT_NEW_DEVICE

If no message loop exists, a USB camera will not be automatically detected after reconnecting.

4.8 Lists

- [Complete list of all return values](#)
- [Error codes of AVI functions](#)
- [Linux: not supported functions](#)

4.8.1 Complete List of All Return Values

No	Error	Description
-1	IS_NO_SUCCESS	General error message
0	IS_SUCCESS	Function executed successfully
1	IS_INVALID_CAMERA_HANDLE	Invalid camera handle Most of the uc480 SDK functions expect the camera handle as the first parameter.
2	IS_IO_REQUEST_FAILED	An IO request from the uc480 driver failed. Possibly the versions of the uc480.dll (API) and the driver file (uc480_usb.sys) do not match.
3	IS_CANT_OPEN_DEVICE	An attempt to initialize or select the camera failed (no camera connected or initialization error).
11	IS_CANT_OPEN_REGISTRY	Error opening a Windows registry key
12	IS_CANT_READ_REGISTRY	Error reading settings from the Windows registry
15	IS_NO_IMAGE_MEM_ALLOCATED	The driver could not allocate memory.
16	IS_CANT_CLEANUP_MEMORY	The driver could not release the allocated memory.
17	IS_CANT_COMMUNICATE_WITH_DRIVER	Communication with the driver failed because no driver has been loaded.
18	IS_FUNCTION_NOT_SUPPORTED_YET	The function is not supported yet.
32	IS_INVALID_CAPTURE_MODE	The function can not be executed in the current camera operating mode (free run, trigger or standby).
49	IS_INVALID_MEMORY_POINTER	Invalid pointer or invalid memory ID
50	IS_FILE_WRITE_OPEN_ERROR	File cannot be opened for writing or reading.
51	IS_FILEREAD_OPEN_ERROR	The file cannot be opened.
52	IS_FILE_READ_INVALID_BMP_ID	The specified file is not a valid bitmap file.
53	IS_FILE_READ_INVALID_BMP_SIZE	The bitmap size is not correct (bitmap too large).
108	IS_NO_ACTIVE_IMG_MEM	No active image memory available. You must set the memory to active using the <code>is_SetImageMem()</code> function or create a sequence using the <code>is_AddToSequence()</code> function.
112	IS_SEQUENCE_LIST_EMPTY	The sequence list is empty and cannot be deleted.

No	Error	Description
113	IS_CANT_ADD_TO_SEQUENCE	The image memory is already included in the sequence and cannot be added again.
117	IS_SEQUENCE_BUF_ALREADY_LOCKED	The memory could not be locked. The pointer to the buffer is invalid.
118	IS_INVALID_DEVICE_ID	The device ID is invalid. Valid IDs start from 1 for USB cameras, and from 1001 for GigE cameras.
119	IS_INVALID_BOARD_ID	The board ID is invalid. Valid IDs range from 1 through 255.
120	IS_ALL_DEVICES_BUSY	All cameras are in use.
122	IS_TIMED_OUT	A timeout occurred. An image capturing process could not be terminated within the allowable period.
123	IS_NULL_POINTER	Invalid array
125	IS_INVALID_PARAMETER	One of the submitted parameters is outside the valid range or is not supported for this sensor or is not available in this mode.
127	IS_OUT_OF_MEMORY	No memory could be allocated.
129	IS_ACCESS_VIOLATION	An internal error has occurred.
139	IS_NO_USB20	The camera is connected to a port which does not support the USB 2.0 high-speed standard. Cameras without a memory board cannot be operated on a USB 1.1 port.
140	IS_CAPTURE_RUNNING	A capturing operation is in progress and must be terminated before you can start another one.
145	IS_IMAGE_NOT_PRESENT	The requested image is not available in the camera memory or is no longer valid.
148	IS_TRIGGER_ACTIVATED	The function cannot be used because the camera is waiting for a trigger signal.
151	IS_CRC_ERROR	A CRC error-correction problem occurred while reading the settings.
152	IS_NOT_YET_RELEASED	This function has not been enabled yet in this version.
153	IS_NOT_CALIBRATED	The camera does not contain any calibration data.
154	IS_WAITING_FOR_KERNEL	The system is waiting for the kernel driver to respond.
155	IS_NOT_SUPPORTED	The camera model used here does not support this function or setting.
156	IS_TRIGGER_NOT_ACTIVATED	The function is not possible as trigger is disabled.
157	IS_OPERATION_ABORTED	The dialog was canceled without a selection so that no file could be saved.

No	Error	Description
158	IS_BAD_STRUCTURE_SIZE	An internal structure has an incorrect size.
159	IS_INVALID_BUFFER_SIZE	The image memory has an inappropriate size to store the image in the desired format.
160	IS_INVALID_PIXEL_CLOCK	This setting is not available for the currently set pixel clock frequency.
161	IS_INVALID_EXPOSURE_TIME	This setting is not available for the currently set exposure time.
162	IS_AUTO_EXPOSURE_RUNNING	This setting cannot be changed while automatic exposure time control is enabled.
163	IS_CANNOT_CREATE_BB_SURF	The BackBuffer surface cannot be created.
164	IS_CANNOT_CREATE_BB_MIX	The BackBuffer mix surface cannot be created.
165	IS_BB_OVLMEM_NULL	The BackBuffer overlay memory cannot be locked.
166	IS_CANNOT_CREATE_BB_OVL	The BackBuffer overlay memory cannot be created.
167	IS_NOT_SUPP_IN_OVL_SURF_MODE	Not supported in BackBuffer Overlay mode.
168	IS_INVALID_SURFACE	Back buffer surface invalid.
169	IS_SURFACE_LOST	Back buffer surface not found.
170	IS_RELEASE_BB_OVL_DC	Error releasing the overlay device context.
171	IS_BB_TIMER_NOT_CREATED	The back buffer timer could not be created.
172	IS_BB_OVL_NOT_EN	The back buffer overlay was not enabled.
173	IS_ONLY_IN_BB_MODE	Only possible in BackBuffer mode.
174	IS_INVALID_COLOR_FORMAT	Invalid color format
175	IS_INVALID_WB_BINNING_MODE	Mono binning/mono sub-sampling do not support automatic white balance.
176	IS_INVALID_I2C_DEVICE_ADDRESS	Invalid I ² C device address
177	IS_COULD_NOT_CONVERT	The current image could not be processed.
178	IS_TRANSFER_ERROR	Transfer error. Frequent transfer errors can mostly be avoided by reducing the pixel rate.
179	IS_PARAMETER_SET_NOT_PRESENT	Parameter set is not present.
180	IS_INVALID_CAMERA_TYPE	The camera type defined in the .ini file does not match the current camera model.
181	IS_INVALID_HOST_IP_HIBYTE	Invalid HIBYTE of host address
182	IS_CM_NOT_SUPP_IN_CURR_DISPLAYMODE	The color mode is not supported in the current display mode.
183	IS_NO_IR_FILTER	No IR filter available
184	IS_STARTER_FW_UPLOAD_NEEDED	The camera's starter firmware is not compatible with the driver and needs to be updated.
185	IS_DR_LIBRARY_NOT_FOUND	The DirectRenderer library could not be found.
186	IS_DR_DEVICE_OUT_OF_MEMORY	Not enough graphics memory available.

No	Error	Description
187	IS_DR_CANNOT_CREATE_SURFACE	The image surface or overlay surface could not be created.
188	IS_DR_CANNOT_CREATE_VERTEX_BUFFER	The vertex buffer could not be created.
189	IS_DR_CANNOT_CREATE_TEXTURE	The texture could not be created.
190	IS_DR_CANNOT_LOCK_OVERLAY_SURFACE	The overlay surface could not be locked.
191	IS_DR_CANNOT_UNLOCK_OVERLAY_SURFACE	The overlay surface could not be unlocked.
192	IS_DR_CANNOT_GET_OVERLAY_DC	Could not get the device context handle for the overlay.
193	IS_DR_CANNOT_RELEASE_OVERLAY_DC	Could not release the device context handle for the overlay.
194	IS_DR_DEVICE_CAPS_INSUFFICIENT	Function is not supported by the graphics hardware.
195	IS_INCOMPATIBLE_SETTING	Because of other incompatible settings the function is not possible.
196	IS_DR_NOT_ALLOWED_WHILE_DC_IS_ACTIVE	A device context handle is still open in the application.
197	IS_DEVICE_ALREADY_PAIED	The device is already paired.
198	IS_SUBNETMASK_MISMATCH	The subnet mask of the camera and PC network card are different.
199	IS_SUBNET_MISMATCH	The subnet of the camera and PC network card are different.
200	IS_INVALID_IP_CONFIGURATION	The configuration of the IP address is invalid.
201	IS_DEVICE_NOT_COMPATIBLE	The device is not compatible to the drivers.
202	IS_NETWORK_FRAME_SIZE_INCOMPATIBLE	The settings for the image size of the camera are not compatible to the PC network card.
203	IS_NETWORK_CONFIGURATION_INVALID	The configuration of the network card is invalid.
204	IS_ERROR_CPU_IDLE_STATES_CONFIGURATION	The configuration of the CPU idle has failed.
205	IS_DEVICE_BUSY	The camera is busy and cannot transfer the requested image.

4.8.2 Error Codes of AVI Functions

No.	Error	Description
300	IS_AVI_NO_ERR	Function executed successfully.
301	IS_AVI_ERR_INVALID_FILE	The file has no valid AVI format.
302	IS_AVI_ERR_NEW_FAILED	No memory could be allocated for the AVI file.
303	IS_AVI_ERR_CREATESTREAM	No AVI stream could be created.
304	IS_AVI_ERR_PARAMETER	One of the submitted parameters is outside the valid range.
305	IS_AVI_ERR_NO_CODEC_AVAIL	The maximum number of instances allowed in this system has been reached. It is not possible to create another instance.
306	IS_AVI_ERR_INVALID_ID	The indicated AVI instance could not be found. Either the AVI ID is invalid or the instance has already been deleted using isavi_ExitAVI() .
307	IS_AVI_ERR_COMPRESS	The last image compression failed.
309	IS_AVI_ERR_CAPTURE_RUNNING	Another capturing operation is in progress or an AVI file is still open.
310	IS_AVI_ERR_CAPTURE_NOT_RUNNING	No capturing operation is running or no AVI file is opened.
311	IS_AVI_ERR_PLAY_RUNNING	A playback is already running.
312	IS_AVI_ERR_PLAY_NOT_RUNNING	No playback is running.
313	IS_AVI_ERR_WRITE_INFO	The AVI file could not be modified.
314	IS_AVI_ERR_INVALID_VALUE	The q parameter is outside the range of 1...100.
315	IS_AVI_ERR_ALLOC_MEMORY	No memory could be allocated.
316	IS_AVI_ERR_INVALID_CM	The submitted color mode is not supported for AVI capturing.
317	IS_AVI_ERR_COMPRESSION_RUN	The current image could not be processed since compression is still in progress.
318	IS_AVI_ERR_INVALID_SIZE	The submitted size is invalid.
319	IS_AVI_ERR_INVALID_POSITION	The submitted position is invalid.
320	IS_AVI_ERR_INVALID_UEYE	No DCx camera was found.
321	IS_AVI_ERR_EVENT_FAILED	The event could not be initialized. The Windows <code>SetEvent</code> function failed.

4.8.3 Linux: Not Supported Functions

The following uc480 API functions are not supported by the Linux driver version 4.82 and older:

is_DirectRenderer()
is_GetColorDepth()
<code>is_GetDC()</code>
<code>is_ReleaseDC()</code>

[is_RenderBitmap\(\)](#)

[is_SetDisplayMode\(\)](#)

[is_SetDisplayPos\(\)](#)

[is_SetHwnd\(\)](#)

[is_SetOptimalCameraTiming\(\)](#)

[is_UpdateDisplay\(\)](#)

5 Specifications

This chapter lists the specifications of the available DCx camera models.

- [Model comparison](#)
- [Camera and sensor data](#)
- [Mechanical specifications](#)
- [Electrical specifications](#)
- [Accessories](#)

5.1 Model Comparison

The following table outlines the key features of DCx camera series for direct comparison (see also the [DCx Camera Family](#) chapter).

	DCC1545M DCC1645C	DCU223x DCU224x	DC1240x	DC3240x	DC3260x
Sensor Type	CMOS	CCD	CMOS	CMOS	CMOS
Objective Mount	CS (C*)	C	C	C	C
EMC compliance	CE B FCC B	CE A FCC A	CE B FCC A	CE A FCC A	CE A FCC A
WxHxD [mm] (CCD size)	44 x 44 x 25.6	32 x 34 x 30.2 (37.2)	32 x 34 x 30.2	29 x 29 x 29	29 x 29 x 29
Mounting holes					
bottom	1	2	2	2	2
top	0	3	3	3	3
side	0	3	3	3	3
Thread diameter	1/4" (#8-32, M4)*)	#8-32, M4	#8-32, M4	1/4" (#8-32, M4)*)	1/4" (#8-32, M4)*)
Adjustable flange back distance	+	-	-	-	-
IP protection class	30	30	30	30	30
Interface	USB 2.0	USB 2.0	USB 2.0	USB 3.0 USB 2.0	USB 3.0 USB 2.0
Power supply	USB	USB	USB	USB	USB
Lockable connector	-	+	+	+	+
I/O connector	10-pin connector	9-pin micro D- SUB	9-pin micro D- SUB	8-pin HR25	8-pin HR25
Optocoupler for I/O	-/-	1/1	1/1	1/1	1/1
Optocoupler speed	-	100 µs	100 µs	High (trigger)	High (trigger)
Max. cable length (m)	5	5	5	5	5
Dig. I/O/GPIO	0/0/0	1/1/0	1/1/0	1/1/2	1/1/2
PWM	-	-	-	+	+
RS-232	-	-	-	+ (GPIO)	+ (GPIO)
PC bus	-	-	-	-	-
Hot pixel correction	Software	Software	Software	Hardware	Hardware
Color calculation	Software	Software	Software	Hardware	Hardware
Hardware: Dig. Subsampling/ Binning	-	-	-	+	+
Bit depth: Internal/ transferred	8/8	8/8	8/8	16/16	12/12

	DCC1545M DCC1645C	DCU223x DCU224x	DC1240x	DC3240x	DC3260x
LUT: Internal/ transferred	-	-	-	12/12	12/12
Max. pixel clock (MHz) at full resolution	43	43	43	480	118
Max. data rate (MBytes/s)	38	38	38	300	300

*) with supplied adapter

5.2 Model Naming Conventions

In the table below, you will find an overview of all the models in the DCx camera series.

		DC	X	##	#	X
Short for "Digital Camera"						
C	CMOS					
U	CCD					
Sensor and shutter type						
12	CMOS Global Shutter					
15/16	CMOS Rolling Shutter					
22	CCD Progressive Scan					
32	CMOS Global Shutter, USB3					
Resolution						
3	XGA					
4, 40, 45	SXGA					
6						
Color / Mono						
M	Monochrome sensor					
C	Color sensor					
N	Near Infrared sensor					

5.3 Camera and Sensor Data

In this chapter the technical properties of the sensors are listed. You can look up e.g. binning factors. Also you will find the parameters for the different interfaces.

Cameras with CMOS sensors

- [DCC1240x / DCC3240x](#)
- [DCC1545M](#)
- [DCC1645C](#)

Cameras with CCD sensors

- [DCU223x](#)
- [DCU224x](#)

Note:

The diagrams shown in the sensor specifications section indicate the **relative** sensitivities of the DCx Cameras in the spectral range. Therefore, the characteristic curves cannot be compared to each other.

5.3.1 DCC3260x

Sensor specification	
Sensor type	CMOS
Shutter system	Global shutter
Characteristic	Linear
Readout mode	Progressive scan
Resolution class	2 MP
Resolution	2.35 MPix
Resolution (h x v)	1936 x 1216 pixels
Aspect ratio	16:10
ADC (analog-to-digital converter)	12 bits
Optical sensor class	1/1,2 inch
Optical area	11.340 x 7.130 mm
Optical sensor diagonal	13.395 mm (1/1.19 inch)
Pixel size	5.86 μm , square
Micro lens shift	0°
Manufacturer	Sony
Sensor name, monochrome	IMX249LLJ-C
Sensor name, color	IMX249LQJ-C
Special features	<ul style="list-style-type: none"> • Supports overlap trigger mode for high frame rates in trigger mode • Using the extended pixel clock
Gain	
Monochrome model (master gain)	24x
Color model (master gain/RGB)	24x/4x
Analog gain boost	-

Camera timing		USB 3 uc480 CP Rev. 2
Pixel clock range ^{*1}	MHz	30 – 118
Frame rate freerun mode ^{*2}	fps	41
Frame rate trigger mode (continuous) ^{*2}	fps	41
Frame rate trigger mode (maximum) ^{*2}	fps	41
Exposure time (min ^{*2} , max ^{*3})	ms	0.033 – 30 000
Maximum long exposure	ms	30 000
AOI		
Mode		Horizontal + vertical ^{*4}
AOI image width / step width	pixels	96 – 1936 / 8
AOI image height / step width	pixels	2 – 1216 / 2
AOI position grid / horizontal, vertical	pixels	8 / 2
Binning		
Mode		-

Subsampling		
Mode		Horizontal + vertical
Method		Monochrome/color: automatic
Factor		2x, 3x, 4x, 5x, 6x, 8x, 16x
Hardware trigger		
Mode		Asynchronous
Trigger delay with rising edge	μs	3 ±0.25* ⁶
Trigger delay with falling edge	μs	21 ±0.25* ⁶
Adjustable trigger delay	μs	15 – 4 000 000
Power consumption*⁵		USB 3 uc480 CP Rev. 2
	W	1.4 – 2.4

** Not yet defined.

*1 The maximum possible pixel clock frequency depends on the PC hardware used.

*2 Requires maximum pixel clock frequency.

*3 Requires minimum pixel clock frequency.

*4 Use of this function increases the frame rate.

*5 The power consumption depends on the sensor model and the pixel clock setting.

*6 Not yet confirmed



Please see also the DCC3260 application notes chapter.

Relative sensor sensitivity




Fig. 1: Sensor sensitivity DCC3260/ (monochrome)




Fig. 2: Sensor sensitivity DCC3260/ (color)

The cut-off wavelength of the IR cut filter is at 650 nm (if the IR cut filter is used in the camera model).

5.3.2 DCC1240x / DCC3240x

Sensor specification	
Sensor type	CMOS
Shutter system	Electronic global and rolling shutter
Characteristic	Linear
Readout mode	Progressive scan
Resolution class	SXGA
Resolution	1280 x 1024 pixels (1.3 Megapixel)
Aspect ratio	5:4
Bit depth	10 bits *8
Optical sensor class	1/1.8 inch
Exact sensitive area	6.784 mm x 5.427 mm
Exact optical sensor diagonal	8.69 mm (1/1.84 inch)
Pixel size	5.30 μ m, square
Micro lens shift	12°
Sensor name, monochrome	e2v EV76C560ABT
Sensor name, color	e2v EV76C560ACT
Sensor name, NIR *7	e2v EV76C661ABT
Special features	<ul style="list-style-type: none"> • Automatic hotpixel correction in the sensor, see is_HotPixel() • Multi AOI with 2 or 4 AOI, see Camera basics: AOI • Sequence AOI • Sensor internal image scaler, downscaling by factor 1...4, see is_SetSensorScaler() and uc480 Viewer: Size)

		<ul style="list-style-type: none"> Allows to switch between global and rolling shutter readout, see is_DeviceFeature() 	
Gain			
Monochrome model (master gain)		4x	
Color model (master / RGB)		4x / 3.96x	
Gain boost		2x	
Camera timing		DCC1240x	DCC3240x
Pixel clock range (allowed/recommended)	MHz	7 to 35 / 35 ^{*1}	5 to 85 / 85 ^{*1}
Max. pixel clock with subsampling/binning	MHz	85 ^{*1}	85 ^{*1}
Frame rate (freerun mode)	fps	25.8 ^{*2}	60.0 ^{*2}
Frame rate (trigger mode, 1 ms exposure)	fps	24.7 ^{*2}	56.9 ^{*2}
Exposure time in freerun mode	ms	0.009 ^{*2} to 2000 ^{*3}	0.009 ^{*2} to 2000 ^{*3}
Exposure time in trigger mode	ms	0.009 ^{*2} to 2000 ^{*3}	0.009 ^{*2} to 2000 ^{*3}
AOI			
Mode		Horizontal + Vertical ^{*4}	
AOI image width, step width	Pixels	16 to 1280, 4	16 to 1280, 4
AOI image height, step width	Pixels	4 to 1024, 2	4 to 1024, 2
AOI position grid horizontal, vertical	Pixels	2, 2	2, 2
AOI frame rate, 640 x 480 pixels (VGA)	fps	52.0	123.0
AOI frame rate, 320 x 240 pixels (CIF)	fps	98.0	229.0
Binning			
Mode		Horizontal + Vertical ^{*4}	
Method		H + V combined, mono/color binning, H: additive. V: averaging	
Factors		2x	
Frame rate with 2x binning, 640 x 480 pixels (VGA)	fps	60.0	64.0
Subsampling			
Mode		Scaler	
Hardware trigger			
Mode		Asynchronous	Asynchronous
Trigger delay with rising edge	μs	20 ±0.25	3 ±0.25 ^{*6}
Trigger delay with falling edge	μs	33 ±0.25	21 ±0.25 ^{*6}
Additive trigger delay (optional)	μs	15 μs...4 s	15 μs...4 s
Power consumption ^{*5}			

	W	0.3 to 0.7	1.3 ^{*5}
--	---	------------	-------------------

** Not yet defined.

*1 The maximum possible pixel clock frequency depends on the PC hardware used.

*2 Requires maximum pixel clock frequency.

*3 Requires minimum pixel clock frequency.

*4 Use of this function increases the frame rate.

*5 The power consumption depends on the sensor model and the pixel clock setting.


*6 Not yet confirmed

*7 DCC3240N only


*8 Only for USB3.0 transmission (DCC3240x), with USB2.0 - bit depth is 8 bit.

Please see also the [DCC1240x / DCC3240x Application Notes](#) chapter.

Relative sensor sensitivity



Sensor sensitivity of the DCC1240M (monochrome) / DCC3240x (monochrome and NIR)



Sensor sensitivity of the DCC1240C / DCC3240C

Note

The colored part of above diagrams just indicates the IR wavelength range in order to tell it from the visible.

5.3.3 DCC1545M

Sensor specification		
Sensor type	CMOS	
Shutter system	Electronic rolling shutter	
Readout mode	Progressive scan	
Resolution class	SXGA	
Resolution	1280 x 1024 pixels (1.31 Megapixel)	
Aspect ratio	5:4	
Bit depth	10 bits	
Optical sensor class	1/2 inch	
Exact sensitive area	6.656 x 5.325 mm	
Exact optical sensor diagonal	8.52 mm (1/1.88 inch)	
Pixel size	5.20 μm , square	
Sensor name	Aptina MT9M001 (monochrome)	
Gain		
Monochrome model (master gain)	13x	
Analog gain boost	1.5x	
Camera timing		
Pixel clock range	MHz	5 to 43 ^{*1}
Max. pixel clock with subsampling/binning	MHz	48 ^{*1}
Frame rate (freerun mode)	fps	25.0 ^{*2}
Frame rate (trigger mode, 1 ms exposure)	fps	25.0 ^{*2}
Exposure time in freerun mode	ms	0.037 ^{*2} to 983 ^{*2}
Exposure time in trigger mode	ms	0.037 ^{*2} to 983 ^{*2}
AOI		
Mode		Horizontal ^{*4} + Vertical ^{*4}
AOI image width, step width	Pixels	32 to 1280, 4
AOI image height, step width	Pixels	4 to 1024, 2
AOI position grid horizontal, vertical	Pixels	4, 2
AOI frame rate, 640 x 480 pixels (VGA)	fps	84
Binning		
Mode		none
Subsampling		
Mode		Horizontal ^{*4} + Vertical ^{*4}
Method		H + V: Color subsampling
Factor		2x, 4x, 8x
Frame rate w/ 2x subsampling,	fps	94

640 x 480 pixels		
Frame rate w/ 4x subsampling, 320 x 240 pixels	fps	258
Hardware trigger		
Mode		Asynchronous
Trigger delay with rising edge	μs	22.0 \pm 0.25
Trigger delay with falling edge	μs	40.3 \pm 0.25
Additive trigger delay (optional)	μs	15 μs ...4 s
Power consumption ^{*5}		
	W	0.5 to 1.0

*1 The maximum possible pixel clock frequency depends on the PC hardware used.

*2 Requires maximum pixel clock frequency.


*3 Requires minimum pixel clock frequency.

*4 Use of this function increases the frame rate.

*5 The power consumption depends on the sensor model and the pixel clock setting.

Please see also the [DCC1545M Application Notes](#) chapter.

Relative sensor sensitivity



Sensor sensitivity of the DCC1545M (monochrome)

Note

The colored part of above diagram just indicates the IR wavelength range in order to tell it from the visible.

5.3.4 DCC1645C

Sensor specification		
Sensor type	CMOS	
Shutter system	Electronic rolling shutter	
Characteristic	Linear	
Readout mode	Progressive scan	
Resolution class	SXGA	
Resolution	1280 x 1024 pixels (1.31 Megapixel)	
Aspect ratio	5:4	
Bit depth	10 bits	
Optical sensor class	1/3 inch	
Exact sensitive area	4.608 x 3.686 mm	
Exact optical sensor diagonal	5.90 mm (1/2.71 inch)	
Pixel size	3.60 μm , square	
Micro lens shift	25°	
Sensor name	Aptina MT9M131 (color)	
Gain		
Color model (master/RGB)	4.27x / 3.1x	
Analog gain boost	2.0x	
Camera timing		
Pixel clock range	MHz	5 to 40 ^{*1}
Max. pixel clock with subsampling/ binning	MHz	40 ^{*1}
Frame rate (freerun mode)	fps	25.0 ^{*2}
Frame rate (trigger mode, 1 ms exposure)	fps	24.9 ^{*2}
Exposure time in freerun mode	ms	0.037 ^{*2} to 10122 ^{*3}
Exposure time in trigger mode	ms	0.037 ^{*2} to 10122 ^{*3}
AOI		
Mode		Horizontal ^{*4} + Vertical ^{*4}
AOI image width, step width	Pixels	32 to 1280, 4
AOI image height, step width	Pixels	4 to 1024, 2
AOI position grid horizontal, vertical	Pixels	4, 2
AOI frame rate, 1280 x 720 pixels (HD 720)	fps	34
AOI frame rate, 800 x 600 pixels (SVGA)	fps	61
Binning		
Mode		none
Subsampling		
Mode		Horizontal ^{*4} + Vertical ^{*4}

Method		H + V: Color subsampling
Factor		2x, 4x
Frame rate w/ 2x subsampling, 640 x 480 pixels	fps	89
Frame rate w/ 4x subsampling, 320 x 240 pixels	fps	263
Hardware trigger		
Mode		Asynchronous
Trigger delay with rising edge	μs	180.9 ±0.25
Trigger delay with falling edge	μs	199.3 ±0.25
Additive trigger delay (optional)	μs	15 μs...4 s
Power consumption^{*5}		
	W	0.3 to 0.8

*1 The maximum possible pixel clock frequency depends on the PC hardware used.

*2 Requires maximum pixel clock frequency.


*3 Requires minimum pixel clock frequency.

*4 Use of this function increases the frame rate.

*5 The power consumption depends on the sensor model and the pixel clock setting.

Please see also the [DCC1645C Application Notes](#) chapter.

Relative sensor sensitivity



Sensor sensitivity of the DCC1645C

Note

The colored part of above diagram just indicates the IR wavelength range in order to tell it from the visible.

5.3.5 DCU223x

Sensor specification		
Sensor type	CCD	
Shutter system	Electronic global shutter	
Characteristic	Linear	
Readout mode	Progressive scan	
Resolution class	XGA	
Resolution	1024 x 768 pixels (0.79 Megapixel)	
Aspect ratio	4:3	
Bit depth	12 bits	
Optical sensor class	1/3 inch	
Exact sensitive area	4.762 x 3.571 mm	
Exact optical sensor diagonal	5.95 mm (1/2.69 inch)	
Pixel size	4.65 μm , square	
Sensor name, monochrome	Sony ICX204AL	
Sensor name, color	Sony ICX204AK	
Gain		
Monochrome model (master gain)	10.47x	
Color model (master/RGB)	7.59x/4.0x	
Analog gain boost	2.0x (monochrome model only)	
Camera timing		
Pixel clock range (allowed/recommended)	MHz	5 to 30 / 10 to 20 ^{*1}
Pixel clock for optimal image quality	MHz	15 ^{*1}
Max. pixel clock with subsampling/binning	MHz	30 ^{*1}
Frame rate (freerun mode)	fps	30.0 ^{*2}
Frame rate (trigger mode, 1 ms exposure)	fps	28.7 ^{*2}
Exposure time in freerun mode	ms	0.030 ^{*2} to 773 ^{*3}
Exposure time in trigger mode	ms	0.030 ^{*2} to 10 min ^{*3}
AOI		
Mode		Horizontal + Vertical ^{*4}
AOI image width, step width	Pixels	16 to 1024, 4
Mono: AOI image height, step width	Pixels	120 to 768, 1
Color: AOI image height, step width	Pixels	120 to 768, 2
Mono: AOI position grid horizontal, vertical	Pixels	1, 1

Color: AOI position grid horizontal, vertical	Pixels	2, 2
AOI frame rate, 800 x 600 pixels (SVGA)	fps	37
AOI frame rate, 640 x 480 pixels (VGA)	fps	45
AOI frame rate, 320 x 240 pixels (CIF)	fps	78
Binning		
Mode		Vertical ^{*4}
Method		V: Monochrome binning, additive
Factor		2x, 3x, 4x
Frame rate with 2x binning, 1024 x 384 pixels	fps	53
Frame rate with 3x binning, 1024 x 256 pixels	fps	71
Frame rate with 4x binning, 1024 x 192 pixels	fps	85
Subsampling		
Mode		-
Hardware trigger		
Mode		Asynchronous
Trigger delay with rising edge	μs	39.5 ±2.6
Trigger delay with falling edge	μs	57.9 ±2.6
Additive trigger delay (optional)	μs	15 μs...4 s
Power consumption^{*5}		
	W	1.0 to 1.7

*1 The maximum possible pixel clock frequency depends on the PC hardware used.

*2 Requires maximum pixel clock frequency.


*3 Requires minimum pixel clock frequency.

*4 Use of this function increases the frame rate.

*5 The power consumption depends on the sensor model and the pixel clock setting.

Please see also the [DCU223x Application Notes](#) chapter.

Relative sensor sensitivity

**Note**

The colored part of above diagrams just indicates the IR wavelength range in order to tell it from the visible.

5.3.6 DCU224x

Sensor specification		
Sensor type	CCD	
Shutter system	Electronic global shutter	
Characteristic	Linear	
Readout mode	Progressive scan	
Resolution class	SXGA	
Resolution	1280 x 1024 pixels (1.31 Megapixel)	
Aspect ratio	5:4	
Bit depth	12 bits	
Optical sensor class	1/2 inch	
Exact sensitive area	5.952 x 4.762 mm	
Exact optical sensor diagonal	7.62 mm (1/2.1 inch)	
Pixel size	4.65 μm , square	
Sensor name, monochrome	Sony ICX205AL	
Sensor name, color	Sony ICX205AK	
Gain		
Monochrome model (master gain)	13.66x	
Color model (master/RGB)	8.9/4.0x	
Analog gain boost	2.0x (monochrome model only)	
Camera timing		
Pixel clock range (allowed/recommended)	MHz	5 to 30 / 10 to 20 ^{*1}
Pixel clock for optimal image quality	MHz	14 ^{*1}
Max. pixel clock with subsampling/binning	MHz	30 ^{*1}
Frame rate (freerun mode)	fps	15.0 ^{*2}
Frame rate (trigger mode, 1 ms exposure)	fps	17.0 ^{*2}
Exposure time in freerun mode	ms	0.066 ^{*2} to 1360 ^{*3}
Exposure time in trigger mode	ms	0.066 ^{*2} to 10 min ^{*3}
AOI		
Mode		Horizontal + Vertical ^{*4}
AOI image width, step width	Pixels	16 to 1280, 4
Mono: AOI image height, step width	Pixels	120 to 1024, 1
Color: AOI image height, step width	Pixels	120 to 1024, 2
Mono: AOI position grid horizontal, vertical	Pixels	1, 1
Color: AOI position grid horizontal,	Pixels	2, 2

vertical		
AOI frame rate, 1024 x 768 pixels (XGA)	fps	18
AOI frame rate, 640 x 480 pixels (VGA)	fps	28
AOI frame rate, 320 x 240 pixels (CIF)	fps	38
Binning		
Mode		Vertical ^{*4}
Method		V: Monochrome binning, additive
Factor		2x, 3x, 4x
Frame rate with 2x binning, 1280 x 512 pixels	fps	23
Frame rate with 3x binning, 1280 x 340 pixels	fps	28
Frame rate with 4x binning, 1280 x 256 pixels	fps	31
Subsampling		
Mode		Vertical ^{*4}
Method		V: Color subsampling
Factor		4x
Frame rate w/ 4x subsampling, 1280 x 256 pixels	fps	31
Hardware trigger		
Mode		Asynchronous
Trigger delay with rising edge	µs	39.9 ±2.5
Trigger delay with falling edge	µs	57.7 ±2.5
Additive trigger delay (optional)	µs	15 µs...4 s
Power consumption^{*5}		
	W	1.1 to 2.1

*1 The maximum possible pixel clock frequency depends on the PC hardware used.

*2 Requires maximum pixel clock frequency.



*3 Requires minimum pixel clock frequency.

*4 Use of this function increases the frame rate.

*5 The power consumption depends on the sensor model and the pixel clock setting.

Please see also the [DCU224x Application Notes](#) chapter.

Relative sensor sensitivity



Note

The colored part of above diagrams just indicates the IR wavelength range in order to tell it from the visible.








5.4 Mechanical Specifications

- [DCU223x, DCU224x](#)
- [DCC1545M, DCC1645C](#)
- [DCC1240x, DCC3240x](#)

General



- [Flange back distance](#)
- [Position accuracy of the sensor](#)
- [Filter glasses](#)
- [Ambient conditions](#)

5.4.1 DCU223x, DCU224x





all dimensions are mm
 all threads M3, unless otherwise noted

5.4.2 DCC1240x







DETAIL A
SCALE 2:1


- MICRO D-SUB
- 1: DIGITAL OUT (-)
 - 2: DIGITAL IN (+)
 - 3: SHIELD
 - 4: USB VCC (5V)
 - 5: USB GROUND
 - 6: DIGITAL OUT (+)
 - 7: DIGITAL IN (-)
 - 8: USB DATA (+)
 - 9: USB DATA (-)



5.4.3 DCC3240x




PIN	DESCRIPTION
1	GROUND
2	FLASH OUTPUT, OPTO-DECOUPLED (-)
3	GPIO 1, 3.3V LVCMOS
4	TRIGGER INPUT, OPTO-DECOUPLED (-)
5	FLASH OUTPUT, OPTO-DECOUPLED (+)
6	GPIO 2, 3.3V LVCMOS
7	TRIGGER INPUT, OPTO-DECOUPLED (+)
8	OUTPUT SUPPLY VOLTAGE, 5V (100mA)



DETAIL A
SCALE 2:1
REAR VIEW

5.4.4 DCC1545M, DCC1645C



5.4.5 Flange Back Distance


- [Calculating the flange back distance](#)
- [Maximum immersion depth for lenses](#)

5.4.5.1 Calculating the Flange Back Distance

To correctly determine the flange back distance of a DCx camera, you need to consider the distance between the lens flange and the active area of the sensor and, additionally, the type and thickness of any materials inserted into the optical path.

The "distance in air" between the threaded flange and the active area is 17.526 mm with C-mount lenses and 12.526 mm with CS-mount lenses.

This "mechanical distance" can change due to the material-specific refractive index of the inserted materials. The glass cover of the sensor and all filters inserted into the optical path (see [Filter types](#) table) must be taken into account in the calculation.



Calculating the flange back distance (schematic illustration)

Designation	Description
a	Distance from threaded flange to active sensor area (flange back distance) 17.526 mm*1 for C-mount 12.526 mm*1 for CS-mount
x	Distance from threaded flange to PCB
e	Distance from active sensor area to PCB
d	Thickness of the glass cover of the sensor
f	Filter thickness (optional)
n	Refractive index
h	Maximum sensor height above the PCB

*1 This distance describes the equivalent in air (see introduction above)

Hint

You can use the following formula to calculate the mechanical flange back distance:

$$x = a + \frac{d \times (n_{\text{glass}} - 1)}{n_{\text{glass}}} + \frac{f \times (n_{\text{filter}} - 1)}{n_{\text{filter}}} + e$$

The tolerances for the position accuracy of DCx camera sensors are given in the [Position accuracy](#) chapter.

Calculating the flange back distance for DCx Cameras with C-mount

Camera model	Thickness sensor glass [mm]	Active sensor area to PCB [mm]	Flange to active sensor area without filter glass[mm]	Flange to active sensor area with filter glass [mm]	Sensor height above the PCB [mm]
Designation	d	e	x'	x	h
CMOS					
DCC1240x/ DCC3240x	0.550	1.130	18.840	19.190	2.530
CCD					
DCU223x	0.750	2.810	20.590	20.930	4.780
DCU224x	0.750	2.810	20.590	20.930	4.830

Calculation example: UI-154x-xx with IR-cut filter

(a = 17.526 mm, d = 0.525 mm, nGlass = 1.50, f = 1 mm, nFilter = 1.53; see [Filter types](#) table)

$$x = 17.526 \text{ mm} + \frac{0.525 \text{ mm} \times (1.50 - 1)}{1.50} + \frac{1.00 \text{ mm} \times (1.53 - 1)}{1.53} + 1.27 \text{ mm} = 19.32 \text{ mm}$$

Calculating the flange back distance for DCC1545M and DCC1645C cameras with CS-mount**Note**

For these cameras with CS-mount, the flange back distance is only 12.526 mm.

Camera model	Thickness sensor glass [mm]	Active sensor area to PCB [mm]	Threaded flange to active sensor area	Flange to active sensor area without filter glass[mm]	Flange to active sensor area with filter glass [mm]	Sensor height above the PCB [mm]
Designation	d	e	a	x	x'	h
DCC1545M	0.525	1.270	13.23	15.95	17.50	2.480
DCC1645C	0.550	1.400	13.10	16.08	17.63	2.500

5.4.5.2 Maximum Immersion Depth for Lenses


Some C-mount lenses reach deep into the camera flange. This may cause the lens to push against the back of the filter glass inside the camera or even make it impossible to screw in the lens.

The table below indicates the maximum possible immersion depth for each DCx camera model. The actual immersion depth of a lens is given in the relevant data sheet. As lens parts with a small diameter are allowed to reach deeper into the camera flange, the immersion depths are specified based on the diameter.

Beside the immersion depth also the back focal length has to be considered, that means the distance between the last lens and the sensor (named "bfl" in the image below). The back focal length can be calculated for C-mount with the following formula:

$$bfl = 17,526 - x$$

x stands for the maximum immersion depth (see table below).



Note

Front Panel of the DCU223x / DCU224x housing

Models introduced before 2008 have a different front panel. On these models, the filter glass is held in the C-mount lens connector by two screws. This front panel version is indicated by (V2) in the table below.

Camera	Type	Thread depth	for diameter at lens end	resulting max. immersion depth (mm)		min. required back focal length (mm)	
			(Ø A, B, C in mm)	CMOS	CCD	CMOS	CCD
DCU223X DCU224X DCC1240 X	C-mount	5 mm	24.0	6.9	6.4	6.4	6.9
			20.0	9.4	8.9		
			14.1	11.1	10.6		
DCU223X (V2)	C-mount	5 mm	24.0	-	8.4	-	7.8

DCU224X (V2)			17.1	-	9.2		
			14.1	-	9.7		
DCC1545 M DCC1645 C	CS-mount	4 mm	24.0	6.1	-	4.9	-
			14.6	7.6	-		
	C-mount with extension ring	4 mm* ¹	22.0* ¹	11.1	-		
			14.6* ¹	12.6	-		
DCC3240 x	C-mount	5 mm	24.0	7.2	-	6.8	-
			20.4	9.7	-		
			14.6	10.7	-		

*¹ May vary depending on the inside diameter of the extension ring used

*² Without IR cut filter


Note

The data given in the table include the following tolerances as a safety clearance:


- Immersion depth: 0.2 mm
- Diameter: 0.2 mm

5.4.6 Position Accuracy of the Sensor

The following illustration shows the tolerance margins of the sensor position relative to the camera housing. A maximum error in all directions (rotation, translation) cannot occur simultaneously.



Position accuracy of the sensor (1)



Position accuracy of the sensor (2)

Position accuracy inside the camera housing, in each direction	±0.3	mm
Horizontal/vertical rotation (°)	±1.0	
Translational rotation (°)	±1.0	
Flange back distance	±0.05	mm

Note

C-mount lenses can also be subject to inaccuracies of the flange back distance. The tolerance usually is ±0.05 mm. In some cases, however, the inaccuracies of camera and lens might add up,

resulting in a total error > 0.05 mm.

5.4.7 Filter Glasses

- [Filter types](#)
- [Mounting the filter](#)
- [Cleaning the filter glasses](#)

5.4.7.1 Filter Types


Each DCx camera has a filter glass in the front flange to prevent the entry of dust particles. Color cameras by default use an IR cut filter, which is required to ensure correct color rendering. The default filter glass in monochrome cameras has no filter effect. Every camera model is available with different filter variants such as daylight cut filters (type DL).

The following table shows an overview of the different optical filters used in the DCx Cameras:

Filter type	Name	Refractive index (n_{Filter})	Glass type	Thickness (f)	Cut-off wavelength	Non- reflective
IR cut filter (old)	BG	1.53	BG40	1 mm	650 nm	-
IRcut filter (new)	HQ	1.53	D263	1 mm	650 nm	On one side
Daylight cut filter	DL	1.53	RG665	1 mm	695 nm	-
Glass	GL	1.53	D263	1 mm	380 nm	On both sides


Notes

- All sensors have a D263 type cover glass. This glass is opaque to wavelengths below 330 nm.
- You can tell the filter type from visually by its coloration:
 - Reddish glass: HQ filter
 - Bluish glass: BG filter
 - Opaque glass: DL filter
 - Plain glass: GL filter
- New DCx color cameras use an IR cut filter of the type HQ by default. This filter offers an improved accuracy of the infrared content. HQ filters achieve a higher image brightness and better color rendering compared with the BG filter.
- uc480 drivers of version V3.24 and higher determine automatically which the IR filter is used in a camera. The corresponding color correction is selected automatically.


Infrared cut filter (type BG)**Note**

The colored part of above diagram just indicates the IR wavelength range in order to tell it from the visible.

Infrared cut filter (type HQ)




Daylight cut filter (type DL)



Note

The colored part of above diagrams just indicates the IR wavelength range in order to tell it from the visible.

Plain glass filter (type GL)**Notes**

- The colored part of above diagram just indicates the IR wavelength range in order to tell it from the visible.
- A different scale is used for the Y-axis of the glass filter curve, to show the curves between 400 nm and 1100 nm better.


5.4.7.2 Mounting the Filter

Attention

It is recommended to have the filter changed under clean room conditions, otherwise dust might enter the sensor area and become visible in images.

When completely removing the adjustment ring and the filter glass, the rubber gasket should also be removed. Due to expansion of the rubber gasket during installation, it's difficult to reinstall then the filter glass.

The threaded ring presses the filter glass on a rubber gasket. A properly mounted threaded ring will seal off the sensor. The threaded ring is screwed into the adjusting ring from the front with a torque of 0.2 Nm.



Notes

- A special tool is required for adjusting the threaded ring!
- Some DCx camera models have a different design with a separate filter glass that is secured by two screws.

5.4.7.3 Cleaning the Filter Glasses

When using the DCx camera with its lens removed, the filter glass may become soiled from outside. This might be visible in captured images. In such case, the filter glass needs cleaning.

Note

It is strongly recommended to return the cameras to manufacturer for professional cleaning. The manufacturer is not liable for any damage resulting from cleaning the filter glasses. This even applies if the following instructions have been observed.

Instructions for cleaning DCx camera filter glasses


- The filter glasses can be cleaned only from outside. If remove the filter glasses, the sensor might become soiled. Thorlabs is not liable for any damage to the sensor resulting from removal of the filter glasses.
- First, remove dirt particles on the glass using compressed air. Do not use compressed air from compressors or spray cans since it often contains oil droplets or droplets of other liquids. For best results, use purified nitrogen from nitrogen bottles.
- Only use lint-free wipes or cotton-free swabs for cleaning. Never touch the filter glasses with bare fingers - it's mostly difficult to remove fingerprints completely!
- We recommend to use pure alcohol for cleaning. 100% isopropyl alcohol evaporates without leaving any residues. Only add small quantities of alcohol to the wipe. Never pour alcohol directly onto the camera.

Attention

Never use cleaning agents containing acetone for cleaning the filter glasses! Acetone may damage the filter glass coating and may deteriorate the optical quality of the glasses.

Cameras with fixed filter glass

Use a wipe to remove dust particles in a single sweep over the edge of the filter glass (see figure below).



Cleaning fixed DCx filter glasses

Cameras with replaceable filter glass

Use a wipe to remove dust particles in a circular sweep (see figure below).



Cleaning interchangeable DCx filter glasses

5.4.8 Ambient Conditions

Attention

- Avoid high air humidity levels and rapid temperature changes when using DCx Cameras.
- Temperatures below +4 °C (39 °F) combined with excessive relative air humidity levels can cause icing.
- At ambient temperatures above 45 °C (113 °F), the image quality could be reduced due to increased thermal noise. It is recommended to mount the camera to a heat-dissipating unit when high ambient temperatures prevail.

Note

The temperature values given above refer to the ambient temperature. The internal camera temperature is usually higher than the ambient temperature and must not exceed than 70 °C (158 °F).

	Min.	Max.	
Ambient temperature	-5	50	°C
	23	122	°F
Storage temperature	-20	60	°C
	-4	140	°F
Relative humidity*1	20	80	%

*1 Non-condensing

Note

Non-condensing means that the relative air humidity must be below 100 %. Otherwise, moisture will form on the camera surface. If, for example, air has a relative humidity of 40 % at 35 °C (95 °F), the relative humidity will increase to over 100 % if the air cools down to 19.5 °C (67 °F); condensation begins to form.

Vibration and shock resistance

Vibration and shock resistance of the USB DCx Cameras were tested as specified in DIN EN 60068-2-6(1996-05), DIN EN 60068-2-27(1995-03) and DIN EN 60068-2-29(1995-03). The mechanical shock was at 80 g; the vibration testing was performed with sinusoidal vibration at a frequency between 30 Hz-500 Hz and an amplitude of 10 g.

5.5 Camera Interface

This section of the manual contains information on connecting the cameras and wiring IOs.

For information on a camera's **power consumption**, please refer to the [Camera and sensor data](#) section. This section contains information on all camera models sorted by sensor type.

- [DCU223x, DCU224x, DCC1240x](#)
- [DCC3240x](#)
- [EEPROM Specification](#)

5.5.1 DCU223x, DCU224x, DCC1240x

In this section the additional digital input / output of these cameras is described in detail:


- [I/O Connector - Pin Assignment](#)
- [Digital Input \(Trigger\) Circuit](#)
- [Digital Output \(Flash\) Circuit](#)

5.5.1.1 I/O Connector - Pin Assignment


Attention

USB cables with non-standard connectors must be connected to the camera first and then to the PC. Otherwise the camera might not be recognized correctly.

9-pin micro D-Sub socket

Pin	Description	
1	Digital output (-)	 <p>DCU22xX / DCC1240X Micro D-Sub socket male, camera rear view</p>
2	Digital input (+)	
3	Shielding	
4	USB power supply (VCC) 5 V	
5	USB ground (GND)	
6	Digital output (+)	
7	Digital input (-)	
8	USB data (+)	
9	USB data (-)	

Pin assignment of the CAB-DCU-Tx cable for USB 2.0, trigger and flash

Pin	Description	Cable color	
1	Digital output (-)	green	 <p>DCU22xX / DCC1240X</p>

2	Digital input (+)	white	<i>Micro D-Sub connector female, connecting side view</i>
6	Digital output (+)	yellow	
7	Digital input (-)	brown	

For a comprehensive list of all cables and connectors available for DCU22xX / DCC1240X cameras, please refer to the [DCU22xX / DCC1240X Accessories](#) section.

5.5.1.2 Digital Input (Trigger) Circuit

Digital input specifications


USB board revision *)	1.2		2.0 or higher		
	Min.	Max.	Min.	Max.	
Level low	0	2	0	2	V
Level high	9	24	5	24	V
Voltage range	0	30	0	30	V
Trigger pulse width (edge)	100	-	100	-	μ s
Trigger edge steepness	35		35		V/ms
Breakdown voltage		50		50	V
Input current	10	-	10	-	mA

Note

*) For information on how to determine the USB board revision, please refer to the [DCx Driver Compatibility](#) chapter.

For interpreting the trigger signal, either the negative or the falling edge can be used. The digital input is galvanically isolated using an opto coupler to protect the camera and the PC against surges. Only DC voltages may be applied to the digital input.

Digital input wiring



5.5.1.3 Digital Output (Flash) Circuit

Digital output specifications

USB board revision *)	1.2	2.0 or higher	
	Max.	Max.	
Output current (short-time)	50	500	mA
Output current (permanent)	15	150	mA
Output voltage	30	30	V
Breakdown voltage	50	50	V
Collector power dissipation	100	125	mW

NOTE


*) For information on how to determine the USB board revision, please refer to the [DCx Driver Compatibility](#) chapter.


The digital output is galvanically isolated using an opto coupler to protect the camera and the PC against surges. Only DC voltages may be applied to the digital input.

The output of the opto coupler can be used as an open collector or open emitter output. This means that the output signal can be connected to ground or to the supply voltage. The output signal is active if the collector-emitter switch is closed (software setting: Flash high active, see also the Camera Properties: Input/Output section).


Digital output wiring

The following figures show examples of how the digital output is wired.





Wiring of the digital output as an open collector output (rev. 2.0)



Wiring of the digital output as an open emitter output (rev. 2.0)

5.5.2 DCC3240x

- [I/O connector Pin Assignment](#)
- [GPIO Interface](#)
- [Digital Input \(Trigger\) Circuit](#)
- [Digital Output \(Flash\) Circuit](#)
- [RS-232 Serial Interface](#)

5.5.2.1 I/O Connector Pin Assignment

Attention


The General Purpose IO are not potential-free and have no protective circuit. Faulty wiring (overvoltage, undervoltage or inverting the wiring when used as serial interface) can result in a damage in the electronics.

During operation as serial interface only LVCMOS levels are allowed to the connector pins. To get a serial RS-232 compliant interface, an external level shifter (LVCMOS/RS-232) is required.

Applying RS-323 levels directly to the pins as well as mixing up the signals RxD and TxD can destroy the camera electronics!

Pin assignment of the 8-pin Hirose connector HR25 for trigger, flash and GPIO

Pin	Description	Cable color
1	Ground (GND)	gray
2	Flash output, opto-decoupled (-)	green
3	GPIO 1, 3.3 V LVCMOS	blue
4	Trigger input, opto-decoupled (-)	brown
5	Flash output, opto-decoupled (+)	yellow
6	GPIO 2, 3.3 V LVCMOS	red
7	Trigger input, opto-decoupled (+)	white
8	Output supply voltage, 5 V (100 mA)	pink



Hirose connector male, camera rear view

For a comprehensive list of all cables and connectors available for DCC3240X cameras, please refer to the [DCC3240X Accessories](#) section.

5.5.2.2 GPIO Interface

GPIO specifications

The two GPIOs (General Purpose I/O) can be used as inputs or outputs. This selection is made by software using the corresponding SDK API functions. Please observe the following criteria:

- Input: 3.3 V LVCMOS, max. input voltage 4.0 V
- Output: 3.3 V LVCMOS, max. 10 mA

Attention


The General Purpose IO are not potential-free and have no protective circuit. Faulty wiring (overvoltage, undervoltage or inverting the wiring when used as serial interface) can result in a damage in the electronics.

During operation as serial interface only LVCMOS levels are allowed to the connector pins. To get a serial RS-232 compliant interface, an external level shifter (LVCMOS/RS-232) is required.

Applying RS-232 levels directly to the pins as well as mixing up the signals RxD and TxD can destroy the camera electronics!


GPIO wiring as input

The following figures illustrate GPIO wiring examples.



	Min.	Max.	
Signal level "Low"	0	0.8	V
Signal level "High"	2.0	4.0	V

GPIO wiring as output



	Min.	Max.	
Signal level "Low"	0	0.8	V

Signal level "High"	2.4	3.3	V
Output current	0	8.0	mA


5.5.2.3 Digital Input (Trigger) Circuit

Digital input specifications

	Min.	Max.	
Level low	0	1	V
Level high	5	24	V
Voltage range	0	-	V
Trigger pulse width (edge)	10	-	μs
Trigger edge steepness	35	-	V/ms
Breakdown voltage	-	50	V
Input current	10	-	mA

For interpreting the trigger signal, either the negative or the falling edge can be used. The digital input is galvanically isolated using an optocoupler to protect the camera and the PC against surges. Only DC voltages may be applied to the digital input.

Digital input wiring



5.5.2.4 Digital Output (Flash) Circuit

Digital output specifications


	Max.	
Output current (short-time)	500	mA
Output current (permanent)	150	mA
Output voltage	30	V
Breakdown voltage	50	V
Collector power dissipation	125	mW

The digital output is galvanically isolated using an optocoupler to protect the camera and the PC against surges. Only DC voltages may be applied to the digital output.


The output of the optocoupler can be used as an open collector or open emitter output. This means that the output signal can be connected to ground or to the supply voltage. The output signal is active if the collector-emitter switch is closed (software setting: Flash high active, see also the Camera properties: Input/output section).

Digital output wiring

The following figures show examples of how the digital output is wired.



Wiring of the digital output as an open collector output



Wiring of the digital output as an open emitter output

5.5.2.5 RS-232 Serial Interface

Attention

The General Purpose IO are not potential-free and have no protective circuit. Faulty wiring (overvoltage, undervoltage or inverting the wiring when used as serial interface) can result in a damage in the electronics.

During operation as serial interface only LVCMOS levels are allowed to the connector pins. To get a serial RS-232 compliant interface, an external level shifter (LVCMOS/RS-232) is required.

Applying RS-232 levels directly to the pins as well as mixing up the signals RxD and TxD can destroy the camera electronics!

Serial interface specification


Minimum output voltage		Min.	Max.	
	Signal level "Low"	0	0.8	V
	Signal level "High"	2.0	4.0	V
Maximum input voltage		Min.	Max.	
	Signal level "Low"	0	0.8	V
	Signal level "High"	2.0	4.0	V
Supported Baud rates	9.600 19.200 38.400 57.600 115.200			baud
Transmission mode	Full duplex, 8N1			
Data bits	8			
Stop bits	1			
Parity	None			

Note


With the 8N1 mode, the maximum payload data rate achievable is 80% of the selected baud rate.

Serial interface wiring (UART)

The following figure shows the wiring of the serial interface with GPIO 1 as camera-side output (TxD) and GPIO 2 as camera-side input (RxD). The GPIO must be configured accordingly (see [is_IO\(\)](#)).



Serial interface wiring (RS-232)



5.5.3 Camera EEPROM Specification

DCx Cameras have an EEPROM memory where the camera manufacturer, type, and serial number are stored. A 64 byte memory space can be used freely by the user.

EEPROM Specifications	
Data retention	10 years
Read/write cycles	100,000
Size of user data space	64 bytes

5.6 Accessories for DCx cameras

Lenses

Thorlabs also supplies a wide variety of lenses and [objectives](#) from leading manufacturers. Please contact [Thorlabs](#) for technical support and a detailed quote tailored to your needs.


USB cables

All Thorlabs DCx cameras are shipped with a 1.5m USB2.0, A to Mini B, cable.


5.6.1 Accessories for DCU22xX / DCC1240X

For information on the pin assignment of the cables and connectors see chapter [Pin Assignment I/O Connector](#).

USB Cables with Cables for digital I/Os



CAB-DCU-T1 (angled Micro D SUB camera connector; trigger in / out)



CAB-DCU-T2 (straight Micro D SUB camera connector, trigger in only)


Type	Length	Cable Type	Function	Connector Camera Side	Connector PC Side
CAB-DCU-T1	3 m	USB cable, AWG 28, single shielded,	USB & Trigger in / out	Micro D-Sub for screw-mounting, <u>angled</u>	USB 2.0 Type A
CAB-DCU-T2	3 m	additional cable for <u>digital I/Os</u> , 4-wire, open wires	USB & Trigger in	Micro D-Sub for screw-mounting, <u>straight</u>	

5.6.2 Accessories for DCC1x45X

For USB cables and accessories see also [Accessories for all DCx cameras](#).

DCC1545M and DCC1645C CS mount cameras are delivered with both CS / C objective mount adapter and CS / SM1 1" optics adapters.


5.6.3 Accessories for DCC3240x / DCC3260x



Type	Length	Cable Type	Function	Connector Camera Side
CAB-DCU-T3	2 m	Shielded high-flexible control cable 8 x 0.1mm, Ø 4.9mm	GPIO digital in (trigger) digital out (flash)	Hirose 8 pin (HR25)

Pin assignment of the 8-pin Hirose connector HR25 for trigger, flash and GPIO

Pin	Description	Cable color
1	Ground (GND)	gray
2	Flash output, opto-decoupled (-)	green
3	GPIO 1, 3.3 V LVCMOS	blue
4	Trigger input, opto-decoupled (-)	brown
5	Flash output, opto-decoupled (+)	yellow
6	GPIO 2, 3.3 V LVCMOS	red
7	Trigger input, opto-decoupled (+)	white
8	Output supply voltage, 5 V (100 mA)	pink



*Hirose connector female,
cable tip view*

6 Appendix

- [Troubleshooting/FAQ](#)
- [Status LED on USB DCx Cameras](#)
- [Color and memory formats](#)
- [uc480 parameter file \(ini file\)](#)
- [History of API functions](#)
- [Certifications and Compliances](#)
- [Thorlabs 'End of Life' Policy \(WEEE\)](#)
- [Exclusion of Liability and Copyright](#)
- [Thorlabs Worldwide Contacts](#)

6.1 Troubleshooting/FAQ

Installation and connection

- ☐ Installation of the Thorlabs camera software fails.

You need administrator privileges to install the software. Operating the cameras, however, does not require administrator privileges.

- ☐ The camera is connected to the PC, but cannot be opened in Thorcam

Check the [LED on the camera](#):

- LED is **red**: Camera detection failed. Check whether the TSI driver software has been installed. Disconnect and reconnect the camera to the USB cable. The camera should then be correctly recognized.
If the camera is still not listed in the [uc480 Camera Manager](#), open the Windows Device Manager to check whether the camera has been correctly recognized. If recognition was successful, you will find an entry in the format “uc480.....cameras” under “Universal Serial Bus Controllers.” A question mark or exclamation mark before the entry indicates that camera was not correctly recognized. You can remove the entry using the shortcut menu (right-click). Disconnect and reconnect the camera. The Found New Hardware Wizard will detect it as a new device and install the appropriate drivers. The camera should then be correctly recognized.
- LED is **green**: The camera is fully operational. Check whether the camera has been opened in a different application.
- LED is **off**: No power supply to the camera. Check the cable, the connectors and, if applicable, the power supply to the hubs. In case of a DCU22xX or DC1240X camera, check whether any pins of the [micro D-sub connector](#) have been bent.
- LED **flashes**: A fault has occurred in the camera hardware. Please contact the [Thorlabs](#).

USB DCx camera operation

- ☐ The camera can be opened in the software, but captures images sporadically or not at all.

Check for dropped frames in the settings window of the ThorCam software. If significant percentage of dropped frames are reported, the camera speed settings are too high for the system you are using. Check the following:

- Use only USB 2.0 (USB 3 in case of DCC3240x) certified cables and hubs.
- Do not use any passive extension cables.
- Do not connect the camera to the USB ports on the front of the PC, but to the ones directly on the main board. You will find those USB ports at the back of the PC.

In addition, check the following camera settings in the software:

- Pixel clock frequency: Reduce the [pixel clock](#) if data transfer errors occur. When you are operating more than one USB camera on one port, the pixel clock of all the cameras added together should not exceed about 40 MHz.
- ☐ Is it possible to operate an older USB camera with Windows 7?

Look for the serial number to see if your camera can be operated with Windows 7:

- The support for Windows 7 was introduced with driver version 3.50. This driver can be used with cameras from serial number 400 26 27000 on.

6.1.1 PCs with Energy Saving CPU Technology

This application note is related to all DCx USB cameras connected to PC systems using current CPU models that implement modern energy saving technologies.

Symptoms:

- Low USB bandwidth provided by the PC system
- TransferFailed errors occurring even at moderate pixel clock settings
- Camera operates at low speed only

Summary:

Current CPUs with modern energy saving features can cause bandwidth limitations on USB. The only available approach to this issue is to disable CPU sleep states. Unfortunately this is not possible for all systems.

Detailed explanation:

Modern CPUs like Intel i5 & i7 and others make use of advanced energy saving technologies ensuring a low power consumption and long battery life for mobile

devices. Additionally those CPU implement features for increasing the performance of single cores if there is enough thermal headroom available when other cores have little load.

A basic idea to achieve this is to put a CPU core to sleep while there is nothing to do for it. Various different activity states of CPU cores are available in modern CPUs. These CPU states are referred to as "C-states". C0 is the working state of a core.

Increasing numbers refer to less activity and longer wake up times. Current CPU fall down to variations of the C3 state which are referred to as "Sleep", "Deep Sleep" and similar.

Unfortunately negative effects of the sleep states have shown up. It is observed that the available bandwidth of PC busses drops significantly when part of the CPU enters these states.

The operation of DCx USB cameras is affected by the sleep states because they reduce the speed of the USB system. The available bandwidth on the USB may drop down to around 30% of the maximum bandwidth when the CPU, or one of its cores, enters sleeping states.

One would expect that a CPU core will not fall into a sleep state while it is obviously needed for the operation of the USB. But obviously USB data transfers do not prevent the CPU from falling to sleep. If the code execution load of a CPU core is low enough it will fall asleep and immediately reduce the USB bus speed.

For operation at high frame rates DCx cameras require an adequate USB bandwidth which might not be available when CPU cores are in sleep states.

Advice:

If you seem to be running into this low bandwidth issue please check and try the following. These first hints are general recommendations for issues with the USB

data transfer. You can check the USB performance with the "Optimum" pixel clock settings checkbox in uc480 Demo software. A good USB system should be able to reach a pixel clock

setting near the maximum value.

- Please remove other USB devices from the system (USB keyboard and mouse are fine). Run tests with only one camera connected at once.
- Make sure using a USB port directly on the mainboard. Front panel or other ports are connected to the mainboard with poor cabling quality frequently.
- Make sure to use USB2.0 certified cables to connect the camera.
- If you are using USB hubs or extensions: Run a test without these devices, connect the camera directly to the PC.
- Disable other equipment that is connected via USB. For example WLAN and Bluetooth adapters might use USB to connect.
- If you are using a mobile PC: run it on mains power, not battery.
- Check your energy saving options in the operating system. Disable energy saving features and set the available features to “full performance” or similarly named options.

If you checked the above and still observe low USB performance you might be experiencing the issue with CPU sleep states.

6.2 Status LED on USB DCx Cameras

DCU223x, DCU224x and DCC1240x

The LED on the rear side of the USB DCx camera indicates whether

- the DCx camera is powered on – LED lights up red (only USB board rev. 2.0 or higher).
- the uc480 driver has been loaded and the camera is operational – LED lights up green
- an error has occurred – green LED flashes:
- 2x flash: unknown sensor, please contact our [Thorlabs](#) team.

If the LED does not light up green, please check the following:

- Has the camera been connected correctly?
- Have the driver and the camera been installed properly in the [uc480 Camera Manager](#) on the host PC?
- Does the host PC meet all [system requirements](#)?



Revision 1.2 (green Status LED)



Revision 2.0 (red/green LED)

DCC3240x

The LED on the DCC3240x flashes 2x green if the camera is connected to a USB 2.0 port. If the camera is connected to a USB 3.0 port the LED flashes 3x green.

6.3 Color and Memory Formats

Attention

Obsolete parameters

The following parameters for color formats are obsolete and should no longer be used (see also [is_SetColorMode\(\)](#)):

- IS_SET_CM_RGB32
- IS_SET_CM_RGB24
- IS_SET_CM_RGB16
- IS_SET_CM_RGB15
- IS_CM_UYVY_MONO_PACKED
- IS_CM_UYVY_BAYER_PACKED
- IS_CM_BAYER_RG8 (identical to IS_CM_SENSOR_RAW8)
- IS_CM_BAYER_RG12 (identical to IS_CM_SENSOR_RAW12)
- IS_CM_BAYER_RG16 (identical to IS_CM_SENSOR_RAW16)
- IS_CM_BGR555_PACKED (has been renamed to IS_CM_BGR5_PACKED)

Each color format supported by the DCx camera defines a different memory format. The following table shows the byte arrangement in memory:

API constant	Description	Byte 0	Byte 1	Byte 2	Byte 3	Byte 4	Byte 5	Byte 6	Byte 7
		0 7 8	15 16	23 24	31	32 39 40	47 48	55 56	63
IS_CM_SENSOR_RAW16	Raw Bayer (16)	16 *		16 *		16 *		16 *	
	Odd row	16 *		16 *		16 *		16 *	
IS_CM_SENSOR_RAW12	Raw Bayer (12)	12		12		12		12	
	Odd row	12		12		12		12	
IS_CM_SENSOR_RAW8	Raw Bayer (8)	8		8		8		8	
	Odd row	8		8		8		8	
IS_CM_MONO16	Grey value (16)	16 *		16 *		16 *		16 *	
IS_CM_MONO12	Grey value (12)	12		12		12		12	
IS_CM_MONO8	Grey value (8)	8		8		8		8	
IS_CM_RGBA12_PACKED	RGB48 (12 12 12)	12		12		12		12	
IS_CM_RGB12_PACKED	RGB36 (12 12 12)	12		12		12		12	
IS_CM_RGB10_PACKED	RGB30 (10 10 10)	10		10		10		10	
IS_CM_RGBA8_PACKED	RGB32 (8 8 8)	8		8		8		8	
IS_CM_RGBY8_PACKED	RGBY (8 8 8 8)	8		8		8		8	
IS_CM_RGB8_PACKED	RGB24 (8 8 8)	8		8		8		8	
IS_CM_BGRA12_PACKED	BGR48 (12 12 12)	12		12		12		12	
IS_CM_BGR12_PACKED	BGR36 (12 12 12)	12		12		12		12	
IS_CM_BGR10_PACKED	BGR30 (10 10 10)	10		10		10		10	
IS_CM_BGRA8_PACKED	BGR32 (8 8 8)	8		8		8		8	
IS_CM_BGR8_PACKED	BGR24 (8 8 8)	8		8		8		8	
IS_CM_BGRY8_PACKED	BGRY (8 8 8)	8		8		8		8	
IS_CM_BGR565_PACKED	BGR16 (5 6 5)	5 6 5		5 6 5		5 6 5		5 6 5	
IS_CM_BGR5_PACKED	BGR15 (5 5 5)	5 5 5		5 5 5		5 5 5		5 5 5	
IS_CM_UYVY_PACKED	YUV 4:2:2 (8 8)	8		8		8		8	
IS_CM_UYVY_MONO	YUV 4:2:2 (8 8)	8		8		8		8	
IS_CM_UYVY_BAYER	YUV 4:2:2 (8 8)	8		8		8		8	
IS_CM_CBCRY_PACKED	YCbCr 4:2:2 (8 8)	8		8		8		8	

Color codes	Red channel	Y / grey channel	Cr component
	Green channel	U component	Cb component
	Blue channel	V component	Unused (0)

Colour and memory formats

Note

An asterisk (*) identifies formats which are filled starting with the most significant bit (MSB) but which may have less than the indicated number of payload bits, depending on the camera model.

For the RGB16 and RGB15 data formats, the MSBs of the internal 8-bit R, G and B colors are used.

The first pixel in the first line with the index (0,0) is always a red pixel at color cameras.

The list above does not contain the IS_CM_RGB8_PLANAR color format. In planar RGB the image is

saved as 8 bit RGB. The channels red, green, and blue are stored separately, i.e. first all red information, then all green information and at last all blue information are saved.

6.4 uc480 Parameter File (ini file)

Using the [is_ParameterSet\(\)](#) function, you can save the currently set DCx camera parameters to a file in the ini format (*.ini) or load an ini file.

Attention

Only camera-specific ini files can be loaded.

The ini file you want to load has to match the paired camera model.

When loading an ini file, make sure that the image size (AOI) and color depth parameters in the ini file match those in the allocated memory. Otherwise, display errors may occur.

uc480 parameter files can also be created and edited manually. The following table shows the structure of the parameter file. The entries in square brackets [] indicate sections. If a section does not exist in the ini file, the corresponding camera parameters will not be modified when you load the file.

Hint

You can use wildcards in the ini file:

If you specify * as first character, then the highlighted characters will be interpreted in the camera name only, e.g. ***DCU224C**. Thus this string applies for DCU223C and DCU224C.

Structure of a uc480 parameter file

Parameter	Description	Value range	Example
[Versions]			
uc480.dll	File version of the uc480 API	-	4.00.0000
uc480_eth.sys	(Not applicable to DCx Cameras)	-	4.00.0000
uc480_usb.sys	File version of the USB uc480 driver	-	4.00.0000
uc480_boot.sys	File version of the USB uc480 boot loader	-	4.00.0000
[Sensor]			
Sensor	Full name of the camera model	-	DCC3240C
[Image size] Image size settings			
Start X	Start point (X coordinate) in AOI mode	0...(max. width ^{*1} -width)	100
Start Y	Start point (Y coordinate) in AOI mode	0...(max. height ^{*1} -Height)	100
Start X absolute	Activate absolute AOI positioning in the memory (see is_AOI())	0, 1	1
Start Y absolute	Activate absolute AOI positioning in the memory (see is_AOI())	0, 1	1
Width	Width of the AOI	Sensor-dependent ^{*1}	2460

Parameter	Description	Value range	Example
Height	Height of the AOI	Sensor-dependent ¹	1820
Binning	Activate binning mode and select factor	Sensor-dependent ²	0
Subsampling	Activate subsampling mode and select factor	Sensor-dependent ²	0
[Scaler]	The internal image scaling is only supported by sensors of the DCC1240x / DCC3240x camera series.		
Mode	Enable/disable scaling	0 = Scaling off 1 = Scaling on	0
Factor	Scaling factor		
[Multi AOI]			
Enabled	Activate/deactivate multi AOI	0 = Multi AOI off 1 = Multi AOI on	0
Mode	Mode of multi AOI Currently only IS_AOI_MULTI_MODE_AXES is supported		
x1...x4, y1...y4	Axis for multi AOI mode		
[Shutter]			
Mode	Shutter mode	Sensor-dependent (only supported by DCC1240x / DCC3240x models) 1 = Rolling shutter 2 = Global shutter 4 = Fast linescan 64 = Rolling shutter with global start 128 = Global shutter (alternative timing)	1
Linscan number	Line which is used in the linescan mode. The maximum possible line number depends on the height of the selected AOI.	Sensor-dependent (only supported by DCC1240x / DCC3240x models)	512
[Log Mode]			
Mode	Log mode (only supported by DCC1240x / DCC3240x models)	0 = Factory-default with anti-blooming 1 = Off (no anti-blooming) 2 = Manual Log mode	0
Manual value	Log mode value	Only in combination with manual Log mode	
Manual gain	Log mode gain		
[Timing]	Timing parameter settings		
Pixelclock	Current pixel clock of the camera	Sensor-dependent ¹	103
Framerate	Current frame rate	Depends on Pixelclock	15.104458

Parameter	Description	Value range	Example
		and image geometry	
Exposure	Current exposure time	Depends on Framerate	0.334059
Long exposure	Activates long exposure If the long exposure is active, then the range of exposure changes.	Not supported by DCx cameras	0
[Selected Converter]	Sets the type of Bayer conversion for the specified color format when using color cameras (see is_SetColorConverter()). For a description of all color formats, see the Color and memory formats section.		
IS_SET_CM_RGB32	Color format	0, 1, 2, 4	2
IS_SET_CM_RGB24	Color format	0, 1, 2, 4	2
IS_SET_CM_RGB16	Color format	0, 1, 2, 4	2
IS_SET_CM_RGB15	Color format	0, 1, 2, 4	2
IS_SET_CM_Y8	Color format	0, 1, 2, 4	2
IS_SET_CM_RGB8	Color format	0, 1, 2, 4	2
IS_SET_CM_BAYER	Color format	0, 1, 2, 4	8
IS_SET_CM_UYVY	Color format	0, 1, 2, 4	2
IS_SET_CM_UYVY_MONO	Color format	0, 1, 2, 4	2
IS_SET_CM_UYVY_BAYER	Color format	0, 1, 2, 4	2
IS_CM_CBYCRY_PACKED	Color format	0, 1, 2, 4	8
IS_SET_CM_RGBY	Color format	0, 1, 2, 4	8
IS_SET_CM_RGB30	Color format	0, 1, 2, 4	8
IS_SET_CM_Y12	Color format	0, 1, 2, 4	8
IS_SET_CM_BAYER12	Color format	0, 1, 2, 4	8
IS_SET_CM_Y16	Color format	0, 1, 2, 4	8
IS_SET_CM_BAYER16	Color format	0, 1, 2, 4	8
IS_CM_RGBA8_PACKED	Color format	0, 1, 2, 4	2
IS_CM_RGB8_PACKED	Color format	0, 1, 2, 4	2
IS_CM_RGBY8_PACKED	Color format	0, 1, 2, 4	8
IS_CM_RGB10V2_PACKED	Color format	0, 1, 2, 4	8
[Parameters]	Additional image parameter settings		
Colormode	Sets the current color mode	see Color and memory formats	11
Brightness	Software correction of image brightness ^{*3}	0...255	100
Contrast	Software correction of image contrast ^{*3}	0...511	215
Gamma	Software correction of the gamma value	0.01...10.0	1.000000
Hardware Gamma	Sensor-based hardware correction of the gamma	0, 1	0

Parameter	Description	Value range	Example
	value		
Blacklevel Mode	Mode for black level correction of the sensor	0, 1, 32 ²	1
Blacklevel Offset	Manual offset for black level correction of the sensor	0...255	0
Hotpixel Mode	Mode for hot pixel correction	0, 1, 2, 4 ²	2
Hotpixel Threshold	Not used	-	0
Sensor Hotpixel	Activates the sensor-internal hot pixel correction	Sensor-dependent ¹ 1 = on 0 = off	0
GlobalShutter	Enables the Global Start shutter of the sensor	0, 1 Not supported by DCx models, they return "7" (not supported).	0
[Gain]	Sets the sensor gain control for image brightness		
Master	Master gain	0...100	0
Red	Red gain	0...100	6
Green	Green gain	0...100	0
Blue	Blue gain	0...100	6
GainBoost	Activate gain boost	0, 1	0
[Processing]	Parameters for image pre-processing in the driver		
EdgeEnhancementFactor	Enable edge enhancement	0...2	0
RopEffect	Image geometry change (Rop = raster operation), e.g. mirroring	0, 8, 16, 32, 64 ²	0
Whitebalance	Enable software white balance	0, 1, 2, 4 ²	0
Whitebalance Red	Red factor for software white balance	double value	1.000000
Whitebalance Green	Green factor for software white balance	double value	1.000000
Whitebalance Blue	Blue factor for software white balance	double value	1.000000
Color correction	Enable color correction	0, 1, 2, 4, 80 ²	1
Color_correction_factor	Set the color correction factor	0.0...1.0	1.000000
Color_correction_satU	Saturation-U (see also is_SetSaturation())	0...200 100 = Saturation 1.0 200 = Saturation 2.0	0
Color_correction_satV	Saturation-V (see also is_SetSaturation())	0...200 100 = Saturation 1.0 200 = Saturation 2.0	0

Parameter	Description	Value range	Example
Bayer Conversion	Sets the size of the Bayer conversion mask for the current color format when using color cameras	1, 2 ²	1
[Auto features]	Sets the parameters for automatic image control		
Auto Framerate control	Enable frame rate control	0, 1	0
Brightness exposure control	Enable exposure time control	0, 1	0
Brightness gain control	Enable sensor gain control	0, 1	0
Auto Framerate Sensor control	Enable the sensor-internal control for frame rates (see also is_SetAutoParameter())	0 = off 1 = on	0
Brightness exposure Sensor control	Enable the sensor-internal brightness control	0 = off 1 = on	0
Brightness gain Sensor control	Enable the sensor-internal gain control	0 = off 1 = on	0
Brightness exposure Sensor control photometry	Not supported by DCx cameras		
Brightness gain Sensor control photometry	Not supported by DCx cameras		
Brightness control once		0, 1	0
Brightness reference	Reference value for brightness control	0...255	128
Brightness speed	Brightness control speed	0...100	50
Brightness max gain	Maximum gain for brightness control	0...100	100
Brightness max exposure	Maximum exposure time for brightness control	Depends on <code>Pixelclock</code> and image geometry	66.082816
Brightness Aoi Left	X start point of reference AOI for brightness control	0...(max. width ¹ -Aoi Width)	0
Brightness Aoi Top	Y start point of reference AOI for brightness control	0...(max. height ¹ -Aoi Height)	0
Brightness Aoi Width	Width of reference AOI for brightness control	Sensor-dependent ¹	2560
Brightness Aoi Height	Height of reference AOI for brightness control	Sensor-dependent ¹	1920
Brightness Hysteresis	Hysteresis value for auto exposure and gain (see IS_SET_AUTO_HYSTERESIS)	0...10	2
Brightness Skip Frames	Number of images that will be not analyzed for the control (see IS_SET_AUTO_SKIPFRAMES)	0...1000	4

Parameter	Description	Value range	Example
Auto WB control	Enable white balance control	0, 1	0
Auto WB type	White balance mode		
Auto WB RGB color mode	Color space of white balance (only active, if the mode "Auto (Kelvin)" is set)		
Auto WB offsetR	Red offset for white balance control	0...100	0
Auto WB offsetB	Blue offset for white balance control	0...100	0
Auto WB gainMin	Minimum gain for white balance control	$0 \leq \text{gainMin} \leq \text{gainMax} \leq 100$	0
Auto WB gainMax	Maximum gain for white balance control	$0 \leq \text{gainMin} \leq \text{gainMax} \leq 100$	100
Auto WB speed	White balance control speed	0...100	50
Auto WB Aoi Left	X start point of reference AOI for white balance control	$0 \dots (\text{max. width}^1 - \text{Aoi Width})$	0
Auto WB Aoi Top	Y start point of reference AOI for white balance control	$0 \dots (\text{max. height}^1 - \text{Aoi Height})$	0
Auto WB Aoi Width	Width of reference AOI for white balance control	Sensor-dependent ¹	2560
Auto WB Aoi Height	Height of reference AOI for white balance control	Sensor-dependent ¹	1920
Auto WB Once	Carry out white balance control once	0, 1	0
Auto WB Hysteresis	Hysteresis value for auto white balance (see IS_SET_AUTO_WB_HYSTERESIS)	0...10	2
Auto WB Skip Frames	Number of images that will not be analyzed for control (see IS_SET_AUTO_WB_SKIPFRAMES)	0...1000	4
[Trigger and Flash]	Sets the digital inputs/outputs		
Trigger mode	Trigger mode		
Trigger timeout	Timeout value for triggered image capture in 10 ms steps		
Trigger delay	Delay of triggered image capture in μs	Sensor-dependent ¹	15
Trigger debounce mode	Not supported by DCx		

Parameter	Description	Value range	Example
	cameras		
Trigger debounce delay time	Not supported by DCx cameras		10
Trigger burst size	Not supported by DCx cameras		
Flash strobe	Activate flash output	0...6 ^{*2}	0
Flash delay	Delay of the flash signal in μ s	Depends on sensor setting, can be queried using is_IO()	0
Flash duration	Duration of the flash signal in μ s	Depends on sensor setting, can be queried using is_IO()	200
PWM mode	PWM mode	see is_IO()	
PWM frequency	Frequency of the pulse-width modulation		
PWM dutycycle	Duty cycle of the pulse-width modulation	0.0...1.0 (1.0 corresponds to 100 %)	
GPIO state	State of the GPIO	see is_IO()	
GPIO direction	Direction of the GPIO	see is_IO()	
[Sequence AOI]	The sequence AOI is supported by DCC1240x / DCC3240x models only.		
NumberUsedAOI	Number of used AOIs	see is_AOI()	
StartX1...StartX3	X position of AOI 1-3		
StartY1...StartY3	Y position of AOI 1-3		
Gain1...Gain3	Gain of AOI 1-3		
Exposure1...Exposure3	Exposure of AOI 1-3		
ReadoutCycle1... ReadoutCycle3	Number of read-out cycles of AOI 1-3		
BinningModel1... BinningMode3	Binning mode of AOI 1-3		
SubsamplingModel1... SubsamplingMode3	Subsampling mode of AOI 1-3		
ScalerFactor1... ScalerFactor3	Scaling factor of AOI 1-3		
DetachImageParameter1... DetachImageParameter3	Changes of exposure time and gain are transferred to AOI 1-3	0 = Every change is transferred 1 = Changes are not transferred see also is_AOI()	
[Transfer]	Not supported by DCx cameras		
ImageDelay_us			
PacketInterval_us			

*1 This information is provided in [Camera and sensor data](#) chapter.

*2 For the parameters, please refer to the uc480.h header file provided in the \Develop\include

folder of the uc480 installation directory (see also [Programming notes](#)).

*3 Function obsolete, see chapter [Obsolete functions](#).

6.5 Definition of IP Protection Classes

The housing of the DCx Cameras comply with **IP 30**:

- Protection against the ingress of small particles (diameter > 2.5 mm)
- No protection against water

6.6 History of API functions

New functions in software version 4.50

- [israw_AddFrame\(\)](#)
- [israw_CloseFile\(\)](#)
- [israw_ExitFile\(\)](#)
- [israw_GetFrame\(\)](#)
- [israw_GetImageInfo\(\)](#)
- [israw_GetSize\(\)](#)
- [israw_InitFile\(\)](#)
- [israw_OpenFile\(\)](#)
- [israw_SeekFrame\(\)](#)
- [israw_SetImageInfo\(\)](#)

New functions in software version 4.40

- [is_Gamma\(\)](#)
- [is_LUT\(\)](#)

New functions in software version 4.20

- [is_ImageBuffer\(\)](#)
- [is_Measure\(\)](#)

New functions in software version 4.01

- [is_Blacklevel\(\)](#)

New functions in software version 4.00

- [is_AutoParameter\(\)](#)
- [is_Convert\(\)](#)
- [is_DeviceInfo\(\)](#)
- [is_EdgeEnhancement\(\)](#)
- [is_ImageFile\(\)](#)
- [is_ParameterSet\(\)](#)
- [is_PixelClock\(\)](#)
- [is_Stream\(\)](#)

New functions in software version 3.90

- [is_BootBoost\(\)](#)

-
- [is_CaptureStatus\(\)](#)
 - [is_Exposure\(\)](#)
 - [is_IO\(\)](#)

New functions in software version 3.81

- [is_Configuration\(\)](#)

New functions in software version 3.80

- [is_AOI\(\)](#)
- [is_DeviceFeature\(\)](#)
- [is_HotPixel\(\)](#)

New functions in software version 3.70

- [is_ColorTemperature\(\)](#)
- [is_TriggerDebounce\(\)](#)

New functions in software version 3.61

- `is_ScenePreset()`
- `is_Saturation()`
- `is_Sharpness()`
- `is_Zoom()`

New functions in software version 3.52/3.60

- `is_FaceDetection()`
- `is_Focus()`
- [is_ImageFormat\(\)](#)
- `is_ImageStabilization()`

New functions in software version 3.40

- [is_DirectRenderer\(\)](#)
- [is_GetImageInfo\(\)](#)
- `is_GetDuration()`
- [is_GetSensorScalerInfo\(\)](#)
- [is_SetSensorScaler\(\)](#)

New functions in software version 3.33

- [is_Direct3D\(\)](#)
- [is_GetTimeout\(\)](#)

New functions in software version 3.32

- [is_GetTimeout\(\)](#)
- [is_SetTimeout\(\)](#)
- [is_SetTriggerCounter\(\)](#)

New functions in software version 3.30

- `is_GetCameraLUT()`
- `is_GetCaptureErrorInfo()`
- [is_GetColorConverter\(\)](#)

- `is_GetComportNumber()`
- [is_GetSupportedTestImages\(\)](#)
- [is_GetTestImageValueRange\(\)](#)
- `is_ResetCaptureErrorInfo()`
- `is_SetCameraLUT()`
- [is_SetColorConverter\(\)](#)
- [is_SetSensorTestImage\(\)](#)

New function in software version 3.20

- [is_SetOptimalCameraTiming\(\)](#)

New functions in software version 3.10

- `is_EnableHdr()`
- `is_GetHdrKneepointInfo()`
- `is_GetHdrKneepoints()`
- `is_GetHdrMode()`
- `is_SetHdrKneepoints()`

New functions in software version 3.00

- `is_GetEthDeviceInfo()`
- `is_SetAutoCfgIpSetup()`
- `is_SetPacketFilter()`
- `is_SetPersistentIpCfg()`
- `is_SetStarterFirmware()`

6.7 Certifications and Compliances

Compliance with the directives is demonstrated by meeting the following standards:

Product type	EMC immunity	EMC emission	UL certification
DCC1240x ^{*1}	EN 61000-6-2:2005	EN 61000-6-3:2001 + A11:2004	UL 60950-1, 2nd Edition, 2011-12-19 CSA C22.2 No. 60950-1-07, 2nd Edition, 2011-12
DCU223x ^{*1} DCU224x ^{*1}	EN 61000-6-2:2001	EN 61000-6-4:2001	UL 60950-1, 2nd Edition, 2011-12-19 CSA C22.2 No. 60950-1-07, 2nd Edition, 2011-12
DCC1545M ^{*1} DCC1645C ^{*1}	EN 61000-6-2:2005	EN 61000-6-3:2001 + A11:2004	UL 60950-1, 2nd Edition, 2011-12-19 CSA C22.2 No. 60950-1-07, 2nd Edition, 2011-12
DC3240x ^{*1}	EN 61000-6-2:2005	EN 61000-6-3:2007	UL 60950-1, 2nd Edition, 2011-12-19 CSA C22.2 No. 60950-1-07, 2nd Edition, 2011-12

^{*1} This equipment has been tested and found to comply with part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential

area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

Modifications not expressly approved by the manufacturer could void the user's authority to operated the equipment under FCC rules.

Certificates of Conformity available upon request.

6.8 Thorlabs 'End of Life' Policy (WEEE)

As required by the WEEE (Waste Electrical and Electronic Equipment Directive) of the European Community and the corresponding national laws, Thorlabs Scientific Imaging offers all end users in the EC the possibility to return “end of life” units without incurring disposal charges.

This offer is valid for Thorlabs Scientific Imaging electrical and electronic equipment

- sold after August 13th 2005
- marked correspondingly with the crossed out “wheelie bin” logo (see figure below)
- sold to a company or institute within the EC
- currently owned by a company or institute within the EC
- still complete, not disassembled and not contaminated

As the WEEE directive applies to self contained operational electrical and electronic products, this “end of life” take back service does not refer to other Thorlabs Scientific Imaging products, such as

- pure OEM products, that means assemblies to be built into a unit by the user (e. g. OEM laser driver cards)
- components
- mechanics and optics
- left over parts of units disassembled by the user (PCB's, housings etc.).

Waste treatment on your own responsibility

If you do not return an “end of life” unit to Thorlabs Scientific Imaging, you must hand it to a company specialized in waste recovery. Do not dispose of the unit in a litter bin or at a public waste disposal site.

WEEE Number (Germany) : DE97581288

Ecological background

It is well known that waste treatment pollutes the environment by releasing toxic products during decomposition. The aim of the European RoHS Directive is to reduce the content of toxic substances in electronic products in the future.

The intent of the WEEE Directive is to enforce the recycling of WEEE. A controlled recycling of end-of-life products will thereby avoid negative impacts on the environment.



6.9 Exclusion of Liability and Copyright

Thorlabs Scientific Imaging has taken every possible care in preparing this Operation Manual. We however assume no liability for the content, completeness or quality of the information contained therein. The content of this manual is regularly updated and adapted to reflect the current status of the software. We furthermore do not guarantee that this product will function without errors, even if the stated specifications are adhered to.

Should you require further information about this product or encounter specific problems that are not discussed in sufficient detail in the User Manual, please contact your nearest Thorlabs office.

All rights reserved. This manual may not be reproduced, transmitted or translated to another language, either as a whole or in parts, without the prior written permission of *Thorlabs Scientific Imaging*.

Copyright © Thorlabs Scientific Imaging 2018. All rights reserved.

6.10 Thorlabs Worldwide Contacts

For technical support or sales inquiries, please visit us at www.thorlabs.com/contact for our most up-to-date contact information.

Index

C

C#	426
C/C++	425, 426, 427

Camera	
close	89
EEPROM	77
electrical data	472
Information	88
open	89
parameters	81
camera data	439
Camera EEPROM	109
Camera ID	77
Camera Manager	73
camera list	74
Control center	74
Camera parameters	77
load at start-up	148
Camera status	
get	109
Capture	
freerun synchronisation	93
live mode (freerun)	18
overlap trigger	20
single snap mode	18
trigger mode	20
Capture video	141
Characteristics	
gamma curve	44
linear	44
logarithmic	44
C-mount/CS-mount	460
Color depth	341
Color formats	490
COM port -> Serial interface	57
Contrast adjustment	47
CPU idle state	77
CRA correction	26

D

DCC3240x	
digital input	480
digital output	481
electrical data	477
flash	481
GPIO	478
I/O connector	477

A

ActiveX programming	427
Ambient conditions	471
AOI	
for automatic image control	120
AOI (Area of interest)	38
AOI sequence	38
multi AOI	38
AOI sequence	120
Application notes	
DCC1240x	68
DCC1545M	71
DCC1645C	72
DCC3240x	68
DCC3260	67
DCU223x	72
DCU224x	72
Area of interest -> AOI	38
Area scan camera	37
Area scan sensor	37
Auto exposure shutter (AES)	51
Auto frame rate (AFR)	51
Auto gain (AGC)	51
Auto gain control (AGC)	51
Auto white balance (AWB)	51
Automatic image control	
hysteresis	51
programming	108, 324
AVI functions	
errors	433
AVI recording	395

B

Bad pixel -> hot pixel	260
Bayer filter	
Bayer conversion	29
Binning	38, 43
Bit depth	47, 341
Black level correction	51
Blinking codes	488
Boot boost	77
Burst trigger mode	383

DCC3240x		Event handling	101
pin assignment	477	Exposure	
serial interface	482	dual	208
trigger	480	fine	205
DCC3260		long	207
camera data	439	set	200, 202
DCU223x / DCU224x		Exposure time	50
electrical data	472		
I/O connector	472	F	
Pin assignment	472	Fast line scan	68
DCx		Fill factor	26
driver compatibility	65	Filter glasses	464
DCx camera model	13	clean	464
DCC3240x	13	DL filter	464
DCU223x / DCC224x	14	IR cut filter	464
DCx camera models	13, 435	Firmware	53
comparison table	436	Flange back distance	460
DCx Model Naming	438	calculate	460
Conventions		Flash	20, 33, 286
DCx naming convention	438	Flash timing	
Delphi	427	in freerun mode	56
DIB	22	in trigger mode	54
Digital input		Frame rate	50, 353
trigger	54	Freerun -> Capture: live mode	18
Digital input/output	21, 54	Freerun mode	141
programming	351		
Digital output		G	
flash	54	Gain	358
Direct3D	22, 177, 345	linearity	51
DirectShow	77	Gamma characteristic	44
DirectX	64, 177	Gamma-Kennlinie	216
Display		Global shutter	33
bitmap mode (DIB)	22	GPIO	286
color formats	490	Graphics card	22, 64
DIB mode	92		
Direct3D	22, 64, 92	H	
DirectDraw	422	Histogram	47
OpenGL	22	Hot pixel correction	31, 260
overlay	92	Hotpixel	
overlay display	177	add	82
Display overlay	92	edit	82
Driver version	76		
E		I	
EEPROM of the camera	77, 483	Image capture -> Capture	20
Errors		Image display -> Display	22
error messages	429	Image format	
transmission errors	50	set	269
troubleshooting	486		

Image memory		is_GetActiveImageMem	217
allocate	89	is_GetActSeqBuf	219
INI file -> Parameter file	493	is_GetAutoInfo	221
Interface		is_GetBusSpeed	225
serial	482	is_getCameraInfo	226
IP norm	500	is_GetCameraList	228
IR cut filter	464	is_GetCameraLUT	230
is_AddToSequence	117	is_GetColorConverter	231
is_AllocImageMem	118	is_GetColorDepth	232
is_AOI	120	is_GetDLLVersion	233
is_AutoParameter	130	is_GetError	234
is_Blacklevel	133	is_GetFramesPerSecond	235
is_CameraStatus	136	is_GetFrameTimeRange	236
is_CaptureStatus	138	is_GetImageHistogram	238
is_CaptureVideo	141	is_GetImageMem	245
is_ClearSequence	143	is_GetImageMemPitch	246
is_ColorTemperature	144	is_GetNumberOfCameras	248
is_Configuration	148	is_GetOsVersion	249
is_Convert	155	is_GetSensorInfo	250
is_CopyImageMem	157	is_GetSensorScalerInfo	252
is_CopyImageMemLines	158	is_GetSupportedTestImages	253
is_DeviceFeature	159	is_GetTestImageValueRange	255
AOI merge mode	161	is_GetTimeout	256
image memory	170	is_GetUsedBandwidth	257
level controlled trigger	168	is_GetVsyncCount	258
line scan mode	172	is_HasVideoStarted	259
Log mode	165	is_HotPixel	260
shutter mode	169	is_ImageFile	265
timestamp	173	is_ImageFormat	269
is_DeviceInfo	174	is_InitCamera	276
is_DirectRenderer	177	is_InitEvent	280
is_DisableEvent	187	is_InitImageQueue	282
is_EdgeEnhancement	188	is_InquireImageMem	284
is_EnableAutoExit	191	is_IO	286
is_EnableEvent	192	is_IsVideoFinish	299
is_EnableMessage	195	is_LockSeqBuf	301
is_ExitCamera	197	is_LUT	302
is_ExitEvent	198	is_LUT()	422
is_ExitImageQueue	199	is_Measure	305
is_Exposure	200	is_ParameterSet	309
dual exposure	208	is_PixelClock	312
exposure time	202	is_ReadEEPROM	315
fine increment	205	is_RenderBitmap	317
long exposure	207	is_ResetToDefault	319
is_ForceTrigger	210	IS_SET_EVENT_TRANSFE	192
is_FreeImageMem	212	R_FAILED	
is_FreezeVideo	214	is_SetAllocatedImageMem	321
is_Gamma	216		

is_SetAutoParameter	324	isavi_OpenAVI	408
is_SetBinning	332	isavi_ResetFrameCounters	409
is_SetCameraID	335	isavi_SetFrameRate	410
is_SetColorConverter	337	isavi_SetImageQuality	411
is_SetColorCorrection	339	isavi_SetImageSize	412
is_SetColorMode	341	isavi_StartAVI	414
is_SetDisplayMode	345	isavi_StopAVI	415
is_SetDisplayPos	349	israw_AddFrame	416
is_SetErrorReport	350	israw_CloseFile	416
is_SetExternalTrigger	351	israw_ExitFile	417
is_SetFrameRate	353	israw_GetFrame	417
is_SetGainBoost	355	israw_GetImageInfo	418
is_SetGamma	356	israw_GetSize	419
is_SetHardwareGain	358	israw_InitFile	419
is_SetHWGainFactor	361	israw_OpenFile	420
is_SetImageMem	363	israw_SeekFrame	421
is_SetOptimalCameraTiming	364	israw_SetImageInfo	421
is_SetRopEffect	366		
is_SetSaturation	368	L	
is_SetSensorScaler	369	LED	286
is_SetSensorTestImage	372	LED -> Status LED	
is_SetSubSampling	374	DCx Status LED	488
is_SetTimeout	378	Lens	
is_SetTriggerCounter	380	immersion depth	462
is_SetTriggerDelay	381	Line scan	
is_StopLiveVideo	382	fast line scan	68
is_Trigger	383	Line scan mode	37, 159
is_TriggerDebounce	385	Linux	64
is_UnlockSeqBuf	389	not supported functions	433
is_WaitEvent	390	Live mode	141
is_WaitForNextImage	391	Log mode	159
is_WriteEEPROM	393	Lookup table (LUT)	44, 45
isavi_AddFrame	395		
isavi_CloseAVI	396	M	
isavi_DisableEvent	396	Memory board	422
isavi_EnableEvent	397	Memory formats	490
isavi_ExitAVI	398	Micro lenses	26
isavi_ExitEvent	399	Multi AOI	38, 120
isavi_GetAVIFileName	400	Multi camera systems	
isavi_GetAVIFileNameW	401	system requirements	64
isavi_GetAVISize	402		
isavi_GetnCompressedFrames	403	N	
isavi_GetnLostFrames	404	Network card	64
isavi_InitAVI	405	O	
isavi_InitEvent	406	Obsolete functions	422
isavi_OpenAVI	407		

OpenGL	22, 345	DCC3240x	57
Operating system	64	Setting the area of interest ->	120
Overlap trigger -> Capture	20	AOI	
Overlay -> Display	22	Shock resistance	471
Overlay: Direct3D		Shutter	
display (programming)	92	global	33
Overlay: display	177	rolling	33
		Standby mode of the camera	21
		Starter firmware	53, 276
		Status LED	65
		DCx camera	488
		Subsampling	38, 42
		support	486
		Support file	
		create	77
		System requirements	64
P			
Parameter			
new	53		
Parameter file (ini file)	493		
Pixel clock	50		
Pixel pre-processing	45		
Port	77		
Ports	57		
Processor			
idle states	148		
support for multitasking	148		
Programming languages	425, 426, 427		
Programming notes	425		
PWM	286		
Q			
Quick start			
connection	62		
image capture	62		
programming	86		
R			
Raw Bayer -> Bayer filter	29		
RAW file	416		
Reference AOI			
for automatic image control	120		
Return values	429		
RGB gain	51		
ROI -> AOI	38, 120		
Rolling shutter	33		
Rolling shutter (global start)	33		
RS-232 -> Serial interface	57		
S			
Sensor			
position accuracy	463		
Sensor formats	25		
Serial interface			
		T	
		Temperature range	471
		Thread programming	427
		Trigger	286
		level controlled	168
		Trigger mode	383
		Trigger mode -> Image capture	20
		U	
		UART	482
		uc480 Camera Manager	73
		Control Center	74
		uc480 Hotpixel Editor	82
		uc480 programs	72
		uc480.h	425, 426, 427
		uc480_Api.dll	425, 426, 427
		USB	
		bandwidth	60
		cable	59
		connection	59
		hub	65
		standard	58
		topology	58
		USB 3.0	
		cable	60
		connction	60

V

VB .NET	427
Vibration resistance	471

W

Windows	64
---------	----